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LYMPHOTRACK® ASSAY	LYMPHOTRACK® LOW POSITIVE CONTROL	LYMPHOQUANT® INTERNAL CONTROL
IGHV Leader, IGH FR1 1/2/3	LymphoTrack® B-cell Low Positive Control Catalog #:4-088-0098	LymphoQuant® B-cell Internal Control Catalog #:4-088-0118
TRG, TRB	LymphoTrack® T-cell Low Positive Control Catalog #:4-088-0108	LymphoQuant® T-cell Internal Control Catalog #:4-088-0128

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Note: Total **DNA Input** and **Total Read Depth** are divided across replicates tested, and a replicate is an independent PCR reaction with input DNA from the same subject. Read depth requirements assume 100% target cells. Read depth requirements may vary dependent on replicate numbers used to reach total DNA inputs.

Also, note: A DNA concentration step is often required to achieve the higher levels of DNA input.

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THE RELATIONSHIP BETWEEN DNA INPUT, READ DEPTH, AND CONFIDENCE LEVEL MUST BE CONSIDERED IN MRD STUDIES.

The Project Planner tool in the LymphoTrack MRD Software can be used to create an experimental design that allows a user to optimize confidence in the result using the available sample and read depth. This tool allows for user-defined permutations to DNA input, read depth, and/or replicates; performing calculations based on these permutations to determine confidence levels at various sensitivity levels. These numerical outputs can be used to generate confidence models as shown in **Figure 1**.

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
CONFIDENTLY REPORT TRUE NEGATIVE RESULTS WITH 10^{-6} SENSITIVITY

Increased DNA inputs are required to achieve greater sensitivity in MRD testing due to the probability that a measured specimen will contain a subject's cancer cells. **Figure 2** explains in detail the cell sampling required to reach 10^{-6} sensitivity with 95% confidence.


HOW MANY CELLS MUST BE TESTED TO REACH 1×10^{-6} SENSITIVITY WITH >95% CONFIDENCE THAT IT IS A TRUE NEGATIVE RESULT?

10% Maroon Marble in 90% Teal Marbles

If we pick 1 marble we have a 10% probability of picking a **maroon** marble and 90% probability of picking a **teal** marble.



10% Clonal Cancer Cells in 90% Non-Clonal Cells


This is akin to testing 667 ng of DNA (1/30th of the 20µg total suggested, or ~100k cells). The calculated probability of detecting the cancer cell is only 9.75%.

If we pick 3 marbles we have a 27% probability of picking at least 1 **maroon** marble and 73% probability of picking all **teal** marbles.

This is akin to testing 2000 ng of DNA (3/30th of the 20µg total suggested, or ~310,000 cells). The calculated probability of detecting the cancer cell is 26.5%.

If we pick 30 marbles we have a 95.8% probability of picking at least 1 **maroon** marble and only a 4.2% probability of picking all **teal** marbles.

This is akin to testing 20 µg of DNA (30/30th of the 20µg total suggested, or ~3.1m cells). The calculated probability of detecting the cancer cell is 95.38%.

Figure 2: 20 µg DNA input minimum must be tested to establish a true negative result with 1×10^{-6} sensitivity (>95% confidence).

In this example, the maroon marble(s) found in the teal marble background represents the rare clonal cancer cell(s) in a subject's otherwise non-clonal cell background. The 3 scenarios presented here demonstrate that laboratories must test 20 µg of DNA for 95% confidence that they can find 1 rare clonal cell in 1 million cells (1×10^{-6} sensitivity). It is important to note that the probability is entirely dependent on the assumed proportion of maroon marbles (or clonal cells). If we KNOW there should be 10% maroon marbles, then picking 30 teal marbles tells us we are 95.8% sure that the jar does not have 10% maroon marbles. This is the same logic that is used for confidence in a negative MRD signal using 20 µg of DNA. The actual proportion might be less than 1×10^{-6} sensitivity, but we are 95% sure it is not equal to or higher than 1×10^{-6} .

6. DESIGN OF EXPERIMENT: SAMPLE BATCHING AND ASSAY MULTIPLEXING

LymphoTrack® Assays are available for use with the MiSeq®, Ion S5™ and Ion PGM™ sequencing platforms. The LymphoTrack® Assays are designed to provide the end user efficient and flexible workflow options. Users can design cost-effective MRD runs by multiplexing assays and batching samples, while remaining mindful of the NGS flow cell capacity and the desired MRD sensitivity.

Assay Multiplexing: The full LymphoTrack® Clonality Suite consists of 7 independent assays each targeting one of the following respective loci (*IGHV* Leader, *IGH* FR1/FR2/FR3, *IGK*, *TRG* and *TRB*). These assays were specifically designed so that up to 7 clonality targets can be multiplexed together in a single LymphoTrack® MiSeq® run. Up to 5 independent targets (*IGH* FR1/FR2/FR3, *IGK*, and *TRG*) can be multiplexed together if using the LymphoTrack Assays for Ion S5™ or Ion PGM™.

Sample Batching: LymphoTrack® one-step PCR incorporates indices (molecular barcodes) onto each amplicon. By design, a unique index is used for each specimen. The sample indices are read during sequencing and used by the LymphoTrack® (MRD) Software to de-multiplex samples.

» LymphoTrack® Assay or Panel kits for the MiSeq® are provided with up to 24 different indices (up to 48 with *IGH* FR1). Thus, when tracking clonal sequences, these panels allow 22 different subject samples to be run along with 2 external controls on a single flow cell (or up to 46 different FR1 subject samples).

» LymphoTrack® Assay kits for the Ion S5/PGM™ are provided with 12 different indices, thus allowing 10 different subject samples to be run along with 2 external controls on a single sequencing chip.

7. LYMPHOTRACK® MRD SOFTWARE

The LymphoTrack® MRD Software seamlessly automates control and sample analyses. The powerful bioinformatics software facilitates the longitudinal assessment of up to 5 clonal populations by providing multiple functionalities to the user for project planning, data analyses and reporting. Project planner features assist the user in design of experiment and project management with a built in "project" feature to save and load "projects" as additional time points are tested. When Invivoscribe® LymphoTrack® MRD controls are used, the software automates MRD tracking, reports clonal counts and frequencies, and provides multiple graphs to visualize changes in clonal frequencies over time. **Refer to the LymphoTrack® MRD Software Technical Bulletin for additional details.**

PROJECT PLANNER TOOL

This application provides design of experiment features including a confidence interval calculator. The software embedded calculator considers sample replicate count, resequencing count, read depth, and DNA input amount to determine the confidence level of a true MRD negative sample. This tool also allows the user the ability to customize the desired levels of confidence to ensure the experimental design meets the user-defined specification for combined sensitivity and confidence for MRD monitoring.

AUTOMATED BIOINFORMATICS DATA ANALYSES

The LymphoTrack® MRD software facilitates tracking of up to 5 defined clonal sequences as well as the degree of mismatches (0, 1, 2). It further calculates clonal read frequency, analyzes multiple samples and/or resequencing replicates, and reports the level of confidence when the clonal sequence is not detected at various levels.

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
Figure 3: Estimated MRD Levels

NOTE: Graphical representation only

LYMPHOTRACK® MRD SAMPLE REPORT

The MRD Sample Report provides a summary table for the latest time point for all clonal sequences. Details are included for each sequence analyzed including the Clonal Frequency if DETECTED and the % Confidence if NOT DETECTED. Detailed information applicable to analysis is summarized for the most recent time point as follows:

- » Clonal Sequence(s) Detected or Not Detected
- » Number of PCR Replicates
- » Total Read Number Analyzed
- » Cumulative Target Read Number
- » LymphoQuant® Internal Control Read Number
- » When a clonal sequence of interest is Not Detected, the call confidence is reported for various sensitivities.

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©2021 Invivoscribe, Inc (Invivoscribe®) 無断複写・転載を禁じます。ここに記載されている商標はInvivoscribe®およびその関連会社、またはここに記載されている他社の商標については、それぞれその所有者に帰属します。ILLUMINA®およびMiSeq®はIllumina, Inc.の登録商標です。Ion S5TMおよびIon PGM™はThermo Fisher Scientificの商標です。
Invivoscribe, Inc.