

Instructions for Use




## IdentiClone® *TCRB* + *TCRG* T-cell Clonality Assay

For identification of clonal T-cell receptor beta chain and gamma chain gene rearrangements.

**IVD** For *In Vitro* Diagnostic Use



 Storage Conditions: **-85°C to -65°C**

(DNA controls may be separated from assay kits and stored at 2°C to 8°C)

Catalog#	Products	Quantity
<b>REF</b> 92000011	IdentiClone <i>TCRB</i> + <i>TCRG</i> T-cell Clonality Assay – ABI Fluorescence Detection	33 Reactions
<b>REF</b> 92000021	IdentiClone <i>TCRB</i> + <i>TCRG</i> T-cell Clonality MegaKit – ABI Fluorescence Detection	330 Reactions

## Table of Contents

<b>1.</b>	<b>INTENDED USE .....</b>	<b>3</b>
<b>2.</b>	<b>SUMMARY AND EXPLANATION OF THE TEST .....</b>	<b>3</b>
2.1.	Background.....	3
2.2.	Summary .....	3
<b>3.</b>	<b>PRINCIPLES OF THE PROCEDURE .....</b>	<b>4</b>
3.1.	Polymerase Chain Reaction (PCR).....	4
3.2.	Differential Fluorescence Detection .....	4
<b>4.</b>	<b>REAGENTS .....</b>	<b>5</b>
4.1.	Reagent Components .....	5
4.2.	Warnings and Precautions .....	6
4.3.	Storage and Handling.....	6
<b>5.</b>	<b>INSTRUMENTS .....</b>	<b>7</b>
5.1.	Thermal cycler .....	7
5.2.	ABI Capillary Electrophoresis Instruments .....	7
<b>6.</b>	<b>SPECIMEN COLLECTION AND PREPARATION .....</b>	<b>8</b>
6.1.	Precautions .....	8
6.2.	Interfering Substances .....	8
6.3.	Specimen Requirements and Handling.....	8
6.4.	Sample Preparation .....	8
6.5.	Sample Storage.....	8
<b>7.</b>	<b>ASSAY PROCEDURE .....</b>	<b>9</b>
7.1.	Materials Provided.....	9
7.2.	Materials Required (not provided) .....	9
7.3.	Reagent Preparation .....	10
7.4.	Amplification.....	11
7.5.	ABI Fluorescence Detection.....	11
7.6.	Quality Control .....	12
7.7.	Recommended Positive Controls .....	12
<b>8.</b>	<b>INTERPRETATION OF RESULTS.....</b>	<b>13</b>
8.1.	Analysis.....	13
8.2.	Sample Interpretation .....	14
<b>9.</b>	<b>LIMITATIONS OF PROCEDURE .....</b>	<b>14</b>
<b>10.</b>	<b>EXPECTED VALUES.....</b>	<b>14</b>
10.1.	Expected Size of Amplified Products.....	14
10.2.	Sample Data.....	15
<b>11.</b>	<b>PERFORMANCE CHARACTERISTICS .....</b>	<b>16</b>
<b>12.</b>	<b>TECHNICAL AND CUSTOMER SERVICE .....</b>	<b>17</b>
<b>13.</b>	<b>REFERENCES.....</b>	<b>17</b>
<b>14.</b>	<b>SYMBOLS.....</b>	<b>17</b>
<b>15.</b>	<b>LEGAL NOTICE.....</b>	<b>18</b>
15.1.	Warranty and Liability .....	18
15.2.	Patents and Trademarks .....	18

# 1. Intended Use

The IdentiClone *TCRB* + *TCRG* T-cell Clonality Assay is an *in vitro* diagnostic product intended for PCR-based detection of clonal T-cell receptor beta chain and gamma chain gene rearrangements in patients with suspect lymphoproliferations. Specifically, the *TCRB* + *TCRG* T-Cell Clonality Assay can be used to:

- Identify clonality in suspect lymphoproliferations
- Support a differential diagnosis between reactive lesions and T-cell and some immature B-cell malignancies
- Determine lineage involvement in mature lymphoproliferative disorders
- Monitor and evaluate disease recurrence

# 2. Summary and Explanation of the Test

## 2.1. Background

Rearrangements of the antigen receptor genes occur during ontogeny in B and T lymphocytes. These gene rearrangements generate products that are unique in length and sequence for each cell. Therefore, polymerase chain reaction (PCR) assays can be used to identify lymphocyte populations derived from a single cell by detecting the unique V-J gene rearrangements present within these antigen receptor loci.<sup>1</sup> This IdentiClone PCR assay employs multiple consensus DNA primers that target conserved genetic regions within the T cell receptor beta chain and gamma chain genes. This test is used to detect the vast majority of clonal T-cell malignancies from DNA. Test products can be analyzed using a variety of detection formats, including gel and capillary electrophoresis.

Gene rearrangement analysis can also be performed by Southern Blot (SB)-based techniques. Although SB analysis is very reliable, it is increasingly replaced by PCR techniques because of the greater efficiency and sensitivity of PCR. Moreover, PCR is relatively easy, less labor intensive and requires much lower quantities of high molecular weight DNA than SB tests. In addition, PCR can often be performed on DNA isolated from paraffin-embedded tissue, whereas SB cannot be performed because the DNA is often degraded. Therefore, there is a strong need to replace SB analysis with reliable PCR techniques.

## 2.2. Summary

Invivoscribe's IdentiClone assays represent a new approach to PCR-based clonality testing. These standardized assays were carefully optimized testing positive and negative control samples using multiplex master mixes. Assay development was followed by extensive validation including the testing of more than 400 clinical samples using Revised European/American Lymphoma (REAL) Classification. Testing was performed at more than thirty prominent independent testing centers throughout Europe in a collaborative study known as the BIOMED-2 Concerted Action. Results from this BIOMED-2 study appear in *Leukemia*, a leading peer-reviewed journal<sup>2</sup>. In a *Leukemia* 2007 article, testing for both *TCRB* and *TCRG* gene rearrangements led to 94% sensitivity, compared to 91% for *TCRB* and 89% for *TCRG* when tested alone. It may also increase the reliability of the tests as it is more likely that the clonal products will be detected in more than one tube.<sup>4</sup>

The ABI detection based assays cannot reliably detect clonal populations comprising less than 1% of the total lymphocyte cell population. Always interpret the results of molecular clonality tests in the context of clinical, histological and immunophenotypic data.

This test kit includes 6 master mixes. *TCRB* Tubes A and B target framework regions within the variable region and the joining region of the *TCR* beta chain locus. *TCRB* Tube C targets the diversity and joining regions of the *TCR* beta chain locus. *TCRG* Tube A contains primers that target the *Vγ1-8* + *Vγ10* genes and *Jγ1.1*, *Jγ1.3*, *Jγ2.1* and *Jγ2.3* genes (also known as *JγP1*, *Jγ1*, *JγP2* and *Jγ2* respectively). *TCRG* Tube B contains primers that target the *Vγ9* + *Vγ11* genes and *Jγ1.1*, *Jγ1.3*, *Jγ2.1* and *Jγ2.3* genes. Lastly, the Specimen Control Size Ladder master mix, targets multiple genes and generates a series of amplicons of approximately 96, 197, 297, 397 and 602 base pairs to ensure that the quality and quantity of input DNA is adequate to yield a valid result. A single thermal cycler program and similar detection methodologies are used with all of our Gene Clonality Assays. This improves consistency and facilitates cross training on a broad range of different assays.


This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.



### 3. Principles of the Procedure

#### 3.1. Polymerase Chain Reaction (PCR)

PCR assays are routinely used for the identification of clonal T-cell populations. These tests amplify the DNA between primers that target the conserved variable (V) regions and the conserved joining (J) regions (*TCRB* Tubes A & B and *TCRG* Tubes A & B), as well as the diversity (D) and joining regions (*TCRB* Tube C). These conserved regions lie on either side of an area within the V-J region where programmed genetic rearrangements occur during maturation of all B and T lymphocytes. The antigen receptor genes that undergo rearrangement are the immunoglobulin heavy chain and light chains in B-cells and the T cell receptor genes in T-cells. Each B- and T-cell has a single productive V-J rearrangement that is unique in both length and sequence. Therefore, when DNA from a normal or polyclonal population is amplified using DNA primers that flank the V-J region, a bell-shaped curve (Gaussian distribution) of amplicon products within an expected size range is produced. This Gaussian distribution reflects the heterogeneous population of V-J rearrangements. (In certain cases, where lymphocyte DNA is not present, no product is seen.) For DNA from samples containing a clonal population, the yield is one or two prominent amplified products (amplicons) within a diminished polyclonal background.



**Figure 1.** This is a simplified diagram of a representative rearranged T cell receptor beta gene and the T cell receptor gamma gene showing the approximate placement of the upstream and downstream DNA primers. The numbers of primers and their specificity are listed for master mix *TCRB* Tubes A, B and C and *TCRG* Tubes A and B. (The V $\gamma$ 1f primer is a consensus primer that targets V $\gamma$ 1 through V $\gamma$ 8)..

Since the antigen receptor genes are polymorphic (consisting of a heterogeneous population of related DNA sequences), it is difficult to employ a single set of DNA primer sequences to target all of the conserved flanking regions around the V-J rearrangement. N-region diversity and somatic mutation further scramble the DNA sequences in these regions. Therefore, multiplex master mixes, which target several FR regions, are required to identify the majority of clonal rearrangements. As indicated, clonal rearrangements are identified as prominent, single-sized products within the background of different-sized amplicon products that form a Gaussian distribution around a statistically favored, average-sized rearrangement.

#### 3.2. Differential Fluorescence Detection

Differential fluorescence detection is commonly used to resolve the different-sized amplicon products using a capillary electrophoresis instrument. Primers can be conjugated with several different fluorescent dyes (fluorophores) so that they can produce different emission spectra upon excitation by a laser in the capillary electrophoresis instrument. In this manner, different fluorescent dyes can correspond to different targeted regions. This detection system results in unsurpassed sensitivity, single nucleotide resolution, differential product detection and relative quantification. In addition, the use of agarose and polyacrylamide gels, as well as the use of carcinogens such as ethidium bromide, can virtually be eliminated. Further, differential detection allows accurate, reproducible and objective interpretation of

primer-specific products and automatic archiving of data. Inter-assay and intra-assay reproducibility in size determination using capillary electrophoresis is approximately 1 to 2 base pair. This reproducibility and sensitivity coupled with the automatic archiving of specimen data allows for the monitoring, tracking and comparison of data from individual patients over time.



## 4. Reagents

### 4.1. Reagent Components

Table 1. Available Kits

Catalog #	Product	Quantity
<b>REF</b> 92000011	IdentiClone <i>TCRB</i> + <i>TCRG</i> Gene Clonality Assay – ABI Fluorescence Detection	33 Reactions
<b>REF</b> 92000021	IdentiClone <i>TCRB</i> + <i>TCRG</i> Gene Clonality Assay MegaKit – ABI Fluorescence Detection	330 Reactions

Table 2. Reagent Components

Reagent	Catalog #	Reagent Components (active ingredients)	Unit Quantity	92000011 # of Units	92000021 # of Units	Storage Temp.
Master Mixes	22050011CE	<b><i>TCRB</i> Tube A – 6FAM &amp; HEX</b> Multiple oligonucleotides targeting the V $\beta$ + J $\beta$ 1 + J $\beta$ 2 regions of the T cell receptor beta gene in a buffered salt solution.	1500 $\mu$ L	1	10	
	22050021CE	<b><i>TCRB</i> Tube B – 6FAM</b> Multiple oligonucleotides targeting the V $\beta$ + J $\beta$ 2 regions of the T cell receptor beta gene in a buffered salt solution.	1500 $\mu$ L	1	10	
	22050031CE	<b><i>TCRB</i> Tube C – 6FAM &amp; HEX</b> Multiple oligonucleotides targeting the D $\beta$ + J $\beta$ 1 + J $\beta$ 2 regions of the T cell receptor beta gene in a buffered salt solution.	1500 $\mu$ L	1	10	
	22070031CE	<b><i>TCRG</i> Tube A – 6FAM &amp; HEX</b> Multiple oligonucleotides targeting the V $\gamma$ 1-8, V $\gamma$ 10 + multiple J $\gamma$ regions of the T cell receptor gamma gene in a buffered salt solution.	1500 $\mu$ L	1	10	
	22070041CE	<b><i>TCRG</i> Tube B – 6FAM &amp; HEX</b> Multiple oligonucleotides targeting the V $\gamma$ 9, V $\gamma$ 11 + multiple J $\gamma$ regions of the T cell receptor gamma gene in a buffered salt solution.	1500 $\mu$ L	1	10	
Template Amplification Control Master Mix	20960021	<b>Specimen Control Size Ladder – 6FAM</b> Multiple oligonucleotides targeting housekeeping genes.	1500 $\mu$ L	1	10	
Positive Control DNAs	40881210	<b>IVS-0021 Clonal Control DNA</b> 200 $\mu$ g/mL of DNA in 1/10 <sup>th</sup> TE solution	100 $\mu$ L	1	5	
	40880490	<b>IVS-0009 Clonal Control DNA</b> 200 $\mu$ g/mL of DNA in 1/10 <sup>th</sup> TE solution	100 $\mu$ L	1	5	
	40880190	<b>IVS-0004 Clonal Control DNA</b> 200 $\mu$ g/mL of DNA in 1/10 <sup>th</sup> TE solution	100 $\mu$ L	1	5	
Negative (Normal) Control DNA	40920010	<b>IVS-0000 Polyclonal Control DNA</b> 200 $\mu$ g/mL of DNA in 1/10 <sup>th</sup> TE solution	100 $\mu$ L	1	5	

**Note:** There are no preservatives used in the manufacture of this kit.

## 4.2. Warnings and Precautions

- **IVD** This product is For *In Vitro* Diagnostic Use
- Use the assay kit as a system; do not substitute other manufacturer's reagents. Dilution, reducing amplification reaction volumes or other deviation in this protocol may affect the performance of this test and/or nullify any limited sublicense that comes with the purchase of this testing kit.
- Materials are stable until the labeled expiration date when stored and handled as directed. Do not use kits beyond their expiration date.
- Close adherence to the protocol will assure optimal performance and reproducibility. Ensure correct thermal cycler program is used, as other programs may provide inaccurate/faulty data, such as false positive and false negative results.
- Do not mix or combine reagents from kits with different lot numbers.
- Laboratory personnel are reminded to wear appropriate personal protective equipment and follow good laboratory practices and universal precautions when working with specimens. Handle specimens in approved biological safety containment facilities and open only in certified biological safety cabinets. Use molecular biology grade water with the preparation of specimen DNA.
- Due to the analytical sensitivity of this test, use extreme care to avoid the contamination of reagents or amplification mixtures with samples, controls or amplified materials. Closely monitor all reagents for signs of contamination (e.g., negative controls giving positive signals). Discard reagents suspected of contamination.
- To minimize contamination, wear clean gloves when handling samples and reagents and routinely clean work areas and pipettes prior to doing PCR.
- Autoclaving does not eliminate DNA contamination. Follow uni-directional work flow between separate areas of the PCR laboratory: begin with Master Mix Preparation, move to the Specimen Preparation, then to the Amplification and finally to Detection. Do not bring amplified DNA into the areas designated for master mix or specimen preparation.
- Dedicate all pipettes, pipette tips and any equipment used in a particular area to that area of the laboratory.
- Use sterile, disposable plastic ware whenever possible to avoid RNase, DNase or cross-contamination.

## 4.3. Storage and Handling

- For any duration other than immediate use, **store assay kits at -85°C to -65°C.**
- The optimum storage temperature for DNA controls is 2°C to 8°C, but DNA controls can be stored at -85°C to -65°C.
- All reagents and controls must be thawed and vortexed or mixed thoroughly prior to use to ensure that they are resuspended completely. Excessive vortexing may shear DNA and cause labeled primers to lose their fluorophores.
- Materials are stable until the labeled expiration date when stored and handled as directed. Do not use kits beyond their expiration date.
- Due to high salt concentrations, PCR master mixes are sensitive to freeze/thaw cycles. Aliquot master mixes into sterile o-ring screw-cap tubes if necessary.

## 5. Instruments

### 5.1. Thermal cycler

- Use or Function: Amplification of DNA samples
- Suggested Instrument: Veriti™ Thermal Cycler or equivalent
- Performance Characteristics and Specification:
  - Minimum Thermal Range: 15°C to 96°C
  - Minimum Ramping Speed: 0.8°C/sec
- Follow manufacturer's installation, operation, calibration and maintenance procedures.
- See section 7.4 *Amplification* for thermal cycler program.

### 5.2. ABI Capillary Electrophoresis Instruments

- Use or Function: Fragment detection and analysis
- Performance Characteristics and Specification:
  - The following capillary electrophoresis instruments will meet the performance needs for this assay:
    - ABI 310 Genetic Analyzer (1-capillary)
    - ABI 3100 Avant Genetic Analyzer (4-capillaries)
    - ABI 3100 Genetic Analyzer (16-capillaries)
    - ABI 3130 Genetic Analyzer (4-capillaries)
    - ABI 3130xL Genetic Analyzer (16-capillaries)
    - ABI 3500 Genetic Analyzer (8-capillaries)
    - ABI 3500xL Genetic Analyzer (24-capillaries)
- Follow manufacturer's installation, operation, calibration and maintenance procedures.
- The ABI instrument used must be calibrated with appropriate Matrix Standards as outlined in section 7.2 *Materials Required (not provided)*.
- Use the default settings for your polymer and capillary type.
- See section 7.5 ABI Fluorescence Detection.

## 6. Specimen Collection and Preparation

### 6.1. Precautions

Biological specimens from humans may contain potentially infectious materials. Handle all specimens in accordance with the OSHA Standard on Bloodborne Pathogens or Biosafety Level 2.

### 6.2. Interfering Substances

The following substances are known to interfere with PCR:

- Divalent cation chelators
- Low retention pipette tips
- EDTA (not significant at low concentrations)
- Heparin

### 6.3. Specimen Requirements and Handling

This assay tests **genomic DNA** from the following sources:

- 5 cc of peripheral blood, bone marrow biopsy or bone marrow aspirate anti-coagulated with heparin or EDTA (stored at 2°C to 8°C and shipped at ambient temperature)
- Minimum 5 mm cube of tissue (stored and shipped frozen; or stored and shipped in RPMI 1640 at ambient temperature or on ice)
- 3 µg of genomic DNA (stored at 2°C to 8°C and shipped at ambient temperature)
- Formalin-fixed paraffin embedded tissue or slides (stored and shipped at ambient temperature)

### 6.4. Sample Preparation

Extract the genomic DNA from patient specimens as soon as possible. Resuspend DNA to a final concentration of 100 µg to 400 µg per mL in 1/10<sup>th</sup> TE (1 mM Tris-HCl, pH 8.0; 0.1 mM EDTA) or in molecular biology grade or USP water. This is a robust assay system. A wide range of DNA concentrations will generate a valid result. Therefore, quantifying and adjusting DNA concentrations is generally not necessary. Testing sample DNA with the Specimen Control Size Ladder master mix will ensure that DNA of sufficient quality and quantity was present to yield a valid result.

### 6.5. Sample Storage











Store genomic DNA at 2°C to 8°C or at -85°C to -65°C until use.



## 7. Assay Procedure

### 7.1. Materials Provided

Table 3. Materials Provided

Catalog #	Description
 22050011CE	TCRB Tube A – 6FAM & HEX
 22050021CE	TCRB Tube B – 6FAM
 22050031CE	TCRB Tube C – 6FAM & HEX
 22070031CE	TCRG Tube A – 6FAM & HEX
 22070041CE	TCRG Tube B – 6FAM & HEX
 20960021	Specimen Control Size Ladder – 6FAM
 40881210	IVS-0021 Clonal Control DNA
 40880490	IVS-0009 Clonal Control DNA
 40880190	IVS-0004 Clonal Control DNA
 40920010	IVS-0000 Polyclonal Control DNA

### 7.2. Materials Required (not provided)

Table 4. Materials Required (not provided)

Reagent/Material	Recommended Reagents/Materials and Suppliers	Catalog #	Notes
DNA Polymerase	Roche: <ul style="list-style-type: none"> <li>EagleTaq DNA Polymerase</li> </ul> Invivoscribe: <ul style="list-style-type: none"> <li>FalconTaq DNA Polymerase or equivalent</li> </ul>	05206944190  60970130	N/A
Glass Distilled De-ionized Molecular Biology Grade or USP Water	N/A	N/A	DNase / RNase free
Calibrated Pipettes	Rainin: <ul style="list-style-type: none"> <li>P-2, P-20, P-200 and P-1000 pipettes</li> <li>Or SL-2, SL-20, SL-200 and SL-1000 pipettes</li> </ul>	N/A	Must be able to accurately measure volumes between 1 µL and 1000 µL.
Thermal cycler	Thermo Fisher Scientific: <ul style="list-style-type: none"> <li>Veriti Dx Thermal Cycler</li> </ul> Bio-Rad: <ul style="list-style-type: none"> <li>MJ Research PTC-100 or PTC-200, PTC-220, PTC-240</li> </ul> Perkin-Elmer <ul style="list-style-type: none"> <li>PE 9600 or PE 9700</li> </ul>	N/A	N/A
Vortex Mixer	N/A	N/A	N/A
PCR plates or tubes	N/A	N/A	Sterile
Filter barrier pipette tips	N/A	N/A	Sterile, RNase/DNase/Pyrogen-free
Microcentrifuge tubes	N/A	N/A	Sterile
ABI Capillary Electrophoresis Instrument	Thermo Fisher Scientific: <ul style="list-style-type: none"> <li>ABI 310, 3100 or 3500 series</li> </ul>	N/A	N/A
Hi-Di Formamide	Thermo Fisher Scientific: <ul style="list-style-type: none"> <li>Hi-Di™ Formamide</li> </ul>	4311320	N/A
Size Standards	Invivoscribe: <ul style="list-style-type: none"> <li>Hi-Di Formamide w/ROX size standards for ABI 3100</li> </ul> Thermo Fisher Scientific: <ul style="list-style-type: none"> <li>For ABI 3100 or 3130 instruments:               <ul style="list-style-type: none"> <li>GeneScan™ - 400HD [ROX]™</li> </ul> </li> <li>For ABI 3500 instruments:               <ul style="list-style-type: none"> <li>GeneScan - 600 [LIZ]™ v2.0</li> </ul> </li> </ul>	60980061  402985  4408399	N/A

Table 4. Materials Required (not provided)

Reagent/Material	Recommended Reagents/Materials and Suppliers	Catalog #	Notes
Spectral Calibration Dye Sets	Thermo Fisher Scientific: <ul style="list-style-type: none"> <li>For ABI 3100 and 3130 instruments: <ul style="list-style-type: none"> <li>DS-30 Matrix Standard Kit (Dye Set D)</li> </ul> </li> <li>For ABI 310 instruments: <ul style="list-style-type: none"> <li>NED Matrix Standard</li> <li>And Fluorescent Amidite Matrix Standards [6FAM, TET, HEX, TAMRA, ROX]</li> </ul> </li> <li>For ABI 3500 instruments: <ul style="list-style-type: none"> <li>DS-33 Matrix Standard Kit (Dye Set G5)</li> </ul> </li> </ul>	4345827  402996 401546  4345833	N/A
Polymer	Thermo Fisher Scientific: <ul style="list-style-type: none"> <li>POP-4™ Polymer: <ul style="list-style-type: none"> <li>POP-4 for 310 Genetic Analyzers</li> <li>POP-4 for 3100/3100-Avant Genetic Analyzers</li> <li>POP-4 for 3130/3130xL Genetic Analyzers</li> </ul> </li> <li>POP-7™ Polymer: <ul style="list-style-type: none"> <li>POP-7 for 3130/3130xL Genetic Analyzers</li> <li>POP-7 for 3500/3500xL Genetic Analyzers</li> </ul> </li> </ul>	402838 4316355 4352755  4352759 4393714	N/A
Buffer	Thermo Fisher Scientific: <ul style="list-style-type: none"> <li>10X Genetic Analyzer Buffer with EDTA</li> </ul>	402824	Dilute 1:10 in sterile water before use

### 7.3. Reagent Preparation

- Test all unknown samples using the Specimen Control Size Ladder master mix. This is to ensure that no inhibitors of amplification are present and there is DNA of sufficient quality and quantity to generate a valid result.
- Singlicate test results are valid; however, we recommend **duplicate** testing when possible. If duplicate testing provides inconsistent results, re-testing or re-evaluation of the sample is necessary.
- Test positive, negative and no template** controls for each of the master mixes.
- It is recommended to batch multiple samples in a run to avoid running out of the negative control (IVS-0000 Polyclonal Control DNA). If batching samples is not practical in your laboratory setting, IVS-0000 Polyclonal Control DNA is also available for purchase separately.

7.3.1. Using gloved hands, remove the master mixes from the freezer. Allow the tubes to thaw completely; then gently vortex to mix.

7.3.2. In containment hood or dead air box remove an appropriate aliquot from each master mix to individual clean, sterile microcentrifuge tubes.

- Aliquot volumes are 45 µL for each reaction.
- Add an additional reaction for every 15 reactions to correct for pipetting errors.
- For each master mix (except for the Specimen Control Size Ladder), the number of reactions (**n**) is:

<b>n = 2 × # of samples</b>	(run each sample in duplicate)
+ 1	positive control DNA (refer to Table 6)
+ 1	negative control DNA (IVS-0000 Polyclonal Control DNA)
+ 1	no template control (water)
+ 1	to correct for pipetting errors

---

**n = 2 × # of samples + 4      Total**

- Therefore, the total aliquot volume for each master mix is **n × 45 µL**.
- For the Specimen Control Size Ladder master mix, the number of reactions (**m**) is:

<b>m = # of samples</b>	(run each sample in singlicate)
+ 1	negative control DNA (IVS-0000 Polyclonal Control DNA)
+ 1	no template control (water)
+ 1	to correct for pipetting errors

---

**m = # of samples + 3      Total**

- Therefore, the total aliquot volume for the Specimen Control Size Ladder master mix is **m × 45 µL**.

- 7.3.3. **For *TCRB* Tubes A and B:** Add 2.25 U (or 0.45  $\mu\text{L}$  @5 U/ $\mu\text{L}$ ) of Taq DNA polymerase per reaction to each master mix.
- The total Taq DNA polymerase added to each master mix is  $n \times 0.45 \mu\text{L}$ .
  - Gently vortex to mix.
- 7.3.4. For *TCRB* Tube C and *TCRG* Tubes A and B and Specimen Control Size Ladder: Add 1.25 U (or 0.25  $\mu\text{L}$  @5 U/ $\mu\text{L}$ ) of Taq DNA polymerase per reaction.
- The total Taq DNA polymerase added to each master mix is  $n \times 0.25 \mu\text{L}$  for the *TCRB* Tube C master mix and  $m \times 0.25 \mu\text{L}$  for the Specimen Control Size Ladder master mix.
  - Gently vortex to mix.
- 7.3.5. For each reaction, aliquot 45  $\mu\text{L}$  of the appropriate master mix + DNA polymerase solution into individual wells in a PCR plate or tube.
- 7.3.6. Add 5  $\mu\text{L}$  of appropriate template (sample DNA, positive control DNA, negative control DNA or water) to the individual wells containing the respective master mix solutions. Pipette up and down several times to mix.
- 7.3.7. Cap or cover the PCR plate.
- Samples are now ready to be amplified on a thermal cycler.
  - If amplification cannot be performed immediately following reagent preparation, the PCR plate or tubes can be stored at 2°C to 8°C for up to 24 hours.

#### Quick Guide

For each master mix and n reactions, mix:

$n \times 45 \mu\text{L}$

Master Mix

$n \times 0.25 \mu\text{L}$  or  $0.45 \mu\text{L}^*$

Taq DNA polymerase

Vortex gently to mix.

Aliquot 45  $\mu\text{L}$  of master mix + Taq DNA polymerase solution into each reaction well.

Add 5  $\mu\text{L}$  of appropriate Template to each well.

**Total reaction volume = 50  $\mu\text{L}$**

**\*Note:** Use 0.45  $\mu\text{L}$  of Taq DNA polymerase for *TCRB* Tubes A and B and 0.25  $\mu\text{L}$  of Taq DNA polymerase for *TCRB* Tube C and Specimen Control Size Ladder.

## 7.4. Amplification

- 7.4.1. Amplify the samples using the following PCR program:

- Use the **calculated** option for temperature measurement with the BioRad MJ PTC thermal cyclers.

**Table 5.** Thermal cycling conditions

Step	Temperature	Duration	Cycles
1	95°C	7 minutes	1
2	95°C	45 seconds	35
3	60°C	45 seconds	
4	72°C	90 seconds	
5	72°C	10 minutes	1
6	15°C	$\infty$	1

- 7.4.2. Remove the amplification plate or tubes from the thermal cycler.

- Although amplified DNA is stable at room temperature for extended periods of time, store PCR products at 2°C to 8°C until detection. Detection must be within 30 days of amplification.

## 7.5. ABI Fluorescence Detection

Please note that for ABI fluorescence detection a preceding peak is often seen and is an artifact due to the detection method the ABI platforms use. Preceding peaks are sometimes skewed and have bases that slope on the right side towards the real peak. This is especially evident in the Specimen Control Size Ladder master mix where the 96-base pair peak has a preceding peak that shows up at 84 base pair.

- Do not multiplex PCR products from different master mixes; this will result in overall reduced sensitivity of the assay.

### ABI 310, 3100 or 3130 Platforms

- 7.5.1. In a new microcentrifuge tube, mix an appropriate amount (for a total of 10 µL per reaction) of Hi-Di Formamide with ROX Size Standards. Vortex well.
- 7.5.2. In a new 96-well PCR plate, add 10 µL of Hi-Di Formamide with ROX size standards to individual wells for each PCR.
- 7.5.3. Transfer 1 µL of each reaction to the wells containing Hi-Di Formamide with ROX size standards. Add only one sample per well. Pipette up and down to mix.
- 7.5.4. Cap or cover the PCR plate or tubes.
- 7.5.5. Heat denature the samples at 95°C for 2 minutes then snap chill on ice for 5 minutes.
- 7.5.6. Prepare a **sample sheet** and **injection list** for the samples.
- 7.5.7. Run the samples on an ABI capillary electrophoresis instrument according to the user manual.
  - Data are automatically displayed as size and color specific peaks.
- 7.5.8. Review profile and controls, report results. (See sections 8 *Interpretation of Results* and 10 *Expected Values*)

### ABI 3500 Platforms:

**Note:** Due to instrument to instrument variation in the performance of the ABI 3500 platform, the amount of formamide, sample and size standard listed in the protocol is intended to be a starting point. The protocol may need to be optimized for specific ABI 3500 Platforms.

- 7.5.9. In a new microcentrifuge tube, mix an appropriate amount (9.5 µL per reaction) of Hi-Di Formamide with LIZ Size Standards. Vortex well.
- 7.5.10. In a new 96-well PCR plate, add 9.5 µL of Hi-Di Formamide with LIZ size standards to individual wells for each PCR reaction.
- 7.5.11. Transfer 0.5 µL of each reaction to the wells containing Hi-Di Formamide with LIZ size standards. Add only one sample per well. Pipette up and down to mix.
- 7.5.12. Cap or cover the PCR plate.
- 7.5.13. Heat denature the samples at 95°C for 3 minutes, then snap chill on ice for 5 minutes.
- 7.5.14. Prepare a sample sheet and injection list for the samples.
- 7.5.15. Run the samples on an ABI 3500 capillary electrophoresis instrument according to the user manual.
- 7.5.16. Data are automatically displayed as size and color specific peaks. Review profile and controls, report results. (see sections 8 *Interpretation of Results* and 10 *Expected Values*)

### 7.6. Quality Control

Positive and negative (or normal) controls are furnished with the kit and are run in singlicate each time the assay is performed to ensure proper performance of the assay. In addition, a no template control (e.g. water) is also included to test for contamination of the master mix or cross-contamination of PCR reactions due to improper sterile technique. A buffer control may also be added to ensure that no contamination of the buffer used to resuspend the samples has occurred. The values for the positive controls are provided under section 10.1 *Expected Size of Amplified Products*. Additional controls and sensitivity controls (dilutions of positive controls into our negative control) are available from Invivoscribe.

### 7.7. Recommended Positive Controls

The amplicon sizes listed were determined using an ABI platform. Amplicon sizes seen on your specific capillary electrophoresis instrument may differ 1 to 4 base pair (bp) from those listed depending on the platform of detection and the version of the analysis software used. Once identified, the amplicon size as determined on your specific platform will be consistent from run to run. This reproducibility is extremely useful when monitoring disease recurrence.

**Note:** “Color” indicates the color of products generated with the master mix when using the default color assignment on ABI fluorescence detection systems. The amplicon sizes listed were determined using an ABI platform.

Table 6. Recommended Positive Controls

Master Mix	Target	Color	Control DNA	Catalog #	Product Size (bp)
<i>TCRB</i> Tube A	Vβ + Jβ1/2	Blue (Jβ2.X) + Green (Jβ1.X)	Valid Size Range IVS-0009 Clonal Control DNA	--- 40880490	240 - 285 264
<i>TCRB</i> Tube B	Vβ + Jβ2	Blue (Jβ2.X)	Valid Size Range IVS-0004 Clonal Control DNA IVS-0021 Clonal Control DNA	--- 40880190 40881210	240 - 285 253 267
<i>TCRB</i> Tube C	Dβ + Jβ1/2	Blue (Jβ2.X) + Green (Jβ1.X)	Valid Size Range IVS-0009 Clonal Control DNA	--- 40880490	170 - 210 (Dβ2), 285 - 325 (Dβ1) 309
<i>TCRG</i> Tube A	Vγ1-8 + Vγ10 + multiple Jγ regions	Blue (Jγ 1.1/2.1) + Green (Jγ 1.3/2.3)	Valid Size Range IVS-0021 Clonal Control DNA	--- 40881210	145 - 255 211 (Vγ1-8 + Jγ1.3/2.3)
<i>TCRG</i> Tube B	Vγ9 Vγ11 + multiple Jγ regions	Blue (Jγ 1.1/2.1) + Green (Jγ 1.3/2.3)	Valid Size Range IVS-0021 Clonal Control DNA	--- 40881210	80 - 220 167 (Vγ9 + Jγ1.3/2.3)
Specimen Control Size Ladder	Multiple Genes	Blue	Valid Size Range IVS-0000 Polyclonal Control DNA	--- 40920010	96, 197, 297, 397, 602 <sup>a</sup> 96, 197, 297, 397, 602 <sup>a</sup>

<sup>a</sup>Note: Because smaller PCR fragments are preferentially amplified, it is not unusual for the 602 bp fragment to have a diminished signal or to be missing entirely.

## 8. Interpretation of Results

Although positive results are highly suggestive of malignancy, interpret both positive and negative results in the context of all clinical information and laboratory test results. The size range for each of the master mixes has been determined by testing positive and negative control samples. For accurate and meaningful interpretation it is important to ignore peaks that occur outside of the valid size range for each of the master mixes.

### 8.1. Analysis

- 8.1.1. Report samples that fail to amplify following repeat testing as “A result cannot be reported on this specimen because there was DNA of insufficient quantity or quality for analysis.”
- 8.1.2. Repeat samples that test negative if the positive control reaction failed.
- 8.1.3. If samples run in duplicate yield differing results, re-test and/or re-evaluate the samples for sample switching.
- 8.1.4. All assay controls must be examined prior to interpretation of sample results.

Table 7. The following describes the analysis of each of the controls and the decisions necessary based upon the results.

Type of Control	Expected Result	Aberrant Result
No Template Control	No amplification present, continue with analysis	Amplification present. Repeat the assay.
Polyclonal Control	Product size is consistent with expected size listed in section 10.1 <i>Expected Size of Amplified Products</i> . No clonal rearrangements are present. Continue with analysis.	Clonal rearrangements are present. Repeat the assay
Positive Control (This can also be an extraction control if positive control material is taken through extraction processes)	Product size is consistent with expected size listed in section 10.1 <i>Expected Size of Amplified Products</i> . Continue with analysis.	Product not within the valid size range. Repeat the assay.
Specimen Control Size Ladder (This amplification control is essential for samples of unknown quantity and quality.)	If all of the 96, 197, 297, 397 and 602 bp peaks are seen, continue with analysis. Because smaller PCR fragments are preferentially amplified, it is not unusual for the 600 bp fragment to have a diminished signal or to be missing entirely. Continue with analysis.	If no bands are seen, repeat the assay <u>unless specimen tests positive</u> . If only 1, 2, or 3 bands are seen, re-evaluate sample for DNA degradation <u>unless specimen tests positive</u> .

## 8.2. Sample Interpretation

Given that the controls produce expected results, interpret the clinical samples as follows:

- One or two prominent positive bands<sup>a</sup> within the valid size range are reported as:

**“Positive for the detection of clonal T cell receptor beta chain or gamma chain gene rearrangement(s) consistent with the presence of a clonal cell population.”** In the context of overall diagnostic criteria, clonal cell populations can indicate the presence of hematologic malignancy.”

- An absence of positive bands<sup>a</sup> within the valid size range is reported as:

**“Negative for the detection of clonal T cell receptor beta chain or gamma chain gene rearrangement(s).”**

<sup>a</sup>**Note:** Criteria for defining a positive band are as follows:

- Products generated from **diagnostic samples** that fall within the valid size range and are at least three times the amplitude of the third largest peak in the polyclonal background are consistent with a positive peak.
- Products generated from **samples collected after initial diagnosis** that fall within the valid size range and are either: 1) at least three times the amplitude of the third largest peak; or, 2) exceed the amplitude of adjacent neighboring peaks and are identical in size to clonal amplicon products previously generated from the same patient using the same master mix, are consistent with a positive peak. Products generated from samples that fall within the valid size range and produce a discrete band(s) distinct from any background smear are consistent with a positive band.

## 9. Limitations of Procedure

- This assay does not identify 100% of clonal cell populations.
- This assay cannot reliably detect less than 1 positive cell per 100 normal cells.
- Always interpret the results of molecular clonality tests in the context of clinical, histological and immunophenotypic data.
- PCR-based assays are subject to interference by degradation of DNA or to inhibition of PCR due to EDTA, heparin and other agents.

## 10. Expected Values

### 10.1. Expected Size of Amplified Products

The amplicon sizes listed were determined using an ABI platform. Amplicon sizes seen on your specific capillary electrophoresis instrument may differ 1 to 4 base pair (bp) from those listed depending on the platform of detection and the version of the analysis software used. Once identified, the amplicon size as determined on your specific platform will be consistent from run to run. This reproducibility is extremely useful when monitoring disease recurrence.

**Note:** “Color” indicates the color of products generated with the master mix when using the default color assignment on ABI fluorescence detection systems. The amplicon sizes listed were determined using an ABI platform.

**Table 8.** Expected Size of Amplified Products

Master Mix	Target	Color	Control DNA	Catalog #	Product Size (bp)
<i>TCRB</i> Tube A	V $\beta$ + J $\beta$ 1/2	Blue (J $\beta$ 2.X) + Green (J $\beta$ 1.X)	<b>Valid Size Range</b>	---	<b>240 - 285</b>
			IVS-0000 Polyclonal Control DNA	40920010	240 - 285, 271 <sup>a</sup>
			IVS-0009 Clonal Control DNA	40880490	264
			IVS-0004 Clonal Control DNA	40880190	295
<i>TCRB</i> Tube B	V $\beta$ + J $\beta$ 2	Blue (J $\beta$ 2.X)	<b>Valid Size Range</b>	---	<b>240 - 285</b>
			IVS-0000 Polyclonal Control DNA	40920010	240 - 285, 221 <sup>b</sup>
			IVS-0009 Clonal Control DNA	40880490	---
			IVS-0004 Clonal Control DNA	40880190	253
<i>TCRB</i> Tube C	D $\beta$ + J $\beta$ 1/2	Blue (J $\beta$ 2.X) + Green (J $\beta$ 1.X)	<b>Valid Size Range</b>	---	<b>170 - 210 (D<math>\beta</math>2), 285 - 325 (D<math>\beta</math>1)</b>
			IVS-0000 Polyclonal Control DNA	40920010	128 <sup>b</sup> , 170 - 210, 285 - 325, 337 <sup>b</sup>
			IVS-0009 Clonal Control DNA	40880490	309
			IVS-0004 Clonal Control DNA	40880190	295

Table 8. Expected Size of Amplified Products

Master Mix	Target	Color	Control DNA	Catalog #	Product Size (bp)
<b>TCRG Tube A</b>	Vγ1-8 + Vγ10 + multiple Jγ regions	Blue (Jγ 1.1/2.1) + Green (Jγ 1.3/2.3)	<b>Valid Size Range</b>	---	<b>145 - 255</b>
			IVS-0000 Polyclonal Control DNA	40920010	<b>230 - 255</b> (Vγ1-8 + Jγ 1.1/2.1), <b>195 - 230</b> (Vγ1-8 + Jγ 1.3/2.3), <b>175 - 195</b> (Vγ10 + Jγ 1.1/2.1), <b>145 - 175</b> (Vγ10 + Jγ 1.3/2.3)
			IVS-0009 Clonal Control DNA	40880490	212 (Vγ1-8 + Jγ 1.3/2.3)
			IVS-0021 Clonal Control DNA	40881210	211 (Vγ1-8 + Jγ 1.3/2.3)
<b>TCRG Tube B</b>	Vγ9 Vγ11 + multiple Jγ regions	Blue (Jγ 1.1/2.1) + Green (Jγ 1.3/2.3)	<b>Valid Size Range</b>	---	<b>80 - 220</b>
			IVS-0000 Polyclonal Control DNA	40920010	<b>195 - 220</b> (Vγ9 + Jγ 1.1/2.1), <b>160 - 195<sup>c</sup></b> (Vγ9 + Jγ 1.3/2.3), <b>110 - 140<sup>d</sup></b> (Vγ11 + Jγ 1.1/2.1), <b>80 - 110<sup>d</sup></b> (Vγ11 + Jγ 1.3/2.3)
			IVS-0009 Clonal Control DNA	40880490	115 <sup>e</sup> (Vγ11 + Jγ 1.3/2.3)
			IVS-0021 Clonal Control DNA	40881210	167 (Vγ9 + Jγ 1.3/2.3)
<b>Specimen Control Size Ladder</b>	Multiple Genes	Blue	<b>Valid Size Range</b> IVS-0000 Polyclonal Control DNA	---	<b>96, 197, 297, 397, 602<sup>f</sup></b> 96, 197, 297, 397, 602 <sup>f</sup>

<sup>a</sup>Note: The 271 bp band is particularly seen in samples with low numbers of contaminating lymphoid cells.

<sup>b</sup>Note: Under sub-optimal conditions nonspecific products of 221 bp (in Tube B) and 128 and 337 bp (in Tube C) can be detected. If present, these bands will normally be faint.

<sup>c</sup>Note: Because smaller PCR fragments are preferentially amplified, it is not unusual for the 602 bp fragment to have a diminished signal or to be missing entirely.

<sup>d</sup>Note: Amplicon product is often not seen in this size range. This is an extremely restricted repertoire.

<sup>e</sup>Note: This may be seen as a weak amplicon.

<sup>f</sup>Note: Because smaller PCR fragments are preferentially amplified, it is not unusual for the 600 bp fragment to have a diminished signal or to be missing entirely.

## 10.2. Sample Data

The data shown in Figures 2 - 7 were generated using the master mixes indicated. Amplified products were run on an ABI instrument.




Figure 2. TCRB Tube A – 6FAM &amp; HEX




Figure 3. TCRB Tube B – 6FAM




Figure 4 TCRB Tube C – 6FAM &amp; HEX




Figure 5 TCRG Tube A – 6FAM &amp; HEX




Figure 6 TCRG Tube B – 6FAM &amp; HEX




Figure 7 Specimen Control Size Ladder – 6FAM master mix

## 11. Performance Characteristics

This IdentiClone *TCRB* + *TCRG* T-Cell Clonality PCR test is a rapid and reliable procedure that is far more sensitive than Southern Blot (SB) analysis in detecting clonality in suspect lymphoproliferations. The final clinico-histopathological diagnosis correlates well with PCR results in a higher number of patients in comparison with SB results. This is evidenced by two notable papers, one published in 2003 in *Leukemia* by van Dongen et al. and one published in 2005 in the *Journal of Molecular Diagnostics* (JMD) by Sandberg et al.

Table 9. Comparison of PCR and Southern Blot detection

PCR/SB concordance ( <i>Leukemia</i> ): <sup>2</sup>		PCR/SB concordance ( <i>JMD</i> ): <sup>3</sup>	
<i>IGH</i> :	93% sensitivity/ 92% specificity	<i>IGH</i> + <i>IGK</i> :	85% sensitivity
<i>IGK</i> :	90% sensitivity/ 90% specificity		
<i>IGL</i> :	86% sensitivity/ 92% specificity		
<i>TCRB</i> :	86% sensitivity/ 98% specificity	<i>TCRB</i> :	85% sensitivity
<i>TCRG</i> :	89% sensitivity/ 94% specificity		
<i>TCRD</i> :	83% sensitivity/ 95% specificity		



Table 10. PCR vs. SB analysis relative to histopathology and final diagnosis

	PCR/SB concordance:	PCR sensitivity:	SB sensitivity
IGH + IGK:	85%	98%	39%
TCRB:	85%	96%	35%

The study by Sandberg *et al.* was an independent study of 300 patient samples from a variety of sample types. In cases where both PCR and SB analyses were done and results could be correlated with histopathology and a final diagnosis, the diagnostic accuracy of selected IdentiClone tests was determined to be at least 96%. This was far more accurate than SB analysis, which in this study missed 23 clear cases of malignancy and 7 probable malignancies. There were no clear false positive results generated using the IdentiClone tests and there was a high level of precision.<sup>3</sup> In addition, a clear benefit of this assay was that clonal results generated allowed for subsequent detection of patient- and tumor-specific gene rearrangements for minimal residual disease detection.

## 12. Technical and Customer Service

Technical and Customer Service Representatives are available Monday through Friday to answer phone, e-mail or website inquiries.

### Contact Information



Invivoscribe, Inc

10222 Barnes Canyon Road | Building 1 | San Diego | California 92121-2711 | USA

Phone: +1 858 224-6600 | Fax: +1 858 224-6601 | Business Hours: 7:00AM - 5:00 PM PST/PDT

Technical Service: [support@invivoscribe.com](mailto:support@invivoscribe.com) | Customer Service: [sales@invivoscribe.com](mailto:sales@invivoscribe.com) | Website: [www.invivoscribe.com](http://www.invivoscribe.com)

## 13. References

1. Miller, JE, Wilson, SS, Jaye, DJ, Kronenberg, M. (1999). An automated semiquantitative B and T cell clonality assay. *Molecular Diagnostics* 4, 101-117.
2. Van Dongen, JJM *et al.* (2003). Design and standardization of PCR primers and protocols for detection of clonal immunoglobulin and T-cell receptor gene recombinations in suspect lymphoproliferations: Report of the BIOMED-2 Concerted Action BMH4-CT98-3936. *Leukemia*. 17, 2257-2317.
3. Sandberg, Y, *et al.* (2005). BIOMED-2 multiplex immunoglobulin/T-cell receptor polymerase chain reaction protocols can reliably replace Southern Blot analysis in routine clonality diagnostics. *J. Molecular Diagnostic* 7, 495-503.
4. van Krieken, JHJM *et al.* (2007). Improved reliability of lymphoma diagnostics via PCR-based clonality testing: – Report of the BIOMED-2 Concerted Action BHM4-CT98-3936. *Leukemia* 21, 201-206.

## 14. Symbols

The following symbols are now used in labeling for Invivoscribe diagnostic products.

	Catalog Number		Expiration Date
	Reagent Volume		Authorized Representative in the European Community
	Lot Number		Consult Instructions for Use
	Storage Conditions		For <i>In Vitro</i> Diagnostic Use
	Unique Device Identifier		Manufacturer
	UK Conformity Assessed		UK Responsible Person
	Swiss Authorized Representative		European Conformity

## 15. Legal Notice

### 15.1. Warranty and Liability

Invivoscribe, Inc. (Invivoscribe®) is committed to providing the highest quality products. Invivoscribe® warrants that the products meet or exceed the performance standards described in the Instructions For Use, as to products with such an insert. If a product is covered by product specifications and does not perform as specified, our policy is to replace the product or credit the full purchase price. No other warranties of any kind, expressed or implied, are provided by Invivoscribe®. Invivoscribe® liability shall not exceed the purchase price of the product. Invivoscribe shall have no liability for direct, indirect, consequential or incidental damages arising from the use, results of use, or inability to use its products; product efficacy under purchaser controlled conditions in purchaser's laboratory must be established and continually monitored through purchaser defined and controlled processes including but not limited to testing of positive, negative, and blank controls every time a sample is tested. Ordering, acceptance, and use of product constitutes purchaser acceptance of sole responsibility for assuring product efficacy and purchaser agreement to the limitation of liability set forth in this paragraph.

This product is an in vitro diagnostic product is not available for sale or use within North America.

### 15.2. Patents and Trademarks

This product is covered by one or more of the following: European Patent Number 1549764, European Patent Number 2418287, European Patent Number 2460889, Japanese Patent Number 4708029, United States Patent 8859748, and related pending and future applications. All of these patents and applications are licensed exclusively to Invivoscribe®. Additional patents licensed to Invivoscribe covering some of these products apply elsewhere. Many of these products require nucleic acid amplification methods such as Polymerase Chain Reaction (PCR). No license under these patents to use amplification processes or enzymes is conveyed expressly or by implication to the purchaser by the purchase of this product.

Identiclone® is a registered trademark of Invivoscribe®

©2023 Invivoscribe, Inc. All rights reserved. The trademarks mentioned herein are the property of Invivoscribe, Inc. and/or its affiliates or (as to the trademarks of others used herein) their respective owners.