

Instructions for Use

TCRD Gene Clonality Assay


RUO

For identification of clonal T-cell receptor delta chain gene rearrangements.

RUO For RESEARCH USE ONLY. Not for use in diagnostic procedures.

Manufactured in U.S.A.



 Storage Conditions: **-85°C to -65°C**
(DNA controls may be separated from assay kits and stored at 2°C to 8°C)

Catalog#	Products	Quantity
REF 12060010	<i>TCRD</i> Gene Clonality Assay for Gel Detection	33 Reactions
REF 12060011	<i>TCRD</i> Gene Clonality Assay for ABI Fluorescence Detection	33 Reactions
REF 12060020	<i>TCRD</i> Gene Clonality Assay MegaKit for Gel Detection	330 Reactions
REF 12060021	<i>TCRD</i> Gene Clonality Assay MegaKit for ABI Fluorescence Detection	330 Reactions

TABLE OF CONTENTS

1.	ASSAY USE	3
2.	SUMMARY AND EXPLANATION OF THE TEST.....	3
2.1.	Background	3
2.2.	Summary.....	3
3.	PRINCIPLES OF THE PROCEDURE.....	3
3.1.	Polymerase Chain Reaction (PCR).....	3
3.2.	Gel Detection.....	4
3.3.	Differential Fluorescence detection.....	4
4.	REAGENTS.....	5
4.1.	Reagent Components.....	5
4.1.	Warnings and Precautions	6
4.2.	Storage and Handling.....	6
5.	INSTRUMENTS.....	7
5.1.	Thermal cycler	7
5.2.	Electrophoresis Unit.....	7
5.3.	UV Illumination Unit.....	7
5.4.	ABI Capillary Electrophoresis Instruments.....	7
6.	SPECIMEN COLLECTION AND PREPARATION	8
6.1.	Precautions.....	8
6.2.	Interfering Substances	8
6.3.	Specimen Requirements and Handling	8
6.4.	Sample Preparation.....	8
6.5.	Sample Storage	8
7.	ASSAY PROCEDURE.....	8
7.1.	Materials Provided.....	8
7.2.	Materials Required (not provided)	9
7.3.	Reagent Preparation.....	10
7.4.	Amplification	11
7.5.	Detection.....	11
7.6.	Template Amplification Control	13
7.7.	Recommended Positive Controls.....	13
8.	INTERPRETATION AND REPORTING	13
8.1.	Sample Interpretation	13
8.2.	Control Interpretation.....	13
8.3.	Sample Reporting	13
9.	LIMITATIONS OF PROCEDURE.....	14
10.	EXPECTED VALUES.....	14
10.1.	Approximate product sizes with <i>TCRD</i> primers.....	14
10.2.	Expected Size of Amplified Products.....	14
10.3.	Sample Data.....	15
11.	TECHNICAL AND CUSTOMER SERVICE	15
12.	REFERENCES	16
13.	SYMBOLS	16
14.	LEGAL NOTICE	16
14.1.	Warranty and Liability	16
14.2.	Associated Patents	16
15.	<i>TCRD</i> GENE CLONALITY ASSAY: SINGLE PAGE GUIDE.....	17

1. Assay Use

This Research Use Only assay identifies T-cell receptor delta chain clonality and is useful for the study of identifying clonal T-cell populations and evaluation of new research and methods in malignancy studies.

2. Summary and Explanation of the Test

2.1. Background

Rearrangements of the antigen receptor genes occur during ontogeny in B and T lymphocytes. These gene rearrangements generate products that are unique in length and sequence for each cell. Therefore, polymerase chain reaction (PCR) assays can be used to identify lymphocyte populations derived from a single cell by detecting the unique V-J gene rearrangements present within these antigen receptor loci.¹ This PCR assay employs multiple consensus DNA primers that target conserved genetic regions within the T-cell receptor delta chain gene. The human T-cell receptor delta (*TCRD*) locus comprises a cluster of 10 genes located on chromosome 14 at 14q11.2. The cluster spans 60 kb and is localized between the T-cell receptor alpha (*TCRA*) variable and joining region genes.

This test is used to detect the vast majority of clonal T-cell malignancies from DNA. Test products can be analyzed using a variety of detection formats, including gel and capillary electrophoresis. Invivoscribe's gene clonality assays represent a simple approach to PCR-based clonality testing. These standardized assays were carefully optimized testing positive and negative control samples using multiplex master mixes.

2.2. Summary

This test kit includes two (2) master mixes. The *TCRD* Tube multiplex master mix contains 12 individual primers that target conserved regions within the variable (V), diversity (D), and the joining (J) regions that flank the unique hypervariable antigen-binding region 3 (CDR3). This standardized multiplex PCR assay detects the vast majority of clonal *TCR* delta gene rearrangements using a single multiplex master mix (Figure 1). The Specimen Control Size Ladder master mix, targets multiple genes and generates a series of amplicons of approximately 100, 200, 300, 400, and 600 base pairs (bp) to ensure that the quality and quantity of input DNA is adequate to yield a valid result. A single thermal cycler program and similar detection methodologies are used with all of our gene clonality assays which improves consistency and facilitates cross training on a broad range of different assays.

This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.



3. Principles of the Procedure


3.1. Polymerase Chain Reaction (PCR)

Polymerase chain reaction (PCR) assays are routinely used for the identification of clonal B-cell populations. These tests amplify the DNA between primers that target the variable (V) and joining (J) regions. The conserved regions lie on either side of an area within the V-J region where programmed genetic rearrangements occur during maturation of all B and T lymphocytes. The antigen receptor genes that undergo rearrangement are the immunoglobulin heavy chain and light chains genes in B-cells, and the T-cell receptor genes in T-cells. Each B- and T-cell has a single productive V-J rearrangement that is unique in both length and sequence. Therefore, when this region is amplified using DNA primers that flank this region, a clonal population of cells yields one or two prominent amplified products (amplicons) within the expected size range.

Two products are generated in cases when the initial rearrangement was non-productive and was followed by rearrangement of the other homologous chromosome. In contrast, DNA from a normal or polyclonal (many clones) population produces a bell-shaped curve of amplicon products (or Gaussian distribution) that reflect the heterogeneous population of V-J region rearrangements.

3.2. Gel Detection

Gel electrophoresis, such as agarose gel electrophoresis or non-denaturing polyacrylamide gel electrophoresis (PAGE), is commonly used to resolve the different amplicon products based on their size, charge, and conformation. Since DNA is negatively charged, when an electrical potential (voltage) is applied across the gel containing PCR products, the electrical field causes the amplicons to migrate through the gel. Smaller DNA fragments are able to easily migrate through the gel matrix, whereas larger DNA fragments migrate more slowly. This causes a separation of the amplicon products based on size. Ethidium bromide or other DNA intercalating dyes can then be used to stain and detect these products in the gel.



***TCRD* Tube: 6 $V\delta$ and 1 $D\delta 2$ primers + 4 $J\delta$ and 1 $D\delta 3$ primers**

Figure 1. This is a simplified diagram of a representative rearranged T-cell receptor delta gene on chromosome 14 (14q11.2) showing the approximate placement of the upstream and downstream DNA primers. The numbers of primers and their specificity are listed for the *TCRD* Tube master mix.

A heteroduplex analysis can also be performed and run on a polyacrylamide gel to differentiate clonal and non-clonal PCR products. A heteroduplex analysis involves denaturing the PCR products at a high temperature, then quickly re-annealing the DNA strands by suddenly reducing the temperature. This causes a large portion of DNA strands to incorrectly bind to other non-homologous strands creating loops in the DNA. These loops cause a significant reduction in the ability of the DNA to migrate through a polyacrylamide gel. However, if the majority of the PCR products are clonal, when a heteroduplex analysis is performed, most of these PCR products will correctly re-anneal with a homologous strand. These PCR products will run normally through the polyacrylamide gel. Therefore, in a clonal sample with a polyclonal background, a heteroduplex analysis will cause most of the polyclonal product to run much slower through the polyacrylamide gel, thereby increasing their separation and the ability to identify the clonal band(s).

3.3. Differential Fluorescence detection

Gel electrophoresis is commonly used to resolve the different-sized amplicon products and ethidium bromide or other DNA intercalating dyes to stain and detect these products. A powerful alternative method is use of differential fluorescence detection with primers conjugated with fluorescent dyes that correspond to different targeted regions. Reaction products from several different master mixes can be pooled, fractionated using capillary electrophoresis, and detected simultaneously. This detection system results in unsurpassed sensitivity, single base resolution, differential product detection and quantification. In addition, the laboratory can eliminate the use of agarose and polyacrylamide gels, as well as the use of carcinogens such as ethidium bromide. Further, differential detection allows accurate, reproducible and objective interpretation of primer-specific products and automatic archiving of data.

4. Reagents

4.1. Reagent Components

Table 1. Available Kits











Catalog #	Description	Quantity
 12060010	<i>TCRD</i> Gene Clonality Assay for Gel Detection	33 Reactions
 12060020	<i>TCRD</i> Gene Clonality Assay MegaKit for Gel Detection	330 Reactions
 12060011	<i>TCRD</i> Gene Clonality Assay for ABI Fluorescence Detection	33 Reactions
 12060021	<i>TCRD</i> Gene Clonality Assay MegaKit for ABI Fluorescence Detection	330 Reactions

Table 2. Reagent Components

Reagent	Catalog # 	Reagent Components (active ingredients)	Unit of Quantity	Assay Kit # of Units	MegaKit # of Units	Storage Temp.
Positive Control DNA	40881210	IVS-0021 Clonal Control DNA 200 µg/mL of DNA in 1/10 th TE	100 µL	1	5	
Negative (Normal) Control DNA	40920010	IVS-0000 Polyclonal Control DNA 200 µg/mL of DNA in 1/10 th TE	100 µL	1	5	 or 
<i>Gel Detection</i>						
Master Mixes	22060010	<i>TCRD</i> Tube – Unlabeled Multiple oligonucleotides targeting the Vδ + Dδ + Jδ regions of the T-cell receptor delta gene in a buffered salt solution.	1500 µL	1	10	
Template Amplification Master Mix	20960020	Specimen Control Size Ladder – Unlabeled Multiple oligonucleotides targeting housekeeping genes.	1500 µL	1	10	
<i>Differential Fluorescence Detection</i>						
Master Mixes	22060011	<i>TCRD</i> Tube – 6FAM & HEX Multiple oligonucleotides targeting the Vδ + Dδ + Jδ regions of the T-cell receptor delta gene in a buffered salt solution.	1500 µL	1	10	
Template Amplification Master Mix	20960021	Specimen Control Size Ladder – 6FAM Multiple oligonucleotides targeting housekeeping genes.	1500 µL	1	10	

Note: There are no preservatives used in the manufacture of this kit.

4.1. Warnings and Precautions

- **RUO** For RESEARCH USE ONLY. Not for use in diagnostic procedures
- Use this assay kit as a system. Do not substitute other manufacturer's reagents. Dilution, reducing amplification reaction volumes, or other deviation in this protocol may affect the performance of this test and/or nullify any limited sublicense that comes with the purchase of this testing kit.
- Materials are stable until the labeled expiration date when stored and handled as directed. Do not use kits beyond their expiration date.
- Close adherence to the protocol will assure optimal performance and reproducibility. Ensure use of correct thermal cycler program, as other programs may provide inaccurate/faulty data, such as false positive and false negative results.
- Do not mix or combine reagents from kits with different lot numbers.
- Wear appropriate personal protective equipment and follow good laboratory practices and universal precautions when working with specimens. Handle specimens in approved biological safety containment facilities and open only in certified biological safety cabinets.
- Due to the analytical sensitivity of this test, use extreme care to avoid the contamination of reagents or amplification mixtures with samples, controls or amplified materials. Closely monitor all reagents for signs of contamination (e.g., negative controls giving positive signals). Discard reagents suspected of contamination.
- To minimize contamination, wear clean gloves when handling samples and reagents and routinely clean work areas and pipettes prior to doing PCR.
- Autoclaving does not eliminate DNA contamination. Follow uni-directional workflow in the PCR laboratory; begin with master mix preparation, move to specimen preparation, then to amplification, and finally to detection. Do not bring amplified DNA into the areas designated for master mix or specimen preparation.
- Dedicate all pipettes, pipette tips, and any equipment used in a particular area to that area of the laboratory.
- Use sterile, disposable plastic ware whenever possible to avoid RNase, DNase, or cross-contamination.

4.2. Storage and Handling

- For any duration other than immediate use, **store assay kits at -85°C to -65°C.**
- The optimum storage temperature for DNA controls is 2°C to 8°C, but for long term storage DNA controls can be stored at -85°C to -65°C.
- All reagents and controls must be thawed and vortexed or mixed thoroughly prior to use to ensure that they are resuspended completely. Excessive vortexing may shear DNA and cause labeled primers to lose their fluorophors.
- Materials are stable until the labeled expiration date when stored and handled as directed. Do not use kits beyond their expiration date.
- Due to high salt concentrations, PCR master mixes are sensitive to freeze/thaw cycles. Aliquot master mixes into sterile o-ring screw-cap tubes if necessary.

5. Instruments

5.1. Thermal cycler

- Use or function: Amplification of DNA samples
- Performance characteristics and specification:
 - Minimum thermal range: 15°C to 96°C
 - Minimum ramping speed: 0.8°C/sec
- Follow manufacturer's installation, operation, calibration, and maintenance procedures.
- See section 7.4 *Amplification* for thermal cycler program.

5.2. Electrophoresis Unit

- Use or function: DNA fragment separation
- Performance characteristics and specification:
 - Capable of running at 35V to 135V for extended times
- Follow manufacturer's installation, operation, calibration, and maintenance procedures.

5.3. UV Illumination Unit

- Use or function: DNA detection
- Performance characteristics and specification:
 - Capable of emitting light at a wavelength of ~302 nm
- Follow manufacturer's installation, operation, calibration, and maintenance procedures.

5.4. ABI Capillary Electrophoresis Instruments

- Use or function: Fragment detection and analysis
- Performance characteristics and specification:
 - The following capillary electrophoresis instruments will meet the performance needs for this assay:
 - ABI 3100 Avant Genetic Analyzer (4-capillaries)
 - ABI 3100 Genetic Analyzer (16-capillaries)
 - ABI 3130 Genetic Analyzer (4-capillaries)
 - ABI 3130xL Genetic Analyzer (16-capillaries)
 - ABI 3500 Genetic Analyzer (8-capillaries)
 - ABI 3500xL Genetic Analyzer (24-capillaries)
- Follow manufacturer's installation, operation, calibration, and maintenance procedures.
- The ABI instrument used must be calibrated with appropriate Matrix Standards as outlined in section 7.2 *Materials Required But Not Provided*.
- Use the default settings for your polymer and capillary type.
- See section 7.5 *ABI Fluorescence Detection* for further details.

6. Specimen Collection and Preparation

6.1. Precautions

Biological specimens from humans may contain potentially infectious materials. Handle all specimens in accordance with the OSHA Standard on Bloodborne Pathogens or Biosafety Level 2.

6.2. Interfering Substances

The following substances are known to interfere with PCR:

- Divalent cation chelators
- Low retention pipette tips
- EDTA (not significant at low concentrations)
- Heparin

6.3. Specimen Requirements and Handling

This assay tests extracted and purified genomic DNA (gDNA). Common sources of gDNA include:

- 5 cc of peripheral blood, bone marrow biopsy or bone marrow aspirate anti-coagulated with heparin or EDTA. Ship at ambient temperature; OR
- Formalin-fixed paraffin embedded tissue or slides

6.4. Sample Preparation

Extract the gDNA from specimens as soon as possible. Resuspend DNA to a final concentration of 100 µg to 400 µg per ml in 1/10th TE (1 mM Tris-HCl, pH 8.0; 0.1 mM EDTA) or in molecular biology grade or USP water. This is a robust assay system. A wide range of DNA concentrations will generate a valid result. Therefore, quantifying and adjusting DNA concentrations is generally not necessary. Testing sample DNA with the Specimen Control Size Ladder master mix will ensure that DNA of sufficient quality and quantity was present to yield a valid result.

6.5. Sample Storage

Store gDNA at 2°C to 8°C or long term at -85°C to -65°C.

7. Assay Procedure

7.1. Materials Provided

Table 3: Kit components

<i>Gel Detection</i>	
Catalog #	Description
REF 22060010	<i>TCRD</i> Tube – Unlabeled
REF 20960020	Specimen Control Size Ladder – Unlabeled
REF 40881210	IVS-0021 Clonal Control DNA
REF 40920010	IVS-0000 Polyclonal Control DNA

<i>Differential Fluorescence Detection</i>	
Catalog #	Description
REF 22060011	<i>TCRD</i> Tube – 6FAM & HEX
REF 20960021	Specimen Control Size Ladder – 6FAM
REF 40881210	IVS-0021 Clonal Control DNA
REF 40920010	IVS-0000 Polyclonal Control DNA

7.2. Materials Required (not provided)

Table 4: Materials Required (not provided)

Reagent/Material	Recommended Reagents/Materials and Suppliers	Catalog # (REF)	Notes
DNA Polymerase	Roche: <ul style="list-style-type: none"> EagleTaq DNA Polymerase or equivalent 	05206944190	N/A
Glass Distilled De-ionized Molecular Biology Grade or USP Water	N/A	N/A	Sterile and free of DNase and RNase.
Calibrated Pipettes	N/A	N/A	Must be able to accurately measure volumes between 1 μ L and 1000 μ L.
Thermal cycler	Thermo Fisher Scientific: <ul style="list-style-type: none"> Veriti Dx Thermal Cycler Bio-Rad: <ul style="list-style-type: none"> MJ Research PTC-100 or PTC-200, PTC-220, PTC-240 Perkin-Elmer <ul style="list-style-type: none"> PE 9600 or PE 9700 	N/A	N/A
Vortex Mixer	N/A	N/A	N/A
PCR plates or tubes	N/A	N/A	Sterile
Filter barrier pipette tips	N/A	N/A	Sterile, RNase/DNase/Pyrogen-free
Microcentrifuge tubes	N/A	N/A	Sterile
Gel Electrophoresis Unit	N/A		For polyacrylamide gels
<i>Gel Detection</i>			
Ethidium Bromide	Thermo Fisher Scientific®: <ul style="list-style-type: none"> UltraPure® 10 mg/mL Ethidium Bromide 	15585-011	N/A
6% Polyacrylamide Gels	Thermo Fisher Scientific: <ul style="list-style-type: none"> Novex® TBE Gels (6%, 12 well) 	EC62652Box	N/A
TBE Running Buffer	Invitrogen: <ul style="list-style-type: none"> Novex TBE Running Buffer (5X) 	LC6675	Dilute 1:5 prior to use.
Gel Loading Buffer	Thermo Fisher Scientific: <ul style="list-style-type: none"> 10X BlueJuice™ Gel Loading Buffer Novex Hi-Density TBE Sample Buffer (5X) 	10816-015 LC6678	N/A
100 bp DNA Ladder	Invitrogen: <ul style="list-style-type: none"> TrackIt™ 100 bp DNA Ladder 	10488-058	N/A
<i>Differential Fluorescence Detection</i>			
ABI Capillary Electrophoresis Instrument	Applied Biosystems: <ul style="list-style-type: none"> ABI 310, 3100, or 3500 series 	S.O.	N/A
Hi-Di Formamide	Applied Biosystems: <ul style="list-style-type: none"> Hi-Di™ Formamide 	4311320	N/A

Table 4: Materials Required (not provided)

Reagent/Material	Recommended Reagents/Materials and Suppliers	Catalog # (REF)	Notes
Size Standards	Invivoscribe, Inc.:		
	• Hi-Di Formamide w/ROX size standards for ABI 310	60980051	
	• Hi-Di Formamide w/ROX size standards for ABI 3100	60980061	
	Applied Biosystems:		N/A
	• For ABI 3100 or 3130 instruments:		
	○ GeneScan™ - 400HD [ROX]™	402985	
	• For ABI 3500 instruments:		
	○ GeneScan - 600 [LIZ]™ v2.0	4408399	
Spectral Calibration Dye Sets	Applied Biosystems:		
	• For ABI 3100 and 3130 instruments:		
	○ DS-30 Matrix Standard Kit (Dye Set D)	4345827	
	• For ABI 310 instruments:		
	○ NED Matrix Standard	402996	N/A
○ And Fluorescent Amidite Matrix Standards [6FAM, TET, HEX, TAMRA, ROX]	401546		
• For ABI 3500 instruments:			
○ DS-33 Matrix Standard Kit (Dye Set G5)	4345833		
Polymer	Applied Biosystems:		
	• POP-4™ Polymer:		
	○ POP-4 for 310 Genetic Analyzers	402838	
	○ POP-4 for 3100/3100-Avant Genetic Analyzers	4316355	
	○ POP-4 for 3130/3130xL Genetic Analyzers	4352755	N/A
	• POP-7™ Polymer:		
○ POP-7 for 3130/3130xL Genetic Analyzers	4352759		
○ POP-7 for 3500/3500xL Genetic Analyzers	4393714		
Buffer	Applied Biosystems:		
	• 10X Genetic Analyzer Buffer with EDTA	402824	Dilute 1:10 in sterile water before use

7.3. Reagent Preparation

- Test all unknown samples using the template amplification control (**Specimen Control Size Ladder**) master mix to ensure that no inhibitors of amplification are present, and there is DNA of sufficient quality and quantity to generate a valid result.
 - Test all samples in **singlicate**.
 - **Test positive, negative and no template** controls for each master mix.
- 7.3.1. Using gloved hands, remove the master mixes from the freezer. Allow the tubes to thaw; then gently vortex to mix.
- 7.3.2. In containment hood or dead air box, remove an appropriate aliquot to clean, sterile microfuge tube (one tube for each master mix).
- Aliquot volumes are 45 µL for each sample + 135 µL (3 x 45 µL) for the positive, negative and no template controls.
 - Include an additional 20 µL to correct for pipetting errors.
- 7.3.3. Add the appropriate amount of Taq DNA polymerase (0.25 µL @5 U/µL per 50 µL total PCR volume) to the master mix and gently mix by inverting several times or gently vortexing.
- The master mixes are now ready for distribution to reaction tubes or plate, and addition of sample.

7.4. Amplification

- 7.4.1. Aliquot 45 μL of the master mix/enzyme solutions into individual PCR wells or tubes.
- 7.4.2. Add 5 μL of sample or control DNA to the individual tubes or wells containing the respective master mix reactions.
 - Pipette up and down several times to mix.
- 7.4.3. Amplify the reactions using the following PCR program
 - Use the **calculated** option for temperature measurement with the PTC instruments.

Table 5: Thermal cycling conditions

Step	Temperature	Duration	Cycles
1	95°C	7 minutes	1
2	95°C	45 seconds	35
3	60°C	45 seconds	
4	72°C	90 seconds	
5	72°C	10 minutes	1
6	15°C	∞	1

- 7.4.4. Remove the amplification plate from the thermal cycler.

7.5. Detection

Gel Detection – Agarose TBE Gels

- 7.5.1. Prepare a 2% MetaPhor or NuSieve 3:1 agarose/TBE gel.
- 7.5.2. Individually mix 20 μL from each amplification reaction mixed with 4 μL of 6X gel loading buffer.
- 7.5.3. Load 20 μL of this mixture into separate wells of the gel, flanked by DNA size standards.
- 7.5.4. Detect products using ethidium bromide or an equivalent dye.
- 7.5.5. Photograph the gel and interpret the data.

Gel Detection – Polyacrylamide TBE Gels

- 7.5.6. Assemble electrophoresis unit using a 6% non-denaturing polyacrylamide TBE gel and 0.5X TBE running buffer.
- 7.5.7. Add 5 μL ice-cold non-denaturing bromophenol blue loading buffer to samples.
- 7.5.8. Load 20 μL mixture into wells of the gel.
- 7.5.9. Run gel at 110V for 2-3 hours or 40-50V overnight.
 - Voltage and electrophoresis time depend on the PCR amplicon size, acrylamide gel thickness, and type of PCR equipment.
 - Voltage and run time can be adapted accordingly.
- 7.5.10. Stain gels in 0.5 $\mu\text{g}/\text{mL}$ ethidium bromide (in water or 0.5X TBE Buffer) for 5-10 minutes.
- 7.5.11. Destain gels 2X in water for 5-10 minutes.
- 7.5.12. Visualize the gel using UV illumination.
- 7.5.13. Photograph the gel and interpret the data.

Gel Detection – Heteroduplex Analysis (RECOMMENDED)

- 7.5.14. Denature 20 µL PCR products at 94°C for 5 minutes.
- 7.5.15. Re-anneal PCR products at 4°C for 60 minutes.
- 7.5.16. Assemble electrophoresis unit using a 6% non-denaturing polyacrylamide TBE gel and 0.5X TBE running buffer.
- 7.5.17. Add 5 µL ice-cold non-denaturing bromophenol blue loading buffer to samples.
- 7.5.18. Load 20 µL of mixture into wells of the gel.
- 7.5.19. Run gel at 110V for 2-3 hours or 40-50V overnight.
 - Voltage and electrophoresis time depend on the PCR amplicon size, acrylamide gel thickness, and type of PCR equipment.
 - Voltage and run time can be adapted accordingly.
- 7.5.20. Stain gels in 0.5 µg/mL ethidium bromide (in water or 0.5X TBE Buffer) for 5-10 minutes.
- 7.5.21. Destain gels 2X in water for 5-10 minutes.
- 7.5.22. Visualize the gel using UV illumination.
- 7.5.23. Photograph the gel and interpret the data.

ABI Fluorescence Detection with ABI 310, 3100, & 3130 instruments

- 7.5.24. In a new microcentrifuge tube, mix an appropriate amount (10 µL per reaction) of Hi-Di Formamide with ROX Size Standards. Vortex well.
- 7.5.25. In a new 96-well PCR plate, add 10 µL of Hi-Di Formamide with ROX size standards to individual wells for each reaction.
- 7.5.26. Transfer 1 µL of each reaction to the wells containing Hi-Di Formamide with ROX size standards.
 - Add only one sample per well.
 - Pipette up and down to mix.
- 7.5.27. Cap or cover the PCR plate.
- 7.5.28. Heat denature the samples at 95 °C for 2 minutes, then snap chill on ice for 5 minutes.
- 7.5.29. Prepare a **sample sheet** and **injection list** for the samples.
- 7.5.30. Run the samples on an ABI 310/3100/3130 capillary electrophoresis instrument according to its user manual.
 - Data are automatically displayed as size and color specific peaks.
- 7.5.31. Review profile and controls, report results.

ABI Fluorescence Detection with ABI 3500 instruments

- Note:** Due to instrument to instrument variation in the performance of the ABI 3500 platform, the amount of formamide, sample and size standard listed in the protocol is intended to be a starting point. The protocol may need to be optimized for specific ABI 3500 Platforms.
- 7.5.32. In a new microcentrifuge tube, mix an appropriate amount (9.5 µL per reaction) of Hi-Di Formamide with LIZ Size Standards. Vortex well.
 - 7.5.33. In a new 96-well PCR plate, add 9.5 µL Hi-Di Formamide with LIZ size standards to individual wells for each reaction.
 - 7.5.34. Transfer 0.5 µL of each reaction to the wells containing Hi-Di Formamide with LIZ size standards.
 - Add only one sample per well.
 - Pipette up and down to mix.
 - 7.5.35. Cap or cover the PCR plate.
 - 7.5.36. Heat denature the samples at 95°C for 3 minutes, then snap chill on ice for 5 minutes.
 - 7.5.37. Prepare a **sample sheet** and **injection list** for the samples.
 - 7.5.38. Run the samples on an ABI 3500 capillary electrophoresis instrument according to its user manual.
 - Data are automatically displayed as size and color specific peaks.
 - 7.5.39. Review profile and controls, report results.

7.6. Template Amplification Control

The **Specimen Control Size Ladder** master mix primers are available unlabeled for gel detection and labeled with a fluorescent dye (6-FAM) for differential fluorescent detection. This label is detected as **BLUE** using the differential fluorescence software. The amplicons produced with this master mix are at ~100, 200, 300, 400, and 600 bp.

- The ~100 bp band is comprised of 84 bp and 96 bp bands. Both of these bands co-migrate on a gel.
- Run the products of this master mix separately.

7.7. Recommended Positive Controls

Table 6: Recommended Positive Control Amplified Product Sizes

Master Mix	Target	Color	Control DNA	Catalog #	Product Size (bp)
TCRD Tube	V δ + J δ + D δ	Blue (Jδ) & Green (Dδ3)	Valid Size Range IVS-0021 Clonal Control DNA	--- 40881210	120 - 280¹ 158 (minor), 181
Specimen Control Size Ladder	Multiple Genes	Blue	Valid Size Range IVS-0000 Polyclonal Control DNA	--- 40920010	84, 96, 200, 300, 400, 600 84, 96, 200, 300, 400, 600

Note: The amplicon sizes listed above were determined using an ABI platform. Amplicon sizes observed on each specific CE instrument may differ 1-4 bp from those listed above depending on the platform of detection (ABI) and the version of the analysis software used. Once identified, the amplicon size as determined on each specific instrument will be consistent from run to run.

¹**Note:** A 90 bp peak may be observed in normal (polyclonal) samples.

8. Interpretation and Reporting

This assay is for research use only; not intended for diagnostic purposes. PCR based testing does not identify 100% of clonal cell populations; therefore, repeat testing by Next Generation Sequencing (NGS) may be advisable to rule out clonality.

8.1. Sample Interpretation

Results can be reported as “**Positive**” or “**Negative**” for “**Detection of clonal T-cell receptor delta chain gene rearrangement consistent with the presence of a clonal cell population**”.

- 8.1.1. Report samples that fail to amplify following repeat testing as “**A result cannot be reported on this specimen because there was DNA of insufficient quantity or quality for analysis**”.
- 8.1.2. All assay controls must be examined prior to interpretation of sample results. If the controls do not yield the correct results, the assay is not valid and the samples cannot be interpreted.

8.2. Control Interpretation

The following describes the analysis of each of the controls and the decisions necessary based upon the results.

- **Negative Control:** (Polyclonal control, water or no template blank). If the negative control is:
 - **Positive:** Possible contamination of all PCR amplification reactions. Do not continue with the interpretation of results. Prepare fresh master mix and repeat amplification.
 - **Negative:** Continue with the analysis.
- **Positive Control:** (This can also be an extraction control if positive control material is taken through extraction processes) If the positive control is:
 - **Positive:** Continue with analysis.
 - **Negative:** Repeat the assay.
- **Specimen Control Size Ladder:** (This is run on unknown samples only) If the specimen control size ladder:
 - **Generates 84, 96, 300, 400, 500, and 600 bp products:** Continue with analysis.
 - **Does not generate specified product:** Repeat assay and/or re-extract the specimen.

8.3. Sample Reporting

Following the acceptance of the controls, the samples are interpreted as follows: One or two prominent bands within the valid size range are reported as “**Detection of clonal T-cell receptor delta chain gene rearrangement consistent with the presence of a clonal cell population.**”

9. Limitations of Procedure

The assay is subject to interference by degradation of DNA or inhibition of PCR due to heparin or other agents.

10. Expected Values

The size range for each master mix has been determined testing positive control samples. For accurate and meaningful interpretation, ignore peaks that occur outside of the valid size range for each master mix. Peaks that are outside of the range cannot be assumed to be valid.

10.1. Approximate product sizes with *TCRD* primers

Table 7: Approximate product sizes with TCRD primers

	J _D 1	J _D 2	J _D 3	J _D 4	DD3-3'
V _D 1	170 bp - 210 bp	170 bp - 210 bp	180 bp - 220 bp	160 bp - 200 bp	~220 bp
V _D 2new3	200 bp - 240 bp	200 bp - 240 bp	200 bp - 240 bp	190 bp - 230 bp	~240 bp
V _D 3	230 bp - 270 bp	230 bp - 270 bp	240 bp - 280 bp	220 bp - 260 bp	~280 bp
V _D 4	180 bp - 220 bp	180 bp - 220 bp	190 bp - 230 bp	170 bp - 210 bp	~230 bp
V _D 5	220 bp - 260 bp	220 bp - 260 bp	230 bp - 270 bp	210 bp - 250 bp	~270 bp
V _D 6	200 bp - 240 bp	200 bp - 240 bp	210 bp - 250 bp	190 bp - 230 bp	~250 bp
DD2-5'	~130 bp	~140 bp	~150 bp	~130 bp	~190 bp

*These cell lines give these APPROXIMATE product sizes:

<u>Loucy:</u>	D2J1 150 bp,	V6J2 210 bp
<u>Nalm-16:</u>	D2D3 170 bp,	V2D3 230 bp
<u>REH:</u>	V2D3 240 bp	

*Courtesy of Dr. Louise Lavender, Molecular Pathology Unit, Southampton University Hospitals

10.2. Expected Size of Amplified Products

Note: “Color” indicates the color of products generated with the master mix when using differential fluorescence detection format (e.g., ABI instruments).

Table 8: Expected Size of Amplified Products

Master Mix	Target	Color	Control DNA	Catalog #	Product Size (bp)
<i>TCRD</i> Tube	V δ + J δ + D δ	Blue (J δ) & Green (D δ 3)	Valid Size Range IVS-0000 Polyclonal Control DNA IVS-0021 Clonal Control DNA	--- 40920010 40881210	120 - 280¹ 130 - 280, 160 - 280² 158 (minor), 181
Specimen Control Size Ladder	Multiple Genes	Blue	Valid Size Range Any human DNA	--- 40920010	84, 96, 200, 300, 400, 600 84, 96, 200, 300, 400, 600

Note: The amplicon sizes listed above were determined using an ABI platform. Amplicon sizes observed on each specific CE instrument may differ 1-4 bp from those listed above depending on the platform of detection (ABI) and the version of the analysis software used. Once identified, the amplicon size as determined on each specific instrument will be consistent from run to run.

¹**Note:** A 90 bp peak may be observed in normal (polyclonal) samples.

²**Note:** A nonspecific 163 bp peak may be observed in normal (polyclonal) samples.

10.3. Sample Data

Figure 2. Gel Detection

The image shown at right were generated using the *TCRD* Tube master mix. Amplified products were run on a 6% polyacrylamide gel.

TCRD Tube

Lane 1 = 100% IVS-0014
Lane 2 = 100% IVS-0021
Lane 3 = 10% IVS-0021
Lane 4 = 100% IVS-0000

Valid Size Range =
120-280 bp





Figure 3. The data shown at left were generated using the *TCRD* Tube master mix and run on an ABI instrument.

11. Technical and Customer Service

Thank you for purchasing our ***TCRD* Gene Clonality Assay**. We appreciate your business. We are happy to assist you in the validation of this assay and will provide ongoing technical assistance to keep the assays performing efficiently in your laboratory. Technical assistance is most rapidly obtained using our Internet site: <http://www.inivoscribe.com> or by sending an email inquiry to: support@inivoscribe.com.

Contact Information

Inivoscribe, Inc
10222 Barnes Canyon Road, Bldg. 1
San Diego, California 92121-2711
USA

Phone: +1 858 224-6600
Fax: +1 858 224-6601
Technical Service: support@inivoscribe.com
Customer Service: sales@inivoscribe.com
Website: www.inivoscribe.com
Business Hours: 7:00AM – 5:00PM PST/PDT

Authorized Representative and EU Technical Assistance

EC REP Inivoscribe Technologies, SARL
Le Forum – Bât B
515 Avenue de la Tramontane
ZI Athélia IV
13600 La Ciotat, France










Phone: +33 (0)4 42 01 78 10
Fax: +33 (0)4 88 56 22 89
Technical Service: support@inivoscribe.com
Customer Service: sales-eu@inivoscribe.com
Website: www.inivoscribe.com
Business Hours: 9:00AM – 5:00PM CET/CEST

12. References

1. Miller, JE, Wilson, SS, Jaye, DJ, Kronenberg, M. An automated semiquantitative B and T cell clonality assay. *Molecular Diagnostics*. 1999, **4(2)**:101-117.
2. van Dongen, JJM *et al.* Design and standardization of PCR primers and protocols for detection of clonal immunoglobulin and T-cell receptor gene recombinations in suspect lymphoproliferations: Report of the BIOMED-2 Concerted Action BMH4-CT98-3936. *Leukemia*. 2003, **17(12)**:2257-2317.
3. van Krieken JH, Langerak AW, Macintyre EA, Kneba M, Hodges E, Sanz RG, Morgan GJ, Parreira A, Molina TJ, Cabeçadas J, Gaulard P, Jasani B, Garcia JF, Ott M, Hannsmann ML, Berger F, Hummel M, Davi F, Brüggemann M, Lavender FL, Schuurin E, Evans PA, White H, Salles G, Groenen PJ, Gameiro P, Pott Ch, van Dongen JJM. Improved reliability of lymphoma diagnostics via PCR-based clonality testing: report of the BIOMED-2 Concerted Action BHM4-CT98-3936. *Leukemia*. 2007; **21(2)**:201-6.

13. Symbols

The following symbols are used in labeling for Invivoscribe products.

	Catalog Number		Expiration Date
	Reagent Volume		Authorized Representative in the European Community
	Lot Number		Manufacturer
	Storage Conditions		Consult Instructions for Use
	Research Use Only		

14. Legal Notice

14.1. Warranty and Liability

Invivoscribe, Inc. (Invivoscribe®) is committed to providing the highest quality products. Invivoscribe® warrants that the products meet or exceed the performance standards described in the Instructions For Use, as to products with such an insert. If a product is covered by product specifications and does not perform as specified, our policy is to replace the product or credit the full purchase price. No other warranties of any kind, expressed or implied, are provided by Invivoscribe®. Invivoscribe® liability shall not exceed the purchase price of the product. Invivoscribe® shall have no liability for direct, indirect, consequential or incidental damages arising from the use, results of use, or inability to use its products; product efficacy under purchaser controlled conditions in purchaser's laboratory must be established and continually monitored through purchaser defined and controlled processes including but not limited to testing of positive, negative, and blank controls every time a sample is tested. Ordering, acceptance, and use of product constitutes purchaser acceptance of sole responsibility for assuring product efficacy and purchaser agreement to the limitation of liability set forth in this paragraph.

This product is for Research Use Only; not for use in diagnostic procedures.

14.2. Associated Patents

This product is covered by one or more of the following: European Patent Number 1549764, European Patent Number 2418287, European Patent Number 2460889, Japanese Patent Number 4708029, United States Patent 8859748, and related pending and future applications. All of these patents and applications are licensed exclusively to Invivoscribe®. Additional patents licensed to Invivoscribe® covered in some of these products apply elsewhere. Many of these products require nucleic acid amplification methods such as Polymerase Chain Reaction (PCR). No license under these patents to use amplification processes or enzymes is conveyed expressly or by implication to the purchaser by the purchase of this product.

©2020 Invivoscribe, Inc. All rights reserved. The trademarks mentioned herein are the property of Invivoscribe, Inc. and/or its affiliates, or (as to the trademarks of others used herein) their respective owners.

15. *TCRD* Gene Clonality Assay: Single Page Guide

- 15.1. Using gloved hands, remove the master mixes from the freezer. Allow the tubes to thaw; then gently vortex to mix.
- 15.2. In a containment hood or dead air box remove an appropriate aliquot to clean, sterile microfuge tube (one tube for each master mix).
 - Aliquot volumes are 45 μ L for each sample + 135 μ L for the positive, negative and no template controls.
 - Include an additional 20 μ L to correct for pipetting errors.
- 15.3. Add the appropriate amount of Taq DNA polymerase (0.25 μ L @5 U/ μ L per 50 μ L total reaction volume) to each master mix and gently mix by inverting several times or gently vortexing.
- 15.4. Aliquot 45 μ L of master mix to individual wells of a PCR plate.
- 15.5. Add 5 μ L of DNA from the unknown and control samples to individual tubes or wells containing the respective master mix reactions and pipette up and down several times to mix.
- 15.6. Amplify target DNA using the universal thermal cycler program.

Gel Detection – Heteroduplex Analysis

- 15.7. Denature 20 μ L PCR products at 94°C for 5 minutes.
- 15.8. Re-anneal PCR products at 4°C for 60 minutes.
- 15.9. Assemble electrophoresis unit using a 6% non-denaturing polyacrylamide TBE gel and 0.5X TBE running buffer.
- 15.10. Add 5 μ L of ice-cold non-denaturing bromophenol blue loading buffer to samples
- 15.11. Load 20 μ L of mixture into wells of the gel.
- 15.12. Run gel at 110V for 2-3 hours or 40-50V overnight.
- 15.13. Gels are stained in 0.5 μ g/mL ethidium bromide (in water or 0.5X TBE Buffer) for 5-10 minutes.
- 15.14. Gels are destained 2X in water for 5-10 minutes.
- 15.15. UV illumination is used for visualization.
- 15.16. Gel is photographed and data are interpreted.

ABI Fluorescence Detection with ABI 310, 3100 & 3130 instruments

- 15.17. In a new microcentrifuge tube, mix an appropriate amount (10 μ L per reaction) of Hi-Di Formamide with ROX Size Standards. Vortex well.
- 15.18. In a new 96-well PCR plate, add 10 μ L of Hi-Di Formamide with ROX size standards to individual wells for each reaction.
- 15.19. Transfer 1 μ L of each reaction to the wells containing Hi-Di Formamide with ROX size standards. Add only one sample per well. Pipette up and down to mix.
- 15.20. Cap or cover the PCR plate.
- 15.21. Heat denature the samples at 95 °C for 2 minutes, then snap chill on ice for 5 minutes.
- 15.22. Prepare a **sample sheet** and **injection list** for the samples.
- 15.23. Run the samples on an ABI 310/3100/3130 capillary electrophoresis instrument according to its user manual.
- 15.24. Review profile and controls, report results.

ABI Fluorescence Detection with ABI 3500 instruments

- 15.25. In a new microcentrifuge tube, mix an appropriate amount (9.5 μ L per reaction) of Hi-Di Formamide with LIZ Size Standards. Vortex well.
- 15.26. In a new 96-well PCR plate, add 9.5 μ L of Hi-Di Formamide with LIZ size standards to individual wells for each reaction.
- 15.27. Transfer 0.5 μ L of each reaction to the wells containing Hi-Di Formamide with LIZ size standards. Add only one sample per well. Pipette up and down to mix.
- 15.28. Cap or cover the PCR plate.
- 15.29. Heat denature the samples at 95°C for 3 minutes, then snap chill on ice for 5 minutes.
- 15.30. Prepare a sample sheet and injection list for the samples.
- 15.31. Run the samples on an ABI 3500 capillary electrophoresis instrument according to the user manual.
- 15.32. Review profile and controls, report results.