

Instructions for Use

General Purpose Reagents: Controls

For use as qualitative PCR controls.

| Catalog Number | Product Name | Quantity |
|---------------------|-------------------------------------|------------------|
| REF 40880008 | <i>IGH</i> SHM Pos (+) | 1 Tube - 45 µL |
| REF 40880010 | IVS-0001 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40880190 | IVS-0004 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40880370 | IVS-0007 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40880430 | IVS-0008 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40880490 | IVS-0009 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40880550 | IVS-0010 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40880730 | IVS-0013 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40881090 | IVS-0019 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40881210 | IVS-0021 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40881390 | IVS-0024 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40881690 | IVS-0029 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40881750 | IVS-0030 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40881810 | IVS-0031 Clonal Control DNA | 1 Tube - 100 µL |
| REF 40883320 | 5% <i>TCRG</i> Positive Control DNA | 1 Tube - 50 µL |
| REF 40920020 | <i>TCRG</i> Negative Control DNA | 1 Tube - 50 µL |
| REF 40920010 | IVS-0000 Polyclonal Control DNA | 1 Tube - 100 µL |
| REF R0880230 | <i>FLT3</i> ITD Positive DNA | 1 Tube - 100 µL |
| REF R0880240 | <i>FLT3</i> TKD Positive DNA | 1 Tube - 100 µL |
| REF R0880250 | <i>FLT3</i> Extraction Control | 1 Tube - 1800 µL |
| REF 40890100 | IVS-0002 Clonal Control RNA | 1 Tube - 100 µL |
| REF 40890190 | IVS-0003 Clonal Control RNA | 1 Tube - 100 µL |
| REF 40890910 | IVS-0011 Clonal Control RNA | 1 Tube - 100 µL |
| REF 40891270 | IVS-0015 Clonal Control RNA | 1 Tube - 100 µL |
| REF 40891720 | IVS-0020 Clonal Control RNA | 1 Tube - 100 µL |
| REF 40892800 | IVS-0032 Clonal Control RNA | 1 Tube - 100 µL |
| REF 40893070 | IVS-0035 Clonal Control RNA | 1 Tube - 100 µL |

Table of Contents

| | |
|--|----------|
| 1. PROPRIETARY NAME | 3 |
| 2. INTENDED USE | 3 |
| 3. REAGENTS | 4 |
| 3.1. Tissue DNA Controls..... | 4 |
| 3.2. Cell Line DNA Controls | 4 |
| 3.3. RNA Controls | 5 |
| 3.4. Warnings and Precautions..... | 6 |
| 3.5. Storage and Stability..... | 6 |
| 4. PROCEDURE | 6 |
| 5. INTERPRETATION OF RESULTS | 6 |
| 6. SYMBOLS | 6 |
| 7. TECHNICAL AND CUSTOMER SERVICE | 7 |
| 8. LEGAL NOTICE | 7 |

1. Proprietary Name

Table 1. Proprietary Product Names and Unique Device Identification (UDI)

| Catalog Number | Product Name | UDI | Volume |
|---------------------|-------------------------------------|--------------|---------|
| REF 40880008 | <i>IGH</i> SHM Pos (+) | 810022731642 | 45 µL |
| REF 40880010 | IVS-0001 Clonal Control DNA | 810022731666 | 100 µL |
| REF 40880190 | IVS-0004 Clonal Control DNA | 810022731710 | 100 µL |
| REF 40880370 | IVS-0007 Clonal Control DNA | 810022731727 | 100 µL |
| REF 40880430 | IVS-0008 Clonal Control DNA | 810022731734 | 100 µL |
| REF 40880490 | IVS-0009 Clonal Control DNA | 810022731741 | 100 µL |
| REF 40880550 | IVS-0010 Clonal Control DNA | 810022731758 | 100 µL |
| REF 40880730 | IVS-0013 Clonal Control DNA | 810022731765 | 100 µL |
| REF 40881090 | IVS-0019 Clonal Control DNA | 810022731789 | 100 µL |
| REF 40881210 | IVS-0021 Clonal Control DNA | 810022731796 | 100 µL |
| REF 40881390 | IVS-0024 Clonal Control DNA | 810022731802 | 100 µL |
| REF 40881690 | IVS-0029 Clonal Control DNA | 810022731819 | 100 µL |
| REF 40881750 | IVS-0030 Clonal Control DNA | 810022731826 | 100 µL |
| REF 40881810 | IVS-0031 Clonal Control DNA | 810022731833 | 100 µL |
| REF 40883320 | 5% <i>TCRG</i> Positive Control DNA | 810022731857 | 50 µL |
| REF 40920020 | <i>TCRG</i> Negative Control DNA | 810022732076 | 50 µL |
| REF 40920010 | IVS-0000 Polyclonal Control DNA | 810022732052 | 100 µL |
| REF R0880230 | <i>FLT3</i> ITD Positive DNA | 810022732373 | 100 µL |
| REF R0880240 | <i>FLT3</i> TKD Positive DNA | 810022732380 | 100 µL |
| REF R0880250 | <i>FLT3</i> Extraction Control | 810022732472 | 1800 µL |
| REF 40890100 | IVS-0002 Clonal Control RNA | 810022731949 | 100 µL |
| REF 40890190 | IVS-0003 Clonal Control RNA | 810022731956 | 100 µL |
| REF 40890910 | IVS-0011 Clonal Control RNA | 810022731963 | 100 µL |
| REF 40891270 | IVS-0015 Clonal Control RNA | 810022731970 | 100 µL |
| REF 40891720 | IVS-0020 Clonal Control RNA | 810022731987 | 100 µL |
| REF 40892800 | IVS-0032 Clonal Control RNA | 810022731994 | 100 µL |
| REF 40893070 | IVS-0035 Clonal Control RNA | 810022731642 | 100 µL |

2. Intended Use

Inivoscribe's General Purpose Reagents (GPR) are intended for validation by the end user for use in molecular assays. These controls are specific to gene rearrangements, mutations, and/or translocations from human DNA or RNA indicated in section 3: *Reagents*.

These products are General Purpose Reagents and are for Laboratory Use.


3. Reagents

3.1. Tissue DNA Controls

Tissue DNA controls are extracted from normal, disease-free tissue and are tested extensively to ensure quality and reproducibility of your test results.

- Tissue DNA Controls are adjusted to the final concentration specified in Table 2 with 1/10 TE (1 mM Tris-HCl (pH 8.0), 0.1 mM EDTA), which provides sufficient buffering capacity and EDTA to protect the DNA controls without interfering with the Mg^{2+} concentrations required for robust amplification reactions.

Table 2. Tissue DNA Controls

| Catalog Number | Description | Concentration | Characteristics | Storage Conditions |
|---------------------|---------------------------------|---------------|---|---|
| REF 40920010 | IVS-0000 Polyclonal Control DNA | 200 µg/mL | <ul style="list-style-type: none"> genomic DNA isolated from the tissue of normal human tonsils negative control for gene rearrangements, chromosome translocations, and mutation tests |  8°C * |

*Minimize the number of freeze-thaw cycles.

3.2. Cell Line DNA Controls

Cell Line DNA Controls are extracted from established cell lines grown under standard cell culture conditions and extensively tested to ensure quality and reproducibility.

- DNA Clonal Controls are adjusted to the final concentration specified in Table 3 with 1/10th TE (1 mM Tris-HCl (pH 8.0), 0.1 mM EDTA), which provides sufficient buffering capacity and EDTA to protect the DNA controls without interfering with the Mg^{2+} concentrations required for robust amplification reactions.

Table 3. Cell Line DNA Controls











| Catalog Number | Description | Concentration | Characteristics | | | Storage Conditions |
|---------------------|--|---------------|---|---------------------------------------|---------------|---|
| | | | Gene Rearrangement | Chromosome Translocation | Gene Mutation | |
| REF 40880008 | <i>IGH</i> SHM Pos (+) [‡] | 15 µg/mL | <i>IGH, IGK</i> | <i>IGH-BCL2</i> t(14;18) Mbr | n/c |  8°C * |
| REF 40880010 | IVS-0001 Clonal Control DNA | 200 µg/mL | n/c | <i>IGH-BCL2</i> t(14;18) mcr | n/c | |
| REF 40880190 | IVS-0004 Clonal Control DNA | 200 µg/mL | <i>TRB, TRG</i> | n/c | n/c | |
| REF 40880370 | IVS-0007 Clonal Control DNA | 200 µg/mL | <i>IGH, IGK, IGL</i> | <i>IGH-BCL2</i> t(14;18) Mbr | n/c | |
| REF 40880430 | IVS-0008 Clonal Control DNA [†] | 200 µg/mL | <i>IGH</i> D _H -J _H , <i>TRB, TRG</i> | n/c | n/c | |
| REF 40880490 | IVS-0009 Clonal Control DNA | 200 µg/mL | <i>TRB, TRG</i> | n/c | n/c | |
| REF 40880550 | IVS-0010 Clonal Control DNA | 200 µg/mL | <i>IGH, IGK, IGL</i> | <i>IGH-BCL1</i> t(11;14) [§] | n/c | |
| REF 40880730 | IVS-0013 Clonal Control DNA | 200 µg/mL | <i>IGH, IGK, IGL</i> | n/c | n/c | |
| REF 40881090 | IVS-0019 Clonal Control DNA | 200 µg/mL | <i>IGH, IGK</i> | n/c | n/c | |
| REF 40881210 | IVS-0021 Clonal Control DNA | 200 µg/mL | <i>TRG, TRB, TRD</i> | n/c | n/c | |
| REF 40881390 | IVS-0024 Clonal Control DNA | 200 µg/mL | <i>IGH, IGK</i> | n/c | n/c | |
| REF 40881690 | IVS-0029 Clonal Control DNA | 200 µg/mL | <i>IGH, IGK, IGL</i> | n/c | n/c | |

Table 3. Cell Line DNA Controls

| Catalog Number | Description | Concentration | Characteristics | | | Storage Conditions |
|--|--|---------------|---|------------------------------|----------------------|---|
| | | | Gene Rearrangement | Chromosome Translocation | Gene Mutation | |
|  40881750 | IVS-0030 Clonal Control DNA [†] | 200 µg/mL | <i>IGH, IGK</i> | <i>IGH-BCL2</i> t(14;18) Mbr | n/c |  |
|  40881810 | IVS-0031 Clonal Control DNA | 200 µg/mL | <i>IGH, IGK</i> | <i>IGH-BCL2</i> t(14;18) mcr | n/c | |
|  40883320 | 5% <i>TCRG</i> Positive Control DNA | 50 µg/mL | <i>TRG</i> | n/c | n/c | |
|  40920020 | <i>TCRG</i> Negative Control DNA | 50 µg/mL | negative control for <i>TRG</i> gene rearrangements | n/c | n/c | |
|  R0880230 | <i>FLT3</i> ITD Positive DNA | 10 µg/mL | n/c | n/c | <i>FLT3</i> ITD |  |
|  R0880240 | <i>FLT3</i> TKD Positive DNA | 10 µg/mL | n/c | n/c | <i>FLT3</i> TKD | |
|  R0880250 | <i>FLT3</i> Extraction Control | 50 µg/mL | n/c | n/c | <i>FLT3</i> negative | |

Note: n/c is used to indicate that the control has not been fully characterized; there may be additional rearrangements, translocations or mutations associated with the control.

*Minimize the number of freeze-thaw cycles.

[†]This control does not contain a complete *IGH* V_H-J_H rearrangement and may only be suitable for *IGH* D_H-J_H rearrangements.

[§]*IGH-BCL1* was previously referred to as *BCL1*/J_H









[‡]These controls can be used as SHM positive controls with ≥2% mutational rates compared to the germline sequence.

3.3. RNA Controls

Inivoscribe RNA Controls are extracted from established cell lines grown under standard cell culture conditions and extensively tested to ensure quality and reproducibility. RNA

- Controls are adjusted to the final concentration specified in Table 4. with molecular grade water.

Table 4. RNA Controls

| Catalog Number | Description | Concentration | Target | Storage Conditions |
|--|-----------------------------|---------------|---|---|
|  40890100 | IVS-0002 Clonal Control RNA | 400 µg/mL | <i>E2A-PBX1</i> t(1;19)(q23;p13) |  |
|  40890190 | IVS-0003 Clonal Control RNA | 400 µg/mL | <i>BCR-ABL1</i> t(9;22) p210 e13a2 (b2a2) | |
|  40890910 | IVS-0011 Clonal Control RNA | 400 µg/mL | <i>BCR-ABL1</i> t(9;22) p210 e14a2 (b3a2) | |
|  40891270 | IVS-0015 Clonal Control RNA | 400 µg/mL | <i>CBFB-MYH11</i> inv(16) | |
|  40891720 | IVS-0020 Clonal Control RNA | 400 µg/mL | <i>PML-RARα</i> t(15;17)(q22;q11) | |
|  40892800 | IVS-0032 Clonal Control RNA | 400 µg/mL | <i>BCR-ABL1</i> t(9;22) p190 e1a2 | |
|  40893070 | IVS-0035 Clonal Control RNA | 400 µg/mL | negative for <i>BCR-ABL</i> and <i>PML-RARα</i> chromosome translocations | |

*Minimize the number of freeze-thaw cycles.

3.4. Warnings and Precautions

- Invivoscribe's GPRs are for Laboratory Use.
- Invivoscribe's GPR Controls are not intended for use as a substitute for the internal controls provided by *in vitro* diagnostic kit manufacturers.
- Establish standard operating procedures and instructions for using the Invivoscribe GPR DNA and RNA controls in molecular assays.
- Reagents are stable until the labeled expiration date when stored and handled as directed. Do not use reagents beyond their expiration date.
- Perform all quality control requirements in conformance with local, state and/or federal regulations or accreditation requirements.
- Wear appropriate personal protective equipment and follow good laboratory practices and universal precautions when working with specimens.
- Handle specimens in approved biological safety containment facilities and open only in certified biological safety cabinets.
- Use extreme care to avoid the contamination of reagents with samples, controls or amplified materials. Closely monitor all reagents for signs of contamination (*e.g.*, negative controls giving positive signals). Discard reagents suspected of contamination.
- To minimize contamination, wear clean gloves when handling samples and reagents and routinely clean work areas and pipettes prior to doing PCR.
- Autoclaving does not eliminate DNA contamination.
- Follow uni-directional workflow in the PCR laboratory; begin with master mix preparation, move to specimen preparation, then to amplification, and finally to detection. Do not bring amplified DNA into the areas designated for master mix or specimen preparation.
- Dedicate all pipettes, pipette tips, and any equipment used in a particular area to that area of the laboratory.
- Use sterile, disposable plastic ware whenever possible to avoid RNase, DNase, or cross-contamination.

3.5. Storage and Stability

- When stored at the intended storage conditions, Invivoscribe's GPR controls are stable until the expiration date indicated on the vial label.

4. Procedure

- 4.1. Allow the GPR control to equilibrate to room temperature.
- 4.2. Vortex gently, then pulse-spin in a centrifuge (4 to 6 seconds) to collect the contents at the bottom of the tube.
- 4.3. Introduce the GPR control as an independent sample at the template addition step in the workflow.
- 4.4. Handle the GPR control similarly to nucleic acids extracted from routine samples and run in parallel with routine samples.

5. Interpretation of Results

- 5.1. Results generated by Invivoscribe's GPR controls may differ according to the molecular test method.
- 5.2. To establish a baseline performance, incorporate results from multiple runs under different conditions (*e.g.*, operator, run, day) to determine a valid size range specific to the assay used.
- 5.3. Once the validated size range is ascertained, the expected size range can be used to verify each subsequent run result of the GPR control.

6. Symbols

The following symbols are used in labeling for Invivoscribe products.

| | | | |
|---|--------------------|---|---|
|  | Storage Conditions |  | Expiration Date |
|  | Catalog Number |  | Authorized Representative in the European Community |
|  | Reagent Volume |  | Manufacturer |
|  | Lot Number |  | Consult Instructions for Use |

7. Technical and Customer Service

Technical and Customer Service Representatives are available Monday through Friday to answer phone, e-mail or website inquiries. Please do not hesitate to contact sales@invivoscribe.com for assistance evaluating controls to suit your testing needs.

Contact Information

Invivoscribe, Inc.
10222 Barnes Canyon Road, Bldg. 1
San Diego, CA 92121-2777
USA

Phone: +1 858 224-6600
Fax: +1 858 224-6601
Technical Service: support@invivoscribe.com
Customer Service: sales@invivoscribe.com
Website: www.invivoscribe.com
Business Hours: 7:00AM – 5:00PM PST/PDT

Authorized Representative and EU Technical Assistance

EC|REP Invivoscribe Technologies, SARL
Le Forum – Bât B
515 Avenue de la Tramontane
ZI Athélia IV
13600 La Ciotat, France

Phone: +33 (0)4 42 01 78 10
Fax: +33 (0)4 88 56 22 89
Technical Service: support@invivoscribe.com
Customer Service: sales-eu@invivoscribe.com
Website: www.invivoscribe.com
Business Hours: 9:00AM – 5:00PM CET/CEST

8. Legal Notice

Invivoscribe, Inc. (Invivoscribe®) is committed to providing the highest quality products. Invivoscribe® warrants that the products meet or exceed the performance standards described in the Instructions For Use, as to products with such an insert. If a product is covered by product specifications and does not perform as specified, our policy is to replace the product or credit the full purchase price. No other warranties of any kind, expressed or implied, are provided by Invivoscribe®. Invivoscribe® liability shall not exceed the purchase price of the product. Invivoscribe shall have no liability for direct, indirect, consequential or incidental damages arising from the use, results of use, or inability to use its products; product efficacy under purchaser controlled conditions in purchaser's laboratory must be established and continually monitored through purchaser defined and controlled processes including but not limited to testing of internally validated positive, negative, and blank controls every time a sample is tested. Ordering, acceptance and use of product constitutes purchaser acceptance of sole responsibility for assuring product efficacy and purchaser agreement to the limitation of liability set forth in this paragraph.

This product is a General Purpose Reagent (GPR).

©2020 Invivoscribe, Inc. All rights reserved. The trademarks mentioned herein are the property of Invivoscribe, Inc. and/or its affiliates or (as to the trademarks of others used herein) their respective owners.