

# Validate an NGS Based Measurable Residual Disease (MRD) Assay for B-Cell Lymphoproliferative Diseases


P. Shah<sup>1</sup>, L. Hebert<sup>1</sup>, O. Diaz II<sup>1</sup>, A.A. Street<sup>1</sup>, J. Li<sup>1</sup>, Ekaterina Rudenko<sup>1</sup>, A.M. Zlotnicki<sup>1</sup>, C. Morrison<sup>2</sup>, S. Glenn<sup>2</sup>, R. Bomben<sup>3</sup>, A. Medina Herrera<sup>4</sup>, A. Wong<sup>5</sup>, and Y. Huang<sup>1</sup> | <sup>1</sup>Invivoscribe, Inc., San Diego, United States; <sup>2</sup>Roswell Park Comprehensive Cancer Center, New-York, United States; <sup>3</sup>Reference Oncologic Center, National Cancer Institute, Aviano, Italy; <sup>4</sup>Salamanca University Hospital, Salamanca, Spain; <sup>5</sup>Lahey Hospital & Medical Center, Burlington, United States.

## Background

MRD assessment in B-cell lymphoproliferative diseases has proven utility in assessing response to therapy, refining treatment and predicting clinical outcome in patients previously treated with a combination of treatments (Garcia-Marco et al. Haematologica 2019; Pui et al. Leukemia 2017). Next-generation sequencing (NGS) based approaches for MRD detection have provided high sensitivity and specificity to identify and track measurable residual disease using clonal immunoglobulin heavy chain (*IGH*) gene rearrangements (Rawstron, Leukemia, 2016). Here, we report a sensitive NGS-based B-cell MRD Assay to detect MRD in B-cell malignancies such as acute lymphoblastic leukemia (ALL), chronic lymphocytic leukemia (CLL), multiple myeloma (MM) and mantle cell lymphoma (MCL) with *IGH* FR1 and FR3.

## Methods

### Sample Prep




- **Sample Input volume:** ~1-2 mL
- **Specimen type:** Bone Marrow Aspirates
- **Two methods for DNA extraction**
  - Manual using QIAmp DNA Blood Midi Kit
  - Automated using the Kingfisher Flex

### Assay Workflow



- **Sample DNA Input:** 0.5 µg to 20 µg
- **MRD Positive & Negative Controls:** Run controls
- **MRD Internal Controls:** Quantitative Controls (Spiked into samples and Positive Control)
- **B-cell MRD Software:** Designed to track clonal sequences and analyze data in follow-up samples tested with the *IGH* FR1/FR3 MRD assay on an Illumina NextSeqDx® instrument, and generate reports for use by the clinical lab to report on the assessed MRD signal of each sample.

## Results: Cell vs DNA Contriving



Cell Line  
● LP-1  
● U-266

Transformed Fit to Log	
Slope	Intercept
1.09	- 0.32

- Study was conducted with two Multiple Myeloma cell lines – U-266 and LP-1 contrived in healthy bone marrow background that covered dynamic range from 2.0E-01 to 3.0E-07.
- MRD detection of cell and DNA contrived samples was observed to be equivalent with R<sup>2</sup> of 0.87 and slope of 1.09 in regards to calculated (calc.) MRD Frequency (i.e Clonal MRD Frequency)

## Results: LOB, LOD, and LOQ

- **LOB Results:** No clonal sequence was detected in 10 healthy bone marrow samples tested at 0.5 µg and 20 µg (N= 180 per master mix lot) with both *IGH* FR1 and FR3 gene targets. Thus, LOB was zero.

### IGH FR1- LOD/LOQ Results

DNA Input (µg)	Clonal cell estimates	Clonal MRD Frequency	Clonal cell estimates	Calculated (Clonal) MRD Frequency
	LOD [95% CI]	LOD [95% CI]	LOQ [95% CI]	LOQ [95% CI]
0.5	7.0 [4.9 - 12.1]	9.1E-05 [6.4E-05 - 1.6E-04]	5.4 [3.8 - 8.5]	7.0E-05 [4.9E-05 - 1.1E-04]
2	4.7 [3.5 - 7.7]	1.5E-05 [1.1E-05 - 2.5E-05]	4.6 [2.6 - 7.9]	1.5E-05 [8.5E-06 - 2.6E-05]
20	5.4 [4.0 - 8.5]	1.8E-06 [1.3E-06 - 2.7E-06]	7.3 [5.4 - 9.2]	2.4E-06 [1.7E-06 - 3.0E-06]

- LOD is highest value at 20 µg : 7.0 cells which is equivalent to 2.3E-6
- LOQ is highest value at 20 µg: 7.3 cells which is equivalent to 2.4E-6


### IGH FR3 - LOD/LOQ Results

DNA Input (µg)	Clonal cell estimates	Clonal MRD Frequency	Clonal cell estimates	Calculated (Clonal) MRD Frequency
	LOD [95% CI]	LOD [95% CI]	LOQ [95% CI]	LOQ [95% CI]
0.5	3.9 [3.0-7.2]	5.11E-05 [3.9E-05 – 9.4E-05]	6.4 [5.1-8.0]	8.3E-05[6.6E-05-1.0E-04]
2	4.7 [3.7-7.4]	1.54E-05[1.2E-05 - 2.4E-05]	12.1 [5.2-28.2]	3.9E-05[1.7E-05-9.2E-05]
20	5.7 [4.2-10.1]	1.84E-06[1.4-06 – 3.3E-06]	51.4*[6.7-396.9]	1.7E-05[2.2E-06-1.3E-04]

- LOD is highest value at 20 µg : 5.7 cells which is equivalent to 1.8E-6
- LOQ is highest value at 20 µg: 51.4 cells which is equivalent to 1.7E-5

\*LOQ is 14.4 cells which is equivalent 4.7E-6 when removing data for one sample that was an outlier

## Results: Linearity



DNA Input (µg)	Log-transformed Fit		
	Slope	Intercept	R <sup>2</sup>
0.5	1.06	-0.005	0.91
2	1.06	0.021	0.88
20	1.01	-0.542	0.85

DNA Input (µg)	Log-transformed Fit		
	Slope	Intercept	R <sup>2</sup>
0.5	1.05	-0.017	0.80
2	1.06	-0.020	0.82
20	1.05	-0.588	0.74*

\*R<sup>2</sup> = 0.85 when removing data for one sample that was an outlier

- MRD Detection was linear across the entire detection range of the B-cell MRD Assay for both *IGH* FR1 and FR3 (8.0E-03 to 2.6E-06)

## Results: Clinical Accuracy

- 10 Healthy bone marrow and 20 contrived MRD clinical samples (CLL, MM, and MCL) were tested with B-Cell MRD Assay for the *IGH* FR1 and FR3.
- The clonal sequences detected during baseline testing with the reference method (clonality assay) were used to track MRD in contrived clinical samples.
- Positive Percent agreement (PPA) and Negative Percent Agreement (NPA) was 100% for both *IGH* FR1 and FR3 gene targets.

N =220		IGH FR1	
		Detected	Not Detected
FR1 Clonality Assay	Clonal	20	0
	Non-Clonal	0	200

N =220		IGH FR3	
		Detected	Not Detected
FR3 Clonality Assay	Clonal	20	0
	Non-Clonal	0	200

## Results: Carryover and Crossover Contamination


Study	Read Cutoff (MRD Detection)	Minimum Threshold Reads	Rate 95% CI ± 2 sd	Average Rate
Crossover	≥ 2 reads	≥165,160	1.21E-05	1.11E-06
Carryover	≥ 2 reads	≥ 19,730	1.01E-04	7.91E-06

- The clonal sequence of a sample is searched in a previous run for carryover analysis and within the same run for crossover contamination. No two samples should have the same clonal sequence in a previous or within the same run.
- If target reads of a sample are greater than minimum threshold reads, a Z-score is calculated to determine the probability (p-value) for a sequence to be detected.
  - If p-value < 0.05 then the sample is considered detected, it's a true positive
  - If p-value ≥ 0.05 then the sample is considered not detected, and signal is a false positive.

## Results: Precision


### Sadler's Precision Profile

#### IGH FR1



N (Total data points)	%CV per Malignant Cell Level											
	216	216	216	216	216	216	216	216	216	216	216	216
Malignant Cell Input	615.4	215.4	61.5	21.5	8.0	6.2	4.0	2.2	1.5	0.6		
DNA Input	12.9	10.0	14.3	10.4	0.0	5.3	13.3	0.0	22.6	24.1		
Operator	3.4	0.0	0.0	0.0	5.1	0.0	0.0	0.0	0.0	3.8		
MM Lot/Instrument /Run/Day	0.0	0.5	0.0	0.0	7.7	2.0	9.4	0.0	4.6	0.0		
Random Error	35.4	34.3	37.2	52.9	59.0	63.6	71.4	97.2	102.3	173.8		
Precision (Total Variation)	37.8	35.8	39.8	53.9	59.2	63.8	72.7	97.2	104.7	175.5		

#### IGH FR3



N (Total data points)	%CV per Malignant Cell Level											
	216	216	216	216	216	216	216	216	216	216	216	216
Malignant Cell Input	615.4	215.4	61.5	21.5	8.0	6.2	4.0	2.2	1.5	0.6		
DNA Input	25.8	30.5	30.4	19.5	30.8	26.5	26.4	24.6	25.9	30.2		
Operator	13.3	14.7	0.0	7.0	13.1	16.3	0.0	0.0	6.5	0.0		
MM Lot/Instrument /Run/Day	10.7	11.1	9.4	8.6	14.6	0.0	0.0	0.0	0.0	35.1		
Random Error	41.9	44.2	46.2	60.3	64.0	72.5	78.3	100.6	110.4	178.4		
Precision (Total Variation)	52.1	56.8	56.1	64.3	73.6	78.9	82.6	103.6	113.6	184.4		

- Study was conducted with two master mix lots, 2 operators, 2 instruments, 3 DNA inputs with 9 clinical samples (CLL, MM and MCL).
- CV% decreases with increase in Clonal Malignant cell level

## Conclusions

- *IGH* FR1 exhibited a LoB at 0, LoD at 2.3E-06, LoQ at 2.4E-06, good linearity with R<sup>2</sup>>0.88 across 4.5 logs and precision with %CV ranging from 27.6% to 63.5% at the LoD.
- *IGH* FR3 exhibited a LoB at 0, LoD at 1.8E-06, LoQ at 1.7E-05, good linearity with R<sup>2</sup>>0.74 across 4.5 logs and precision with %CV ranging from 46.6% to 89.4% at the LoD.
- The B-cell MRD Assay demonstrated excellent clinical performance; yielding 100% agreement in an accuracy study that included 20 clinical specimens (CLL, MM and MCL) and 10 healthy donors.