* invivoscribe

Ô

Improving Lives with Precision Diagnostics*

Product Catalog

invivoscribe.com

Table of Contents

Introduction

- **03** Table of Contents
- 05 Global Invivoscribe Locations
- 07 Board of Directors
- 09 Executive Leadership
- 11 Leadership Team
- 13 Best-in-Class Assays and Reagents
- 16 LabPMM Clinical Lab Services
- 17 Streamlined CDx Partnerships
- **19** Global Distribution Network

Clonality Overview

21 Background23 Test Algorithm for Suspect Lymphoproliferations

Hematological Test Menu

25 Hematological Test Menu

Next Generation Sequencing (NGS)

27 LymphoTrack® Dx (CE-IVD) Assays

- 29 LymphoTrack Dx IGHV Leader Somatic Hypermutation Assay
- 31 LymphoTrack Dx IGH FR1/2/3 Assays
- **35** LymphoTrack Dx *IGK* Assays
- 37 LymphoTrack Dx TRG Assays
- **39** LymphoTrack Dx *TRB* Assays
- 41 LymphoTrack Dx Bioinformatics Software

43 LymphoTrack® (RUO) Assays

- 45 LymphoTrack IGHV Leader Somatic Hypermutation Assay
- 47 LymphoTrack *IGH* FR1/FR2/FR3 Assays
- **51** LymphoTrack *IGK* Assays
- 53 LymphoTrack TRG Assays
- **55** LymphoTrack *TRB* Assays
- 57 LymphoTrack Software S5/PGM™
- 58 LymphoTrack Software MiSeq™
- **59** LymphoTrack Software Analysis Outputs
- 60 LymphoTrack Flex (RUO) Assays
- 62 LymphoTrack Flex TRG Assay
- 62 LymphoTrack Flex *IGH* FR1 Assay
- 62 LymphoTrack Flex IGHV Leader Assay
- 64 LymphoTrack Enterprise Software

66 MRD Solutions (RUO)

68 LymphoTrack MRD Software 70 FLT3 ITD MRD Assay 72 NPM1 MRD Assay

Capillary (CE-IVD) Assays

74 IdentiClone® (CE-IVD) Assays

- 76 IdentiClone IGH + IGK B-Cell Clonality Assays
- **78** IdentiClone *IGH* Gene Clonality Assays **80** IdentiClone *IGK* Gene Clonality Assays
- 82 IdentiClone *IGL* Gene Clonality Assays
- 84 IdentiClone *TCRB* Gene Clonality Assays
- 86 IdentiClone T-Cell Receptor Gamma Gene Rearrangement Assay 2.0
- 88 IdentiClone TCRD Gene Clonality Assays

90 LeukoStrat[®] CDx Mutation Assays

- 92 LeukoStrat CDx FLT3 Mutation Assay (IVD) USA
- 94 LeukoStrat CDx FLT3 Mutation Assay (CE 2797 IVD)
- 96 LeukoStrat CDx FLT3 Mutation Assay Japan

Gel and Capillary (RUO) Assays

- 100 IGH + IGK B-Cell Clonality Assays
 102 IGH Gene Rearrangement Assays
 104 IGH Gene Clonality Assays
 106 IGK Gene Clonality Assays
 108 IGL Gene Clonality Assays
 110 TCRB Gene Clonality Assays
 112 T-Cell Receptor Gamma Gene Rearrangement Assay 2.0
 114 T-Cell Receptor Gamma Gene Rearrangement Assays
 116 TCRD Gene Clonality Assays
 118 BCL1/J_H Translocation Assay
 120 BCL2/J_H t(14;18) Translocation Assay
 124 BCR/ABL t(9;22) Translocation Assays
 126 PML/RARa t(15;17) Translocation Assays
- 128 IGH Somatic Hypermutation Assays v2.0

FLT3 Mutation (RUO) Assays

130 FLT3 Mutation Assays

Analyte Specific Reagents

132 Analyte Specific Reagents

Controls, Reagents, and Enzymes

- 136 DNA Controls
- 139 LymphoTrack[®] Low Positive Controls
- 139 LymphoQuant® Internal Controls
- 140 RNA Controls
- 143 RNA Sensitivity Panels
- 144 Master Mix Controls
- 144 Enzymes

Custom Products

145 Custom Products

Reference

- 147 Clonality Testing Workflow Methods
- 149 Illumina® MiSeq[™] Workflow Summary
- **150** Thermo Fisher Scientific® Ion S5^m PGM^m Workflow Summary
- 151 LymphoTrack Flex Workflow Summary
- 152 Next-Generation Sequencing Menu
- 153 Gel and Capillary Electrophoresis Menu
- 154 Common Technical Support Questions

Product List

156 Product List by Catalog Number

For Legal Notices related to these products, visit: https://invivoscribe.com/legalnotices/

Improving Lives with Precision Diagnostics®

A technology-driven international company focused on precision diagnostics with best-in-class products and services that optimize patient care and accelerate drug approvals worldwide.

LEARN MORE: invivoscribe.com/

Global Invivoscribe® Locations

A long and proud history of working with scientists, laboratorians and clinicians from around the world. We offer product support in over 160 countries. Our distributor partners are in locations around the globe.

LEARN MORE:

invivoscribe.com/global-distributors/

Munich, Germany



Shanghai, China

Tokyo, Japan



United States

Invivoscribe®, Inc.

San Diego, United States

ivoscriba

10222 Barnes Canyon Rd, Bldg. 1, San Diego, CA 92121, United States Toll Free: +1 866 623 8105 | Phone: +1 858 224 6600 Fax: +1 858 224 6601 | Email: customerservice@invivoscribe.com

LabPMM[®], LLC

10222 Barnes Canyon Rd, Bldg. 1, San Diego, CA 92121, United States Phone: +1 858 224 6650 | Fax: +1 858 224 6655 Email: support@labpmm.com

Germany

LabPMM[®], GmbH

Zeppelinstr. 1, 85399 Halbergmoos (Munich Area), Germany Phone: +49 (0) 89 904 299 800 Email: info@labpmm.de

China

Invivoscribe[®] Diagnostic Technologies (Shanghai) Co., Ltd.

Unit 6302, #6 Building, 338 Jia Li Lue Road, Zhangjiang Hi-Tech Park, Pu Dong District 201203 Shanghai, China Phone: +86 21 6162 9669 | Email: sales@invivoscribe.com

Japan LabPMM®, GK

Life Science and Environment Research Center (LiSE) 4th Floor, 3-25-13 Tonomachi, Kawasaki-ku, Kawasaki-shi, Kanagawa 210-0821, Japan Phone: +81 (0) 44 281 1500 | Fax: +81 (0) 3 6745 9346 Email: services@labpmm.co.jp or info@labpmm.de

France

Invivoscribe® Technologies, SARL

c/o Ficorec Domiciliation Services, 327, Boulevard Michelet, 13009 Marseille, France Phone: +33 (0) 4 42 01 8 10 Email: sales-eu@invivoscribe.com

%invivoscribe® Board of Directors

LEARN MORE: invivoscribe.com/company/board-of-directors/



Jeffrey Edward Miller, Ph.D.

Dr. Miller is a scientist, inventor, and entrepreneur focused on Improving Lives with Precision Diagnostics[®] – coupling drug trials and therapeutic treatment regimens with optimized clinically-actionable diagnostic methods in order to select the correct patients and then monitor and track their response throughout the course of their disease. He received his undergraduate degree in Biochemistry from UCLA and a combined Ph.D. in Biochemistry & Molecular Biology from UCSB. Prior to starting Invivoscribe, Dr. Miller had more than twenty years of combined experience in protein biochemistry, cellular and molecular immunology, cardiac physiology, virology, and molecular biology, experience he had developed working in laboratories at the Medical School, Department of Chemistry, Molecular Biology Institute at UCLA, while earning his Ph.D. at UCSB, and as a postdoctoral scientist at Applied Molecular Evolution. He also spent several years at Quest Diagnostics Nichols Institute, setting up the molecular oncology laboratory and developing and launching PCR-based molecular assays for infectious disease and hematopathology.



James Isaacs, Jr., JD.

James B. Isaacs, Jr. has practiced law since 1983. He currently serves as Licensing and Contracts Counsel at Invivoscribe. Jim attended Stanford University and Yale Law School, then began his career at the Los Angeles law firm of O'Melveny & Myers. As a trial lawyer and later in-house counsel with a focus on intellectual property disputes, Jim has successfully defended a myriad of businesses and individuals; in plaintiffs' actions he has obtained and collected multi-million dollar judgments in the United States and abroad. As a businessman and co-founder of Invivoscribe, Jim has been active in the legal and commercial affairs of the company since 1995.



Gary Clouse, JD.

Gary Clouse has practiced as a litigator and business attorney in Southern California for more than three decades. He currently serves as Corporate Secretary and Legal Counsel for Special Projects at Invivoscribe. Clouse is a graduate of Indiana State University and Northwestern University School of Law. Following law school, he clerked for the federal Seventh Circuit Court of Appeals in Chicago. He began his legal practice at the law firm O'Melveny & Myers in Los Angeles. Clouse is one of the founders of invivoscribe.





Mitchell Kronenberg, Ph.D.

Dr. Kronenberg received a B.A. from Columbia University, a Ph.D. from the California Institute of Technology, and served on the faculty of the UCLA School of Medicine from 1986-1997. He joined the La Jolla Institute for Immunology in 1997, and currently is a Professor, the President Emeritus and a member of the institute Board of Directors. Dr. Kronenberg was the President of the institute from 2003 to 2021 and Chief Scientific Officer from 2021-2023. The Institute has grown in accomplishment and reputation under his leadership. Dr. Kronenberg's research interests include natural killer T cells, other innate lymphocytes such as MAIT cells and ILC, regulation of mucosal immunology and the pathogenesis of inflammatory bowel disease. He has co-authored more than 400 publications, and is a fellow of the American Association for the Advancement of Science (AAAS), a Distinguished Fellow of the American Association of Immunologists, recipient of an NIH MERIT award and is an Institute for Scientific Information (ISI) Highly Cited Scientist. In 2016, he was named the most admired CEO (large nonprofit organization category) by the San Diego Business Journal. He is an advisor to a number of for-profit and nonprofit organizations, including service as a member of the Board of Scientific Counselors for Basic Science, National Cancer Institute. He has been involved with Invivoscribe since its founding.

Stephen Wilson, Ph.D.

Dr. Wilson is an immunologist and biotech entrepreneur with more than 30 years of experience in biomedical research, development and executive management. He earned his Ph.D. from the University of Arizona's College of Medicine and joined the La Jolla Institute for Immunology (LJI) in 1997 where he was a fellow of the National Multiple Sclerosis Society and National Institutes of Health. Dr. Wilson entered management at LJI, eventually serving as LJI's Chief Operating Officer. During his tenure the Institute grew substantially, becoming a global leader in immunology, conducting nearly \$1B in total research operations and ranking as the #1 place to work in the world by The Scientist magazine. His original research has been published in high impact journals. Dr. Wilson has served as principal or co-principal investigator on more than \$75M in competitive grants and awards, including the development of the 2021 X-Prize winning rapid CoVID-19 diagnostic test. Currently, Dr. Wilson is serving as the Chief Executive Officer of Foundation Immune Engineering's Botnar Institute.



Hardwick Simmons, MBA

Hardwick 'Wick' Simmons retired as the Chairman and CEO of The Nasdaq Stock Market, Inc. in May of 2003. Prior to Nasdaq, he served as President and CEO of Prudential Securities Inc., a major investment management and securities brokerage firm. Simmons is a former chairman of the Securities Industry Association, a former director of the Chicago Board Options Exchange and a former president of the New York Bond Club. He is currently a director of Lionsgate Entertainment Corp., president of Stonetex Oil Corp., and a trustee of Woods Hole Oceanographic Institution. He is a graduate of Harvard College and Harvard Business School and served in the U.S. Marine Corps Reserve from 1960 until 1966.

*invivoscribe Executive Leadership

LEARN MORE: invivoscribe.com/company/leadership-team/







Jeffrey Edward Miller, Ph.D.

Founder, Chief Scientific Officer, and Chief Executive Officer

Dr. Miller is a scientist, inventor, and entrepreneur focused on Improving Lives with Precision Diagnostics® – coupling drug trials and therapeutic treatment regimens with optimized clinically-actionable diagnostic methods in order to select the correct patients and then monitor and track their response throughout the course of their disease. He received his undergraduate degree in Biochemistry from UCLA and a combined Ph.D. in Biochemistry & Molecular Biology from UCSB. Prior to starting Invivoscribe, Dr. Miller had more than twenty years of combined experience in protein biochemistry, cellular and molecular immunology, cardiac physiology, virology, and molecular biology, experience he had developed working in laboratories at the Medical School, Department of Chemistry, Molecular Biology Institute at UCLA, while earning his Ph.D. at UCSB, and as a postdoctoral scientist at Applied Molecular Evolution. He also spent several years at Quest Diagnostics Nichols Institute, setting up the molecular ancelular and eveloping and launching PCR-based molecular as for infectious disease and hematopathology.

Meghna Bhatnagar, MBA Chief Financial Officer

Meghna Bhatnagar joined Invivoscribe in 2010 as Chief Financial Officer. In this role she is responsible for leading the Invivoscribe global finance organization, along with human resources and information technology. Since her arrival, Ms. Bhatnagar has played an integral role in directing all aspects of company strategy, planning and operations. Ms. Bhatnagar has over 20 years of experience building and leading finance and IT teams in global companies. Prior to joining IVS, she served as COO of Radiant Technologies, a technology company focused on providing business management solutions to small and medium sized companies where she was responsible for leadership and development of an entire project delivery team with full P&L responsibility. She played a key role in guiding overall strategy of the company and at the same time provided leadership for operational improvements.

Dr. Bradley Patay, M.D. Chief Medical Officer

Dr. Patay is dedicated to improving health by integrating genomic knowledge into medical care. He has authored numerous articles in this field, presented at multiple conferences and has been featured in Bloomberg Business Week. He has been head of the internal medicine section at Scripps Torrey Pines, worked as an Assistant Professor at STSI, and has been a founding member of the Board, and Vice President of the College of Genomic Medicine, which was established in 2010 to educate physicians and other health care professionals about genomic medicine. His diverse clinical experience prior to joining Scripps Clinic in 2005 includes four years as an internist and pediatricina at Neighborhood Healthcare, a private, nonprofit community healthcare practice, and at Palomar Hospital. At these institutions, Brad cared for a wide range of patients, from neonates to the elderly, in both intensive care and the general wards. Through his service on several committees, he helped improve health care institutions' systems.







Tania Obranovich, BSc (Hons) PhD LLB

Chief IP and Legal Strategy Officer

Dr. Obranovich has more than 25 years' experience developing and managing global IP portfolios in the biotechnology sector, in particular in relation t o molecular, cellular and genomic technologies as they relate to diagnostics and therapeutics. She has acquired this experience both through private practice, having spent more than a decade as a partner in a top tier IP firm, and as legal counsel for a US based diagnostics company where she ran patent and trade mark prosecution and drafting, in addition to overseeing all legal, compliance and corporate governance matters. Tania has practiced throughout her career as both a lawyer and patent attorney. Prior to these roles she completed a PhD and post-doctoral fellowship in immunology. In addition to her legal practice, Tania has lectured extensively in patent law at The University of Melbourne and Monash University and also speaks widely in relation to IP issues. She was heavily involved with navigating the Australian and US legal and policy issues that arcse as a result of the Myriad gene patent litigation in both these jurisdictions. Tania neceived a Federal Government appointment to the Professional Standards Board for Patent and TM Attorneys in 2012 and has held several other Board positions. Tania has served as Chief IP and Legal Strategy Officer to Invivoscribe since 2019.

Dr. Meindert Niemeijer, Ph.D., M.Sc.

Chief Information Officer

Dr. Niemeijer came to IVS from Digital Diagnostics where he led the teams that designed, built, documented and put into production the first autonomous diagnostic Artificial Intelligence (AI), cleared by the FDA. He is an experienced technology leader and deep learning expert passionate about developing AI based medical devices. During his career, he has built a deep expertise in medical device (software) development processes and the way in which these interact with quality management and other regulatory requirements. He loves building and developing technology, processes and teams. Dr. Niemeijer obtained a Ph.D. in Medical Image Analysis and Machine Learning from Utrecht University in the Netherlands. In addition, he has an M.Sc. in Medical Computer Science, also from Utrecht University.

Andreia de Albuquerque, Ph.D.

Chief Commercial Officer

Dr. Andreia de Albuquerque joined Invivoscribe in 2015 and currently serves as the Chief Commercial Officer. In this role, she is responsible for leading Invivoscribe's global partnerships with Pharma and Biotech's, focusing on companion diagnostics development, custom assay development, clinical trial testing, corporate partnerships and other strategic relationships. Prior to joining Invivoscribe, Dr. de Albuquerque has held roles in precision molecular diagnostics, focusing on development, implementation and validation of PCR and NGS assays for diagnostic and research purposes. Her work has resulted in several peer-reviewed publications and participations in international conferences. During this time, she was also responsible for overseeing laboratory diagnostic teams and operations. Her extensive expertise in precision diagnostics and understanding of assay development and clinical testing, paved the way into Business Development, where she focus on accelerating drug approvals, by building successful partnerships, providing technical expertise and working closely with biomarker leaders, sample operations, clinical science and regulatory to implement CDx, IVD and R&D assays in clinical trials.

Dr. de Albuquerque earned her degree in Biology from the University of Coimbra (Portugal) and a combined Ph.D. in Oncology and Molecular Biology from the University of Dresden (Germany).

* invivoscribe Leadership Team

LEARN MORE: invivoscribe.com/company/leadership-team/



Jordan Thornes

V.P., Global Clinical Laboratory Operations

Jordan Thornes joined Invivoscribe in 2004 and currently serves as the Global Director of Clinical Labs. In this role, he is responsible for leading Invivoscribe's global network of clinical labs, The Laboratory for Personalized Molecular Medicine®, (LabPMM®). There are currently four such LabPMM facilities which are strategically located globally in San Diego, USA; Hallbergmoos (Munich), Germany; Shanghai, China and Kawasaki (Tokyo), Japan. LabPMM aids patients and pharmaceutical companies worldwide with globally standardized and clinically-actionable diagnostic methods. LabPMM utilizes proprietary biomarkers to select the correct patients and monitor and track their critically important response throughout the course of their disease.



Ying Huang, Ph.D. V.P., Product Development

Dr. Huang joined Invivoscribe in 2010 and currently serves as the Sr. Director of Product Development. In this role, she is responsible for leading Invivoscribe's product development activities to address unmet clinical needs in the fields of molecular diagnostic and precision medicine. Dr. Huang has over 25 years of research and industry experience with a broad knowledge in technology and asay development. She is passionate about developing novel methods using cutting edge technologies to improve lives with precision diagnostics. She has successfully led development teams through entire product development life cycle, from feasibility to RUO, CLIA and IVD grade product launches and post market support under respective regulatory standards. Prior to Invivoscribe, Dr. Huang worked at Illumina, ACEA Biosciences, Nanogen and MD Anderson Cancer Center. In those roles, she developed a variety of techniques ranging from isolating CTC using microfluidic devices to detecting biowarfare agents using electronic microarray. Dr. Huang was the PI for research grants from DARPA and NIH. She is the co-author of more than 60 publications and meeting abstracts and is the co-inventor of 9 patents. Dr. Huang received her Ph.D. in Electrical Engineering from Bangor University, UK and B.E. in Electronic Engineering from Xidian University, China.



Jason Gerhold

V.P., Global Quality, Regulatory, and Clinical Affairs

Jason Gerhold joined Invivoscribe in 2012 and serves as the Global Director of Regulatory, Quality, and Clinical Affairs focusing on developing high quality diagnostics that meet both international pathology needs and support pharma partner's drug development programs. He directly manages employees in the US, Japan, and China, and is responsible for worldwide compliance to relevant laws, regulations, and standards. Mr. Gerhold has over 20 years of experience in the biotech industry with biologic therapeutics and diagnostics, including process development, regulatory, quality, and clinical roles. Mr. Gerhold has worked in small, medium and large companies; successfully registering IVD and CDx products in several international markets, developing, improving, and maintaining quality management systems while serving as the management representative to regulatory authorities, and leading clinical teams managing international registrational studies. He is known for building strong rapport with regulators and working well with research and pharma teams resulting in highly successful partnerships. He has completed trainings and/or accreditations for ISO 13485 Lead Auditor, Six Sigma Greenbelt, Clinical Trials Administration, Project Management, and CTD writing.

Best-in-Class Assays and Reagents

Invivoscribe[®] provides a full range of standardized CE-marked *in vitro* diagnostic cGMP products for hematology-oncology, as well as RUO assays, analyte specific reagents (ASRs), and DNA & RNA controls.

LEARN MORE: invivoscribe.com/products/

Next-Generation Sequencing (NGS)

NGS is a powerful, high-throughput DNA sequencing technology that allows for massively parallel sequencing of millions of DNA fragments in a single sequencing run. Our NGS assays provide comprehensive solutions from diagnosis through Measurable Residual Disease tracking.

ABI Fluorescence Detection and Gel Detection

We exclusively offer a comprehensive selection of PCR-based assays for ABI fluorescence detection and gel detection, including targeted *FLT3* ITD and TKD mutation assays, B- and T-cell clonality assays (based on EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936), and translocation assays.

Controls and Reagents

Invivoscribe[®] offers an extensive range of General Purpose Reagents (GPRs) and Research Use Only (RUO) reagents. Our high quality, reliable DNA and RNA controls, master mixes, and enzymes, are all manufactured under cGMP conditions to guarantee reliability.

Companion Diagnostics (CDx)

Streamline drug approvals leveraging our Companion Diagnostic (CDx) assay development, manufacturing and regulatory expertise to accelerate drug approvals. Custom assays are designed and manufactured under an ISO certified quality management system.

LabPMM® Clinical Lab Services

Invivoscribe's international network of clinical labs, Laboratory for Personalized Molecular Medicine® (LabPMM®), provides globallystandardized testing of novel and proprietary biomarkers that are critical for patient care.

LEARN MORE: invivoscribe.com/clinical-lab-services/

NGS Cancer Panels

Identification of chromosome abnormalities and actionable driver mutations aid in the stratification of patient risk, and provide critical insight for patients with leukemia or lymphoma.

NGS MRD Tests

Track disease and trend treatment response by monitoring clonal B- (Ig) and T-cell (TCR) gene rearrangements as well as *FLT3* ITD and *NPM1* mutations with unprecedented sensitivity.

Clonality Tests

Malignant lymphoproliferative disorders collectively represent a large portion of hematological cancers and are characterized by reduced population diversity of Ig and TCR rearrangements. Confirm diagnostic suspicions and identify hematological neoplasms with LabPMM's NGS Clonality services, available for Ig, TCR, and somatic hypermutation (SHM).

Companion Diagnostic (CDx) Tests

CDx enable the identification of biomarkers that indicate the patients most likely to benefit (or alternately not benefit) from specific or combinatorial targeted therapies. Their utilization enhances patient care and the potential for improved outcomes. The fms related tyrosine kinase 3 (*FLT3*) is one of the most commonly mutated genes in acute myeloid leukemia (AML), occurring in approximately 30% of patients at the time of diagnosis.¹ Although generally associated with normal cytogenetics where patients have standard risk of relapse, *FLT3* mutations have also been identified in sub-groups of patients with chromosomal abnormalities that are associated with high risk of disease relapse.²⁻³ The LeukoStrat CDx *FLT3* Mutation Assay is a globally standardized test validated for the detection of internal tandem duplications (ITD) and tyrosine kinase domain (TKD) mutations of the *FLT3* gene.

Custom Assays

Invivoscribe's experienced regulatory team has over 50 product registrations and multiple companion diagnostic approvals, allowing for optimal strategy, streamlined development and milestone achievement. With an ever-changing regulatory landscape, why take a chance with anyone else?

Targeted Genes

FLT3 and NPM1 are offered to detect targeted mutations. FLT3 and NPM1 mutations can be used for risk stratification and prognosis of AML.

Multiparametric Flow Cytometry

Complement NGS testing with globally-standardized, 12-color flow cytometry testing. Our comprehensive assay designs allow greater specificity and sensitivity using smaller sample volumes as compared to 4- or 6-fluorochrome tests. Some of the many MFC applications include leukemia and lymphoma diagnosis and screening and rare event analyses such as MRD.

References

- 1. Acute Myeloid Leukemia, *Clinical Practice Guidelines in Oncology*, (v.2.2014) National Comprehensive Cancer Network.
- 2. Lowenberg, B. et al. (1999) "Acute myeloid leukemia." N Engl J Med 341(14):1051-62.
- Thiede, C. et al. (2002) "Analysis of FLT3-activating mutations in 979 patients with acute myelogenous leukemia: association with FAB and identification of subgroups with poor prognosis." *Blood* 99(12): 4326-35.

Streamlined CDx Partnerships

Companion Diagnostics represents a critical milestone in precision medicine and plays a pivotal role in the validation of targeted therapies. Use of a single partner for CDx development minimizes complexity and risk which could affect drug approvals. Our Streamlined CDx[™] program supports all stages of development, validation, regulatory submission, and commercialization. Let us help you accelerate your drug approval.

LEARN MORE: invivoscribe.com/partnerships/

Product Development

Design Controls: All biomarker assays & software developed under design controls.

Bioinformatics Software: Comprehensive solutions to prevent errors, streamline workflow and reduce turn-around times.

Identify and Track: NGS panels and targeted assays to identify clonal architecture and track disease burden.

Custom Development: Custom biomarker assay development.

Complementary MRD Assays: MRD assays may be used as surrogate endpoints for P3 clinical trials.

Manufacturing

cGMP Compliant: US FDA/CDRH registered, EN ISO 13485:2016 certified manufacturing facility based in San Diego.

CE-marked IVDs and CDx Manufacturing: CE-Marked, PMDA and FDA approved CDx, and >50 CE-IVDs including gel, capillary electrophoresis and NGS assays.

Assay Development: IUO & RUO assays, CE-marked IVDs, & CDx.

Controls and Reagents: DNA / RNA controls, ASRs, GPRs, MRD controls & proficiency panels.

Clinical Lab Services

Clinical Lab Experience: More than a dozen years of clinical reference lab experience.

Internationally Standardized: Internationally standardized CDx, flow cytometry and biomarker testing with internationally accredited labs serving the US, Europe, and Asia.

Comprehensive Panels and Targeted Assays: Comprehensive NGS and MFC-based panels to screen, identify and track clinically actionable biomarkers.

MOA Profiling of Novel Compounds: Using primary AML patient samples – IC50 determination, mRNA analysis, genetic, and phenotypic characterization.

Worldwide Enrollment: Testing services have supported hundreds of enrollment sites worldwide.

Global Regulatory Expertise

Accredited and Proven: EN ISO 13485:2016 accredited. Experienced staff & proven Quality Management System. Registered Medical Device Establishment: Registered Medical Device Establishment with the US FDA, KFDA, Saudi Arabia, and the MHLW/PMDA. Multiple CDx Approvals: Multiple companion diagnostics that are CE 2797 IVD certified, FDA and PMDA/MHLW approved. CE-marked IVDs: More than 50 CE-marked IVDs available in the EU and select ROW markets; >60 tests included in the ARTG in Australia. Marketing Authorization Holder (MAH): Marketing Authorization Holder (MAH) and national reimbursement for CDx in Japan.

Global Distribution Network

A long and proud history of working with scientists, laboratorians and clinicians from around the world. We offer product support in over 160 countries. Our distributor partners are in locations around the globe.

LEARN MORE: invivoscribe.com/global-distributors/

North America

Canada Mexico United States

South America

Argentina Brazil Chile Colombia Peru Uruguay

Central America

Caribbean Costa Rica El Salvador Guatemala Honduras Nicaragua Panama

Europe

Austria

Belarus

Belgium

Bosnia &

Bulgaria

Croatia

Cyprus

Denmark

Estonia

Finland

France

Germany

Greece

Hungary

Iceland

Ireland

Kazakhstan

Italy

Latvia

Malta

Moldova

Monaco

Norway Poland

Portugal Romania Russia Serbia Slovakia Slovenia Spain Sweden Switzerland Ukraine

United Kingdom

Netherlands

Lithuania

Luxembourg

Macedonia

Herzegovina

Czech Republic

Middle East

Bahrain Iran Iraq Israel Jordan Kuwait Lebanon Oman Pakistan Palestine Qatar Saudi Arabia Turkey United Arab Emirates

Africa

Algeria Angola Burkina Faso Burundi Cameroon Chad Côte d'Ivoire Egypt Equitorial Guinea Eswatini Ethiopia Gabon Ghana Kenya Madagascar Mali Morocco Mozambique Namibia Nigeria Republic of the Congo Senegal Sierra Leone South Africa Tanzania Tunisia Uganda Zambia Zimbabwe

Asia

Afghanistan Armenia Azerbaijan Bangladesh China Hong Kong India Indonesia Japan Kyrgyzstan Malaysia Nepal Philippines Singapore South Korea Sri Lanka Taiwan Tajikistan Thailand Turkmenistan Uzbekistan Vietnam

Australia and New Zealand

Australia New Zealand

Background

PCR and NGS-Based Assessment of Clonality in Hematologic Malignancies

For nearly 30 years, Invivoscribe® has developed, manufactured, and commercialized the gold-standard molecular hematopathology assays and reagents for gel and capillary electrophoresis detection, and most recently, next-generation sequencing instruments. These standardized, cGMP manufactured assays and reagents were developed and validated using standardized workflow and optimized primer sets, reagents and controls.

A number of our products were developed in collaboration with studies conducted by the EuroClonality BIOMED-2 concerted action group; these capillary based products have provided reliable methods for clonality detection that have withstood the test of time.

We have never accepted the status quo, so our comprehensive menu of clonality assays continues to evolve. All of our NGS-based clonality assays were developed in-house together with accompanying bioinformatics software by our Invivoscribe R&D team. Developed under full ISO 13485 design control, these assays and bioinformatics software were designed to run on several next-generation sequencing platforms. These NGS-based assays are several generations ahead of capillary-based products.

Our comprehensive bioinformatics software not only provides critical information on the presence of clonality, but also identifies the sequence information required to track clones in subsequent samples.

The unique process of gene rearrangement that occurs within the immunoglobulin (Ig) and T-cell receptor (TCR) gene loci during immune cell development and maturation generates a vast pool of genetically distinct cells. The resulting diverse population of lymphocytes displays an astonishing number of diverse antigen receptors, each coded in the DNA by a unique sequence, and each displayed on the cell surface, or as antibodies in the blood unique to a given cell.¹² This diversity allows the adaptive immune system to carry out its role in protecting the human body by recognizing the infinite number of pathogens it might encounter during a lifetime.

In sum, lymphoid malignancies are characterized by size- and sequence-specific rearrangements within these loci, which result from the transformation and subsequent proliferation from a single cell. The associated cellular population typically shares one or more cell-specific or "clonal" antigen-receptor gene rearrangements. The detection of these clonal cells forms the basis for clonality assessment in leukemia, lymphoma, and hematologic disease. These methods can also be used to assess somatic hypermutation (SHM) and to study measurable residual disease (MRD).

Malignant cells that remain in the bone marrow following treatment are a major cause of disease relapse. MRD testing by NGS offers enhanced sensitivity and specificity (compared to MRD testing by flow cytometry), and allows residual cells to be identified at very low levels and monitored throughout the different stages of disease.

Invivoscribe[®] can provide you with the necessary tools to accommodate your needs. From gel detection to NGS, we can help you accurately identify and track hematologic biomarkers. For more information on the detection methods available and the biomarkers offered, please refer to the respective product sections of this catalog.

Immunoglobulin and T-Cell Receptor Gene Rearrangements and the Principle and Method of Clonality Testing

The adaptive vertebrate immune system produces a repertoire of immunoglobulin and T-cell receptor molecules using a relatively limited number of heritable germline gene segments. Somatic gene rearrangement is the fundamental mechanism used to generate different immunoglobulin and T-cell receptor molecules, each with unique binding specificity. Lymphocytes undergo gene rearrangements to assemble CDR3 coding regions that are unique in both size and DNA sequence. Since leukemias and lymphomas arise from the malignant transformation of a single cell, they share clonal rearrangement(s) of the antigen receptor genes. This is the basis for clonality testing.³

Test Algorithm for Suspect Lymphoproliferations

Developed in concert with the EuroClonality/BIOMED-2 group for PCR-based clonality assessment of suspected B- and T-cell lymphoproliferative disorders.



If a definitive diagnosis is not possible following:



IGH $(D_{H}-J_{H})$ preferably with IGL –

No evidence of clonality

Clonal



*Previously known as BCL1/J_H. Results should be considered in the context of all available clinical, histological and immunophenotypic data.

TRD

Hematological Test Menu

References

Kakizuka, A et al. (1991) *Cell*, 66:675-684.
 Rimokh, R et al. (1994) *Blood*, 83:1871-1875.
 Gribben, JG et al. (1994) *Blood*, 83:1871-1875.
 Gribben, JG et al. (1994) *Leukemia*, 8:1688-1895.
 Gallagher, RE et al. (1995) *Blood*, 86:1540-1547.
 De Boer, CJ et al. (1995) *Blood*, 86:1540-1547.
 De Boer, CJ et al. (1995) *Blood*, 86:343-56.
 Schlenk RF et al. (2008) *N Engl J Med.*, 356:1909-191
 Gonzalez, D et al. (2003) *Leukemia*, 17:1051-1057.
 Gonzalez, D et al. (2003) *Leukemia*, 17:1398-1403.
 Ghia, P et al. (2007) *Leukemia*, 21:201-61.
 Evans, PAS et al. (2007) *Leukemia*, 21:215-221.
 Davi, F et al. (2008) *Leukemia*, 22:212-214.
 Döhner, H et al. (2000) *Blood*, 82:1264-1269.
 Hughes, T et al. (2006) *Blood*, 80:2983-2990.
 Ralph, QM et al. (1993) *Blood*, 82:202-20.
 Juérsz-Velázquez, MR et al. (2013) *Clinical*

Epidemiology of Acute Lymphoblastic Leukemic - From the Molecules to the Clinic, 193-235.

- Raimondi, S. (2007) Atlas Genet Cytogenet Oncol Haematol., 11(4):328-339.
- Vardi, A et al. (2015) *Clin Cancer Res.*, 22(1):167-174.
 Tan, BT et al. (2008) *J Mol Diag.*, 10(6):502–512.
 González, D et al. (2007) *Blood*, 110:3112-3121.
- 26. Langerak AW et al. (2007) *Exp Opin Med Diag.* 1(4):451-461.
- 27. Arcila, ME et al. (2019) *J Mol Diag.*, 2: 330-341.
- Campo, E et al. (2011) Blood, 117(19): 5019–5032.
 Smirnova, SY et al. (2016) Acta naturae, 8(4):100–109.
 Szczepanski, T et al. (2001) Leukemia, 15:1415–1423.
 Szczepanski, T et al. (1994) Leukemia, 12:1081–1088.
 Beishuizen, A et al. (1994) Leukemia, 15:121–127.
 Tumkaya T et al. (2001) Leukemia, 15:121–127.
 Felix, CA et al. (1990) J Clin Oncol., 8(3):431-42.
 Rosenquist R et al. (2017) Leukemia, 13:1477–1461.
 Arons, E et al. (2006) Modern Pathology, 19:1498–1505.
 Srour SA et al. (2017) Br J Haematol, 177(2): 331.
 Ohshima, K et al. (1990) Hematol Oncol., 8(2):111-8.
 Van Ongen, IIM et al. (2003) Leukemia, 17, 225–2317.

	Diseases	Gene Rearrangement						Translocation			Mutatior				
		<i>IGH</i> V _H – J _H	IGH D _H - J _H	IGK	IGL	ighv Shm	TRB	TRD	TRG	IGH- BCL1 (CCND1)	IGH- BCL2	BCR- ABL1*	PML RARa*	FLT3	NPM1
	Marginal Zone Lymphoma (MZL), extranodal ^{12,13,27}	88%	58%	84%	29%		23%	10%	16%						
	Marginal Zone Lymphoma (MZL), nodal ¹³	100%	30%	80%	30%		10%	20%	10%						
	Mantle Cell Lymphoma (MCL) ^{2,6,7,12,13,27,37}	100%	11%	100%	44%	*	9%	4%	11%	75%					
	Follicular Lymphoma (FL) ^{3,7,12,13,27,28}	84%	19%	84%	21%		6%	5%	2%		90%				
	Diffuse Large B-cell Lymphoma (DLBCL) ^{3,12,13,27}	80%	30%	80%	28%		21%	14%	15%		30%				
	Multiple Myeloma (MM) and other Plasma Cell Neoplasms (PCN) ^{2,9,10,20,25}	84%	60%	57%	97%					20%					
omo	Chronic Lymphocytic Leukemia (CLL) ^{11,12,13,15,23,27,35}	100%	43%	100%	30%	*	25%	12%	18%						
hqm	B-cell Acute Lymphoblastic Leukemia (B-ALL) ^{4,12,14,19,21,22,27,29,30,31,32,33,34}	96%	57%	95%	20%		81%	86%	75%			30%			
	Suspect B-cell Proliferations ^{12,26,27,33}	93%	93%	90%	40%		20%		20%						
hoid	Peripheral T-cell Lymphoma (PTCL) ^{12,13,14,24}	35%	4%		2%		98%		94%						
ymp	T-cell Acute Lymphoblastic Leukemia (T-ALL) ^{12,14,21,22,29,31}	24%	25%	4%			92%	68%	95%						
	Angioimmunoblastic T-cell Lymphoma (AILT) ^{12,13,14}	19%	11%	30%	5%		99%	35%	92%						
	Adult T-cell Leukemia/Lymphoma ³⁹						97%		96%						
	Anaplastic Large-Cell Lymphoma (ALCL) ^{12,13,14}						74%	12%	74%						
	T-cell Prolymphocytic Leukemia (T-PLL) ^{12,13,14}	3%	3%	3%	3%		100%	6%	94%						
	T-cell Large Granular Lymphocytic Leukemia (T-LGL Leukemia) ^{12,13,14}			4%	4%		97%	29%	96%						
	Suspect T-Cell Proliferations ^{12,26,40}	10%		10%			90%	11%	90%						
	Acute Myeloid Leukemia (AML) ^{8,16}													33%	64%
loid	Acute Promyelocytic Leukemia (APL) ^{1,5,16,17}												90%		
Mye	Chronic Myeloid Leukemia (CML) ^{7,16,19,21,38}											87%			
	Myeloproliferative Neoplasms (MPNs) ³⁸											10%			
Note: Perce	The percentage of samples within a given dise ntages indicate the highest referenced value.	ase catego	ory were de	etected us	ing each	gene tar	get.			Recomm Primary	nended Test		Recomm Seconda	ended ry Test	

*Some translocation assays, LymphoTrack®, and MRD applications are currently Research Use Only; not for use in diagnostic procedures. Illumina® and MiSeq™ are registered trademarks of Illumina®, Inc. Ion S5™ and Ion PGM™ are trademarks of Thermo Fisher Scientific®.

..... LymphoTrack[®] Dx

Next-Generation Squencing (NGS) CE-IVD Assays

LymphoTrack Dx Assay kits are designed for the identification of gene rearrangements in hematologic samples utilizing next-generation sequencing (NGS) technologies.

LEARN MORE: catalog.invivoscribe.com

Introduction

These assays take advantage of the wealth and depth of data generated by the Illumina[®] MiSeq[™], Thermo Fisher Scientific[®] Ion PGM[™] and Ion S5[™] platforms.

The Invivoscribe NGS assays offer significant improvements over conventional fragment analysis of B- and T-cell gene rearrangements by providing detailed information regarding the DNA sequences, sequence frequency, and mutational status (*IGHV* Leader and *IGH* FR1 only) of each clonotype.

LymphoTrack Dx Assay kits are a complete solution. Kits contain ready-to-use indexed amplification master mixes, necessary controls, and bioinformatics software. As primers are designed with barcoded indices and adapters, sequencing libraries can be generated with a single PCR, streamlining the overall workflow, eliminating the need for a post-PCR ligation step, and reducing the potential for sample cross contamination.

The per sample testing costs can be reduced by pooling samples from different LymphoTrack Dx Assays into a single sequencing run. The bioinformatics software will sort the complex NGS data for easy analysis and visualization of individual samples.

Detailed instructions for use are provided with all kits and the Invivoscribe technical support team is always available to answer your questions.

Key Benefits

- » One-step PCR for amplicon and library generation
- » Identify and assess mutation status of B- and T-cell gene rearrangements
- » Sequence amplicons from any LymphoTrack kit together
- » Bioinformatics software is available in desktop or Docker images for analysis and interpretation
- » Same reagents for clonality, somatic hypermutation (SHM), minimal residual disease (MRD) testing, and tracking/ monitoring of immunotherapy constructs

Chapter Contents

NGS Assays

- 29 LymphoTrack Dx IGHV Leader Somatic Hypermutation Assay
- 31 LymphoTrack Dx IGH FR1/2/3 Assays
- 35 LymphoTrack Dx IGK Assays
- **37** LymphoTrack Dx *TRG* Assays
- **39** LymphoTrack Dx *TRB* Assay
- 41 LymphoTrack Dx Bioinformatics Software

The LymphoTrack Dx Assays are *in vitro* diagnostic products and are not available for sale or use within North America. For more information regarding the research use only reagents, please see the Next-Generation Sequencing – LymphoTrack RUO section.



LymphoTrack[®] Dx *IGHV* Somatic Hypermutation Assay

Assay Description

The LymphoTrack Dx *IGHV* Leader Somatic Hypermutation Assay for the Illumina[®] MiSeq[™] is an *in vitro* diagnostic product intended for nextgeneration sequencing (NGS) based determination of the frequency distribution of *IGH* gene rearrangements, as well as the degree of somatic hypermutation (SHM) of rearranged genes in patients suspected of having lymphoproliferative disease. This assay aids in the identification of lymphoproliferative disorders as well as providing an aid in determining disease prognosis. If you would like to test for *IGHV* somatic hypermutation using the Thermo Fisher[®] Ion PGM[™] or Ion S5[™] platform, please refer to the LymphoTrack Dx *IGH* FR1 Assay (Catalog #: 91210007).

Summary and Explanation of the Test

The NGS LymphoTrack Dx *IGHV* Leader Somatic Hypermutation Assay for the Illumina® MiSeq[™] represents a significant improvement over clonality assays using fragment analysis as it efficiently detects the majority of *IGH* gene rearrangements using a single multiplex master mix, identifies the DNA sequence specific for each clonal gene rearrangement, and assesses the somatic hypermutation rate of clonal samples in the same workflow.

The master mixes target the Leader (V_{HL}) and the joining (J_{H}) gene regions of the *IGH* locus and are designed with Illumina[®] adapters and indices (8 included in Kit A and 24 included in the Panel). This allows for a one-step PCR and pooling of amplicons from several different samples and targets into a single Illumina[®] MiSeq[™] run. No post-PCR ligation step is required.

Positive (clonal positive, SHM negative), negative (polyclonal) and SHM (clonal positive, SHM positive) DNA controls are included in the kits.

Background

The human immunoglobulin heavy chain (*IGH*) gene locus on chromosome 14 (14q32.3) includes 46–52 functional and 30 nonfunctional variable (V_H) gene segments, 27 functional diversity (D_H) gene segments, and 6 functional joining (J_H) gene segments spread over 1,250 kilobases.

During B-cell development, genes encoding the *IGH* protein are assembled from multiple polymorphic gene segments that undergo rearrangements and selection, generating cell specific V_H - D_H - J_H rearrangements that are unique in both length and sequence.¹

Since leukemias and lymphomas originate from the malignant transformation of individual lymphoid cells, all leukemias and lymphomas generally share one or more cell-specific or "clonal" antigen receptor gene rearrangements. Therefore, tests that detect *IGH* clonal rearrangements can be useful in the study of B- cell malignancies. An additional level of diversity is further generated in the antigen receptors by introducing point mutations in the variable regions, also named SHM. In instances where there is a high degree of SHM, there is the risk that primers located within the variable region will not be able to bind and clonal products will not amplify. In these cases, the leader primers located upstream of the variable region can be beneficial for the detection of clonal products, due to the conserved nature of the VHL region. In addition, the SHM rate of the entire variable gene can be determined using the V_HL primers.

Determining the immunoglobulin variable heavy chain gene (*IGHV*) hypermutation rate is considered a gold standard for determining the prognosis of patients with chronic lymphocytic leukemia (CLL)² and small lymphocytic lymphoma (SLL). In addition, NGS methods can improve disease stratification.

Specimen Requirement

50 ng of high-quality genomic DNA.

References

- 1. Miller et al., Molecular Genetic Pathology (2nd ed.). Springer Science & Business Media. 2013: 302.2.7.13 and 30.2.7.18.
- 2. Ghia et al., *Leukemia* 21: 2-3 (2007).



Simplified representation of the immunoglobulin heavy chain (IGH) gene locus on chromosome 14. Depicted are the variable (V_H) and downstream consensus joining (J_H) region genes involved in rearrangements. Upstream of the variable gene segments, the leader sequence (V_HL) is also depicted. Diversity region genes are not depicted.

LymphoTrack®**Dx**

29 Next-Generation Sequencing (NGS) | CE-IVD Assays

Reagents - MiSeq[™] Detection

Kit A Components

Panel Components (includes all master mixes from Kit A plus the items below)

Master Mix Name	Index #	Master Mix Name	Index #	Master Mix Name	Index #
IGH Leader MiSeq 01	A001	IGH Leader MiSeq 09	A009	IGH Leader MiSeq 18	A018
<i>IGH</i> Leader MiSeq 02	A002	IGH Leader MiSeq 10	A010	IGH Leader MiSeq 19	A019
IGH Leader MiSeq 03	A003	IGH Leader MiSeq 11	A011	IGH Leader MiSeq 20	A020
<i>IGH</i> Leader MiSeq 04	A004	IGH Leader MiSeq 12	A012	IGH Leader MiSeq 21	A021
IGH Leader MiSeq 05	A005	IGH Leader MiSeq 13	A013	IGH Leader MiSeq 22	A022
<i>IGH</i> Leader MiSeq 06	A006	IGH Leader MiSeq 14	A014	IGH Leader MiSeq 23	A023
IGH Leader MiSeq 07	A007	IGH Leader MiSeq 15	A015	IGH Leader MiSeq 25	A025
IGH Leader MiSea 08	A008	IGH Leader MiSeq 16	A016	IGH Leader MiSeq 27	A027

Controls

IGH SHM POS (+) Qty. 1 IGH POS (+) Qty. 1 NGS NEG (-) Qty. 1

Controls

IGH SHM POS (+) Qty. 3 IGH POS (+) Qty. 3 NGS NEG (-) Qty. 3

Rank	Sequence	Length	Merge Count	V-gene	J-gene	% Total reads	Cumulative %	Mutation Rate Partial V-gene (%)	In-frame (Y/N)	No stop codon (Y/N)	V-coverage
1	TTCTCGTGGTG	455	29603	IGHV4-59_08	IGHJ4_02	9.93	9.93	11.26	Y	Υ	98.63
2	CTCGCCCTCCT	463	205	IGHV5-51_01	IGHJ4_02	0.07	9.99	0.00	Y	Y	99.66
3	GGTTTTCCTTG	484	201	IGHV3-7_01	IGHJ4_02	0.07	10.06	7.77	Y	Y	100.00
4	CTCGCCCTCCT	463	185	IGHV5-51_01	IGHJ5_02	0.06	10.12	6.08	Y	Y	99.32
5	CTCGCCCTCCT	469	170	IGHV5-51_01	IGHJ4_02	0.06	10.18	0.00	Y	Y	99.32
6	CTCGCCCTCCT	466	160	IGHV5-51_01	IGHJ4_02	0.05	10.23	0.00	Y	Y	99.66
7	CTGCTGCTGAC	460	159	IGHV2-5_10	IGHJ5_02	0.05	10.29	8.08	Y	Y	97.64
8	GGTTTTCCTTG	493	156	IGHV3-48_02	IGHJ6_02	0.05	10.34	3.72	Y	Y	98.99
9	CTCGCCCTCCT	334	153	IGHV5-51_02	IGHJ2_01	0.05	10.39	3.72	Y	Ν	27.70
10	CTCGCCCTCCT	334	152	IGHV5-51_02	IGHJ2_01	0.05	10.44	3.38	Y	Ν	26.01

Example Data. The top 10 sequences from a read summary generated by the LymphoTrack Dx Software – MiSeq[™] with the SHM mutation rate and predictions pertaining to whether a sequence is in-frame or contains a premature stop codon are depicted. To learn more about the LymphoTrack Dx software, please refer to the LymphoTrack Dx Bioinformatics Software section.

Ordering Information

Catalog #	Products
9-121-0059	LymphoTrack $^{\circ}$ Dx <i>IGHV</i> Leader Somatic Hypermutation Assay Kit A – MiSeq $^{ m imes}$
9-121-0069	LymphoTrack $^{\otimes}$ Dx <i>IGHV</i> Leader Somatic Hypermutation Assay Panel – MiSeq $^{\scriptscriptstyle \mathbb{M}}$
9-500-0009	LymphoTrack® Dx Software – MiSeq™

Quantity Components

8 indices – 5 sequencing reactions each 24 indices – 5 sequencing reactions each 1 Software Package

LymphoTrack[®] Dx *IGH* FR1/2/3 Assays

Assay Description

LymphoTrack Dx IGH FR1 Assays

The LymphoTrack Dx *IGH* FR1 Assay for the Illumina MiSeqTM or Thermo Fisher Scientific[®] Ion S5TM and Ion PGMTM is an *in vitro* diagnostic product intended for next-generation sequencing (NGS) based determination of the frequency distribution of *IGH* gene rearrangements as well as the degree of somatic hypermutation of rearranged genes in patients suspected with having lymphoproliferative disease. This assay aids in the identification of lymphoproliferative disorders as well as providing an aid in determining disease prognosis.

LymphoTrack Dx IGH FR1/2/3 Assays

The LymphoTrack Dx *IGH* FR1 Assay for the Illumina MiSeq[™] or Thermo Fisher Scientific[®] Ion S5[™] and Ion PGM[™] is an *in vitro* diagnostic product intended for next-generation sequencing (NGS) based determination of the frequency distribution of *IGH* gene rearrangements as well as the degree of somatic hypermutation of rearranged genes in patients suspected with having lymphoproliferative disease. This assay aids in the identification of lymphoproliferative disorders as well as providing an aid in determining disease prognosis.

This LymphoTrack Dx *IGH* FR2 Assay is an *in vitro* diagnostic product intended for next-generation sequencing (NGS) for the Illumina MiSeqTM or Thermo Fisher Scientific[®] Ion S5TM and Ion PGMTM instruments. The assay will determine the frequency distribution of *IGH* V_H-J_H gene rearrangements in patients suspected with having lymphoproliferative disease. This assay aids in the identification of lymphoproliferative disorders.

The LymphoTrack Dx *IGH* FR3 Assay is an *in vitro* diagnostic product intended for next-generation sequencing (NGS) for the Illumina MiSeqTM or Thermo Fisher Scientific[®] Ion S5TM and Ion PGMTM instruments. The assay will determine the frequency distribution of *IGH* V_H-J_H gene rearrangements in patients suspected with having lymphoproliferative disease. This assay aids in the identification of lymphoproliferative disorders.

Summary and Explanation of the Test

The LymphoTrack Dx *IGH* Assays represent a significant improvement over conventional clonality assessment methods utilizing fragment analysis by providing four important and complementary uses in a single workflow:

- 1. Detection of clonal populations.
- 2. Identification of sequence information and gene segment utilization.
- 3. The LymphoTrack Dx IGH framework 1 FR1 master mixes provide the degree of SHM in the immunoglobulin variable heavy chain (IGHV) gene locus.
- 4. The ability to track sequences in subsequent samples with the Invivoscribe LymphoTrack MRD* Software.

These assays utilize a single multiplex master mix to target each conserved *IGH* Framework Region (FR1, FR2, and FR3) within the V_{μ} and the J_{μ} regions described in lymphoid malignancies. Each master mix targets one of the conserved *IGH* framework regions (FR1, FR2, or FR3) within the V_{μ} and the J_{μ} regions described in lymphoid malignancies. Targeting all three framework regions significantly reduces the risk of not being able to detect the presence of clonality, as somatic hypermutations in the primer binding sites of the involved V_{μ} gene segments can impede DNA amplification.¹ The included primers are designed with Illumina® or Thermo Fisher Scientific adapters and indices (8-24 and 12, respectively). This allows up to 24 samples on MiSeqTM and 12 samples on lon PGM[™] and lon S5[™] to be sequenced at the same time with any of the individual FRs.

In addition, amplicons generated with different FR master mixes or Invivoscribe LymphoTrack Dx kits (such as *IGK* or *TRG*) can be pooled together in the same sequencing library to reduce testing costs. Positive clonal (SHM negative) and negative polyclonal DNA controls are included in kits. A clonal SHM positive control can be purchased separately (cat#: 40880008).

Background

The human immunoglobulin heavy chain (*IGH*) gene locus on chromosome 14 (14q32.3) includes 46-52 functional and 30 non-functional variable (V_{μ}) , 27 functional diversity (D_{μ}) , and 6 functional joining (J_{μ}) gene segments. The VH gene segments can be further broken down into three conserved frameworks (FR) and three variable complementarity-determining regions (CDRs).

During development of lymphoid cells, antigen receptor genes undergo somatic gene rearrangements.² Specifically during B-cell development, *IGH* molecules are assembled from multiple polymorphic gene segments that undergo rearrangements generating $V_{H}-D_{H}-J_{H}$ combinations unique in both length and sequence.³ Since leukemias and lymphomas originate from the malignant transformation of individual lymphoid cells, all leukemias and lymphomas generally share one or more cell-specific or "clonal" antigen receptor gene rearrangements.

In addition, the *IGHV* hypermutation status obtained with the LymphoTrack Dx *IGH* FR1 master mixes, provides important prognostic information for patients with chronic lymphocytic leukemia (CLL) and small lymphocytic lymphoma (SLL). The SHM rate has been shown to have clinical relevance for CLL, as there is a clear distinction in the median survival of patients with and without SHM.⁴

*MRD Software can be used to track sequences generated by either LymphoTrack Assays – MiSeq[™] or Ion S5/PGM™. MRD applications are for Research Use Only.



31 Next-Generation Sequencing (NGS) | CE-IVD Assays

Specimen Requirement

50 ng of high-quality genomic DNA.

References

- 1. Evans, P. A. et al., (2007). Leukemia 21, 207-14.
- 2. Tonegawa, S. (1983). Nature 302, 575-581.
- 3. Miller JE. (2013) Molecular Genetic Pathology (2nd Edition., sections 30.2.7.13 and 30.2.7.18).
- 4. Ghia et al., *Blood* 105:1678-1685 (2005).

Simplified Representation of the IGH Gene



Simplified depiction of variable (V_H) and downstream consensus joining (J_H) region genes involved in gene rearrangements.

LymphoTrack Dx IGH FR3 Components

Reagents - MiSeq[™] Detection

The LymphoTrack Dx IGH FR1/2/3 Assays contain components from respective individual FR kit A's or panels.

LymphoTrack Dx *IGH* FR1 Components

LymphoTrack Dx IGH FR2 Components

Master Mix Name	Index		Master Mix Name	Index		Master Mix Name	Index	
IGH FR1 MiSeq 01	A001		IGH FR2 MiSeq 01	A001		IGH FR3 MiSeq 01	A001	
IGH FR1 MiSeq 02	A002		IGH FR2 MiSeq 02	A002		IGH FR3 MiSeq 02	A002	
IGH FR1 MiSeq 03	A003		IGH FR2 MiSeq 03	A003		IGH FR3 MiSeq 03	A003	
IGH FR1 MiSeq 04	A004 🧹		IGH FR2 MiSeq 04	A004 🧹		IGH FR3 MiSeq 04	A004 <	٢
IGH FR1 MiSeq 05	A005 呈		IGH FR2 MiSeq 05	A005		IGH FR3 MiSeq 05	A005	2
IGH FR1 MiSeq 06	A006		IGH FR2 MiSeq 06	A006		IGH FR3 MiSeq 06	A006	
IGH FR1 MiSeq 07	A007		IGH FR2 MiSeq 07	A007		IGH FR3 MiSeq 07	A007	
	A008			A008			A008	
IGH FR1 MiSeq 09	A009		IGH FR2 MiSeq 09	A009		IGH FR3 MiSeq 09	A009	
IGH FR1 MiSeq 10	A010		IGH FR2 MiSeq 10	A010		IGH FR3 MiSeq 10	A010	
IGH FR1 MiSeq 11	A011		IGH FR2 MiSeq 11	A011		IGH FR3 MiSeq 11	A011	
IGH FR1 MiSeq 12	A012	Ę	IGH FR2 MiSeq 12	A012	AEL	IGH FR3 MiSeq 12	A012	VEL
IGH FR1 MiSeq 13	A013	PAI	IGH FR2 MiSeq 13	A013	PAI	IGH FR3 MiSeq 13	A013	PAI
IGH FR1 MiSeq 14	A014		IGH FR2 MiSeq 14	A014		IGH FR3 MiSeq 14	A014	
IGH FR1 MiSeq 15	A015		IGH FR2 MiSeq 15	A015		IGH FR3 MiSeq 15	A015	
IGH FR1 MiSeq 16	A016		IGH FR2 MiSeq 16	A016		IGH FR3 MiSeq 16	A016	
IGH FR1 MiSeq 18	A018		IGH FR2 MiSeq 18	A018		IGH FR3 MiSeq 18	A018	
IGH FR1 MiSeq 19	A019		IGH FR2 MiSeq 19	A019		IGH FR3 MiSeq 19	A019	
IGH FR1 MiSeq 20	A020		IGH FR2 MiSeq 20	A020		IGH FR3 MiSeq 20	A020	
IGH FR1 MiSeq 21	A021		IGH FR2 MiSeq 21	A021		IGH FR3 MiSeq 21	A021	
IGH FR1 MiSeq 22	A022		IGH FR2 MiSeq 22	A022		IGH FR3 MiSeq 22	A022	
IGH FR1 MiSeq 23	A023		IGH FR2 MiSeq 23	A023		IGH FR3 MiSeq 23	A023	
IGH FR1 MiSeq 25	A025		IGH FR2 MiSeq 25	A025		IGH FR3 MiSeq 25	A025	
IGH FR1 MiSeq 27	A027		IGH FR2 MiSeq 27	A027		IGH FR3 MiSeq 27	A027	

Kit A's contain indices IGH FRX A001 to A008. Panels contain all master mixes listed above.

Controls in Individual
FR (1, 2, or 3) Kit A'sControls in Individual
FR (1, 2, or 3) PanelsControls in Combo
FR 1/2/3 Kit AIGH POS (+) Qty. 1IGH POS (+) Qty 3IGH POS (+) Qty. 2NGS NEG (-) Qty. 1NGS NEG (-) Qty 3NGS NEG (-) Qty. 2

*MRD Software can be used to track sequences generated by either LymphoTrack Assays – $MiSeq^{TM}$ or Ion S5/PGM™. MRD applications are for Research Use Only.

Controls in Combo

FR 1/2/3 Panel

IGH POS (+) Qty. 6 NGS NEG (-) Qty. 6

Reagents - Ion S5[™] / PGM[™] Detection

The LymphoTrack Dx IGH FR1/2/3 Assays contain components from respective individual FR Assays.

LymphoTrack Dx IGH FR2 Components LymphoTrack Dx IGH FR3 Component LymphoTrack Dx IGH FR1 Components **Master Mix Name** Index # **Master Mix Name** Index # **Master Mix Name** Index # IGH FR1 Ion S5/PGM 01 lonXpress_001 IGH FR2 Ion S5/PGM 01 lonXpress_001 IGH FR3 Ion S5/PGM 01 lonXpress_001 lonXpress_002 IGH FR2 Ion S5/PGM 02 lonXpress_002 IGH FR1 Ion S5/PGM 02 IGH FR3 Ion S5/PGM 02 lonXpress_002 IGH FR1 Ion S5/PGM 03 lonXpress_003 IGH FR2 Ion S5/PGM 03 lonXpress_003 IGH FR3 Ion S5/PGM 03 lonXpress_003 IGH FR1 Ion S5/PGM 04 lonXpress_004 IGH FR2 Ion S5/PGM 04 lonXpress_004 IGH FR3 Ion S5/PGM 04 lonXpress_004 IGH FR1 Ion S5/PGM 07 lonXpress_007 IGH FR2 Ion S5/PGM 07 lonXpress_007 IGH FR3 Ion S5/PGM 07 lonXpress_007 IGH FR2 Ion S5/PGM 08 IGH FR1 Ion S5/PGM 08 lonXpress_008 lonXpress_008 IGH FR3 Ion S5/PGM 08 lonXpress_008 IGH FR1 Ion S5/PGM 09 lonXpress_009 IGH FR2 Ion S5/PGM 09 lonXpress_009 IGH FR3 Ion S5/PGM 09 lonXpress_009 IGH FR1 Ion S5/PGM 10 lonXpress_010 IGH FR2 Ion S5/PGM 10 lonXpress_010 IGH FR3 Ion S5/PGM 10 lonXpress_010 IGH FR1 Ion S5/PGM 11 lonXpress_011 IGH FR2 Ion S5/PGM 11 lonXpress_011 IGH FR3 Ion S5/PGM 11 lonXpress_011 IGH FR1 Ion S5/PGM 12 lonXpress_012 IGH FR2 Ion S5/PGM 12 lonXpress_012 IGH FR3 Ion S5/PGM 12 lonXpress_012 IGH FR1 Ion S5/PGM 13 lonXpress_013 IGH FR2 Ion S5/PGM 13 lonXpress_013 IGH FR3 Ion S5/PGM 13 lonXpress_013 IGH FR1 Ion S5/PGM 14 lonXpress_014 lonXpress_014 IGH FR2 Ion S5/PGM 14 IGH FR3 Ion S5/PGM 14 lonXpress_014 Controls in Individual FR (1,2, or 3) Kits

IGH POS (+) Qty. 2 NGS NEG (-) Qty. 2

Controls in FR 1/2/3 Kit

IGH POS (+) Qty. 4 NGS NEG (-) Qty. 4



V-I Sequence Frequencies. The LymphoTrack Dx Software provides a stacked bar graph depicting the relative frequencies of the 200 most prevalent $V_{\rm H}$ -J_{\rm H} rearrangements identified in a sample. To learn more about the LymphoTrack Dx software, please refer to the LymphoTrack Dx Bioinformatics Software section.



Comparison of SHM Analysis Methods. The SHM rate of 51 CLL samples was determined by the LymphoTrack Dx IGH FR1 Assay – MiSeq™ and analyzed with both the LymphoTrack Dx Software - MiSeq™ and IMGT analysis.



Ordering Information for Reagents - MiSeq[™] Detection

Catalog #	Products
9-121-0129	LymphoTrack $^{\circ}$ Dx <i>IGH</i> FR1/2/3 Assay Kit A - MiSeq TM
9-121-0139	LymphoTrack $^{\circ}$ Dx <i>IGH</i> FR1/2/3 Assay Panel – MiSeq TM
9-121-0009	LymphoTrack® Dx <i>IGH</i> FR1 Assay Kit A – MiSeq™
9-121-0039	LymphoTrack® Dx <i>IGH</i> FR1 Assay Panel – MiSeq™
9-121-0089	LymphoTrack [®] Dx <i>IGH</i> FR2 Assay Kit A − MiSeq [™]
9-121-0099	LymphoTrack® Dx <i>IGH</i> FR2 Assay Panel - MiSeq™
9-121-0109	LymphoTrack® Dx <i>IGH</i> FR3 Assay Kit A − MiSeq™
9-121-0119	LymphoTrack® Dx <i>IGH</i> FR3 Assay Panel - MiSeq™
9-500-0009	LymphoTrack® Dx Software – MiSeq™

Quantity Components

8 + 8 + 8 indices - 5 sequencing reactions each
24 + 24 + 24 indices - 5 sequencing reactions each
8 indices - 5 sequencing reactions each
24 indices - 5 sequencing reactions each
24 indices - 5 sequencing reactions each
8 indices - 5 sequencing reactions each
24 indices - 5 sequencing reactions each
24 indices - 5 sequencing reactions each
1 Software Package

Ordering Information for Reagents – Ion S5[™] / PGM[™] Detection

Culuidy #	
9-121-0057	
9-121-0007	
9-121-0037	
9-121-0047	
9-500-0007	

Catalog #

Products

LymphoTrack® Dx I/GH FR1/2/3 Assay – S5/PGMTM LymphoTrack® Dx IGH FR1 Assay – S5/PGMTM LymphoTrack® Dx IGH FR2 Assay – S5/PGMTM LymphoTrack® Dx IGH FR3 Assay – S5/PGMTM LymphoTrack® Dx Software – S5/PGMTM

Quantity Components

12 + 12 + 12 indices - 5 sequencing reactions each 1 Software Package

CEIVD These products are CE-IVD assays for *in vitro* diagnostic use and available in regions that accept CE-IVD products.



LymphoTrack[®] Dx *IGK* Assay

Assay Description

The LymphoTrack Dx *IGK* Assays for the Illumina[®] MiSeq[®], or Thermo Fisher Scientific[®] Ion PGM[®] and Ion S5[®] instruments are *in vitro* diagnostic products intended for next-generation sequencing (NGS) based determination of the frequency distribution of *IGK* gene rearrangements in patients suspected of having lymphoproliferative disease. These assays aid in the identification of lymphoproliferative disorders.

Summary and Explanation of the Test

In contrast to the IdentiClone[®] fragment analysis assays for clonality that utilize two master mixes, these NGS assays contain a single multiplex master mix to target conserved regions of the *IGK* gene locus described in lymphoid malignancies, thereby reducing sample DNA requirements and simplifying the testing workflow. The LymphoTrack Dx *IGK* master mix primers are also designed with Illumina[®] or Thermo Fisher Scientific adapters and up to 24 different indices. This allows amplicons generated from different indexed *IGK* master mixes to be pooled into a single library for loading onto one MiSeq[®] flow cell, Ion PGM[®] or Ion S5[®] chips.

Positive clonal and negative polyclonal DNA controls are included in kits.

LymphoTrack Dx Software is capable of sorting complex NGS data generated by LymphoTrack Dx Assays by gene target, providing users the ability to reduce per sample testing costs by sequencing amplicons from any LymphoTrack Dx Assay (*e.g. IGH, IGK, TRB, TRG*) at the same time. In addition, the LymphoTrack Dx Software provides an easy and streamlined method for visualization of data and guidelines provided in the instructions for use allow samples to be interpreted for evidence or no evidence of clonality.

Background

The LymphoTrack Dx *IGK* Assays represent a significant improvement over existing fragment analysis clonality assays by providing two important and complementary uses:

1. Detection of clonal populations.

- 2. Identification of sequence information and gene segment utilization.
- 3. Ability to track sequences in subsequent samples with the use of the Invivoscribe LymphoTrack MRD Software.

The human immunoglobulin kappa (*IGK*) gene locus on chromosome 2 (2p11.2) includes 76 V (variable) region genes spanning 7 subgroups and 5 J (joining) region gene segments upstream of the (Ck) region.

During lymphoid cell development, antigen receptor genes undergo somatic gene rearrangements.¹ Specifically, during B-cell development, genes encoding *IGK* molecules are assembled from multiple polymorphic gene segments that generate V-J combinations unique in both length and sequence.²

In addition, the kappa deleting element (K_{de}), approximately 24 kb downstream of the Jk-Ck region can also rearrange with Vk gene segments and the isolated recombination signal sequence in the Jk-Ck intron (Jk-Ck INTR).

Since leukemias and lymphomas originate from the malignant transformation of individual lymphoid cells, all leukemias and lymphomas generally share one or more cell-specific or "clonal" antigen receptor gene rearrangements. Therefore, tests that detect *IGK* clonal rearrangements can be useful in the study of B- cell malignancies and complement *IGH* testing, as the *IGK* receptor is less susceptible to somatic mutations.

Specimen Requirement

50 ng of high-quality genomic DNA.

References

- 1. Tonegawa, S. (1983). Nature 302, 575-581.
- 2. Miller JE. (2013) Molecular Genetic Pathology (2nd Edition., sections 30.2.7.13 and 30.2.7.18).

Simplified Representation of the IGK Gene



Depicted are the variable region (Vk) or variable intragenic J $_{k-Ck}$ intron (J $_{k-Ck}$ INTR) genes involved in *IGK* gene rearrangements in addition to the downstream consensus joining region genes (J $_{k}$) or kappa deleting element (K $_{de}$).

35 Next-Generation Sequencing (NGS) | CE-IVD Assays

Reagents - MiSeq[™] Detection

Kit A Components		Panel Components	Panel Components (includes all master mixes from Kit A plus the items below)						
Master Mix Name	Index #	Master Mix Name	Index #	Master Mix Name	Index #				
IGK MiSeq 01	A001	IGK MiSeq 09	A009	IGK MiSeq 18	A018				
IGK MiSeq 02	A002	IGK MiSeq 10	A010	IGK MiSeq 19	A019				
IGK MiSeq 03	A003	IGK MiSeq 11	A011	IGK MiSeq 20	A020				
IGK MiSeq 04	A004	IGK MiSeq 12	A012	IGK MiSeq 21	A021				
IGK MiSeq 05	A005	IGK MiSeq 13	A013	IGK MiSeq 22	A022				
IGK MiSeq 06	A006	IGK MiSeq 14	A014	IGK MiSeq 23	A023				
IGK MiSeq 07	A007	IGK MiSeq 15	A015	IGK MiSeq 25	A025				
IGK MiSeq 08	A008	IGK MiSeq 16	A016	IGK MiSeq 27	A027				
Controls		Controls							
IGK POS (+) Qty. 1	NGS NEG (-) Qty. 1	IGK POS (+) Qty. 3		NGS NEG (-) Qty. 3					

Reagents – Ion S5[™]/PGM[™] Detection

Assay Components

Master Mix Name	Index #	Master Mix Name	Index #
IGK S5/PGM 01	lonXpress_001	IGK \$5/PGM 11	lonXpress_011
IGK \$5/PGM 02	lonXpress_002	IGK S5/PGM 12	lonXpress_012
IGK \$5/PGM 04	lonXpress_004	<i>IGK</i> S5/PGM 13	lonXpress_013
IGK \$5/PGM 08	lonXpress_008	<i>IGK</i> S5/PGM 14	lonXpress_014
IGK \$5/PGM 09	lonXpress_009	IGK S5/PGM 16	lonXpress_016
IGK \$5/PGM 010	lonXpress_010	IGK \$5/PGM 17	lonXpress_017

Controls

IGK POS (+) Qty. 2



V-J Sequence Frequencies. The LymphoTrack Dx Software provides a stacked bar graph depicting the relative frequencies for the most prevalent rearrangements identified in a sample. To learn more about the LymphoTrack Dx software, please refer to

NGS NEG (-) Qty. 2

the LymphoTrack Dx Bioinformatics Software section.

Ordering Information

Products
LymphoTrack® Dx <i>IGK</i> Assay Kit A – MiSeq™
LymphoTrack® Dx <i>IGK</i> Assay Panel - MiSeq™
LymphoTrack® Dx Software – MiSeq™
LymphoTrack® Dx <i>IGK</i> Assay - S5/PGM™
LymphoTrack® Dx Software – S5/PGM™

Quantity Components

8 indices - 5 sequencing reactions each 24 indices - 5 sequencing reactions each 1 Software Package 12 indices - 5 sequencing reactions each 1 Software Package

CEIVD These products are CE-IVD assays for *in vitro* diagnostic use and available in regions that accept CE-IVD products.

Next-Generation Sequencing (NGS) | CE-IVD Assays 36
LymphoTrack[®] Dx *TRG* Assay

Assay Description

The LymphoTrack Dx TRG Assays for the Illumina[®] MiSeq[®], Thermo Fisher Scientific[®] Ion PGM[®] or Ion S5[®] instruments are *in vitro* diagnostic products intended for next-generation sequencing (NGS) based determination of the frequency distribution of TRG gene rearrangements in patients suspected with having lymphoproliferative disease. These assays aid in the identification of lymphoproliferative disorders.

Summary and Explanation of the Test

This assay utilizes a single multiplex master mix to target conserved V and J regions of the human *TRG* gene locus described in lymphoid malignancies, thereby reducing sample DNA requirements and simplifying the testing workflow. Primers are designed with Illumina® or Thermo Fisher Scientific adapters and up to 24 different indices; thereby allowing amplicons generated from different *TRG* master mixes to be pooled together for sequencing on a single MiSeq[™] flow cell, or Ion PGM[™] or Ion S5[™] chip.

Positive clonal and negative polyclonal DNA controls are included in kits.

LymphoTrack Dx Software sorts complex NGS data generated by LymphoTrack Dx Assays by gene target,, providing users the ability to reduce per sample testing costs by sequencing amplicons generated with any LymphoTrack Dx Assay (*e.g. IGH, IGK, TRB, TRG*) at the same time. In addition, the LymphoTrack Dx Software provides an easy and streamlined method for data visualization and guidelines provided in the instructions for use allow samples to be interpreted for evidence or no evidence of clonality.

Background

The LymphoTrack Dx TRG Assays represent a significant improvement over existing fragment analysis clonality assays by providing two important and complementary uses:

- 1. Detection of initial clonal populations.
- 2. Identification of sequence information required to track clonal rearrangements in subsequent samples.

The TRG gene locus on chromosome 7 (7q14) includes 14 V (variable region) genes (Group I, II, III, and IV), 5 J (joining region) gene segments, and 2 C (constant region) genes spread over 200 kilobases.

During development of lymphoid cells, antigen receptor genes undergo somatic gene rearrangements.¹ Specifically during T-cell development, genes encoding *TRG* molecules are assembled from multiple polymorphic gene segments that undergo rearrangement generating V-J combinations unique in both length and sequence.² Since leukemias and lymphomas originate from the malignant transformation of individual lymphoid cells, all leukemias and lymphomas generally share one or more cell-specific or "clonal" antigen receptor gene rearrangement. Therefore, tests that detect *TRG* clonal rearrangements can be useful in the study of B- and T-cell malignancies.

Note: For a more thorough explanation of the locus and the targeted deep sequencing strategy, please refer to Principle of Immunoglobulin and T-Cell Receptor Gene Rearrangement.²

Specimen Requirement

50 ng of high-quality genomic DNA.

References

Tonegawa, S. (1983). Nature 302, 575-581.
 Miller JE. (2013) Molecular Genetic Pathology (2nd Edition., sections 30.2.7.13 and 30.2.7.18).

Simplified Representation of the TRG



Depicted are the variable region (V) genes and downstream consensus joining region genes (J) that are involved in TRG gene rearrangements.

LymphoTrack®Dx
 Next-Generation Sequencing (NGS) | CE-IVD Assays

Reagents - MiSeq[™] Detection

Kit A Components		Panel Components (includes all master mixes from Kit A plus the items below)					
Master Mix Name	Index #	Master Mix Name	Index #	Master Mix Name	Index #		
TRG MiSeq 01	A001	TRG MiSeq 09	A009	TRG MiSeq 18	A018		
TRG MiSeq 02	A002	TRG MiSeq 10	A010	TRG MiSeq 19	A019		
TRG MiSeq 03	A003	TRG MiSeq 11	A011	TRG MiSeq 20	A020		
TRG MiSeq04	A004	TRG MiSeq 12	A012	TRG MiSeq 21	A021		
TRG MiSeq 05	A005	TRG MiSeq 13	A013	TRG MiSeq 22	A022		
TRG MiSeq 06	A006	TRG MiSeq 14	A014	TRG MiSeq 23	A023		
TRG MiSeq 07	A007	TRG MiSeq 15	A015	TRG MiSeq 25	A025		
TRG MiSeq 08	A008	TRG MiSeq 16	A016	TRG MiSeq 27	A027		
Controls <i>TRG</i> POS (+) Qty. 1	NGS NEG (-) Qty. 1	Controls TRG POS (+) Qty. 3		NGS NEG (-) Qty. 3			

Reagents - Ion S5[™]/PGM[™] Detection

Master Mix Name	Index #	Master Mix Name	Index #	Master Mix Name	Index #
TRG S5/PGM 01	lonXpress_001	TRG S5/PGM 07	lonXpress_007	TRG S5/PGM 11	lonXpress_011
TRG \$5/PGM 02	lonXpress_002	TRG S5/PGM 08	lonXpress_008	TRG S5/PGM 12	lonXpress_012
TRG \$5/PGM 03	lonXpress_003	TRG \$5/PGM 09	lonXpress_009	TRG S5/PGM 13	lonXpress_013
TRG \$5/PGM 04	lonXpress_004	TRG \$5/PGM 10	lonXpress_010	TRG S5/PGM 14	lonXpress_014
<i>TRG</i> S5/PGM 04	lonXpress_004	<i>TRG</i> \$5/PGM 10	lonXpress_010	<i>TRG</i> \$5/PGM 14	lonXpress_014

NGS NEG (-) Qty. 2

Controls

TRG POS (+) Qty. 2



V-J Sequence Frequencies. The LymphoTrack Dx bioinformatics software provides PDF reports which include Top 10 Merged Read Summary as well as a stacked bar graph depicting the relative frequencies of the V-J rearrangements identified in the sample. To learn more about the LymphoTrack Dx software, please refer to the LymphoTrack Dx Bioinformatics Software section.

Ordering Information

Catalog #	Products
9-227-0019	LymphoTrack® Dx <i>TRG</i> Assay Kit A – MiSeq™
9-227-0009	LymphoTrack® Dx <i>TRG</i> Assay Panel – MiSeq™
9-500-0009	LymphoTrack® Dx Software – MiSeq™
9-227-0007	LymphoTrack [®] Dx <i>TRG</i> Assay - S5/PGM [™]
9-500-0007	LymphoTrack [®] Dx Software – S5/PGM [™]

Quantity Components

8 indices - 5 sequencing reactions each 24 indices - 5 sequencing reactions each 1 Software Package 12 indices - 5 sequencing reactions each 1 Software Package

CEIVD These products are CE-IVD assays for *in vitro* diagnostic use and available in regions that accept CE-IVD products.

Next-Generation Sequencing (NGS) | CE-IVD Assays 38

LymphoTrack[®] Dx *TRB* Assay

Assay Description

The LymphoTrack[®] Dx *TRB* Assay for the Illumina MiSeq[™] is an *in vitro* diagnostic product intended for next-generation sequencing (NGS) based determination of the frequency distribution of *TRB* gene rearrangements in patients suspected of having lymphoproliferative disease. This assay aids in the identification of lymphoproliferative disorders.

Summary and Explanation of the Test

This assay utilizes a single multiplex master mix to target conserved V and J regions of the human *TRB* gene locus described in lymphoid malignancies, thereby reducing sample DNA requirements and simplifying the testing workflow. The LymphoTrack Dx *TRB* master mix primers are designed with Illumina® adapters and 8 (Kit A) or 24 (Panel) different indices. This allows amplicons generated from different indexed *TRB* master mixes to be pooled into a single library for loading onto one MiSeq[™] flow cell.

Positive clonal and negative polyclonal DNA controls are included in kits.

LymphoTrack Dx Software is capable of sorting complex NGS data generated by LymphoTrack Dx Assays – MiSeq[™] by gene target, offering a second layer of multiplexing. This provides users the ability to reduce per sample testing costs by sequencing amplicons from any LymphoTrack Dx Assay from the MiSeq instrument (e.g. *TRB*, *TRG*, *IGH*, *IGK*) at the same time. In addition, the LymphoTrack Dx Software provides an easy and streamlined method for visualization of data. Guidelines to interpret samples for evidence or no evidence of clonality are included in the instructions for use provided with each kit.

Background

The LymphoTrack Dx TRB Assays represent a significant improvement over fragment analysis methods for clonality testing by providing two important and complementary uses:

- 1. Detection of initial clonal populations.
- 2. Identification of sequence information required to track clonal rearrangements in subsequent samples.

Analysis of the rearranged *TRB* locus increases the probability of identifying T cell receptor gene rearrangements, as compared to testing for *TRG* gene rearrangements only. As a result, combining the analysis of *TRB* and *TRG* loci increases the sensitivity of clonality detection.

The human T-cell receptor beta (*TRB*) gene locus on chromosome 7 (7q34) includes 65 V β (variable) gene segments, followed by two separate clusters of genes each containing a D β (diversity) gene, several J β (joining) genes, and a C β (constant) region spread over 685 kilobases. The 2 C β genes, *TRBC1* and *TRBC2*, encode highly homologous products with no functional difference.

During lymphoid cell development, antigen receptor genes undergo somatic gene rearrangements.¹ Specifically, during T-cell development genes encoding *TRB* molecules are assembled from multiple polymorphic gene segments that generate $V\beta$ - $D\beta$ - $J\beta$ combinations unique in both length and sequence.²

Since leukemias and lymphomas originate from the malignant transformation of individual lymphoid cells, all leukemias and lymphomas generally share one or more cell-specific or "clonal" antigen receptor gene rearrangements. Therefore, tests that detect *TRB* clonal rearrangements can be useful in the study of B- and T-cell malignancies.

Note: For a more thorough explanation of the locus and the targeted deep sequencing strategy, please refer to Principle of Immunoglobulin and T-Cell Receptor Gene Rearrangement.²

Specimen Requirement

50 ng of high-quality genomic DNA.

References

- 1. Tonegawa, S. (1983). *Nature* 302, 575-581.
- 2. Miller JE. (2013) Molecular Genetic Pathology (2nd Edition., sections 30.2.7.13 and 30.2.7.18).
- 3. JE Miller et al., Molecular Genetic Pathology (2nd ed.). Springer Science & Business Media. 2013: 30.2.7.13.



Simplified Representation of the TRB Gene

Dβ1

Vβ

Depicted are the variable (VB), diversity (DB), and joining (JB) gene regions involved in TRB gene rearrangements, in addition to the downstream consensus (CB) gene regions.

Cβ1

Dβ2

Cβ2

Reagents - MiSeq[™] Detection

Kit A Components		Panel Components (includes all master mixes from Kit A plus the items below)					
Master Mix Name	Index #	Master Mix Name	Index #	Master Mix Name	Index #		
TRB MiSeq 01	A001	TRB MiSeq 09	A009	TRB MiSeq 18	A018		
TRB MiSeq 02	A002	TRB MiSeq 10	A010	TRB MiSeq 19	A019		
TRB MiSeq 03	A003	TRB MiSeq 11	A011	TRB MiSeq 20	A020		
TRB MiSeq 04	A004	TRB MiSeq 12	A012	TRB MiSeq 21	A021		
TRB MiSeq 05	A005	TRB MiSeq 13	A013	TRB MiSeq 22	A022		
TRB MiSeq 06	A006	TRB MiSeq 14	A014	TRB MiSeq 23	A023		
TRB MiSeq 07	A007	TRB MiSeq 15	A015	TRB MiSeq 25	A025		
TRB MiSeq 08	A008	TRB MiSeq 16	A016	TRB MiSeq 27	A027		
Controls		Controls					
TRB POS (+) Qty. 1	NGS NEG (-) Qty. 1	TRB POS (+) Qty. 3		NGS NEG (-) Qty. 3			



V-J Sequence Frequencies. The LymphoTrack Dx bioinformatics software provides a stacked bar graph depicting the relative frequencies for the 200 most prevalent rearrangements sequenced and identified in the sample. To learn more about the LymphoTrack Dx software, please refer to the LymphoTrack Dx Bioinformatics Software section.

Ordering Information

Cata	log	#
------	-----	---

9-225-0009 9-225-0019 9-500-0009

Products

LymphoTrack[®] Dx *TRB* Assay Kit A - MiSeq[™] LymphoTrack[®] Dx *TRB* Assay Panel - MiSeq[™] LymphoTrack[®] Dx Software - MiSeq[™]

Quantity Components

8 indices - 5 sequencing reactions each 24 indices - 5 sequencing reactions each 1 Software Package

CEIVD These products are CE-IVD assays for *in vitro* diagnostic use and available in regions that accept CE-IVD products.



LymphoTrack[®] Dx Bioinformatics Software

Software Use

The LymphoTrack Dx Bioinformatics Software package is designed to analyze raw FASTQ files generated by LymphoTrack Dx Assays for clonality analysis of single or multiple target data sets (*IGHV* Leader, *IGH* FR1, *IGH* FR2, *IGH* FR3, *IGK*, *TRB*, or *TRG*). For data generated with LymphoTrack Dx *IGHV* Leader or *IGH* FR1 Assays, the software provides additional information, including the rate of somatic hypermutation and whether a clone would be functional based upon the presence of a premature stop codon. The software can also predict whether an open reading frame would be in- or out-of-frame.

The LymphoTrack Dx software is composed of two distinct parts:

- 1. A bioinformatics Data Analysis Application
- 2. Microsoft Excel® Data Visualization spreadsheets and automated Sample-to-PDF Reports for streamlined data analysis for VJ usage and VJ sequence frequency graphs.

j_ymphoTrack® Dx MiSeq Data Analysis version 2.4.3 Copyright 2018 Invivoscribe Inc.	
LymphoTrack [®] Software	
Choose analysis mode(s): Leader 77.6 77.6 IGH FR1 IGH FR2 IGH FR Choose decimal mark: Choose decimal mark:	23 📄 IGK
*," decimal comma, e.g., "1.000,00" Choose an Input Folder for analysis:	Browse Launch Program

Sequence Frequency Graph



The stacked bar graph depicts the top 200 sequencing reads for a sample. Each individual colored bar represents a unique population of cells. Different colors stacked at the same point on the x-axis represent unique sequences that utilize the same V and J gene families. The amplicons of these products vary in sequence and may also vary in product size.

LymphoTrack[®]**Dx**

41 Next-Generation Sequencing (NGS) | CE-IVD Assays

Sequencing Summary

Using the merged read summary, along with the easy to follow flow charts in the LymphoTrack Dx Assay instructions for use (IFU), interpretation is quick and easy.

5	Sample Name										
Total reads = 32,458											
	2			3		4					
Rank	Sequence	Length	Merge count	V-gene	J-gene	% Total reads	Cumulative %	Mutation rate partial V-gene (%)	In-frame (Y/N)	No stop codon (Y/N)	V-coverage
1	TTCTCGTGGTG	455	29603	IGHV4-59_08	IGHJ4_02	9.93	9.93	11.26	Y	Y	98.63
2	CTCGCCCTCCT	463	205	IGHV5-51_01	IGHJ4_02	0.07	9.99	0.00	Y	Y	99.66
3	GGTTTTCCTTG	484	201	IGHV3-7_01	IGHJ4_02	0.07	10.06	7.77	Y	Y	100.00
4	CTCGCCCTCCT	463	185	IGHV5-51_01	IGHJ5_02	0.06	10.12	6.08	Y	Y	99.32
5	CTCGCCCTCCT	469	170	IGHV5-51_01	IGHJ4_02	0.06	10.18	0.00	Y	Y	99.32
6	CTCGCCCTCCT	466	160	IGHV5-51_01	IGHJ4_02	0.05	10.23	0.00	Y	Y	99.66
7	CTGCTGCTGAC	460	159	IGHV2-5_10	IGHJ5_02	0.05	10.29	8.08	Y	Y	97.64
8	GGTTTTCCTTG	493	156	IGHV3-48_02	IGHJ6_02	0.05	10.34	3.72	Y	Y	98.99
9	CTCGCCCTCCT	334	153	IGHV5-51_02	IGHJ2_01	0.05	10.39	3.72	Y	Ν	27.70
10	CTCGCCCTCCT	334	152	IGHV5-51_02	IGHJ2_01	0.05	10.44	3.38	Y	Ν	26.01

1 The sample name is clearly identified and the total number of reads (= Read Depth) generated for the sample is provided. Following the IFU, it is easy to determine whether the data generated for a sample can be assessed for the presence or absence of clonality.

2 The sequence of clonal populations is provided and populations are ranked from most abundant to least prevalent. Sequences that differ by 1-2 bp are automatically merged to account for possible sequencing errors and to improve the accuracy and ease of sample interpretation.

3 Sequences are aligned with reference genes to allow for easy identification of specific types of gene rearrangements such as *IGHV*3-21, which is characteristic of some CLL cases and correlates with a poor prognosis.

The percentage that a unique sequence contributes to the total number of reads for a sample is calculated. Following the guidelines in the IFU, samples can be interpreted for the evidence indicating the presence or absence of clonality.

³ The software for LymphoTrack Dx *IGHV* Leader and *IGH* FR1 Assays automatically calculates the somatic hypermutation status of a sequence by comparing the identified sequence with a germline reference. In addition, predictions on whether the sequence would be functional can be drawn by the provided information regarding the presence of a premature stop codon or an open reading frame that is out-of-frame.

Ordering Information

Catalog #

9-500-0009 9-500-0007

Products

LymphoTrack® Dx Software – MiSeq[™] LymphoTrack® Dx Software – S5/PGM[™]

Quantity Components

1 Software Package 1 Software Package

E IVD These products are CE-IVD assays for *in vitro* diagnostic use and available in regions that accept CE-IVD products.



LymphoTrack[®]

Next-Generation Sequencing (NGS) Research Use Only (RUO) Assays

LymphoTrack Assay kits are designed for the identification of gene rearrangements in hematologic samples utilizing next-generation sequencing (NGS) technologies.

LEARN MORE: catalog.invivoscribe.com

Introduction

These assays take advantage of the wealth and depth of NGS data generated with either the Illumina[®] MiSeq[™] or the Thermo Fisher Scientific[®] Ion PGM[™] and Ion S5[™] platforms. They offer significant improvements over conventional clonality testing methods, by providing the distribution of gene rearrangements, DNA sequences, the mutational status (*IGH* only), and the ability to track specific gene rearrangements all with the same workflow.

Primers included in the master mixes are designed with Illumina® adapters and indices (8 or 24 indices configurations for a total possible 24 or *48 unique indices [*FR1 only] per framework kits) or Thermo Fisher adapters and indices (12 indices per framework kits). By offering multiple kit configurations (8-or-24 indices for MiSeq[™], 12 for Ion S5/ PGM[™]), Invivoscribe provides laboratories the ability to choose the optimal kit for their sample throughput and read-depth requirements. Testing costs can be reduced by sequencing in a single run, with the possibility to combine: a) samples with up to 48 different indices and b) amplicons from other LymphoTrack Assays.

Key Benefits

- » One-step PCR for amplicon and library generation
- » Identify and assess mutation status of B- and T-cell gene rearrangements
- » Sequence amplicons from any LymphoTrack kit together
- » Bioinformatics software for analysis and interpretation
- » Same reagents for clonality, somatic hypermutation (SHM), measurable residual disease (MRD) testing, and tracking/ monitoring of immunotherapy constructs

Chapter Contents

NGS Assays

- $\textbf{45} \hspace{0.1 cm} \texttt{LymphoTrack} \hspace{0.1 cm} \textit{IGHV} \hspace{0.1 cm} \texttt{Leader} \hspace{0.1 cm} \texttt{Somatic} \hspace{0.1 cm} \texttt{Hypermutation} \hspace{0.1 cm} \texttt{Assay}$
- 47 LymphoTrack *IGH* FR1/2/3 Assays
- **51** LymphoTrack *IGK* Assays
- 53 LymphoTrack TRG Assays
- 55 LymphoTrack TRB Assay
- 57 LymphoTrack Software S5/PGM™
 58 LymphoTrack Software MiSeq™
- 58 LymphoTrack Software Analysis Outputs

Research Use Only (RUO) assays are not for sale in Europe and other global markets where equivalent CE-IVD assays are available and registered with the appropriate regulatory agencies. For more information regarding the CE-marked LymphoTrack® Dx Assays, please refer to the preceding pages.

LymphoTrack[®] *IGHV* Leader Somatic Hypermutation Assays

Assay Uses

This research use only (RUO) assay for next-generation sequencing (NGS), identifies clonal *IGH* V_{H} -J_H rearrangements, the associated V_{H} -J_H DNA sequences, and the frequency distribution of V_{H} region and J_H region segment utilization. The assay also uses the Illumina[®] MiSeq[®] platform to define the extent of somatic hypermutation (SHM) present in the *IGHV* gene of analyzed samples. If you would like to test for *IGHV* somatic hypermutation using the Thermo Fisher[®] Ion PGM[®] or Ion S5[®] platform please refer to the LymphoTrack *IGH* FR1 Assay (71210007).

Summary and Explanation of the Test

The LymphoTrack *IGHV* Leader Somatic Hypermutation Assay for NGS provides significant improvements over clonality assays using fragment analysis and Sanger sequencing. The assay efficiently detects the majority of *IGH* gene rearrangements using a single multiplex master mix, identifies the DNA sequence specific for each clonal gene rearrangement, and calculates the degree of SHM for each sample.

The master mixes included in this assay target the Leader (V_HL) and the joining (J_H) gene regions of *IGH* and are designed with Illumina® adapters and indices (8 included in Kit A and 24 included in Panels). Including adapters and indices in the primer design allows for a one-step PCR approach to generate sequence-ready amplicons, followed by direct pooling of samples for sequencing on a Illumina® MiSeq[™] flow cell.

Positive (clonal positive, SHM negative), negative (polyclonal), and SHM positive (clonal positive, SHM positive) controls are included in the kit.

LymphoTrack Software is designed to sort complex NGS data generated from LymphoTrack Assays by gene target, providing easy bioinformatics analyses and data visualization. Once the clonal sequence(s) are identified, they can be monitored/tracked using the LymphoTrack MRD Software in subsequent samples. Please see the LymphoTrack MRD software section to learn how the LymphoTrack Assays can be applied to MRD studies, or email <u>marketing@invivoscribe.com</u>.

Background

The human immunoglobulin heavy chain (*IGH*) gene locus on chromosome 14 (14q32.3) includes 46 – 52 functional and 30 nonfunctional variable (V_H) gene segments, 27 functional diversity (D_H) gene segments, and 6 functional joining (J_H) gene segments spread over 1250 kilobases.

During B-cell development, genes encoding the *IGH* molecules are assembled from multiple polymorphic gene segments that undergo rearrangements and selection, generating V_{H} - D_{H} - J_{H} combinations that are unique in both length and sequence.¹

Since leukemias and lymphomas originate from the malignant transformation of individual lymphoid cells, all leukemias and lymphomas generally share one or more cell-specific or "clonal" antigen receptor gene rearrangements. Therefore, tests that detect *IGH* clonal rearrangements can be useful in the study of B-cell malignancies. An additional level of diversity is generated in the antigen receptors by somatic point mutations in the variable regions and this mutation status provides important prognostic information for chronic lymphocytic leukemia (CLL)² and small lymphocytic lymphoma (SLL). In addition, NGS methods can improve disease stratification and elucidate subclone gene profiles.

Specimen Requirement

50 ng of high-quality genomic DNA.

References

- 1. Miller et al., Molecular Genetic Pathology (2nd ed.). Springer Science & Business Media. 2013: 302.2.7.13 and 30.2.7.18.
- 2. Ghia et al., Leukemia 21: 2-3 (2007).



Simple representation of the organization of the immunoglobulin heavy chain (*IGH*) gene on chromosome 14. Depicted are the variable region (V_H) genes and downstream consensus joining region genes (J_H) that are involved in rearrangements. Upstream of the variable gene segments, the leader sequence (V_HL) is also depicted. Diversity region genes are not depicted.



Reagents – MiSeq[™] Detection

Kit A Components

Panel Components (includes all master mixes from Kit A plus the items below)

Master Mix Name	Index #	Master Mix Name	Index #	Master Mix Name	Index #
IGH Leader MiSeq 01	A001	IGH Leader MiSeq 09	A009	IGH Leader MiSeq 18	A018
IGH Leader MiSeq 02	A002	<i>IGH</i> Leader MiSeq 10	A010	IGH Leader MiSeq 19	A019
IGH Leader MiSeq 03	A003	IGH Leader MiSeq 11	A011	IGH Leader MiSeq 20	A020
<i>IGH</i> Leader MiSeq 04	A004	<i>IGH</i> Leader MiSeq 12	A012	IGH Leader MiSeq 21	A021
IGH Leader MiSeq 05	A005	IGH Leader MiSeq 13	A013	IGH Leader MiSeq 22	A022
IGH Leader MiSeq 06	A006	<i>IGH</i> Leader MiSeq 14	A014	IGH Leader MiSeq 23	A023
IGH Leader MiSeq 07	A007	<i>IGH</i> Leader MiSeq 15	A015	IGH Leader MiSeq 25	A025
IGH Leader MiSeq 08	A008	<i>IGH</i> Leader MiSeq 16	A016	IGH Leader MiSeq 27	A027

Controls

IGH SHM POS (+) Qty. 1 IGH POS (+) Qty. 1 NGS NEG (-) Qty. 1

Controls

IGH SHM POS (+) Qty. 3 IGH POS (+) Qty. 3 NGS NEG (-) Qty. 3

Rank	Sequence	Length	Merge count	V-gene	J-gene	% Total reads	Cumulative %	Mutation rate partial V-gene (%)	In-frame (Y/N)	No stop codon (Y/N)	V-coverage
1	TTCTCGTGGTG	455	29603	IGHV4-59_08	IGHJ4_02	9.93	9.93	11.26	Y	Y	98.63
2	CTCGCCCTCCT	463	205	IGHV5-51_01	IGHJ4_02	0.07	9.99	0.00	Y	Y	99.66
3	GGTTTTCCTTG	484	201	IGHV3-7_01	IGHJ4_02	0.07	10.06	7.77	Y	Y	100.00
4	CTCGCCCTCCT	463	185	IGHV5-51_01	IGHJ5_02	0.06	10.12	6.08	Y	Y	99.32
5	CTCGCCCTCCT	469	170	IGHV5-51_01	IGHJ4_02	0.06	10.18	0.00	Y	Y	99.32
6	CTCGCCCTCCT	466	160	IGHV5-51_01	IGHJ4_02	0.05	10.23	0.00	Y	Y	99.66
7	CTGCTGCTGAC	460	159	IGHV2-5_10	IGHJ5_02	0.05	10.29	8.08	Y	Y	97.64
8	GGTTTTCCTTG	493	156	IGHV3-48_02	IGHJ6_02	0.05	10.34	3.72	Y	Y	98.99
9	CTCGCCCTCCT	334	153	IGHV5-51_02	IGHJ2_01	0.05	10.39	3.72	Y	Ν	27.70
10	CTCGCCCTCCT	334	152	IGHV5-51_02	IGHJ2_01	0.05	10.44	3.38	Y	Ν	26.01

Example Data. Depicted are the top 10 sequences from a read summary generated by the LymphoTrack Software – MiSeqTM. Highlighted columns represent fields that are unique to SHM analysis and include the SHM mutation rate and predictions pertaining to whether a sequence is in-frame or contains a premature stop codon. To learn more about the LymphoTrack software, please refer to the LymphoTrack Bioinformatics Software section.

Ordering Information

Catalog #	Products	Quanti
7-121-0059	LymphoTrack® <i>IGHV</i> Leader Somatic Hypermutation Assay Kit A - MiSeq™	8 indices
7-121-0069	LymphoTrack $^{\otimes}$ IGHV Leader Somatic Hypermutation Assay Panel – MiSeq $^{\scriptscriptstyle \mathbb{M}}$	24 indice
7-500-0009	LymphoTrack® Software - MiSeq™	1 Softwa

Quantity Components

8 indices – 5 sequencing reactions each

- 24 indices 5 sequencing reactions each
- 1 Software Package

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.



LymphoTrack[®] *IGH* FR1/FR2/FR3 Assays

Assay Uses

These research use only (RUO) assays for next-generation sequencing (NGS), identify clonal *IGH* V_H-J_H rearrangements, the associated V_H-J_H region DNA sequences, the frequency distribution of V_H region and J_H region segment utilization. The LymphoTrack FR1 can also identify the degree of somatic hypermutation (SHM) of rearranged genes using the Illumina[®] MiSeq[™], Thermo Fisher Scientific[®] Ion PGM[™] or Ion S5[™]. The LymphoTrack *IGH* FR1, *IGH* FR2, and *IGH* FR3 Assays contain primers that target the conserved framework 1 (FR1), framework 2 (FR2), and framework 3 (FR3) regions, respectively. The LymphoTrack *IGH* FR1/2/3 Assay kits contain the master mixes of all three frameworks.

Summary and Explanation of the Test

The LymphoTrack *IGH* Assays represent a significant improvement over clonality assays that utilize fragment analysis by providing four important and complementary uses:

- 1. Detection of clonal populations.
- 2. Identification of sequence information and $V_{\text{H}}\text{-}J_{\text{H}}$ segment utilization.
- 3. The LymphoTrack IGH FR1 Assays provide the degree of SHM of the immunoglobulin variable heavy chain (IGHV).
- 4. The ability to track sequences in subsequent samples with the Invivoscribe LymphoTrack MRD Software.

Each single multiplex master mix targets one of the conserved *IGH* framework regions (FR1, FR2, or FR3) within the V_H and the J_H regions described in lymphoid malignancies. **Targeting all three framework regions significantly reduces the risk of not being able to detect the presence of clonality**, as somatic hypermutations in the primer binding sites of the involved V_H gene segments can impede DNA amplification.¹

Primers included in the master mixes are designed with Illumina® adapters and indices (8 included in Kit A, 24 included in the Panel, and an independent 24 included in the Panel B) or Thermo Fisher adaptors and indices (12 indices per framework kits). This allows for a one-step PCR reaction to generate sequence-ready amplicons and pooling of several different samples on the same Illumina® MiSeq[™] cell, Ion S5 or Ion PGM chip.

Positive clonal (SHM negative) and negative polyclonal DNA controls are included in kits. A clonal SHM positive control can be purchased separately (Catalog # 40880008).

LymphoTrack Software is designed to sort complex NGS data generated from LymphoTrack Assays by gene target, providing easy bioinformatics analyses and data visualization. Once the clonal sequence(s) are identified, they can be monitored/tracked using the LymphoTrack MRD Software in subsequent samples.. Please see the LymphoTrack MRD software section to learn how the LymphoTrack Assays can be applied to MRD studies, or email <u>marketing@invivoscribe.com</u>.

Background

The human immunoglobulin heavy chain (*IGH*) gene locus on chromosome 14 (14q32.3) includes 46–52 functional and 30 nonfunctional variable (V_H) gene segments, 27 functional diversity (D_H) gene segments, and 6 functional joining (J_H) gene segments spread over 1250 kilobases.

During development of lymphoid cells, the antigen receptor genes go through somatic gene rearrangements.¹ For example, during B-cell development, genes encoding the *IGH* molecules are assembled from multiple polymorphic gene segments that undergo rearrangements and selection, generating V_H-D_H-J_H combinations. Since leukemias and lymphomas originate from the malignant transformation of individual lymphoid cells, all leukemias and lymphomas generally share one or more cell-specific or "clonal" antigen receptor gene rearrangements. Therefore, tests that detect *IGH* clonal populations can be useful in the study of B- and T-cell malignancies.

Specimen Requirement

50 ng of high-quality genomic DNA.

Reference

1. S Tonegawa. Nature 302: 575-581 (1983).



Simple representation of the organization of the immunoglobulin heavy chain (*IGH*) gene locus on chromosome 14. Depicted are the variable region (V_H) genes and downstream consensus joining region genes segments (J_H) that are involved in rearrangements.





V-J Sequence Frequencies. The LymphoTrack Software provides a stacked bar graph depicting the relative frequencies for the 200 most prevalent VH-JH rearrangements identified in a sample. To learn more about the LymphoTrack software, please refer to the LymphoTrack Bioinformatics Software section.

Ordering Information

Catalog #	Products
7-121-0129	LymphoTrack® <i>IGH</i> FR1/2/3 Assay Kit A – MiSeq™
7-121-0139	LymphoTrack® <i>IGH</i> FR1/2/3 Assay Panel - MiSeq™
7-121-0009	LymphoTrack® <i>IGH</i> FR1 Assay Kit A – MiSeq™
7-121-0039	LymphoTrack [®] <i>IGH</i> FR1 Assay Panel - MiSeq [™]
7-121-0149	LymphoTrack® <i>IGH</i> FR1 Assay Panel B - MiSeq™
7-121-0089	LymphoTrack [®] <i>IGH</i> FR2 Assay Kit A - MiSeq™
7-121-0099	LymphoTrack® <i>IGH</i> FR2 Assay Panel - MiSeq™
7-121-0109	LymphoTrack [®] <i>IGH</i> FR3 Assay Kit A - MiSeq™
7-121-0119	LymphoTrack® <i>IGH</i> FR3 Assay Panel - MiSeq™
7-500-0009	LymphoTrack® Software – MiSeq™
7-121-0057	LymphoTrack® <i>IGH</i> FR1/2/3 Assay – S5/PGM™
7-121-0007	LymphoTrack® <i>IGH</i> FR1 Assay - S5/PGM™
7-121-0037	LymphoTrack® <i>IGH</i> FR2 Assay - S5/PGM™
7-121-0047	LymphoTrack® <i>IGH</i> FR3 Assay - S5/PGM™
7-500-0007	LymphoTrack® Software - S5/PGM™
7-500-0008	LymphoTrack® MRD Software

Quantity Components

Indices 1-8 (5 sequencing reactions each) Indices 1-24 (5 sequencing reactions each) Indices 1-8 (5 sequencing reactions each) Indices 1-24 (5 sequencing reactions each) Indices 25-48 (5 sequencing reaction each) Indices 1-8 (5 sequencing reactions each) Indices 1-24 (5 sequencing reactions each) Indices 1-8 (5 sequencing reactions each) Indices 1-24 (5 sequencing reactions each) 1 Software Package 12 + 12 + 12 indices - 5 sequencing reactions each 12 indices – 5 sequencing reactions each 12 indices – 5 sequencing reactions each 12 indices – 5 sequencing reactions each 1 Software Package 1 Software Package

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.



Reagents - MiSeq[™] Detection

The LymphoTrack IGH FR1/2/3 Assays contain components from respective individual FR Kit A's, Panel's or Panel B.

LymphoTrack IGH FR1 Components		LymphoTrack IGH FR	2 Components	LymphoTrack IGH FR3 Components		
Master Mix Name	Index #	Master Mix Name	Index #	Master Mix Name	Index #	
IGH FR1 MiSeq 01	A001	IGH FR2 MiSeq 01	A001	IGH FR3 MiSeq 01	A001	
IGH FR1 MiSeq 02	A002	IGH FR2 MiSeq 02	A002	IGH FR3 MiSeq 02	A002	
IGH FR1 MiSeq 03	A003	IGH FR2 MiSeq 03	A003	IGH FR3 MiSeq 03	A003	
IGH FR1 MiSeq 04	A004 <	IGH FR2 MiSeq 04	A004 <	IGH FR3 MiSeq 04	A004 <	
IGH FR1 MiSeq 05	A005	IGH FR2 MiSeq 05	A005 💆	IGH FR3 MiSeq 05	A005 💆	
IGH FR1 MiSeq 06	A006	IGH FR2 MiSeq 06	A006	IGH FR3 MiSeq 06	A006	
IGH FR1 MiSeq 07	A007	IGH FR2 MiSeq 07	A007	IGH FR3 MiSeq 07	A007	
IGH FR1 MiSeq 08	A008	IGH FR2 MiSeq 08	A008	IGH FR3 MiSeq 08	A008	
IGH FR1 MiSeq 09	A009	IGH FR2 MiSeq 09	A009	IGH FR3 MiSeq 09	A009	
IGH FR1 MiSeq 10	A010	IGH FR2 MiSeq 10	A010	IGH FR3 MiSeq 10	A010	
IGH FR1 MiSeq 11	A011	IGH FR2 MiSeq 11	A011	IGH FR3 MiSeq 11	A011	
IGH FR1 MiSeq 12	A012	IGH FR2 MiSeq 12	A012	IGH FR3 MiSeq 12	A012	
IGH FR1 MiSeq 13	A013	IGH FR2 MiSeq 13	A013 d	IGH FR3 MiSeq 13	A013 4	
IGH FR1 MiSeq 14	A014	IGH FR2 MiSeq 14	A014	IGH FR3 MiSeq 14	A014	
IGH FR1 MiSeq 15	A015	IGH FR2 MiSeq 15	A015	IGH FR3 MiSeq 15	A015	
IGH FR1 MiSeq 16	A016	IGH FR2 MiSeq 16	A016	IGH FR3 MiSeq 16	A016	
IGH FR1 MiSeq 18	A018	IGH FR2 MiSeq 18	A018	IGH FR3 MiSeq 18	A018	
IGH FR1 MiSeq 19	A019	IGH FR2 MiSeq 19	A019	IGH FR3 MiSeq 19	A019	
IGH FR1 MiSeq 20	A020	IGH FR2 MiSeq 20	A020	IGH FR3 MiSeq 20	A020	
IGH FR1 MiSeq 21	A021	IGH FR2 MiSeq 21	A021	IGH FR3 MiSeq 21	A021	
IGH FR1 MiSeq 22	A022	IGH FR2 MiSeq 22	A022	IGH FR3 MiSeq 22	A022	
IGH FR1 MiSeq 23	A023	IGH FR2 MiSeq 23	A023	IGH FR3 MiSeq 23	A023	
IGH FR1 MiSeq 25	A025	IGH FR2 MiSeq 25	A025	IGH FR3 MiSeq 25	A025	
IGH FR1 MiSeq 27	A027	IGH FR2 MiSeq 27	A027	IGH FR3 MiSeq 27	A027	
IGH FR1 MiSeq 17	A017					
IGH FR1 MiSeq 24	A024					
IGH FR1 MiSeq 26	A026					
IGH FR1 MiSeq 28	A028					
IGH FR1 MiSeq 29	A029					
IGH FR1 MiSeq 30	A030 😐					
IGH FR1 MiSeq 31	A031					

A032

A033

A034

A035

A036

A037

IGH FR1 MiSeq 32

IGH FR1 MiSeq 33

IGH FR1 MiSeq 34

IGH FR1 MiSeq 35

IGH FR1 MiSeq 36

IGH FR1 MiSeq 37

Reagents - MiSeq[™] Detection Continued

The LymphoTrack IGH FR1/2/3 Assays contain components from respective individual FR Kit A's, Panel's or Panel B.

LymphoTrack IGH FR1 C	omponents Lymp	hoTrack <i>IGH</i> FR2	Components	LymphoTrack IGH FR3 Component			
Master Mix Name	Index # Maste	er Mix Name	Index #	Master Mix Name	ie Index #		
IGH FR1 MiSeq 38	A038						
IGH FR1 MiSeq 39	A039						
IGH FR1 MiSeq 40	A040						
IGH FR1 MiSeq 41	A041						
IGH FR1 MiSeq 42	A042						
IGH FR1 MiSeq 43	A043						
IGH FR1 MiSeq 44	A044						
IGH FR1 MiSeq 45	A045						
IGH FR1 MiSeq 46	A046						
IGH FR1 MiSeq 47	A047						
IGH FR1 MiSeq 48	A048						
IGH FR1 MiSeq 49	A049						
Controls in Individual FR (1, 2, or 3) Kit A's	Controls in Individual FR (1,2, or 3) Panels	Controls in (FR 1/2/3 Kit	Combo Con A FR 1	trols in Combo /2/3 Panel			
IGH POS (+) Qty. 1	IGH POS (+) Qty. 3	IGH POS (+) Qty	v. 2 IGH P	OS (+) Qty. 6			
NGS NEG (-) Qty. 1	NGS NEG (-) Qty. 3	NGS NEG (-) Qt	y. 2 NGS N	NEG (-) Qty. 6			

Reagents - Ion S5[™]/PGM[™] Detection

The LymphoTrack IGH FR1/2/3 Assays contain components from respective individual FR Assays.

LymphoTrack IGH FR1 Components		LymphoTrack IGH FR2	Components	LymphoTrack IGH FR3 Components		
Master Mix Name	Index #	Master Mix Name	Index #	Master Mix Name	Index #	
IGH FR1 S5/PGM 01	lonXpress_001	IGH FR2 S5/PGM 01	lonXpress_001	IGH FR3 S5/PGM 01	lonXpress_001	
IGH FR1 S5/PGM 02	lonXpress_002	IGH FR2 S5/PGM 02	lonXpress_002	IGH FR3 S5/PGM 02	lonXpress_002	
IGH FR1 S5/PGM 03	lonXpress_003	IGH FR2 S5/PGM 03	lonXpress_003	IGH FR3 S5/PGM 03	lonXpress_003	
IGH FR1 S5/PGM 04	lonXpress_004	IGH FR2 S5/PGM 04	lonXpress_004	IGH FR3 S5/PGM 04	lonXpress_004	
IGH FR1 S5/PGM 07	lonXpress_007	IGH FR2 S5/PGM 07	lonXpress_007	IGH FR3 S5/PGM 07	lonXpress_007	
IGH FR1 S5/PGM 08	lonXpress_008	IGH FR2 S5/PGM 08	lonXpress_008	IGH FR3 S5/PGM 08	lonXpress_008	
IGH FR1 S5/PGM 09	lonXpress_009	IGH FR2 S5/PGM 09	lonXpress_009	IGH FR3 S5/PGM 09	lonXpress_009	
IGH FR1 S5/PGM 10	lonXpress_010	IGH FR2 S5/PGM 10	lonXpress_010	IGH FR3 S5/PGM 10	lonXpress_010	
IGH FR1 S5/PGM 11	lonXpress_011	IGH FR2 S5/PGM 11	lonXpress_011	IGH FR3 S5/PGM 11	lonXpress_011	
IGH FR1 S5/PGM 12	lonXpress_012	IGH FR2 S5/PGM 12	lonXpress_012	IGH FR3 S5/PGM 12	lonXpress_012	
IGH FR1 S5/PGM 13	lonXpress_013	IGH FR2 S5/PGM 13	lonXpress_013	IGH FR3 S5/PGM 13	lonXpress_013	
IGH FR1 S5/PGM 14	lonXpress_014	IGH FR2 S5/PGM 14	lonXpress_014	IGH FR3 S5/PGM 14	lonXpress_014	

Controls in Individual FR (1,2, or 3) Kits

IGH POS (+) Qty. 2 NGS NEG (-) Qty. 2

Controls in FR 1/2/3 Kit

IGH POS (+) Qty. 4 NGS NEG (-) Qty. 4

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.



LymphoTrack[®] IGK Assay

Assay Uses

This research use only (RUO) assay for next-generation sequencing (NGS), identifies clonal *IGK* Vk-Jk, Vk-K_{de}, and intron-K_{de} (INTR-K_{de}) rearrangements, the corresponding DNA sequences, and provides the distribution frequency of Vk-Jk, Vk-K_{de}, and INTR-K_{de} segment utilization using the Illumina[®] MiSeq[™] or Thermo Fisher Scientific[®] Ion PGM[™] and Ion S5[™] platforms.

Summary and Explanation of the Test

The LymphoTrack *IGK* Assay represents a significant improvement over clonality assays that utilize fragment analysis by providing three important and complementary uses:

- 1. Detection of clonal populations.
- 2. Identification of sequence information and gene segment utilization.
- 3. Ability to track sequences in subsequent samples with the use of the LymphoTrack MRD Software.

Unlike conventional fragment analysis assays, this NGS method utilizes a single multiplex master mix to target conserved regions of *IGK* that are described in lymphoid malignancies. Primers are designed with Illumina[®] adapters and indices (8-24) or Thermo Fisher Scientific adaptors and indices (12), thereby allowing for a one-step PCR reaction to generate sequence-ready amplicons. To reduce per sample testing costs, amplicons from different samples (amplified with different indexed master mixes) or LymphoTrack kits can be sequenced together on a single Illumina[®] MiSeq[™] flow cell, Ion S5[™] or PGM[™] chips.

Positive (clonal) and negative (polyclonal) DNA controls are included in kits.

LymphoTrack Software is designed to sort complex NGS data generated from LymphoTrack Assays by gene target, providing easy bioinformatics analyses and data visualization. Once the clonal sequence(s) are identified, they can be monitored/tracked using the LymphoTrack MRD Software in subsequent samples. Please see the LymphoTrack MRD software section to learn how the LymphoTrack Assays can be applied to MRD studies, or email <u>marketing@invivoscribe.com</u>.

Background

The human immunoglobulin kappa (*IGK*) gene locus on chromosome 2 (2p11.2) includes 76 V (variable) region genes spanning 7 subgroups and 5 J (joining) region gene segments upstream of the Ck region. The K_{der} approximately 24 kb downstream of the Jk-Ck region, can also rearrange with Vk gene segments and the isolated recombination signal sequence in the Jk-Ck intronic region.

During development of lymphoid cells, antigen receptor genes undergo somatic gene rearrangements.¹ Specifically during B-cell development, genes encoding *IGK* molecules are assembled from multiple polymorphic gene segments that undergo rearrangements generating V-J combinations unique in both length and sequence.² Since leukemias and lymphomas originate from the malignant transformation of individual lymphoid cells, all leukemias and lymphomas generally share one or more cell-specific or "clonal" antigen receptor gene rearrangements. Therefore, tests that detect *IGK* clonal populations can be useful in the study of B-cell malignancies and complement *IGH* testing, as the *IGK* receptor is less susceptible to somatic mutations.

Specimen Requirement

50 ng of high-quality genomic DNA.

References

- 1. S Tonegawa et al., Nature 302: 575-581 (1983).
- 2. JE Miller et al., Molecular Genetic Pathology (2nd ed.). Springer Science & Business Media. 2013: 30.2.7.13 and 30.2.7.18.

Simplified Representation of the IGK Gene

Depicted are the variable region (VK) genes or variable intragenic JK-CK intron (JK-CK INTR) and downstream consensus joining region genes (JK) or kappa deleting element (Kde) that are involved in *IGK* gene rearrangements.



LymphoTrack[®]

Reagents - MiSeq[™] Detection

VII A	Com		-
NITA	Com	pone	nts

Master Mix Name	Index #	Master Mix Name	Index #	Master Mix Name	Index #
IGK MiSeq 01	A001	IGK MiSeq 09	A009	IGK MiSeq 18	A018
IGK MiSeq 02	A002	IGK MiSeq 10	A010	IGK MiSeq 19	A019
IGK MiSeq 03	A003	IGK MiSeq 11	A011	IGK MiSeq 20	A020
IGK MiSeq 04	A004	IGK MiSeq 12	A012	IGK MiSeq 21	A021
IGK MiSeq 05	A005	IGK MiSeq 13	A013	IGK MiSeq 22	A022
IGK MiSeq 06	A006	IGK MiSeq 14	A014	IGK MiSeq 23	A023
IGK MiSeq 07	A007	IGK MiSeq 15	A015	IGK MiSeq 25	A025
IGK MiSeq 08	A008	IGK MiSeq 16	A016	IGK MiSeq 27	A027
Controls		Controls			
IGK POS (+) Qty. 1	NGS NEG (-) Qty. 1	IGK POS (+) Qty. 3		NGS NEG (-) Qty. 3	

Panel Components (includes all master mixes from Kit A plus the items below)

Reagents - Ion S5/PGM[™] Detection

Assay Components

Master Mix Name	Index #	Master Mix Name	Index #
IGK S5/PGM 01	lonXpress_001	IGK \$5/PGM 11	lonXpress_011
IGK \$5/PGM 02	IonXpress_002	<i>IGK</i> \$5/PGM 12	lonXpress_012
IGK \$5/PGM 04	IonXpress_004	<i>IGK</i> \$5/PGM 13	lonXpress_013
IGK \$5/PGM 08	IonXpress_008	<i>IGK</i> \$5/PGM 14	lonXpress_014
IGK \$5/PGM 09	lonXpress_009	<i>IGK</i> S5/PGM 16	lonXpress_016
IGK \$5/PGM 010	lonXpress_010	IGK \$5/PGM 17	lonXpress_017

Controls

IGK POS (+) Qty. 2



V-J Sequence Frequencies.

NGS NEG (-) Qty. 2

The LymphoTrack bioinformatics software provides a stacked bar graph depicting the relative frequencies for the 200 most prevalent rearrangements sequenced and identified in the sample. To learn more about the LymphoTrack software, please refer to the LymphoTrack Bioinformatics Software section.

Ordering Information

Catalog #	Products
7-122-0009	LymphoTrack® <i>IGK</i> Assay Kit A – MiSeq™
7-122-0019	LymphoTrack [®] <i>IGK</i> Assay Panel - MiSeq™
7-500-0009	LymphoTrack [®] Software – MiSeq™
7-122-0007	LymphoTrack [®] <i>IGK</i> Assay - S5/PGM [™]
7-500-0007	LymphoTrack® Software – S5/PGM™
7-500-0008	LymphoTrack® MRD Software*

*MRD Software can be used to track sequences generated by either LymphoTrack Assays - MiSeq[®] or Ion S5/PGM[®]. These products are for Research Use Only (RUO). Not intended for diagnostic purposes.

Quantity Components

8 indices - 5 sequencing reactions each 24 indices - 5 sequencing reactions each 1 Software Package 12 indices - 5 sequencing reactions each 1 Software Package 1 Software Package



LymphoTrack® TRG Assay

Assay Uses

This research use only (RUO) assay for next-generation sequencing (NGS) identifies clonal *TRG* V-J rearrangements and the associated V-J region DNA sequences. It also provides the frequency distribution of V-J segment utilization using the Illumina[®] MiSeq[™] or Thermo Fisher Scientific[®] Ion S5[™] and PGM[™] instruments.

Summary and Explanation of the Test

The LymphoTrack *TRG* Assay represents a significant improvement over existing clonality assays that utilize fragment analysis by providing three important and complementary uses:

- 1. Detection of clonal populations.
- 2. Identification of sequence information and gene segment utilization.
- 3. Ability to track sequences in subsequent samples with the use of the Invivoscribe MRD Software.

This assay utilizes a single multiplex master mix to target conserved V and J regions of the human *TRG* gene locus described in lymphoid malignancies, thereby reducing sample DNA requirements and simplifying the testing workflow. Primers are designed with Illumina® adapters and indices (8-24) or Thermo Fisher Scientific adaptors and indices (12), thereby allowing for a one-step PCR reaction to generate sequence-ready amplicons. In addition, amplicons from different samples (amplified with different indexed master mixes) or LymphoTrack kits can be sequenced together on a single Illumina® MiSeq[™] flow cell Ion S5[™] or PGM[™] chip to reduce per sample testing costs.

Positive (clonal) and negative (polyclonal) DNA controls are included in kits.

LymphoTrack Software is designed to sort complex NGS data generated from LymphoTrack Assays by gene target, providing easy bioinformatics analyses and data visualization. Once the clonal sequence(s) are identified, they can be monitored/tracked using the LymphoTrack MRD Software in subsequent samples. Please see the LymphoTrack MRD software section to learn how the LymphoTrack Assays can be applied to MRD studies, or email <u>marketing@invivoscribe.com</u>.

Background

The human T-cell receptor gamma (*TRG*) gene locus on chromosome 7 (7q14) includes 14 variable region (Group I, II, III, and IV), 5 joining region (J) gene segments, and 2 constant (C) genes spread over 200 kilobases.

During development of lymphoid cells, the antigen receptor genes (including gene segments within the *TRG* locus), undergo somatic gene rearrangement.¹ These developmentally regulated gene rearrangements generate V-J combinations that are unique for each cell.² Since leukemias and lymphomas originate from the malignant transformation of individual lymphoid cells, all leukemias and lymphomas generally share one or more cell-specific or "clonal" antigen receptor gene rearrangements. Therefore, tests that detect *TRG* clonal populations can be useful in the study of T-cell and certain B-cell malignancies. Since the *TRG* locus rearranges before the *TRB* locus, and all unsuccessful or unproductive rearrangements are retained in the cells, examination of the *TRG* locus can identify both clonal Δ/y as well as clonal α/β T-cells. Clonal α/β T-cells would be expected to have biallelic *TRG* gene rearrangements.

Note | For a more thorough explanation of the locus and the targeted deep sequencing strategy, please refer to Principle of Immunoglobulin and T-Cell Receptor Gene Rearrangement.²

Specimen Requirement

50 ng of high-quality genomic DNA.

References

- 1. Tonegawa, S. (1983). Nature 302, 575-581.
- 2. Miller JE. (2013) Molecular Genetic Pathology (2nd Edition., sections 30.2.7.13 and 30.2.7.18).

Simplified Representation of the TRG Gene



Simple representation of the organization of the T-cell receptor gamma gene on chromosome 7. Depicted are the variable region genes (Vy2–Vy11) and downstream joining region genes (Jy1/2–JyP1/2) that are involved in rearrangements in T-cell lymphomas.

LymphoTrack[®]

Reagents - MiSeq[™] Detection

Kit A Components		Panel Components (includes all master mixes from Kit A plus the items below)				
Master Mix Name	Index #	Master Mix Name	Index #	Master Mix Name	Index #	
TRG MiSeq 01	A001	TRG MiSeq 09	A009	TRG MiSeq 18	A018	
TRG MiSeq 02	A002	TRG MiSeq 10	A010	TRG MiSeq 19	A019	
TRG MiSeq 03	A003	TRG MiSeq 11	A011	TRG MiSeq 20	A020	
TRG MiSeq 04	A004	TRG MiSeq 12	A012	TRG MiSeq 21	A021	
TRG MiSeq 05	A005	TRG MiSeq 13	A013	TRG MiSeq 22	A022	
TRG MiSeq 06	A006	TRG MiSeq 14	A014	TRG MiSeq 23	A023	
TRG MiSeq 07	A007	TRG MiSeq 15	A015	TRG MiSeq 25	A025	
TRG MiSeq 08	A008	TRG MiSeq 16	A016	TRG MiSeq 27	A027	
Controls		Controls				
TRG POS (+) Qty. 1	NGS NEG (-) Qty. 1	TRG POS (+) Qty. 3		NGS NEG (-) Qty. 3		

Reagents – Ion S5/PGM[™] Detection

Assay Components

Master Mix N	lame Index #	Master Mix No	ame Index #
.001 TRG \$5/PGM 07	lonXpress_007	<i>TRG</i> \$5/PGM 11	lonXpress_011
.002 TRG \$5/PGM 08	lonXpress_008	<i>TRG</i> \$5/PGM 12	lonXpress_012
.003 TRG \$5/PGM 09	lonXpress_009	<i>TRG</i> S5/PGM 13	lonXpress_013
.004 TRG S5/PGM 10	lonXpress_010	<i>TRG</i> S5/PGM 14	lonXpress_014
	Master Mix N _001 TRG \$5/PGM 07 _002 TRG \$5/PGM 08 _003 TRG \$5/PGM 09 _004 TRG \$5/PGM 10	Master Mix Name Index # _001 TRG \$5/PGM 07 lonXpress_007 _002 TRG \$5/PGM 08 lonXpress_008 _003 TRG \$5/PGM 09 lonXpress_009 _004 TRG \$5/PGM 10 lonXpress_010	Master Mix Name Index # Master Mix Name _001 TRG \$5/PGM 07 IonXpress_007 TRG \$5/PGM 11 _002 TRG \$5/PGM 08 IonXpress_008 TRG \$5/PGM 12 _003 TRG \$5/PGM 09 IonXpress_009 TRG \$5/PGM 13 _004 TRG \$5/PGM 10 IonXpress_010 TRG \$5/PGM 14

Controls

TRG POS (+) Qty. 2



V - J Sequence Frequencies.

NGS NEG (-) Qty. 2

The LymphoTrack bioinformatics software provides a stacked bar graph depicting the relative frequencies for the V–J rearrangements identified and sequenced in a sample. To learn more about the LymphoTrack software, please refer to the LymphoTrack Bioinformatics Software section.

Ordering Information

Catalog #	Products
7-227-0019	LymphoTrack® <i>TRG</i> Assay Kit A - MiSeq™
7-227-0009	LymphoTrack® <i>TRG</i> Assay Panel – MiSeq™
7-500-0009	LymphoTrack® Software – MiSeq™
7-227-0007	LymphoTrack® <i>TRG</i> Assay - S5/PGM™
7-500-0007	LymphoTrack® Software – S5/PGM™
7-500-0008	LymphoTrack® MRD Software*

Quantity Components

8 indices – 5 sequencing reactions each 24 indices – 5 sequencing reactions each 1 Software Package

- 12 indices 5 sequencing reactions each
- 1 Software Package
- 1 Software Package

*MRD Software can be used to track sequences generated by either LymphoTrack Assays - MiSeq[™] or Ion S5/PGM[™]. These products are for Research Use Only (RUO). Not intended for diagnostic purposes.



LymphoTrack® TRB Assay

Assay Uses

This research use only (RUO) assay for next-generation sequencing (NGS) identifies clonal *TRB* V_{β} –(D_{β} –) J_{β} rearrangements, the associated V_{β} –(D_{β} –) J_{β} region DNA sequences, and provides the frequency distribution of V_{β} , D_{β} , and J_{β} region segment utilization using the Illumina[®] MiSeq[®] platform.

Analysis of the rearranged *TRB* locus increases the probability of identifying T-cell receptor gene rearrangements, as compared to testing for *TRG* gene rearrangements only. As a result, combining the analysis of *TRB* and *TRG* loci increases the sensitivity of clonality detection.

Summary and Explanation of the Test

This assay utilizes a single multiplex master mix to target conserved V, D and J regions of the human *TRB* gene locus described in lymphoid malignancies, thereby reducing sample DNA requirements and simplifying the testing workflow.

The LymphoTrack *TRB* master mix primers are also designed with Illumina[®] adapters and 8 indices (Kit A) or 24 indices (Panel). This allows amplicons generated from different indexed *TRB* master mixes to be pooled into a single sequencing library.

Positive (clonal) and negative (polyclonal) DNA controls are included in kits.

LymphoTrack Software is designed to sort complex NGS data generated from LymphoTrack Assays by gene target, providing easy bioinformatics analyses and data visualization. Once the clonal sequence(s) are identified, they can be monitored/tracked using the LymphoTrack MRD Software in subsequent samples. Please see the LymphoTrack MRD software section to learn how the LymphoTrack Assays can be applied to MRD studies, or email <u>marketing@invivoscribe.com</u>.

Background

The LymphoTrack *TRB* Assay represent a significant improvement over clonality assessment by fragment analysis by providing two important and complementary uses:

1. Detection of clonal populations.

- 2. Identification of sequence information and gene segment utilization.
- 3. Ability to track sequences in subsequent samples with the use of the Invivoscribe MRD Software.

The human T-cell receptor beta (*TRB*) gene locus on chromosome 7 (7q34) includes 65 V β (variable) gene segments, followed by two separate clusters of genes each containing a D β (diversity) gene, several J β (joining) genes, and a C β (constant) region spread over 685 kilobases. The two C β genes, *TRBC1* and *TRBC2*, encode highly homologous products with no functional difference.

During lymphoid cell development, antigen receptor genes undergo somatic gene rearrangements.¹ Specifically, during T-cell development, genes encoding *TRB* molecules are assembled from multiple polymorphic gene segments that generate V_B-D_B-J_B combinations unique in both length and sequence.²

Since leukemias and lymphomas originate from the malignant transformation of individual lymphoid cells, all leukemias and lymphomas generally share one or more cell-specific or "clonal" antigen receptor gene rearrangements. Therefore, tests that detect *TRB* clonal rearrangements can be useful in the study of B- and T-cell malignancies.

Note: For a more thorough explanation of the locus and the targeted deep sequencing strategy, please refer to Principle of Immunoglobulin and T-Cell Receptor Gene Rearrangement.²

Specimen Requirement

50 ng of high-quality genomic DNA.

References

1. Tonegawa, S. (1983). Nature 302, 575-581.

2. Miller JE. (2013) Molecular Genetic Pathology (2nd Edition., sections 30.2.7.13 and 30.2.7.18).



Simplified Representation of the TRB Gene

_ Dβ1 _

Depicted are the variable (V β), diversity (D β), and joining (J β) gene regions involved in *TRB* gene rearrangements, in addition to the downstream consensus (C β) gene regions.

Dβ2

Cβ2

Cβ1

Reagents - MiSeq[™] Detection

٧β

Kit A Components		Panel Components (includes all master mixes from Kit A plus the items below)						
Master Mix Name Index #		Master Mix Name	Index #	Master Mix Name	Index #			
TRB MiSeq 01	A001	TRB MiSeq 09	A009	TRB MiSeq 18	A018			
TRB MiSeq 02	A002	TRB MiSeq 10	A010	TRB MiSeq 19	A019			
TRB MiSeq 03	A003	TRB MiSeq 11	A011	TRB MiSeq 20	A020			
TRB MiSeq 04	A004	TRB MiSeq 12	A012	TRB MiSeq 21	A021			
TRB MiSeq 05	A005	TRB MiSeq 13	A013	TRB MiSeq 22	A022			
TRB MiSeq 06	A006	TRB MiSeq 14	A014	TRB MiSeq 23	A023			
TRB MiSeq 07	A007	TRB MiSeq 15	A015	TRB MiSeq 25	A025			
TRB MiSeq 08	A008	TRB MiSeq 16	A016	TRB MiSeq 27	A027			
Controls		Controls						
TRB POS (+) Qty. 1	NGS NEG (-) Qty. 1	TRB POS (+) Qty. 3		NGS NEG (-) Qty. 3				



V-J Sequence Frequencies.

The LymphoTrack bioinformatics software provides PDF reports which include Top 10 Merged Read Summary as well as a stacked bar graph depicting the relative frequencies of the 200 most prevalent rearrangements sequenced and identified in the sample. To learn more about the LymphoTrack software, please refer to the LymphoTrack Software section.

Ordering Information

Catalog #	Products
7-225-0009	LymphoTrack® <i>TRB</i> Assay Kit A – MiSeq™
7-225-0019	LymphoTrack® <i>TRB</i> Assay Panel - MiSeq™
7-500-0009	LymphoTrack® Software - MiSeq™
7-500-0008	LymphoTrack® MRD Software*

Quantity Components

8 indices – 5 sequencing reactions each 24 indices – 5 sequencing reactions each 1 Software Package 1 Software Package

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.



LymphoTrack[®] Software – S5/PGM™

Software Use

The Software – S5/PGM^M package is used with LymphoTrack Assays – S5/PGM^M to analyze raw FASTQ files for clonality analysis of single or multiple target data sets, *IGH* FR1, *IGH* FR2, *IGH* FR3, *IGK*, *TRG*). Analysis of LymphoTrack IGH FR1 Assay also includes the rate of somatic hypermutation (SHM) and whether a clone will be functional based upon the presence of a premature stop codon. The software can also predict whether an open reading frame would be in- or out-of-frame, so no external data analysis is required for sample interpretation.

The provided software is composed of two distinct parts:

1. A bioinformatics Data Analysis Application

2. Microsoft Excel® Data Visualization spreadsheets for each gene target and automated Sample-to-PDF Reports for streamlined data analysis.

The LymphoTrack Software - S5/PGM Data LymphoTrackSoftware-S5PGM_v2.4.x.RUO Analysis Application provides a simple user interface for selecting data to analyze. 🖾 LymphoTrack[®] Choose analysis mode(s): Leader TRG TRB IGH FR1 IGH FR2 IGH FR3 IGK Choose decimal mark: "." decimal point, e.g., "1,000.00" "," decimal comma, e.g., "1.000,00" Browse Choose an Input Folder for analysis: Launch Program

For Research Use Only. Not for use in Diagnostic Procedure

LymphoTrack Report for assay IGH_FR1

Sample name: id14-7-targets_S14_L001_001_combined

Total Read Count: 93322

	Top 10 Merged Read Summary											
Rank	Sequence	Length	Merge count	V-gene	J-gene	% total reads	Cumulativ e %	Mutation rate to partial V- gene (%)	In-frame (Y/N)	No Stop codon (Y/N)	V- coverage	CDR3 Seq
1	CATCTGGATACAC	295	9172	IGHV1- 46_03	IGHJ4_02	9.83	9.83	0.00	Y	Y	100.00	not found
2	GCCTCTGGATACA	278	120	IGHV3- 30_04	IGHJ5_02	0.13	9.96	13.22	Y	Y	100.00	
3	GCCTCTGGATTCA	278	66	IGHV3- 64_01	IGHJ6_02	0.07	10.03	5.73	Y	Y	100.00	not found
4	GCGTCTGGATTCA	278	58	IGHV3- 33_01	IGHJ6_02	0.06	10.09	0.00	Y	Y	100.00	not found
5	GCCTCTGGATTCA	281	56	IGHV3- 30_18	IGHJ6_02	0.06	10.15	4.85	Y	Y	100.00	not found
6	GCCTCTGGATTCA	260	56	IGHV3- 48_01	IGHJ4_02	0.06	10.21	0.00	Y	Y	98.68	not found
7	GCCTCTGGATTCA	269	56	IGHV3- 30_18	IGHJ4_02	0.06	10.27	0.00	Y	Y	99.12	not found
8	GCCTCTGGATTCA	272	55	IGHV3- 11_05	IGHJ4_02	0.06	10.33	9.25	Y	Y	98.24	not found
9	GCCTCCGAAGTCA	167	55	IGHV3- 43_01	none	0.06	10.39	13.10	n/a	N	64.63	not found
10	GCCTCTGGATTCA	284	54	IGHV3- 48_01	IGHJ6_02	0.06	10.45	0.00	Y	Y	100.00	not found

A PDF report is generated for each sample at the completion of analysis and provides a Top 10 Merged Read Summary, the top 10 reads after merging with the top 500 unique reads that differ by the maximum number of mismatches (2 bp) and the associated V- and J-gene families; a V - J Usage Graph, containing a 3D bar chart showing the V- and I-family usage for the unique sequences included in the analysis (IGH/IGK and TRB) or V- and |gene usage (TRG); a V - J Usage Graph, which is a 3D bar chart showing the Vand J-family usage (IGH/IGK and TRB) or V- and J- gene usage (TRG).

Ordering Information

Catalog # 7-500-0007

Products

LymphoTrack[®] Software – S5/PGM[™]

Quantity Components

1 Software Package

LymphoTrack®

LymphoTrack[®] Software – MiSeq™

Software Use

The research use only (RUO) LymphoTrack Software – MiSeq[™] is a software package designed to provide bioinformatics analysis of FASTQ data files generated from Illumina[®] MiSeq[™] systems using any combination of LymphoTrack Assays – MiSeq[™]. The analysis results include identification of gene rearrangement families and frequencies of single or multiple target data (*IGHV* Leader [or "Leader"], *IGH* FR1, *IGH* FR2, *IGH* FR3, *IGK*, *TRB* or *TRG*) as well as the somatic hypermutation status of samples run with the Leader or *IGH* FR1 assays.

IvmnhoTrack® Software	-	П	×
© Sattings Halp		2	
© Settings neip			
LymphoTrack [®]			
Choose the read length:	2 x 300 bp	,	٣
Choose analysis mode(s):			
IGH FR1 IGH FR2 IGH FR3 IGK			
Leader TRB TRG			
Choose a reference database for analysis:	IMGT/V-Q	UEST	*
Choose the number of unique sequences to analyze:	All		٣
Choose the maximum allowable mismatches between two sequences:	3 bp		*
Choose an Input Folder for analysis:			
	OWSE		
ANALYZE			
RUO For Research Use Only; not for use in diagnostic procedure.			
Copyright @2025 Invivoscribe, Inc. Lym	phoTrack®	Softw	are

LymphoTrack Software - MiSeq™ includes the ability to select different analysis parameters:

- read length
- assay targets
- reference database
- number of unique sequences to analyze
- maximum allowable mismatches

A PDF report is generated for each sample at the completion of analysis and provides a Top 10 Merged Read Summary, the top 10 reads after merging with the selected number of unique reads (top 500 or all) that differ by the selected maximum number of mismatches (0, 1, 2 or 3 bp) and the associated V- and J-gene families; a V - J Usage Graph, containing a 3D bar chart showing the Vand J-family usage for the unique sequences included in the analysis (IGH/IGK and TRB) or Vand J- gene usage for all sequences (TRG); a V - J Usage Graph, which is a 3D bar chart showing the V- and J-family usage for the unique sequences included in the analysis (IGH/IGK and TRB) or Vand J- gene usage for all sequences (TRG); and an area for comments and signatures.

LymphoTrack

IGHV Leader Somatic Hypermutation Assay

	pre reame.	arr berr or				10	tai keau count.	300401	Q30 30016	Collin and all
lanaly	sis R	lead Length:	2 × 250	1 bp			Unique Sequences Ana	slyzed: Al		
aran	ieters R	eference Database	IMGT-	QUEST			Allowable Mismatches	a 31	te	
-										
op 1	10 Merged Re	ad Summary								
nk -	Sequence Lengt	th Merze Count	V-pene	Larene		of Total Reads	Mutation Rate to i	Partial In-frame IV	N) No Stop Coden	V-Coverage (%)
	(bp)				Actual	Cumulative	V-gene (%)		(Y/N)	
_	295	20948	IGHV1-46_03	IGHJ4_02	6.7923	6.7923	0.0000	Y	Y	100.0000
	Sequence:	CATCTGGAT/ AGTTCCAGG GGTTGTATTA	ACACCTTCACCA SCAGAGTCACCA GTACCAGCTGC	GCTACTATATGC ATGACCAGGGAC	ACTOGOTOCOAC ACGTCCACGAGO	AGGCCCCTGGACAA ACAGTCTACATGGA	GOOCTTGAGTGGATG GCTGAGCAGCCTGAG	ATCTGAGGACAC	CCTAGTOGTOGTAGC GCCCGTGTATTACTG	ACAAGCTACGCACAG FGCTAGAGATCTCAC
	CDR3:	TOTOCTADA	DATCTCACAGGT	TOTATTAGTACC	AGCTOCTATCCT	COAACTACTITGAC	TACTOGOCCAGOGA			
_	110	112	IGHV3-30-3_03	IGHJ4 02	0.0363	6.8286	3.6242	n/a	n/a	36.5639

	Sequence:	GCCTCTGGA	TTCACCTTCAGT	ATTGTGTTGTTG	CACTGGGTCCGC	CAGGCTCCAGGCAA	GGGGCTGGAATGGGT	GGCATTTATTATG	TGATGCCTTGACTAC	16666CCA666A6
	Sequence: CDR3:	GCCTCTGGA CT not found	ттелесттелет	AATTGTGTTGTTG	CACTGGGTCCGC	CAGGCTCCAGGCAA	GGGGCTGGAATGGGT	GGCATTTATTATG	TGATGCCTTGACTAC	16666004666466
	Sequence: CDR3: 197	GCCTCTGGA CT not found 104	IGHV3-9_04	IGHJ6_03	0.0337	6.8624	3.0568			28.8210
	Sequence: CDR3: 197 Sequence:	GCCTCTGGA CT not found 104 GCCTCTGAA ACCTTTAGA	IGHV3-9_04	IGHJ6_03 AGAGGCTATGAC TAGAATCAAGCT	0.0337 GGACTCTGTGGAI TAGACTACGACC	6.8624 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC	3.0568 CAAACGAACAGTCTGA DTCT0000CAGAGGG	N CAGCACAGAGGACA		28.8210 GTOCAAGAGACTAAA
	Sequence: CDR3: 197 Sequence: CDR3:	GCCTCTGGA CT not lound 104 GCCTCTGAA ACCTTTAGA TGTGCAAGAG	IGHV3-9_04	IGHU6_03 AGAGOCTATOAC TAGAATCAAGCT TTAGAAGGICTT	0.0337 GGACTCTGTGGA TAGACTACGACC ATTAICTAGAATC	6.8624 SCTCCCCGTATCTAC TCTAATACATGGAC AAGCTTAGACTACG	3.0568 3.0568 SAAACQAACAGTCTQA STCT00000CAGAGOQ ACCTCTAATACATQGA	N CAOCAGAOGACA ACCAC COTCTOOOGCAG	N COOCCTTGTATGAGT	28.8210 GTOCAAGAGACTAAA
_	Sequence: CDR3: 197 Sequence: CDR3: 160	GCCTCTGGA CT not lound 104 GCCTCTGAA ACCTTTAGA TGTGCAAGA 101	IGHV3-9_04 ITCACCGACAGT GGTCTTATTATC DACTAAAAACCT IGHV3-21_06	IGHJ6_03 AGAGOCTATGAC TAGAATGAAGGT TTTAGAAGGTCTT IGHJ5_02	0.0337 DGACTCTGTGGAI TAGACTACGACC ATTATCTAGAATC 0.0327	6.8624 3CTCCC0TATCTAC TCTAATACATOGACT AAGCTTAGACTACG 6.8651	3.0568 CAAACGAACADTCTGA DTCTGOOGCAGAGOG ACCTCTAATACATGGA 3.9548	N CAOCACOAOGACA ACCAC COTCTOOOGCAO n/a	N COOCCTTGTATGAGT AGOG	28.8210 GTOCAAGAGACTAAA 62.9956
	Sequence: CDR3: 197 Sequence: CDR3: 160 Sequence:	GCCTCTGGA CT not found 104 GCCTCTGAA ACCTTTTAGA TGTGCAAGA 101 OCCTCTGGA	IGHV3-9_04 IGHV3-9_04 IGHV3-9_04 GOTOTTATTATO GACTAAAAACCT IGHV3-21_06 ITCACCTTCACT SOCCOATTCACC	IGHJ6_03 AGAGGCTATGAG TAGAATGAAGGTCTT IGHJ5_02 AGGTTAGCATG OTCTCTOTGGC	0.0337 DOACTCT0T0CAC TTACACTACCACC ATTATCTACAACC 0.0327 ACT00CTCC0C CAGOCAACCCT	6.8624 CAGGCTCCAGGCAA CTCCCCCGTATCTAC TCTAATACATGGAC AAGCTTAGACTACG 6.8851 CCAGGCTCCAGGGAC	3.0568 3.0568 CAAACGAACAGTCTGA 31CTG0000CAGAG0G ACCTCTAATACATGGA 3.8548 300000CT00AGT000	N GAGCACACACGAC ACCAC CGTCTGGGGCAG INIB ITCTCATCCATTAC	IN COOCCTTGTATGAGT AGGG INTB CTAGTAGTAGTAGTATGA	28.8210 GTOCAAGAGACTAAA 62.9956 CATATACTACGCAGA
	Sequence: CDR3: 197 Sequence: CDR3: 160 Sequence: CDR3:	GCCTCTGGA CT not found 104 GCCTCTGAA ACCTTTTAGA TGTGCAAGA 101 GCCTCTGGA TGACTGAAGA not found	IGHV3-9_04 IGHV3-9_04 ITCACOGACAGT GGTCTTATTATC DACTAAAAACCT IGHV3-21_06 ITCACCTTCAGT SGCCGATTCACC	IGHJ6_03 AGAGGCTATGAG TAGAATCAAGCT TITTAGAGGTCTT IGHJ5_02 ADGTTTAGCATG OTCTCTGTGGC	0.0337 DOACTCT0T0CAC TTACACTACCACC ATTATCTACAACC 0.0327 IAACTOOGTCC0C CASOGAACCCT	6.9624 6.9624 SCTOCCOGTATCTAC TCTAATACATOGAC AAGCTTAGACTACG 8.9951 SCAGGCTCCAGGGAC	3.0568 SAAACGAACAGTCTGA SICTGGGGCAGAGGG ACCTCTAATACATGGA 3.9548 SOGGGCTGGAGTGGG	N GAGCAGAGGAGA CACCAC CGTCTGGGGCAG NS ITCTCATCCATTAC	IN CODECTTGATGAGT AGOO ING TAGTAGTAGTATGAGTATGA	28.8210 GTOCAAGAGACTAAA 62.9056 CATATACTACOCAGA
	Sequence: CDR3: 197 Sequence: CDR3: 160 Sequence: CDR3: 204	GCCTCTGGA CT not lound 104 OCCTCTGAA ACCTTTACA 101 00CTCTGAA TCACTGAAO not found 90	IGHV3-9_04 IGHV3-9_04 ITCACCDACAGT GGTCTTATTATC GGTCTTATTATC IGHV3-21_08 ITCACCTTCAGT SGCCGATTCACC IGHV3-20_02	IGHJ6_03 ADADOCTATOA TADAATCAAGCT TTTAGAATCAAGCT IGHJ5_02 ADGTTTADCATG GTCTCTGTGGCC Inche	0.0337 DGACTCT0T0GA0 TTAGACTACGACC ATTATCTAGACAATC 0.0327 IAACT0D0TCC00 CAGGGAACCCT 0.0318	6.8024 6.8024 30TCCC00514TCTAC TCTAATACATOGACI AAGCTTAGACTACOG 6.8951 30CAGGCTCCAGGGAA	30568 30568 34ACCAACGTCTCA 31CT0000CK0ACGC 3CCTCTAATACATOCA 39648 30000CTC0AGTC00 0.3472	N CAGCAGAGAGAGACA ACCAC CGTCTGGGGCAG Internet Contraction Internet Contraction	IN COOCCITGACTAC	28.8210 GTGCAAGAGACTAAA 82.9856 CATATACTACOCAGA 78.4722
	Sequence: CDRD: 197 Sequence: CDRD: 160 Sequence: CDRD: 234 Sequence:	GCCTCTGGA CT not found 104 GCCTCTGAA ACCTTTTAGA 101 GCCTCTGAA 101 GCCTCGAAQ not found 90 GCCTCTGAAQ	IGHV3-9_04 IGHV3-9_04 IGHV3-9_04 IGHV3-9_04 IGHV3-21_06 IGHV3-21_06 IGHV3-21_06 IGHV3-20_02 IGHV3-20_02 ITCACCTTGAT	IGHJ6_03 AGAGOCTATGAC TAGAATCAAGCT TITAGAATGAAGCT TITAGAGGTCTT/ IGHJ6_02 AGGTTTAGCATG IGHE SATTATGGCATG ATCTCCAAGAA	0.0337 004CTCTGT00A TradacTAcQAAC ATTATCTAGATC 0.0327 AACT00STCC0C CAG0GAACCCT 0.0318 ACCT00STCC0C CAACGCAAGAA	6.8024 6.8024 30.00001410140 10.000014101404 6.8051 30.400000000000000 6.8059 30.400000000000000000000000000000000000	30568 30568 34ACGAACAGTCTGA 31CT0000CAGAGOG ACCTCTAATACATGGA 39848 300000CT00AGT000 0.3472 400000CT00AGT000 0.3472 400000CT00AGT0000		IN COOCCTTGACTAC	28.8210 GTGCAAGAGACTAAA 62.9856 CATATACTACGCAGAG 78.4722 CACAGGGTTATGCAGA
	Sequence: CDR3: 197 Sequence: CDR3: 160 Sequence: CDR3: 234 Sequence: CDR3: CDR3: 234	GCCTCTGAN CT not found 104 GCCTCTGAN ACCTTTACM 101 GCCTCTGAN GCCTCTGAN Not found 90 GCCTCTGGAN GCCTCTGGAN 105 NG 90 NG	IGHV3-9_04 IGHV3-9_04 IGHV3-9_04 IGACTATATATO GACTAAAAACCT IGHV3-21_06 IGACTAAAAACCT IGHV3-20_02 IGACCTATGAC IGHV3-20_02	INGHUG_03 ADAGOCTATGAC TAGAACCAAGCT TTAGAACAAGCT TTTAGAACAAGCT TTTAGGOCT TTTAGGCATG ADGTTTAGGCATG ATTATGGCATG ATTATGGCATG	0.0337 GGACTCGTGGGA TGACTAGGAC ATTAICTAGAATC AACTOGTCCGC CAGGGAACCCT 0.0318 AGCTGGGTCCGC CAACGCCAAGAAM	6.8624 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC	39566 339568 3400 340 340 340 340 340 340 340 340 34	N GAGCAGAGGACA ACCAC COTE TGGOGCAG NB ITCTCATCCATTAC NB TCTCTCCTCCTATTAA NSCCGAGGMCACC	IN CODECTTGATGAGT AGGO INS TAGTAGTAGTAGTATGA INS TTGGAATGGTGGTAG SOCCTTGTATCACTG	28.8210 GTOCAAGAGACTAAA 62.9856 CATATACTACOCAGA 78.4722 CACAGGTTATGCAGA FOCGAGAGGACGGCC
	Sequence: CDRD: 197 Sequence: CDRD: 190 Sequence: CDRD: 234 Sequence: CDRD: 234 Sequence: CDRD: 234	GCCTCTGGA CT not found 104 GCCTCTGAA GCCTCTGAA GCCTCTGAA TGCCTGGAA TGCTGGAA NG NG SCCTCTGGAA NG NG SCCTCTGGAA NG NG SC SCCTCTGGAA NG SC SCCTCTGGAA NG SC SCCTCTGGAA NG SC SCCTCTGGAA NG NG SC SCCTCTGAA SCCTCTGGAA SCCTCTGGAA SCCTCTGGAA SCCTCTGGAA SCCTCTGGAA SCCTCTGGAA SCCTCTGGAA SCCTCTGGA SCCTCTGGA SCCTCTGAA SCCTCTGGA SCCTCTGAA SCCTCTGGA SCCTCTGAA SCCTCTGGA SCCTCTGGA SCCTCTGAA SCCTCTCTGAA SCCTCTCTCTGAA SCCTCTCTCTGAA SCCTCTCTCTGAA SCCTCTCTCTGAA SCCTCTCTCTGAA SCCTCTCTCTCTGAA SCCTCTCTCTCTCTGAA SCCTCTCTCTGAA SCCTCTCTCTGAA SCCTCTCTCTCTCTCTCTCTGAA SCCTCTCTCTCTCTCTCTCTCTCTCTCTCTCTCTCTCTC	IGHV3-9_04 IGHV3-9_04 ITCACOGACAGI GGTCTIATTATO GGTCTATTATO GGTCTATTATO GGTCTATTATO GGTCTATTATO GGTCCACTCACT IGHV3-20_02 ITCACCTTCATT GGCCGATTCACC IGHV3-74_01	IGHJ6_03 AGAGOCTATGAC TAGAATCAAGCT TITAGAAGGTCTT IGHJ6_02 AGGTTTAGGCATG AGGTTTAGGCATG ATTTAGGCATG ATTTAGGCATG ATTTAGGCATG	I 0.0337 004CTCTGTGGA TAGACTCTGTGGA ATTAICTAGAATC I 0.0327 IAACTODGTCCGC CAGGGAACCCT I 0.0316 CAACGCCAAGAAA I 0.0315	6.803 6.803 6.803 6.803 6.803 6.803 6.809 7.809	30566 330568 3442 53442 5457 554 3567 5000 5454 555 35848 33848 33848 33848 33848 33848 300000 5004 5000 4375 4554 35729	N GAGCAGAGGACA ACCAC COTCTODOGCAG NB TCTCATCCATTAC NB TCTCACCATTAC	IN COOCCTTGATGAGT AGOG Infa TAGTAGTAGTAGTATGAG Infa TAGTAGTAGTAGTAGTAG SOCCTTGTATCACTO Y	28.8210 GTOCAAGAGACTAAA IS2.9856 CATATACTACOCAGA IR6.4722 CACAGGTTATOCAGA ISCOGAGAGGACGGCC I99.5095
	Sequence: CDR3: 197 Sequence: CDR3: 160 Sequence: CDR3: 224 Sequence: CDR3: 224 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 254 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3: 255 Sequence: CDR3:	eccrcroseA cr netitound 104 eccrcroseA deccrcroseA regression for the eccrcroseA for for for for for for for for for for	IGHV3-9_94 IGHV3-9_94 ITCACOGACACT GGTCTTATTATC IGHV3-21_66 ITCACCTTCACT IGHV3-21_66 ITCACCTTCACT IGHV3-20_62 ITCACCTTCACT IGHV3-20_62 ITCACCTCACT IGHV3-20_62 IGHV3-20_62 IGHV3-20_62	IGHJ6_03 AGAGOCTATOA AGAGOCTATOA AGAGOCTATOA TAGAACAAGOC IGHJ5_02 AGGTTAGOCATO AGCTCTGTGGG IGCTCTGTGGG IGCTCTGTGGG IGCTCTGTGGG IGCTCTGGGATG ATCTCCCAGAGA ATCTCCCAGAGA	(a.0337 (a.0337 (a.0337 (a.0337 (a.027 (a	6.8024 6.8024 0crococotArtoraco actorocotArtoraco 6.8051 6.8058 0ccacotocacooa 6.8058 0ccacotocacooa 6.8058 0ccacotocacooa 8.8053 0catoctocotosaa 0catoctocotosaa 0catoctocotosaa	30568 30568 30668 3040004ACARTCT00 306000000404000 0.3472 0000000040040000 0.3472 00000000400000 0.3472 0.57269 000000010010100000 0000001001010000	N CAGCACAGGACA ACCAC COTCTGOOGCAG NS ITCTCATCCATTAC NS TCTCTGCTATTAA NSCCGAGGACAC Y TCTCACCATTAA NSCCGAGGACAC	IN CODECTTGATGAGT AGGG Inte TAGTAGTAGTAGTATGA Inte TAGTAGTAGTAGTATGA Inte TAGTAGTAGTAGTATGA Inte Inte TAGTAGTAGTATAGTG INTE INTE INTE INTE INTE INTE INTE INTE	28.8210 GTGCAAGAGACTAAA 82.9956 CATATACTACOCAGA 78.4722 CACAGOTATACCAGAG 79.65995 TGCAAGAGCATCGTG

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.

Ordering Information

Catalog # 7-500-0009

LymphoTrack® Software- MiSeq™

Products

Quantity Components

1 Software Package

LymphoTrack[®] Software Outputs

Run-, target- and sample-level results are provided as output files, enabling LIMS integration and an automated analytical workflow.

Example Data

Depicted are the top 10 sequences from a read summary generated by the LymphoTrack Enterprise Software. Highlighted columns represent fields that are unique to SHM analysis and include the SHM mutation rate and predictions pertaining to whether a sequence is in-frame or contains a premature stop codon. To learn more about the LymphoTrack Enterprise Software, please refer to the LymphoTrack Enterprise Software section.

Rank	Sequence	Length	Merge count	V-gene	J-gene	% Total reads	Cumulative %	Mutation rate partial V-gene (%)	In-frame (Y/N)	No stop codon (Y/N)	V-coverage
1	TTCTCGTGGTG	455	29603	IGHV4-59_08	IGHJ4_02	9.93	9.93	11.26	Y	Y	98.63
2	CTCGCCCTCCT	463	205	IGHV5-51_01	IGHJ4_02	0.07	9.99	0.00	Y	Y	99.66
3	GGTTTTCCTTG	484	201	IGHV3-7_01	IGHJ4_02	0.07	10.06	7.77	Y	Y	100.00
4	CTCGCCCTCCT	463	185	IGHV5-51_01	IGHJ5_02	0.06	10.12	6.08	Y	Y	99.32
5	CTCGCCCTCCT	469	170	IGHV5-51_01	IGHJ4_02	0.06	10.18	0.00	Y	Y	99.32
6	CTCGCCCTCCT	466	160	IGHV5-51_01	IGHJ4_02	0.05	10.23	0.00	Y	Y	99.66
7	CTGCTGCTGAC	460	159	IGHV2-5_10	IGHJ5_02	0.05	10.29	8.08	Y	Y	97.64
8	GGTTTTCCTTG	493	156	IGHV3-48_02	IGHJ6_02	0.05	10.34	3.72	Y	Y	98.99
9	CTCGCCCTCCT	334	153	IGHV5-51_02	IGHJ2_01	0.05	10.39	3.72	Y	Ν	27.70
10	CTCGCCCTCCT	334	152	IGHV5-51_02	IGHJ2_01	0.05	10.44	3.38	Υ	Ν	26.01

Sequence Frequency Graph



The stacked bar graph depicts the top 200 sequencing reads for a sample. Each individual colored bar represents a unique population of cells. Different colors stacked at the same point on the x-axis represent unique sequences that utilize the same V and J gene families. The amplicons of these products vary in sequence and may also vary in product size.

Ordering Information

Catalog # 7-500-0009 7-500-0007

Products

LymphoTrack® Software- MiSeq[™] LymphoTrack® Software - S5/PGM[™]

Quantity Components

- 1 Software Package
- 1 Software Package

LymphoTrack[®] Flex

Next-Generation Sequencing (NGS) Research Use Only (RUO) Assays

Designed to meet the growing demand for flexible, high-throughput and automation-ready solutions, LymphoTrack Flex Assays identify gene rearrangements and determine somatic hypermutation status.

LEARN MORE: catalog.invivoscribe.com

Introduction

The LymphoTrack[®] Flex Assays for hemato-oncology research represent highly sensitive and scalable NGS reagents that accurately identify gene rearrangement diversity, clonality and somatic hypermutation (for samples tested with IGH FR1 and Leader). These assays offer the ability to multiplex 96 indices per gene target using the highthroughput NextSeq[™]1000 platform. Increased DNA input sample volume enables compatibility with lower concentration samples without the need for additional concentration steps.

Sequencing results are analyzed by LymphoTrack Enterprise Software, a Linux-based bioinformatics application that allows analytical automation, custom visualizations and LIMS integration for downstream reporting. This highthroughput solution was devised from the bottom up to provide input volume flexibility, apply scalability and maximize workflow efficiencies.

Designed to be used with the Illumina® NextSeq™1000 System for high-throughput laboratories, LymphoTrack Flex Assay kits can process a total of 96 samples and controls. The supplementary Linux-based LymphoTrack Enterprise Software analysis generates output files that can be easily transferred into a LIMS for automated reporting.

Key Benefits

- » Scalability using up to 96 indices per target
- » Optimized Design provides flexibility for additional applications
- » Accelerated TAT with results in < 50 hours
- » Standardized Workflow reduces errors and enables automated lab environments
- » Dockerized Software seamlessly integrates with analytical pipelines for LIMS /reporting

Chapter Contents

- LymphoTrack Flex Assays
- 62 LymphoTrack Flex IGHV Leader Assay
- 62 LymphoTrack Flex IGH FR1 Assay
- **62** LymphoTrack Flex *TRG* Assay
- **64** LymphoTrack Enterprise Software

Research Use Only (RUO) assays are not for sale in global markets where equivalent CE-IVD assays are available and registered with the appropriate regulatory agencies. For more information regarding the CE-marked LymphoTrack® Dx Assays, please refer to the preceding pages.

LymphoTrack® Flex Assays

Assay Uses

The research use only (RUO) LymphoTrack Flex Assay portfolio is designed to detect rearranged gene sequences in human genomic DNA (gDNA) from peripheral blood, bone marrow and FFPE specimens. These assays provide reagents to generate sequencer-ready *IGH* FR1, IGHV Leader and/or *TRG* libraries for the Illumina[®] NextSeq[™]1000 System.

Summary and Explanation of the Test

The LymphoTrack Flex Assays for hemato-oncology research represent highly sensitive and scalable NGS library preparation reagents that accurately identify gene rearrangement diversity and clonality. Designed to be used with the Illumina® NextSeq™1000 System to provide flexibility and workflow efficiencies for high-throughput laboratories, as well as additional applications (such as MRD with DNA input volume flexibility). LymphoTrack Flex Assay kits can process a total of 96 samples and controls. The supplementary Linux-based LymphoTrack Enterprise Software analysis generates tabulated output files for sample interpretation that can be easily transferred into a LIMS for automated reporting.

Background

Identifying clonality and determining the frequency of immunoglobulin (Ig) and T-cell receptor (TCR) gene rearrangements are crucial in understanding hematologic malignancies^{1,7} These gene rearrangements serve as markers for clonal lymphocyte populations, aiding in the detection of blood cancers.²

During normal lymphocyte development, Ig and TCR genes undergo rearrangement to generate diverse antigen receptors.¹ Conversely, in lymphoid malignancies, such as B-cell and T-cell lymphomas, a single clone of lymphocytes proliferates uncontrollably, all sharing the same unique Ig or TCR gene rearrangement.^{3,4} Detecting these clonal rearrangements is essential for confirming malignancy.⁵

Traditional methods, such as PCR and Sanger sequencing, have been employed to detect Ig and TCR gene rearrangements for several decades.⁶ However, the introduction of next-generation sequencing (NGS) has enhanced sensitivity and specificity in clonality assessment, allowing a more comprehensive analysis of Ig and TCR gene rearrangements, facilitating the detection of clonal populations even in samples with low tumor cell content. Recent improvements in NGS clonality assays includes providing molecular genetic evidence of clonality in malignant or suspect lymphoproliferative disorders and facilitating minimal residual disease (MRD) assessment.⁷⁹

Specimen Requirement

50 ng of high-quality genomic DNA (200 ng is recommended).

References

- 1. Tonegawa, S. Nature. 1983 Apr;302: 575–581.
- 2. Teramo A et al. Nat Commun. 2022 Jun 8;13(1):3298.
- 3. Logan AC et al. Biol Blood Marrow Transplant. 2014 Sep;20(9):1307-13.
- 4. Crane, GM & Loghavi, S (Eds). Precision molecular pathology of aggressive B-cell lymphomas. Springer International Publishing AG, 2024.
- 5. Rosenquist R et al. Leukemia. 2017 Jul;31(7):1477-1481.
- 6. Langerak, T. (2014 February 11). Complex IG / TCR Clonality Targets: Patterns and Interpretation. 17th Euroclonality Workshop. Nijmegen, Netherlands.
- 7. Groenen PJTA et al. Front Oncol. 2023 Feb 07(13) 2234-2243.
- 8. Montserrat E. J Clin Oncol. 2005 May 1;23(13), pp. 2884–2885.
- 9. Faham M, et al. Blood. 2012 Dec 20;120(26):5173-80.



LymphoTrack[®] Flex IGHV Leader Assay Reagents

Component

IGH SHM (+) Control Polyclonal (-) Control IGHV Leader PCR1 PCR2

Purpose Positive Control

Target PCR

Index PCR

Negative Control

Units in Kit

- $1 \times 45 \ \mu L \ tube$
- 1 x 45 µL tube
- $2 \times 850 \ \mu L \ tubes$
- $2 \times 260 \ \mu L \ tubes$

LymphoTrack[®] Flex IGH FR1 Assay Reagents

Component

IGH (+) Control Polyclonal (-) Control IGH FR1 PCR1 PCR2 Purpose Positive Control Negative Control Target PCR Index PCR

Units in Kit 1 x 45 µL tube

1 x 45 µL tube

 $2 \times 850 \ \mu L \ tubes$

 $2 \times 260 \ \mu L \ tubes$

LymphoTrack[®] Flex TRG Assay Reagents

Component	Purpose	Units in Kit
TRG (+) Control	Positive Control	1 x 45 µL tube
Polyclonal (-) Control	Negative Control	1 x 45 µL tube
TRG PCR1	Target PCR	2 x 850 µL tubes
PCR2	Index PCR	2 x 260 µL tubes

Ordering Information

Catalog #	Products	Quantity Components
7100004	LymphoTrack [®] Flex IGHV Leader Assay	96 reactions
7100003	LymphoTrack [®] Flex <i>IGH</i> FR1 Assay	96 reactions
7100002	LymphoTrack® Flex <i>TRG</i> Assay	96 reactions
S100003	LymphoTrack® Enterprise Software	1 Software Package

These products are for Research Use Only (RUO). Not intended for diagnostic or therapeutic purposes.



LymphoTrack[®] Enterprise Software



Software Use

When used with our comprehensive LymphoTrack® NGS Solutions and LymphoTrack Flex Assays, the LymphoTrack Enterprise Software enables a fully automated pipeline for bioinformatics data analyses using the raw FASTQ files generated from Complete Genomics DNBSEQ-G99 and Illumina® MiSeq™, NextSeq™1000 and NovaSeq™6000 systems with LymphoTrack Assays as well as Illumina NextSeq™1000 with LymphoTrack Flex Assays. This Linux-based software is built for portability, flexibility and can be easily integrated into downstream reporting systems and/or LIMS.

Product Details

The research use only (RUO) LymphoTrack Enterprise Software package was designed to enable pipeline automation for high-throughput sequencing laboratories using LymphoTrack Assays and LymphoTrack MRD Solutions. This Docker image enables clonal analysis and MRD monitoring of single or multiple targets (IGHV Leader, *IGH* FR1, *IGH* FR2, *IGH* FR3, *IGK*, *TRG* and/or *TRB*) and also provides somatic hypermutation (SHM) status of samples run with the Leader and/or *IGH* FR1 assays.

The FASTQ raw data files generated from a sequencing run of samples prepared with LymphoTrack reagents can be automatically processed and output into tab-separated value (TSV) files using the LymphoTrack Enterprise Software. The analysis is initiated by either REST-API or command line interface. The data output can then be visualized in Microsoft® Excel® or imported into custom reports.

Compatible Assays

Assay

LymphoTrack® *IGHV* Leader Somatic Hypermutation Assays – MiSeq[™] LymphoTrack[®] *IGH* FR1/2/3 Assays – MiSeq[™] LymphoTrack[®] *IGK* Assays – MiSeq[™] LymphoTrack[®] *TRG* Assays – MiSeq[™]

LymphoTrack® *TRB* Assays – MiSeq™

Available Configuration(s) 8 indicies / 24 indices

8 indicies / 24 indices / 48 indices 8 indicies / 24 indices 8 indicies / 24 indices 8 indicies / 24 indices

Assay

LymphoTrack® Flex IGHV Assay LymphoTrack® Flex *IGH* FR1 Assay LymphoTrack® Flex *TRG* Assay Available Configuration(s)

96 reactions 96 reactions 96 reactions

This product is for Research Use Only (RUO). Not intended for diagnostic purposes.



LymphoTrack[®] Analysis Outputs

Run-, target- and sample-level results are provided as output files, enabling LIMS integration and an automated analytical workflow.

Sequencing Summary

Using the merged read summary interpretation is quick and easy.

	1 Easy	identifica	tion of spe	ecific types of ge	ene rearrang	ements su	ch as <i>IGHV</i> 3-:	21.			
Rank	Sequence	Length	Merge count	V-gene	J-gene	% Total reads	Cumulative %	Mutation rate partial V-gene (%)	In-frame (Y/N)	No stop codon (Y/N)	V-coverage
1	TTCTCGTGGTG	455	29603	IGHV4-59_08	IGHJ4_02	9.93	9.93	11.26	Y	Y	98.63
2	CTCGCCCTCCT	463	205	IGHV5-51_01	IGHJ4_02	0.07	9.99	0.00	Y	Y	99.66
3	GGTTTTCCTTG	484	201	IGHV3-7_01	IGHJ4_02	0.07	10.06	7.77	Y	Y	100.00
4	CTCGCCCTCCT	463	185	IGHV5-51_01	IGHJ5_02	0.06	10.12	6.08	Υ	Y	99.32
5	CTCGCCCTCCT	469	170	IGHV5-51_01	IGHJ4_02	0.06	10.18	0.00	Y	Y	99.32
6	CTCGCCCTCCT	466	160	IGHV5-51_01	IGHJ4_02	0.05	10.23	0.00	Y	Y	99.66
7	CTGCTGCTGAC	460	159	IGHV2-5_10	IGHJ5_02	0.05	10.29	8.08	Y	Y	97.64
8	GGTTTTCCTTG	493	156	IGHV3-48_02	IGHJ6_02	0.05	10.34	3.72	Y	Y	98.99
9	CTCGCCCTCCT	334	153	IGHV5-51_02	IGHJ2_01	0.05	10.39	3.72	Y	Ν	27.70
10	CTCGCCCTCCT	334	152	IGHV5-51_02	IGHJ2_01	0.05	10.44	3.38	Y 3	Ν	26.01

Identification of clonal sequences for follow up tracking with LymphoTrack MRD Software.

SHM rate and indicators to determine whether a clone is productive. Only provided for *IGHV* Leader and *IGH* FR1.

The read summary provides sequences from a sample ranked from most abundant to least prevalent. The total read count for individual sequences is provided and no independent analysis is required to determine V and J gene families and predictions for SHM when analyzing data from LymphoTrack *IGHV* Leader or *IGH* FR1 Assays. Additionally, the software provides raw and merged data in which reads that differ by 1-2 bp are automatically merged to account for possible sequencing errors and to improve the accuracy and ease of sample interpretation.

Ordering Information

Catalog # S100003 Products

LymphoTrack[®] Enterprise Software- MiSeq™

Quantity Components

1 Software Package

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.



Measurable Residual Disease (MRD) Solutions Research Use Only (RUO)

MRD testing with Next-Generation Sequencing is a proven tool that may be used to develop hematologic malignancy management strategies.

LEARN MORE: catalog.invivoscribe.com



Introduction

A number of investigators have described NGS-based approaches that have demonstrated success in detecting and monitoring MRD in Chronic Lymphocytic Leukemia (CLL), Acute Lymphoblastic Leukemia (ALL) and other lymphoid malignancies.^{1,2} LymphoTrack Assays are NGS-based deep sequencing assays that detect virtually all clonal rearrangements within targeted T-cell receptor (TCR) or immunoglobulin (Ig) antigen receptor loci. Once a specific rearrangement (the clonotype) has been identified, LymphoTrack assays can be used to track these clonotype populations to a sensitivity as low as 10⁻⁶. The LymphoTrack MRD solution includes assays, controls and software to facilitate objective MRD assessments over time.

B- or T-cell bundles may be purchased that include the LymphoTrack Assay, LymphoTrack MRD Bioinformatics Software, a LymphoTrack Low Positive Control and a LymphoQuant Internal Control. Each bundle facilitates the standardization of MRD testing by providing controls suitable for longitudinal MRD tracking with test sensitivity assurance. LymphoTrack MRD Software further simplifies clonal tracking due to rich sequence specific data analyses. This software enables objective monitoring of clonal populations by providing multiple functionalities to the user including project planning features and automated bioinformatics applications.

When monitoring MRD, a highly sensitive detection method such as NGS-based LymphoTrack may aid in the early detection of lymphoproliferative disease relapse. However, MRD test results are dependent on DNA amounts interrogated, as well as the confidence level of the test. Controls tracking MRD test sensitivity are thus necessary when reporting MRD test results. Designed for MRD testing, the LymphoTrack Low Positive Control confirms the sensitivity of respective LymphoTrack MRD runs match or exceed a 10⁻⁴ (or 1 in 10,000) level. Detection of the LymphoTrack Low Positive Control thus lessens false negative reporting concerns at 10⁻⁴, and is further necessary to report MRD negative results with confidence at 10⁻⁴.

Consistent use of a spike-in internal control enables objective monitoring of clonality over time and test standardization. Invivoscribe now offers both a LymphoQuant B-cell and T-cell Internal Control. Addition of a spike-in LymphoQuant B-cell or T-cell Internal Control to samples allows the software to convert percent clonal reads into estimated clonal cell equivalents.

Key Benefits

- » Complete solution for MRD
- » Ensures test sensitivity to enable confidence in reporting
- » Facilitates standardization of clonotype tracking
- » Longitudinal assessment of mutation status of *FLT3* ITD, *NPM1*, B- and T-cell gene rearrangements and somatic hypermutation (SHM)
- » LymphoTrack Assays formatted for both Illumina® and Thermo Fisher® NGS Platforms

Chapter Contents

- RUO MRD Assays

68 LymphoTrack MRD Software 70 *FLT3* ITD MRD Assay 72 *NPM1* MRD Assay

References 1. Leukemia. 27:1659-1665, 2013. 2. Blood. 120:5173-5180, 2012.

LymphoTrack® MRD Software

Software Use

The LymphoTrack MRD Bioinformatics Software package may be purchased to aid in the evaluation of treatment response in many lymphoid malignancies such as Acute Lymphoblastic Leukemia and Multiple Myeloma to minimize the risk of patients relapsing. The exceptional sensitivity and precision of NGS-based MRD tracking can accelerate clinical trials and drug development. This MRD software is intended to detect the presence of clonotype sequences within the output files generated using the Invivoscribe LymphoTrack Assays and accompanying LymphoTrack bioinformatics software; it is not intended to define the significance of these findings. Once a specific rearrangement sequence (the clonotype) has been identified in a primary sample, the MRD software enables streamlined tracking of clonal populations at a sensitivity of 10⁻⁴, or even lower limits provided sufficient DNA is tested. The MRD software can be used to Create, Save and Load projects to objectively track up to 5 clonal sequences to monitor for disease relapse and for use in drug development studies.

The provided software is composed of four distinct parts:

- 1. A Project tool that can be used to plan experiments with sufficient confidence based on read depth, replicate count, and DNA input. Projects can also be Saved and Loaded for use in subsequent time points.
- 2. A bioinformatics data analysis application.
- 3. A PDF Sample Report identifying the clonal sequence(s) if present, a summary of the degree of mismatches, calculations of the read frequency and the degree of confidence if clonal sequence(s) are not detected at various sensitivities.
- 4. A PDF Summary Report that will automatically generate longitudinal graphs of up to 5 clonal frequencies for a Subject and summary tables with estimates of clonal cell equivalents and clonal frequencies for all queried sequences if LymphoQuant Internal Controls is used.

MRD Project Planner

The Project Planner can be used to calculate the confidence of a true negative by adding replicate counts, resequencing counts, sequencing depth, and DNA input amount. The software assumes that the same sequencing depth and DNA input is used for each replicate.

প্ট Project Planner	- 🗆 X					
# of PCR Replicates: 3 # of Resequences:"	Results Resequences × Read Depth = Total Reads Per PCR Replicate: 500000 Total Reads Per PCR Replicate × PCR Replicates = Total Reads: 1500000					
1 Read Depth: 500000 Amount of DNA (co):	Sequence Not Detected % Confidence searched sequence was not detected Confidence at 1E-3 : 100.0% Confidence at 1E-4 : 100.0%					
2000 Calculate Confidence	Confidence at 1E-5 : 68.29% Confidence at 1E-6 : 1.68%					
* Typically a PCR Replicate is only s	equenced once	C				

Research Use Only (RUO) assays are not for sale in Europe and other global markets where equivalent CE-IVD assays are available and registered with the appropriate regulatory agencies.



68 Measurable Residual Disease (MRD) Solution | Research Use Only (RUO)

LymphoTrack MRD Sample Report

This report summarizes the overall call, i.e. if a clonotype was detected or not detected, the number of checked replicates, the total DNA input, total reads analyzed and the location of all output files. The values in the PDF report are also found in the generated text files.

Sequence #3 Details for Subject 1 for Collection/Timepoint: 2020/06/03								
Sequence Name	PCR Replicate(s)	Tot	Total Reads		Gene Target		MRD Result	
Third Clone	1	7629	40		IGH FR1		NOT DETECTED	
GCGTCTGGATTCATTTTCCCTAATGGACAGCCTGAGAGCCGAGGACACGGGTGTGTATAAGTGTGCGAGAAATAGCGTGATGGAAATGCGTGCTTCGTGGTGTCTGGGGGCATAG GAGCCAC								
PCR Replicate Details	Cumulative Target Read Count	Cumulat Total Re	tive % ads	Cumu Read	lative Lympł Count	noQuant	Clonal Frequency	
Exact Match	0	0%		100			0	
1 Mismatch	1	0.0002%		103		0		
2 Mismatch	1	0.0002%	105 0					
							0	
Detection Limit	% Confidence		Detection	Limit		% Confid	lence	
Detection Limit 1E-3	% Confidence 99.999%		Detection 1E-5	Limit		% Confic 53.761%	lence	

MRD Summary Graph



Note: The MRD Report may slightly differ from what is shown.

Ordering Information

Catalog # 7-500-0008

Products

LymphoTrack® MRD Software**

Quantity Components

1 Software Package

** MRD Software can be used to track sequences generated using either LymphoTrack® Assays formatted for either the Illumina® or Thermo Fisher® NGS platforms. MRD applications are for Research Use Only (RUO). Not intended for diagnostic purposes.



Measurable Residual Disease (MRD) Solution | Research Use Only (RUO) 69

FLT3 ITD MRD Assay

Assay Use

This Research Use Only (RUO) assay identifies internal tandem duplications (ITD) of the FLT3 gene.

Summary and Explanation of the Test

FLT3 ITD MRD is an NGS-based assay for use on the Illumina[®] MiSeq[™]. *FLT3* ITD MRD Assay tracks internal tandem duplication mutations of the *FLT3* gene to identify measurable residual disease in subjects following acute myeloid leukemia diagnosis (AML) and treatment. The assay is an amplicon-based approach which uses polymerase chain reaction to amplify the region of interest and next-generation sequencing to detect the region of interest.

The assay kit includes 24 master mixes, along with a ready to use positive and negative controls. Primers included in the master mix are designed with Illumina[®] adapters containing unique dual-indexes allowing for multiplexing at two levels. First, the multiplexing of up to 21 unique individuals in parallel in a single sequencing run, and second, the multiplexing of multiple targets in a single sequencing run.

The *FLT3* ITD MRD Assay is for use with the Illumina[®] MiSeq[™] platform. It was designed to produce sequencing data that can be analyzed using the *FLT3* ITD MRD Software package. The *FLT3* ITD MRD Software is a Dockerized Linux-based package built for portability, flexibility, and can easily be integrated into downstream reporting system and/or LIMS. It provides automated bioinformatic analyses of raw FASTQ files generated by the Illumina[®] MiSeq[™], when used with the *FLT3* ITD MRD Assay. Individual sample and aggregate sample data output is generated in under an hour in the form of .tsv (tab-separated value) files.

This deep-sequencing assay can detect internal tandem duplications up to 252 bp at an RUO validated allelic sensitivity of 5x10⁻⁵.

Background

FLT3 mutations are one of the most common mutations in AML, occurring in roughly ~30% of adult AML subjects.¹ The internal tandem duplication is the most common *FLT3* mutation and is present in about 25% of AML subjects.¹² *FLT3* ITD is considered a driver mutation at diagnosis and in relapsed AML,¹³ conferring a higher relapse rate and poor prognosis.¹ Since *FLT3* ITD mutations are present in 75% of relapsed subjects which were *FLT3* ITD-mutated at diagnosis,³ ITDs can serve as a biomarker to monitor for hematological disease progression and can guide therapeutic decisions.

Specimen Requirements

700 ng DNA isolated from whole blood or bone marrow.

References

- 1. Levis, M. Hematology Am Soc Hematol Educ Program (2013) 2013:220-6.
- 2. Daver, N. et al., *Leukemia* (2019) 33:299-312.
- 3. Krönke, J. et al., Blood (2013) 122(1):100-8.

Figure 1. FLT3 ITD MRD Assay Workflow



Reagents

Dual-indexed Master Mix	Index # (I7 Index ID)	Index # (I5 Index ID)
FLT3 ITD IA13	N701	N502
FLT3 ITD IA14	N702	N504
FLT3 ITD IA15	N703	N506
FLT3 ITD IA16	N704	N508
FLT3 ITD IA17	N705	N501
FLT3 ITD IA18	N706	N503
FLT3 ITD IA19	N707	N505
FLT3 ITD IA20	N708	N507
FLT3 ITD IA21	N709	N502
FLT3 ITD IA22	N710	N504
FLT3 ITD IA23	N711	N506
FLT3 ITD IA24	N712	N508
FLT3 ITD IB01	N702	N501
FLT3 ITD IB02	N703	N503
FLT3 ITD IB03	N704	N505
FLT3 ITD IB04	N705	N507
FLT3 ITD IB05	N706	N502
FLT3 ITD IB06	N707	N504
FLT3 ITD IB07	N708	N506
FLT3 ITD IB08	N709	N508
FLT3 ITD IB09	N710	N501
FLT3 ITD IB10	N711	N503
FLT3 ITD IB11	N712	N505
FLT3 ITD IB12	N701	N507
Controls	Concentration	Units in 96 Reaction Assay

FLT3 ITD MRD POS (+) FLT3 ITD MRD NEG (-) 20 ng/µL 20 ng/µL 2 x 500 µL tube 2 x 500 µL tube

Ordering Information

Catalog

14120019 S100005

Products

FLT3 ITD Assay FLT3 ITD MRD Software

Quantity Components

24 indices - 4 sequencing reactions each 1 Software Package

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.

NPM1 MRD Assay

Assay Use

This Research Use Only (RUO) assay identifies four base pair insertions of the NPM1 gene.

Summary and Explanation of the Test

NPM1 MRD is an NGS-based assay for use on the Illumina[®] MiSeq[™]. *NPM1* MRD Assay tracks four base pair insertions in the *NPM1* gene to identify measurable residual disease in subjects following acute myeloid leukemia diagnosis (AML) and treatment. The assay is an amplicon-based approach which uses polymerase chain reaction to amplify the region of interest and next-generation sequencing to detect the region of interest.

The assay kit includes 24 master mixes, along with a ready to use positive and negative controls. Primers included in the master mix are designed with Illumina[®] adapters containing unique dual-indexes allowing for multiplexing at two levels. First, the multiplexing of up to 21 unique individuals in parallel in a single sequencing run, and second, the multiplexing of multiple targets in a single sequencing run.

The *NPM1* MRD Assay is for use with the Illumina® MiSeq[®] platform was designed to produce sequencing data that can be analyzed using the *NPM1* MRD Software package. The *NPM1* MRD Software is a Dockerized Linux-based package built for portability, flexibility, and can easily be integrated into downstream reporting system and/or LIMS. It provides automated bioinformatic analyses of raw FASTQ files generated by the Illumina® MiSeq[®], when used with the *NPM1* MRD Assay. Individual sample and aggregate sample data output is generated in under an hour in the form of .tsv (tab-separated value) files.

This deep-sequencing assay can detect four base mutation insertions characterized as 'A', 'B', 'D', or 'Other' at an RUO validated allelic sensitivity of 5x10⁻⁵.

Background

One of the most commonly mutated genes in AML, *NPM1* gene mutations occur in about one-third of the cases of primary AML in adults.¹ While risk stratification of *NPM1*-mutated AML is ever evolving,² *NPM1* gene mutations are common,¹ AML specific,³ and its mutant transcripts are expressed at high levels,² making the mutation an ideal target for monitoring measurable residual disease² to guide therapeutic decisions.

Specimen Requirements

700 ng DNA isolated from whole blood or bone marrow.

References

- 1. Falini, B. et al., N Engl J Med. (2005) 352:254-266.
- 2. Falini, B. et al., Leukemia (2021) 35:3113-3126.
- 3. Liso, A. et al., Leukemia (2008) 22:1285-1289.

Figure 1. NPM1 MRD Assay Workflow


Reagents

Dual-indexed Master Mix	Index # (I7 Index ID)	Index # (I5 Index ID)
NPM1 MRD IA01	N701	N501
NPM1 MRD IA02	N702	N503
NPM1 MRD IA03	N703	N505
NPM1 MRD IA04	N704	N507
NPM1 MRD IA06	N706	N504
NPM1 MRD IA07	N707	N506
NPM1 MRD IA08	N708	N508
NPM1 MRD IA09	N709	N501
NPM1 MRD IA10	N710	N503
NPM1 MRD IA11	N711	N505
NPM1 MRD IA12	N712	N507
NPM1 MRD IA25	N705	N502
NPM1 MRD IB13	N702	N502
NPM1 MRD IB14	N703	N504
NPM1 MRD IB15	N704	N506
NPM1 MRD IB16	N705	N508
NPM1 MRD IB17	N706	N501
NPM1 MRD IB18	N707	N503
NPM1 MRD IB19	N708	N505
NPM1 MRD IB21	N710	N502
NPM1 MRD IB22	N711	N504
NPM1 MRD IB23	N712	N506
NPM1 MRD IB24	N701	N508
NPM1 MRD IB26	N709	N507
Controls	Concentration	Units in 96 Reaction Assay
NPM1 MRD POS (+)	20 ng/µL	2 x 500 µL tube
NPM1 MRD NEG (-)	20 ng/µL	2 x 500 µL tube

Ordering Information		
Catalog #	Products	Quantity Components
14160019	NPM1 MRD Assay	24 indices - 4 sequencing reactions each
S100004	NPM1 MRD Software	1 Software Package

For information regarding pricing and availability, please contact Invivoscribe at inquiry@invivoscribe.com

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.

S IdentiClone® Capillary - CE-IVD Assays

IdentiClone Assay kits are CE-marked *in vitro* diagnostic products.*

LEARN MORE: catalog.invivoscribe.com

Introduction

These kits are intended for PCR-based detection of clonal gene rearrangements and translocations in patients with suspected lymphoproliferations, using capillary electrophoresis methods.

The Invivoscribe CE-marked IdentiClone Assays represent a simple, cost-effective approach to PCR-based clonality testing. These standardized assays were carefully optimized testing positive and negative control samples using multiplex master mixes. Assay development was followed by extensive validation including the testing of more than 400 clinical samples using Revised European/American Lymphoma (REAL) Classification. Testing was performed at more than thirty prominent independent testing centers throughout Europe in a collaborative study known as the BIOMED-2 Concerted Action.¹ These PCR-based tests include standardized Instructions For Use (IFUs) with interpretation guidelines describing the use of master mixes and controls. The same thermal cycler program and similar detection methods are used for all IdentiClone kits to improve consistency, reduce human error, and facilitate cross-training.

Chapter Contents

B-cell Assays

- 76 IdentiClone® IGH + IGK B-Cell Clonality Assays
- **78** IdentiClone® *IGH* Gene Clonality Assays
- 80 IdentiClone® IGK Gene Clonality Assays
- 82 IdentiClone® IGL B-Cell Clonality Assays

T-cell Assays

- **84** IdentiClone® *TCRB* Gene Clonality Assays
- 86 IdentiClone® T-Cell Receptor Gama Gene Rearrangement Assay 2.0
- 88 IdentiClone® TCRD Gene Clonality Assays

Reference: 1. JJM van Dongen et al., (2003) *Leukemia* 17:2257-2317.

IdentiClone Assays are *in vitro* diagnostic products and are not available for sale or use within North America. For more information regarding the research use only reagents, please see the Gel & Capillary Research Use Only Assays section.



IdentiClone[®] *IGH* + *IGK* B-Cell Clonality Assays

Assay Description

The IdentiClone *IGH* + *IGK* B-Cell Clonality Assay is an *in vitro* diagnostic product intended for PCR-based detection of clonal immunoglobulin heavy chain and kappa light chain gene rearrangements in patients with suspect lymphoproliferations. Specifically, the *IGH* + *IGK* B-Cell Clonality Assay can be used to:

- Identify clonality in atypical lymphoproliferative disorders
- Support a differential diagnosis between reactive lesions and hematologic malignancies⁴
- Assign presumptive lineage in mature monoclonal lymphoproliferative disorders
- Identify tumor-specific markers (IGH and IGK gene rearrangements) for post-treatment monitoring
- Monitor and evaluate disease recurrence

Summary and Explanation of the Test

These kits include six master mixes to test for rearrangements of both *IGH* and *IGK*. The *IGH* Tube A, B, and C master mixes target the framework 1, 2, 3 regions (respectively) within the variable (V_H) region, and the joining (J_m) region of the immunoglobulin heavy chain locus. The *IGK* Tube A master mix targets the variable (V_K) and the joining (J_K) region. *IGK* Tube B master mix targets kappa deleting element (K_{de}) rearrangements with the variable (V_K) region and the intragenic J_K-C_K region. The resulting V_K-K_{de} and J_K-C_K intron-K_{de} rearrangements are a result of unsuccessful rearrangements retained by the B cell. For best sensitivity, it is recommended to test suspect B-cell malignancies for both *IGH* and *IGK*.³ The included Specimen Control Size Ladder master mix targets multiple genes and generates a series of amplicons of approximately 100, 200, 300, 400, and 600 base pairs to ensure that the quality and quantity of input DNA is adequate to yield a valid result.

Performance Characteristics

Data from two independent studies that tested more than 300 patient samples of varying types suggests the diagnostic accuracy of selected IdentiClone tests to be 96%. In both peer-reviewed studies, there were no clear false-positive results generated using the IdentiClone tests, and there was a high level of precision.² The clinico-histopathological diagnosis correlated well with PCR results in a higher number of patients when compared with Southern Blot (SB) results, as seen below:

PCR/SB concordance:1IGH:93% sensitivity / 92% specificityIGK:90% sensitivity / 90% specificity

PCR vs. SB analysis relative to histopathology and final diagnosis:

PCR/SB con	cordance:2	PCR sensitivity:	SB sensitivity
IGH + IGK:	85%	98%	39%

References

1. JJM van Dongen et al., Leukemia 17:2257-2317 (2003).

- 2. Y Sandberg et al., J. Mol. Diag. 7(4):495-503 (2005).
- 3. Van Krieken, JHJM et al., Leukemia 21:201 206 (2007).

See Euro Clonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.





Figure Legend: Simple representation of the organization of a rearranged immunoglobulin heavy chain (*IGH*) gene on chromosome 14q32.33 and the immunoglobulin kappa light chain gene on chromosome 2p11.2. Black arrows represent the relative positions of primers that target the conserved framework regions (FR1-3) and the downstream consensus J_{H} gene segments for *IGH* and the VK, JK, INTR and Kde primers which are included in the *IGK* master mix tubes.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0030 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0019 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0007 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
<i>IGH</i> Tube A	Framework 1 + J _H	1 x 1500 µL tube	10 x 1500 µL tubes
<i>IGH</i> Tube B	Framework 2 + J _H	1 x 1500 µL tube	10 x 1500 µL tubes
<i>IGH</i> Tube C	Framework 3 + JH	1 x 1500 µL tube	10 x 1500 µL tubes
IGK Tube A	Vƙ-Jƙ	1 x 1500 µL tube	10 x 1500 µL tubes
IGK Tube B	Vƙ-K _{de} , Intron-K _{de}	1 x 1500 µL tube	10 x 1500 µL tubes
Specimen Control Size Ladder	Multiple Genes	1 x 1500 µL tube	10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
9-100-0031	IdentiClone® IGH + IGK B-Cell Clonality Assay - ABI Fluorescence Detection	33 reactions
9-100-0041	IdentiClone® IGH + IGK B-Cell Clonality Assay MegaKit - ABI Fluorescence Detection	330 reactions

These are in vitro diagnostic products and are available in regions that accept CE-IVD products.



IdentiClone® IGH Gene Clonality Assays

Assay Description

The IdentiClone *IGH* Gene Clonality Assay is an *in vitro* diagnostic product intended for PCR-based detection of clonal immunoglobulin heavy chain gene rearrangements in patients with suspect lymphoproliferations.

Specifically, the IGH Gene Clonality Assay can be used to:

- Identify clonality in atypical lymphoproliferative disorders
- Support a differential diagnosis between reactive lesions and hematologic malignancies⁴
- Assign presumptive lineage in mature monoclonal lymphoproliferative disorders
- Identify tumor-specific markers (IGH gene rearrangements) for post-treatment monitoring
- Monitor and evaluate disease recurrence

Summary and Explanation of the Test

These test kits include six master mixes. The *IGH* Tube A, B, and C master mixes target the framework 1, 2, and 3 regions (respectively) within the variable (V_{μ}) region and the joining (J_{μ}) region of the immunoglobulin heavy chain locus. The *IGH* Tube D and E master mixes target the diversity and joining regions. The Specimen Control Size Ladder master mix targets multiple genes and generates a series of amplicons of approximately 100, 200, 300, 400, and 600 base pairs to ensure that the quality and quantity of input DNA is adequate to yield a valid result.

Performance Characteristics

Data from two independent studies that tested more than 300 patient samples of varying types suggests the diagnostic accuracy of selected IdentiClone tests to be 96%. In both peer-reviewed studies, there were no clear false-positive results generated using the IdentiClone tests, and there was a high level of precision.² The clinico-histopathological diagnosis correlated well with PCR results in a higher number of patients when compared with Southern Blot (SB) results, as seen below:

PCR/SB concordance:¹ IGH: 93% sensitivity / 92% specificity

PCR vs. SB analysis relative to histopathology and final diagnosis:

PCR/SB concordance: ²	PCR sensitivity:	SB sensitivity:
IGH + IGK: 85%	98%	39%

References

JJM van Dongen et al., *Leukemia* 17:2257-2317 (2003).
 Y Sandberg et al., *J. Mol. Diag.* 7(4):495-503 (2005).

See Euro Clonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.





Figure Legend: Simple representation of the organization of a rearranged immunoglobulin heavy chain gene on chromosome 14q32.33. Black arrows represent the relative positions of primers that target the conserved framework (FR1-3) and diversity (DH1-7) regions, and the downstream consensus J_{H} gene segments. The amplicon products generated from each of these regions can be differentially detected when fluorescent primer sets are used with capillary electrophoresis instruments that employ differential fluorescence detection.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0030 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0019 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0024 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0008 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
Manada an Milana a	Townsh	Harles In Assess	
Master Mixes	larger	Units in Assay	Units in Assay Megakit
IGH Tube A	Framework 1 + JH	1 x 1500 µL tube	10 x 1500 µL tubes
IGH Tube A IGH Tube B	Framework 1 + J _H Framework 2 + J _H	1 x 1500 μL tube 1 x 1500 μL tube	10 x 1500 μL tubes 10 x 1500 μL tubes
IGH Tube A IGH Tube B IGH Tube C	Framework 1 + JH Framework 2 + JH Framework 3 + JH	1 x 1500 μL tube 1 x 1500 μL tube 1 x 1500 μL tube 1 x 1500 μL tube	10 x 1500 µL tubes 10 x 1500 µL tubes 10 x 1500 µL tubes 10 x 1500 µL tubes
IGH Tube A IGH Tube B IGH Tube C IGH Tube D	Framework 1 + J _H Framework 2 + J _H Framework 3 + J _H D _H 1-6 + J _H	1 x 1500 μL tube 1 x 1500 μL tube 1 x 1500 μL tube 1 x 1500 μL tube 1 x 1500 μL tube	10 x 1500 μL tubes 10 x 1500 μL tubes 10 x 1500 μL tubes 10 x 1500 μL tubes 10 x 1500 μL tubes
IGH Tube A IGH Tube B IGH Tube C IGH Tube D IGH Tube E	Framework 1 + J _H Framework 2 + J _H Framework 3 + J _H D _H 1-6 + J _H D _H 7 + J _H	1 x 1500 μL tube 1 x 1500 μL tube	0 x 1500 μL tubes 10 x 1500 μL tubes
IGH Tube A IGH Tube B IGH Tube C IGH Tube D IGH Tube E Specimen Control Size Ladder	Framework 1 + JH Framework 2 + JH Framework 3 + JH DH1-6 + JH DH7 + JH Multiple Genes	1 x 1500 μL tube 1 x 1500 μL tube	10 x 1500 µL tubes 10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
9–101–0061	IdentiClone® IGH Gene Clonality Assay - ABI Fluorescence Detection	33 reactions
9–101–0081	IdentiClone® IGH Gene Clonality Assay MegaKit - ABI Fluorescence Detection	330 reactions

These are in vitro diagnostic products and are available in regions that accept CE-IVD products.



IdentiClone® IGK Gene Clonality Assays

Assay Description

The IdentiClone *IGK* Gene Clonality Assay is an *in vitro* diagnostic product intended for PCR-based detection of clonal immunoglobulin kappa light chain gene rearrangements in patients with suspect lymphoproliferations.

Specifically, the IGK Gene Clonality Assay can be used to:

- · Identify clonality in atypical lymphoproliferative disorders
- Support a differential diagnosis between reactive lesions andhematologic malignancies
- Assign presumptive lineage in mature monoclonal lymphoproliferative disorders
- Identify tumor-specific markers (IGK and IGK-Kde rearrangements) for post-treatment monitoring
- Monitor and evaluate disease recurrence

Summary and Explanation of the Test

These test kits include three master mixes. The *IGK* Tube A master mix targets the variable (Vk) and the joining (Jk) regions of the immunoglobulin kappa light chain locus, whereas the *IGK* Tube B master mix targets kappa deleting element (K_{de}) rearrangements with the Vk regions and the intragenic Jk-Ck regions. The Vk-K_{de} and Jk-Ck intron-K_{de} rearrangements are a result of unsuccessful rearrangements retained by the B cell. The third master mix, the Specimen Control Size Ladder, targets multiple genes and generates a series of amplicons of approximately 100, 200, 300, 400, and 600 base pairs to ensure that the quality and quantity of input DNA is adequate to yield a valid result.

Performance Characteristics

Data from two independent studies that tested more than 300 patient sample of varying types suggests the diagnostic accuracy of selected IdentiClone tests to be 96%. In both peer-reviewed studies, there were no clear false-positive results generated using the IdentiClone tests, and there was a high level of precision.² The clinico-histopathological diagnosis correlated well with PCR results in a higher number of patients when compared with Southern Blot (SB) results, as seen below:

PCR/SB concordance:¹ IGK: 90% sensitivity / 90% specificity

PCR vs. SB analysis relative to histopathology and final diagnosis:

PCR/SB concordance:2PCR sensitivity:SB sensitivity:IGH + IGK:85%98%39%

References

JJM van Dongen et al., *Leukemia* 17:2257-2317 (2003).
 Y Sandberg et al., *J. Mol. Diag.* 7(4):495-503 (2005).

See Euro Clonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.





IGK tube A: 6 V κ primers + 2 J κ primers



IGK tube B: 6 V κ primers and INTR primer + 1 Kde primer

Figure Legend: Schematic diagram of the immunoglobulin kappa light chain gene complex on chromosome 2p11.2. Shown are the relative positions and orientations for the Vk-Jk, and K_{de} primers, which are included in the *IGK* master mix tubes.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0007 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
Master Mixes IGK Tube A	Target Vk - Jk	Units in Assay 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes
Master Mixes IGK Tube A IGK Tube B	Target Vƙ - Jƙ Vƙ-K _{de} , Intron-K _{de}	Units in Assay 1 x 1500 µL tube 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes 10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
9-102-0021	IdentiClone® IGK Gene Clonality Assay - ABI Fluorescence Detection	33 reactions
9-102-0031	IdentiClone® IGK Gene Clonality Assay MegaKit - ABI Fluorescence Detection	330 reactions

These are in vitro diagnostic products and are available in regions that accept CE-IVD products.



IdentiClone[®] IGL Gene Clonality Assays

Assay Description

The IdentiClone *IGL* Gene Clonality Assay is an *in vitro* diagnostic product intended for PCR-based detection of clonal immunoglobulin lambda light chain gene rearrangements in patients with suspect lymphoproliferations and can be used to:

Specifically, the IdentiClone IGL Gene Clonality Assays can be used to:

- Identify clonality in atypical lymphoproliferative disorders
- Support a differential diagnosis between reactive lesions and hematologic malignancies
- Assign presumptive lineage in mature monoclonal lymphoproliferative disorders
- Identify tumor-specific markers (IGL gene rearrangements) for post-treatment monitoring
- Monitor and evaluate disease recurrence

Summary and Explanation of the Test

These test kits include two master mixes. The *IGL* Tube master mix targets the variable (V λ) region and the joining (J λ) region of the immunoglobulin lambda light chain gene locus (*IGL*). The Specimen Control Size Ladder targets multiple genes and generates a series of amplicons of approximately 100, 200, 300, 400, and 600 base pairs to ensure that the quality and quantity of input DNA is adequate to yield a valid result.

Performance Characteristics

Data from two independent studies that tested more than 300 patient samples of varying types suggests the diagnostic accuracy of selected IdentiClone tests to be 96%. In both peer-reviewed studies, there were no clear false-positive results generated using the IdentiClone tests, and there was a high level of precision.²The clinico-histopathological diagnosis correlated well with PCR results in a higher number of patients when compared with Southern Blot (SB) results, as seen below:

PCR/SB concordance:¹ IGL: 86% sensitivity / 92% specificity

References

1. JJM van Dongen et al., Leukemia 17:2257-2317 (2003).

2. Y Sandberg et al., J. Mol. Diag. 7(4):495-503 (2005).

See Euro Clonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.





IGL tube: 2 V λ primers + 1 J λ primer

Figure Legend: Schematic diagram of the immunoglobulin lambda light chain gene complex on chromosome 22q11.2. Shown are the relative positions and orientations for the Vλ and Jλ primers, which are included in the *IGL* master mix tube. The two Vλ primers only target Vλ1, 2, and 3 because these three V families cover approximately 70% of rearrangeable Vλ gene segments, and approximately 90% of all *IGL* gene rearrangements involve these three families. Similarly, the single Jλ primer only targets Jλ1, 2, and 3 because these three J segments are involved in 98% of all *IGL* gene rearrangements.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0010 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0029 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
<i>IGL</i> Tube	Vλ –Jλ	1 x 1500 µL tube	10 x 1500 µL tubes
Specimen Control Size Ladder	Multiple Genes	1 x 1500 µL tube	10 x 1500 µL tubes

Ordering I	nformation	
Catalog #	Products	Quantity
9-103-0011	IdentiClone® IGL Gene Clonality Assay - ABI Fluorescence Detection	33 reactions

These are in vitro diagnostic products and are available in regions that accept CE-IVD products.



IdentiClone® TCRB Gene Clonality Assays

Assay Description

The IdentiClone *TCRB* Gene Clonality Assay is an *in vitro* diagnostic product intended for PCR-based detection of clonal T-cell receptor beta chain gene rearrangements in patients with suspect lymphoproliferations.

Specifically, the TCRB Gene Clonality Assay can be used to:

- Identify clonality in suspect lymphoproliferations
- Support a differential diagnosis between reactive lesions and T-cell and some immature B-cell malignancies
- Determine lineage involvement in mature lymphoproliferative disorders
- Monitor and evaluate disease recurrence

Summary and Explanation of the Test

These kits include four master mixes. *TCRB* Tubes A and B target framework regions within the variable region, and the joining region (V β) of the *TCR* beta chain locus. *TCRB* Tube C targets the diversity and joining (J β) regions of the *TCR* beta chain locus. The Specimen Control Size Ladder master mix included targets multiple genes and generates a series of amplicons of approximately 100, 200, 300, 400, and 600 base pairs to ensure that the quality and quantity of input DNA is adequate to yield a valid result.

Performance Characteristics

Data from two independent studies that tested more than 300 patient samples of varying types suggests the diagnostic accuracy of selected IdentiClone tests to be 96%. In both peer-reviewed studies, there were no clear false-positive results generated using the IdentiClone tests, and there was a high level of precision.² The clinico-histopathological diagnosis correlated well with PCR results in a higher number of patients when compared with Southern Blot (SB) results, as seen below:

PCR/SB concordance:¹ TRB: 86% sensitivity / 98% specificity

PCR vs. SB analysis relative to histopathology and final diagnosis:

	PCR/SB concordance: ²	PCR sensitivity:	SB sensitivity:
TRB:	85%	96%	35%

References

1. JJM van Dongen et al., *Leukemia* 17:2257-2317 (2003).

2. Y Sandberg et al., J. Mol. Diag. 7(4):495-503 (2005).

See Euro Clonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.





Figure Legend: Simplified diagram of a representative rearranged T-cell receptor beta gene showing the approximate placement of the upstream and downstream DNA primers. The numbers of primers and their specificity are listed for Master Mix Tubes A, B, and C.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0009 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0004 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0021 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
<i>TCRB</i> Tube A	Multiple Vβ + Jβ1/2	1 x 1500 µL tube	10 x 1500 µL tubes
<i>TCRB</i> Tube B	Multiple Vβ + Jβ2	1 x 1500 µL tube	10 x 1500 µL tubes
<i>TCRB</i> Tube C	Multiple Dβ + Jβ1/2	1 x 1500 µL tube	10 x 1500 µL tubes
Specimen Control Size Ladder	Multiple Genes	1 x 1500 µL tube	10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
9-205-0011	IdentiClone® TCRB Clonality Assay - ABI Fluorescence Detection	33 reactions
9-205-0021	IdentiClone® TCRB Clonality Assay MegaKit - ABI Fluorescence Detection	330 reactions

These are in vitro diagnostic products and are available in regions that accept CE-IVD products.



IdentiClone® T-Cell Receptor Gamma Gene Rearrangement Assay 2.0

Assay Description

The IdentiClone T-Cell Receptor Gamma Gene Rearrangement Assay 2.0 is an *in vitro* diagnostic product intended for PCR-based detection of clonal T-cell receptor gamma chain gene rearrangements in patients with suspect lymphoproliferations.

Specifically, the T-Cell Receptor Gamma Gene Rearrangement Assay 2.0 can be used to identify clonality in suspect lymphoproliferations.

Summary and Explanation of the Test

Rearrangements of the antigen receptor genes occur during ontogeny in B- and T-lymphocytes. These gene rearrangements generate products that are unique in length and sequence. Polymerase chain reaction (PCR) assays can be used to identify lymphocyte populations derived from a single cell by detecting the unique V-J gene rearrangements present within these antigen receptor loci.¹

This assay allows for amplification of the *TRG* region with fluorescent labeled primers, yielding products that can be grouped under a single Gaussian distribution when separated by size using capillary electrophoresis. In addition, the product size facilitates increased success when testing FFPE samples. The included analysis algorithm aids in the interpretation of data and identification of significant clonal peaks. Presence or absence of molecular clonality can support the differential diagnosis of reactive lesions and certain B- and T-cell malignancies, provided that the results are interpreted in the context of all available clinical, histological, and immunophenotypic data.

Performance Characteristics

To assess the performance of the *TCRG* 2.0 Assay, testing was performed on cell lines with known clonal rearrangements followed by testing on previously sequenced clinical samples.

When used in combination with the provided *TCRG* Algorithm worksheet, the assay was capable of detecting DNA from 6 control cell lines (200 ng/ μ L) diluted into polyclonal tonsil DNA (200 ng/ μ L) at 5% (v/v).

Furthermore, the performance of the *TCRG* 2.0 Assay was evaluated on clinical samples for which the T-cell receptor gamma gene rearrangement status had been identified by Roche 454 sequencing. For the 7 samples that had been identified as clonal by sequencing, the *TCRG* 2.0 assay had 100% concordance. For the 12 samples that were either negative for a clonal event or were oligoclonal, concordance of the *TCRG* 2.0 assay was 75%. Sample types included peripheral blood, bone marrow, and formalin-fixed, paraffin embedded (FFPE) tissue.

Always interpret the results of molecular clonality tests in the context of clinical, histological and immunophenotypic data.

References

Miller JE, Wilson SS, Jaye DJ, and Kronenberg M. Mol. Diag. 1999, 4(2):101-117.
 Armand, Marine et al. HemaSphere, 2019;3:3.



This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936. The performance of this assay was reviewed and validated by the EuroClonality/BIOMED-2 Group.²





Figure Legend: Simple representation of the organization of the T-cell receptor gamma gene on chromosome 7p14. Black arrows represent the relative positions of primers that target the variable region genes and the downstream joining region gene segments that are involved in rearrangements in T-cell lymphomas. The downstream primers are fluorescently labeled through the incorporation of a 6FAM fluorophore. The amplicon products generated from these rearrangements are detected by capillary electrophoresis.

Reagents

Concentration	Units in Assay	Units in Assay MegaKit
50 μg/mL	1 x 50 µL tube	5 x 50 µL tube
50 μg/mL	1 x 50 µL tube	5 x 50 µL tube
Target	Units in Assay	Units in Assay MegaKit
Vy1-Vy11 + Jy1/Jy2, JyP, JyP1/JyP2	1 x 1500 µL tube	10 x 1500 µL tubes
Multiple Cones	1 x 1500 ul tubo	10 x 1500 ul tubes
	Concentration 50 µg/mL 50 µg/mL Target Vy1-Vy11 + Jy1/Jy2, JyP, JyP1/JyP2	ConcentrationUnits in Assay50 μg/mL1 x 50 μL tube50 μg/mL1 x 50 μL tubeTargetUnits in AssayVy1-Vy11 + Jy1/Jy2, JyP, JyP1/JyP21 x 1500 μL tubeMultiple Genes1 x 1500 μL tube

Ordering Information

Catalog #	Products	Quantity
9-207-0101	IdentiClone® T-Cell Receptor Gamma Gene Rearrangement Assay 2.0 - ABI Fluorescence Detection	33 reactions
9-207-0111	IdentiClone® T-Cell Receptor Gamma Gene Rearrangement Assay MegaKit 2.0 - ABI Fluorescence Detection	330 reactions

These are in vitro diagnostic products and are available in regions that accept CE-IVD products.



IdentiClone® TCRD Gene Clonality Assays

Assay Description

The IdentiClone *TCRD* Gene Clonality Assay is an *in vitro* diagnostic product intended for PCR-based detection of clonal T-cell receptor delta chain gene rearrangements in patients with suspect lymphoproliferations and can be used to:

- Identify clonality in suspect lymphoproliferations
- Support a differential diagnosis between reactive lesions and T-cell and some immature B-cell malignancies
- Determine lineage involvement in mature lymphoproliferative disorders
- Monitor and evaluate disease recurrence

Summary and Explanation of the Test

These test kits include two master mixes. The *TCRD* tube targets the framework regions within the variable region, the diversity region, and the joining region of the T-cell receptor delta chain locus (*TRD*, formerly known as *TCRD*). The Specimen Control Size Ladder master mix targets multiple genes and generates a series of amplicons of approximately 100, 200, 300, 400, and 600 base pairs to ensure that the quality and quantity of input DNA is adequate to yield a valid result.

Performance Characteristics

Data from an independent, peer-reviewed study suggests the diagnostic accuracy of selected IdentiClone tests to be 96%. There were no clear false-positive results generated using the IdentiClone tests, and there was a high level of precision. The clinico- histopathological diagnosis correlates well with PCR results in a higher number of patients when compared with Southern Blot (SB) results, as seen below:

PCR/SB concordance:¹ TRD: 83% sensitivity / 95% specificity

Reference

1. JJM van Dongen et al., Leukemia 17:2257-2317 (2003).

EuroClonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.





TCRD tube: 6 V δ and 1 D δ 2 primers + 4 J δ and 1 D δ 3 primers

Figure Legend: Simplified diagram of a representative rearranged T-cell receptor delta gene showing the approximate placement of the upstream and downstream DNA primers. The numbers of primers and their specificity are listed for the TRD Tube Master Mix tube.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0021 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
Master Mixes TCRD Tube	Target Multiple Võ + Dõ +Jõ	Units in Assay 1 x 1500 µL tube	Units in Assay MegaKit
Master Mixes <i>TCRD</i> Tube Specimen Control Size Ladder	Target Multiple Võ + Dõ +Jõ Multiple Genes	Units in Assay 1 x 1500 µL tube 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes 10 x 1500 µL tubes

Ord	erina	Inform	ation
0.0			

Catalog #	Products	Quantity
9-206-0011	IdentiClone® TCRD Gene Clonality Assay - ABI Fluorescence Detection	33 reactions
9-206-0021	IdentiClone® TCRD Gene Clonality Assay MegaKit - ABI Fluorescence Detection	330 reactions

These are in vitro diagnostic products and are available in regions that accept CE-IVD products.



LeukoStrat[®] CDx FLT3 Mutation Assay

Capillary Detection – CE 2797 IVD Assays and IVD Companion Diagnostic Assays

LeukoStrat CDx *FLT3* Mutation Assay Kits are *in vitro* diagnostic products intended for PCR-based detection of *FLT3* activating mutations (ITD, TKD) in patients with acute myelogenous leukemia (AML) using capillary electrophoresis methods.

LEARN MORE: catalog.invivoscribe.com

For Legal Notices related to this product, visit: https://invivoscribe.com/legal-notices/

Introduction

All kits include PCR master mixes and controls for ITD and TKD detection, along with Instructions For Use (IFUs). Master mixes are composed of a buffered magnesium chloride solution, deoxynucleotides, and primers targeting the gene segments of interest. These assay master mixes are complete other than Taq DNA polymerase.

Companion Diagnostics

The LeukoStrat CDx *FLT3* Mutation Assays are companion diagnostics, and provide a complete solution with software automated PDF reports including mutant identification, mutant to wild type signal ratios, and clinical interpretation.

Chapter Contents

CDx Assays

92 LeukoStrat CDx *FLT3* Mutation Assay (IVD) - United States
94 LeukoStrat CDx *FLT3* Mutation Assay (CE 2797 IVD)
96 LeukoStrat CDx *FLT3* Mutation Assay - Japan

LeukoStrat CDx *FLT3* Mutation Assays are *in vitro* diagnostic products. CE 2797 IVD products are not available for sale or use within North America. IVD products are available for sale or use where registered. Refer to product page for respective detail. For more information regarding the research use only reagents, please see the Gel and Capillary Research Use Only Assays section.



LeukoStrat[®] CDx FLT3 Mutation Assay

FDA Approved Companion Diagnostic

Available in the USA - In Vitro Diagnostic Kit

FDA approved assay for assessment of acute myeloid leukemia (AML) patients eligible for treatment with RYDAPT[®] (midostaurin), XOSPATA[®] (gilteritinib), or VANFLYTA[®] (quizartinib), now available as US distributed kit. This *FLT3* companion diagnostic includes reagents along with software that identifies ITD and TKD mutations, generates mutant/wildtype signal ratios, and predicts response to gilteritinib, midostaurin, and quizartinib.

Clinical Significance of *FLT3* Mutation Status: Each year approximately 21,000 patients in the United States are diagnosed with AML. Of those diagnosed with AML, ~1 out of 3 are expected to have presence of *FLT3* mutations, (*FLT3*mut+). Since *FLT3*mut+ AML is clinically actionable, stratification of AML patients by testing for *FLT3* mutation status has become a standard of care.

Intended Use

The LeukoStrat CDx *FLT3* Mutation Assay is a PCR-based in vitro diagnostic test designed to detect internal tandem duplication (ITD) and tyrosine kinase domain (TKD) mutations D835 and I836 in the *FLT3* gene in genomic DNA extracted from mononuclear cells obtained from peripheral blood or bone marrow aspirates of patients diagnosed with acute myelogenous leukemia (AML).

The LeukoStrat CDx FLT3 Mutation Assay is used as an aid in the assessment of patients with AML for whom RYDAPT® (midostaurin) treatment is being considered.

The LeukoStrat CDx *FLT3* Mutation Assay is used as an aid in the assessment of patients with AML for whom XOSPATA® (gilteritinib) treatment is being considered.

The LeukoStrat CDx *FLT3* Mutation Assay is used as an aid in the assessment of patients with *FLT3*-ITD+ AML for whom VANFLYTA® (quizartinib) treatment is being considered.

The test is for use on the 3500xL Dx Genetic Analyzer.

Summary and Explanation of the Test

Acute myelogenous leukemia (AML) in general has a poor prognosis. Many studies in AML have shown that the presence of *FLT3* activating mutations portends a poor prognosis making it an attractive target for treatment.¹²

The LeukoStrat CDx *FLT3* Mutation Assay targets regions of the *FLT3* gene to identify internal tandem duplication (ITD) mutations and tyrosine kinase domain (TKD) mutations, such as the D835 and I836 mutations.

The LeukoStrat CDx *FLT3* Mutation Assay includes reagents, software and procedures for isolating mononuclear cells and extracting DNA from patient peripheral blood or bone marrow specimens to determine if *FLT3* mutations are present.

DNA is amplified via PCR and the amplicons are detected via capillary electrophoresis. *FLT3* mutation status is determined by the LeukoStrat CDx *FLT3* Software. A *FLT3* ITD and/or TKD mutation is reported as Positive if the mutant:wild-type signal ratio meets or exceeds the clinical cutoff of 0.05.

Method Description

ITD Mutations of FLT3

The LeukoStrat CDx *FLT3* Mutation Assay uses fluorescently labeled primers that are in and around the JM region. Wild-type *FLT3* alleles will amplify and produce a product at 327±1 bp as measured by this assay, while alleles that contain ITD mutations will produce a product that exceeds 327±1 bp (see Figure, right).

TKD Mutations of FLT3

The LeukoStrat CDx *FLT3* Mutation Assay uses primers that lie on either side of the TKD region. The *FLT3* target region is amplified using PCR and then an EcoRV restriction digest is performed. Wild-type alleles of the *FLT3* gene yield digestion products of 79±1 bp whereas mutant alleles yield products of 125±1 bp or 127±1 bp from the original undigested amplicon product of 145±1 bp or 147±1 bp, as measured by this assay (see Figure,



92 Capillary Detection | CE 2797 IVD Assays and IVD Companion Diagnostic Assays

Method Schematic: FLT3 ITD and TKD Mutant Detection



Depicted is a representation of the *FLT3* juxtamembrane (JM) region (TM = transmembrane) and the activating loop of the tyrosine kinase (TK) domain. Black arrows represent the relative positions of primers that target in and around the JM region for ITD or the activating loop of the kinase domain for TKD. Colored dots represent fluorophores on labeled primers. The yellow box has vertical black lines that represent the position of the EcoRV restriction digest sites.

Reagents

Reagent Name	Units in Assay		
FLT3 Extraction Control	1 x 1800 µL tube		
FLT3 ITD Master Mix	1 x 1500 µL tube		
FLT3 TKD Master Mix	1 x 1500 µL tube		
FLT3 ITD Positive Control	1 x 100 µL tube		
FLT3 TKD Positive Control	1 x 100 µL tube		
FLT3 No Template Control	1 x 200 µL tube		
Taq DNA Polymerase	1 x 200 µL tube		
EcoRV Enzyme	1 x 200 µL tube		
NEBuffer r3.1	1 x 1250 µL tube		



Ordering Information

Catalog #	Products	Quantity
K-412-0361	LeukoStrat® CDx FLT3 Mutation Assay	33 reactions
K-412-0371	LeukoStrat® CDx <i>FLT3</i> Software	CD complimentary with purchase
K-412-0401	LeukoStrat® CDx Assay Installer	USB complimentary with purchase

References

- 1. Murphy KM et al., A Clinical PCR/Capillary Electrophoresis Assay for the Detection of Internal Tandem Duplication and Point Mutation of the *FLT3* Gene. J. Mol. Diag. 5:96-102 (2003).
- 2. Yamamoto, Y., et al., Activating mutation of D835 within the activation loop of *FLT3* in human hematologic malignancies. *Blood*, 97(8):2434-9 (2001).



LeukoStrat[®] CDx FLT3 Mutation Assay

CE 2797 IVD Companion Diagnostic

Available in regions that accept CE 2797 IVD products

The only internationally standardized CE 2797 IVD assay for *FLT3* Signal Ratio mutation analysis for assessment of acute myeloid leukemia (AML) patients eligible for treatment with XOSPATA® (gilteritinib fumarate) and/or VANFLYTA® (quizartinib hydrochloride).

Intended Use

The LeukoStrat CDx *FLT3* Mutation Assay is a PCR-based in vitro diagnostic test designed to detect internal tandem duplication (ITD) and tyrosine kinase domain (TKD) mutations D835 and I836 in the *FLT3* gene in genomic DNA extracted from mononuclear cells obtained from peripheral blood or bone marrow aspirates of patients diagnosed with acute myelogenous leukemia (AML). The LeukoStrat CDx *FLT3* Mutation Assay may be used as a companion diagnostic for the following therapeutic:

In regions where XOSPATA® (gilteritinib fumarate) is available, the LeukoStrat CDx *FLT3* Mutation Assay is used as an aid in the assessment of patients with AML for whom XOSPATA (gilteritinib fumarate) treatment is being considered.

In regions where VANFLYTA® (quizartinib hydrochloride) is available, the LeukoStrat CDx *FLT3* Mutation Assay is used as an aid in the assessment of patients with *FLT3*-ITD+ AML for whom VANFLYTA® (quizartinib hydrochloride) treatment is being considered.

The qualitative, non-automated test is for use on the 3500xL or 3500xL Dx Genetic Analyzers.

Summary and Explanation of the Test

AML in general has a poor prognosis. Assessment of the mutation status of the *FLT3* (fms related tyrosine kinase 3) receptor gene in karyotype normal AML is the most important prognostic indicator of disease outcome, which is often substantial, as many studies in AML have shown that the presence of *FLT3* activating mutations portends a poor prognosis.¹²

The LeukoStrat CDx *FLT3* Mutation Assay targets regions of the *FLT3* gene to identify ITD mutations and TKD mutations, such as the D835 and I836 mutations, and has been validated in an international clinical trial.

The LeukoStrat CDx *FLT3* Mutation Assay includes reagents, equipment, software and procedures for isolating mononuclear cells and extracting DNA from patient specimens to determine if *FLT3* mutations are present. DNA is amplified via PCR and the amplicons are detected via capillary electrophoresis. *FLT3* mutation status is determined by the LeukoStrat CDx *FLT3* Software. A *FLT3* ITD and/or TKD mutation is reported as Positive if the mutant:wild-type signal ratio meets or exceeds the clinical cutoff of 0.05.

Method Description

ITD Mutations of FLT3

The LeukoStrat CDx *FLT3* Mutation Assay uses fluorescently labeled primers that are in and around the JM region. Wild-type *FLT3* alleles will amplify and produce a product at 327±1 bp as measured by this assay, while alleles that contain ITD mutations will produce a product that exceeds 327±1 bp (see Figure, right).

TKD Mutations of FLT3

The LeukoStrat CDx *FLT3* Mutation Assay uses primers that lie on either side of the TKD region. The *FLT3* target region is amplified using PCR and then an EcoRV restriction digest is performed. Wild-type alleles of the *FLT3* gene yield digestion products of 79±1 bp whereas mutant alleles yield products of 125±1 bp or 127±1 bp from the original undigested amplicon product of 145±1 bp or 147±1 bp, as measured by this assay (please see Figure, right).

References

1. Murphy KM et al., A Clinical PCR/Capillary Electrophoresis Assay for the Detection of Internal Tandem Duplication and Point Mutation of the *FLT3* Gene. J. Mol. Diag. 5:96-102 (2003).

2. Yamamoto, Y., et al., Activating mutation of D835 within the activation loop of FLT3 in human hematologic malignancies. Blood, 97(8):2434-9 (2001).



94 Capillary Detection | CE 2797 IVD Assays and IVD Companion Diagnostic Assays



Depicted is a representation of the *FLT3* juxtamembrane (JM) region (TM = transmembrane) and the activating loop of the tyrosine kinase (TK) domain. Black arrows represent the relative positions of primers that target in and around the JM region for ITD or the activating loop of the kinase domain for TKD. Colored dots represent fluorophores on labeled primers. The yellow box has vertical black lines that represent the position of the EcoRV restriction digest sites.

Reagents

Units in Assay
1 x 1800 µL tube
1 x 1500 μL tube
1 x 1500 μL tube
1 x 100 µL tube
1 x 100 µL tube
1 x 200 µL tube
1 x 200 µL tube
1 x 200 µL tube
1 x 1250 µL tube

Ordering Information

Cata	log	#
	- 3	

Products

 K-412-0431
 LeukoStrat® CDx FLT3 Mutation Assay

 K-412-0441
 LeukoStrat® CDx FLT3 Software

Quantity

33 reactions CD complimentary with purchase

These are *in vitro* diagnostic products and are available in regions that accept CE-IVD products.



LeukoStrat[®] CDx FLT3 Mutation Assay

PMDA/MHLW Approved Companion Diagnostic

Available in Japan

The only internationally standardized assay for *FLT3* Signal Ratio mutation analysis for assessment of acute myeloid leukemia (AML) patients eligible for treatment with Gilteritinib Fumarate or Quizartinib Hydrochloride.

Intended Use

The LeukoStrat CDx *FLT3* Mutation Assay is a PCR-based, *in vitro* diagnostic test designed to detect internal tandem duplication (ITD) mutations and tyrosine kinase domain (TKD) mutations D835 and I836 in the *FLT3* gene in genomic DNA extracted from mononuclear cells obtained from peripheral blood or bone marrow aspirates of patients diagnosed with acute myelogenous leukemia.

The LeukoStrat CDx *FLT3* Mutation Assay is used as an aid in the assessment of patients with AML for whom Gilteritinib Fumarate treatment is being considered.

The LeukoStrat CDx *FLT3* Mutation Assay is used as an aid in the assessment of patients with AML for whom Quizartinib Hydrochloride treatment is being considered.

Summary and Explanation of the Test

AML in general has a poor prognosis. Assessment of the mutation status of the *FLT3* (fms related tyrosine kinase 3) receptor gene in karyotype normal AML is the most important prognostic indicator of disease outcome, which is often substantial, as many studies in AML have shown that the presence of *FLT3* activating mutations portends a poor prognosis.^{1,2} The LeukoStrat CDx *FLT3* Mutation Assay targets regions of the *FLT3* gene to identify ITD mutations, such as the D835 and I836 mutations, and has been validated in an international clinical trial.

The LeukoStrat CDx *FLT3* Mutation Assay includes reagents, equipment, software and procedures for isolating mononuclear cells and extracting DNA from patient specimens to determine if *FLT3* mutations are present. DNA is amplified via PCR with fluorescently labeled primers, amplicon is enzymatically digested (TKD), and the amplicons are detected via capillary electrophoresis. *FLT3* mutation status is determined by the LeukoStrat CDx *FLT3* Software. A *FLT3* ITD and/or TKD mutation is reported as Positive if the mutant:wild-type signal ratio meets or exceeds the clinical cutoff of 0.05.

Method Description

ITD Mutations of FLT3

The LeukoStrat CDx *FLT3* Mutation Assay uses fluorescently labeled primers that are in and around the JM region. Wild-type *FLT3* alleles will amplify and produce a product at 327±1 bp as measured by this assay, while alleles that contain ITD mutations will produce a product that exceeds 327±1 bp (see Figure, right).

TKD Mutations of FLT3

The LeukoStrat CDx *FLT3* Mutation Assay uses primers that lie on either side of the TKD region. The *FLT3* target region is amplified using PCR and then an EcoRV restriction digest is performed. Wild-type alleles of the *FLT3* gene yield digestion products of 79±1 bp whereas mutant alleles yield products of 125±1 bp or 127±1 bp from the original undigested amplicon product of 145±1 bp or 147±1 bp, as measured by this assay (see Figure, right).

References

- 1. Murphy KM et al., A Clinical PCR/Capillary Electrophoresis Assay for the Detection of Internal Tandem Duplication and Point Mutation of the *FLT3* Gene. J. Mol. Diag. 5:96-102 (2003).
- 2. Yamamoto, Y., et al., Activating mutation of D835 within the activation loop of FLT3 in human hematologic malignancies. Blood, 97(8):2434-9 (2001).





Depicted is a representation of the *FLT3* juxtamembrane (JM) region (TM = transmembrane) and the activating loop of the tyrosine kinase (TK) domain. Black arrows represent the relative positions of primers that target in and around the JM region for ITD or the activating loop of the kinase domain for TKD. Colored dots represent fluorophores on labeled primers. The yellow box has vertical black lines that represent

Reagents

Reagent Name	Units in Assay
FLT3 Extraction Control	1 x 1800 µL tube
FLT3 ITD Master Mix	1 x 1500 µL tube
FLT3 TKD Master Mix	1 x 1500 µL tube
FLT3 ITD Positive Control	1 x 100 µL tube
FLT3 TKD Positive Control	1 x 100 µL tube
FLT3 No Template Control	1 x 200 µL tube
Taq DNA Polymerase Enzyme	1 x 200 µL tube
EcoRV Enzyme	1 x 200 µL tube
NEBuffer r3.1	1 x 1250 µL tube

Ordering Information

Cata	loa	#

Products

K-412-0331LeukoStrat® CDx FLT3 Mutation Assay (Japan)K-412-0341LeukoStrat® CDx FLT3 Software (Japan)

Quantity

33 reactions CD complimentary with purchase

These are in vitro diagnostic products, and are available for sale or use within Japan.



Gel and Capillary Research Use Only (RUO) Assays

Invivoscribe offers an array of assays for B- and T-cell gene clonality/rearrangements, mutations, and chromosome translocations for the study of hematologic malignancies.

LEARN MORE: catalog.invivoscribe.com

Introduction

These (RUO) assays are available for either ABI capillary electrophoresis fluorescence, or PAGE/agarose gel detection, and contain the PCR master mixes, recommended controls, and Instructions For Use.

On the following pages, you will find detailed information on each RUO assay, including: assay use, background information, typical output data, kit contents, and ordering information. These assays are available in regular sizes (30 or 33 reactions). Select assays are also available in high-volume MegaKit formats (300 or 330 reactions).

These pages contain Research Use Only products which are not for use in diagnostic procedures. Research Use Only (RUO) assays are not for sale in Europe and other global markets where equivalent CE-IVD assays are available and registered with the appropriate regulatory agencies. Refer to the preceding pages for information regarding our IdentiClone[®] and LeukoStrat[®] CE 2797 IVD Assays.

For more information, please visit www.invivoscribe.com

Chapter Contents

B-cell Assays

100 IGH + IGK B-cell Clonality Assays
102 IGH Gene Rearrangement Assays
104 IGH Gene Clonality Assays
106 IGK Gene Clonality Assays
108 IGL Gene Clonality Assays

Translocation Assays

118 BCL1/J_H Translocation Assay **120** BCL2/J_H t(14;18) Translocation Assay **122** BCL2/J_H Translocation Assay **124** BCR/ABL t(9;22) Translocation Assays **126** PML/RAR~ t(15;17) Translocation Assays

T-cell Assays

110 TCRB Gene Clonality Assays
112 T-Cell Receptor Gama Gene Rearrangement Assays 2.0
114 T-Cell Receptor Gama Gene Rearrangement Assays
116 TCRD Gene Clonality Assays

SHM Assays

128 IGH Somatic Hypermutation Assays v2.0

FLT3 Mutation Assays

130 FLT3 Mutation Assays

IGH + IGK B-Cell Clonality Assays

Assay Use

Immunoglobulin heavy chain (*IGH*) and Kappa light chain (*IGK*) gene clonality assays are useful for identification of B-cell clonality, studying clonal B-cell populations, and evaluation of new research and methods in malignancy studies.

Summary and Explanation of the Test

Five PCR master mixes are included in these test kits to test for rearrangements of both *IGH* and *IGK*. *IGH* Tubes A, B, and C target the conserved framework 1, 2, and 3 regions (respectively) within the variable (VH) region and the joining (JH) region of the *IGH* locus. *IGK* tubes A and B target the variable (VK), intragenic and joining (JK), and kappa deleting element (Kde) regions of the *IGK* locus.

Positive and negative controls, as well as Specimen Control Size Ladder Master Mix are included. PCR products can be analyzed by capillary electrophoresis or heteroduplex analysis. Clonality is indicated if any one of the master mixes generates clonal products.

Background

The immunoglobulin heavy chain (*IGH*) gene locus on chromosome 14 (14q32.33, formerly 14q32.3) includes 46-52 functional and 30 nonfunctional variable (VH), 27 functional diversity (DH), and 6 functional joining (JH) gene segments spread over 1250 kilobases.^{1,2} The most frequently used VH gene segments in normal and malignant B cells belong to VH3, VH4, and VH1 families, which together cover 75–95% of VH usage. The VH gene segments contain three framework regions (FR) and two complementarity determining regions (CDR). The FRs are characterized by their similarity among the various VH segments, whereas the CDRs are highly different even within the same VH family. The CDRs represent the preferred target sequences for somatic hypermutations; however, somatic mutations can also occur in the FRs. Therefore, family-specific primers in the three different FRs were designed to increase the detection rate of clonal *IGH* B-cell populations and decrease the occurrence of false-negative results due to somatic hypermutation in primer binding sites.¹

The human immunoglobulin kappa (*IGK*) light chain locus on the short arm of chromosome 2 (2p11.2) spans 1820 kb. It is made up of 76 variable (Vk) gene segments belonging to seven subgroups, five joining (Jk) gene segments, and one constant (Ck) gene segment. Productive assembly of the kappa gene is successful in about 60% of human B lymphocytes²; however, even when unsuccessful, clonal B cells generally retain the rearranged kappa genes. The Vk segments encode the first 95 N-terminal amino acids. Positions 96-108 are encoded by one of five joining (Jk) gene segments. The constant (Ck) portion of the kappa light chain (amino acids 109-214) is encoded by a single constant (Ck) region separated from the Jk region by an intron.

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA)

References

- 1. M Hummel et al., Leukemia 17: 2266-2272 (2003).
- 2. AW Langerak et al., Leukemia 17: 2272–2275 (2003).
- 3. EP Rock, PR Sibbald, MM Davis, and YH Chien. J. Exp. Med. 179(1): 323-328 (1994).
- 4. JJM van Dongen et al., *Leukemia* 17: 2257–2317 (2003).

EuroClonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.



INTR primer Figure Legend: Simple representation of the organization of a rearranged immunoglobulin heavy chain (IGH) gene on chromosome 14 and the immunoglobulin kappa

IGK tube B: 6 V κ primers and INTR primer + 1 κ de primer Jк-Cк intron

Kde

Kde primer

light chain gene on chromosome 2p11.2. Black arrows represent the relative positions of primers that target the conserved framework regions (FR1-3) and the downstream consensus JH gene segments for IGH and the VK, JK, INTR and K_{de} primers which are included in the IGK master mix tubes.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0030 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0019 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0007 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
<i>IGH</i> Tube A	Framework 1 + J _H	1 x 1500 µL tube	10 x 1500 µL tubes
<i>IGH</i> Tube B	Framework 2 + J _H	1 x 1500 µL tube	10 x 1500 µL tubes
<i>IGH</i> Tube C	Framework 3 + J _H	1 x 1500 µL tube	10 x 1500 µL tubes
<i>IGK</i> Tube A	Vƙ-Jƙ	1 x 1500 µL tube	10 x 1500 µL tubes
<i>IGK</i> Tube B	Vƙ-K _{de}	1 x 1500 µL tube	10 x 1500 µL tubes
Specimen Control Size Ladder	Multiple Genes	1 x 1500 µL tube	10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
1-100-0010	IGH + IGK B-Cell Clonality Assay - Gel Detection	33 reactions
1-100-0031	IGH + IGK B-Cell Clonality Assay - ABI Fluorescence Detection	33 reactions
1-100-0041	IGH + IGK B-Cell Clonality Assay MegaKit - ABI Fluorescence Detection	330 reactions

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.

IGH Gene Rearrangement Assays

Assay Use

Immunoglobulin heavy chain (*IGH*) gene clonality assays are useful for the study of identifying clonal B-cell populations and evaluation of new research and methods in malignancy studies.

Summary and Explanation of the Test

Genomic DNA is amplified using three PCR master mixes that target the three conserved framework regions (FR1, FR2, and FR3) of the *IGH* gene and the joining (J_{μ}) region. These regions flank the unique, hypervariable, antigen-binding, complementarity determining region 3 (CDR3). All positive and negative DNA controls, as well as an Amplification Control master mix, are included. PCR products can be analyzed by capillary electrophoresis. Clonality is indicated if one or more of the three framework master mixes generates clonal products.

Background

Genes encoding immunoglobulin heavy chain (*IGH*) molecules are assembled from multiple polymorphic gene segments that undergo rearrangement and selection during B-cell development.² Rearrangement of these variable (VH), diversity (DH), and joining (JH) genetic segments result in VDJ products of unique length and sequence.^{1,2} Clonal *IGH* rearrangements can be rapidly identified through analyses of the size distributions of DNA products amplified from conserved sequences that flank this region.² For example, DNA isolated from a normal polyclonal population of B cells produces a Gaussian distribution (bell-shaped size curve) of amplified products; whereas, DNA amplified from a clonal B-cell population generates one or two product(s) of unique size that reflect proliferation of a single rearranged clone.¹

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. JE Miller, SS Wilson, DL Jaye, and M Kronenberg. J. Mol. Diag. 4: 101-117 (1999).
- 2. S Tonegawa. Nature 302: 575-581 (1983).



Figure Legend: Genomic organization of a rearranged immunoglobulin heavy chain gene on chromosome 14. The blue and green arrows represent primers targeting the conserved framework regions within the variable region gene. The relative location, size range of valid products, and colors correspond to the products generated from each of these regions when differential fluorescence detection methods are used.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0030 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0029 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
Master Mixes IGH Framework 1	Target Framework 1 + J _H	Units in Assay 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes
Master Mixes IGH Framework 1 IGH Framework 2	Target Framework 1 + J _H Framework 2 + J _H	Units in Assay 1 x 1500 μL tube 1 x 1500 μL tube	Units in Assay MegaKit 10 x 1500 µL tubes 10 x 1500 µL tubes
Master Mixes IGH Framework 1 IGH Framework 2 IGH Framework 3	Target Framework 1 + J _H Framework 2 + J _H Framework 3 + J _H	Units in Assay 1 x 1500 μL tube 1 x 1500 μL tube 1 x 1500 μL tube	Units in Assay MegaKit 10 x 1500 µL tubes 10 x 1500 µL tubes 10 x 1500 µL tubes
Master Mixes IGH Framework 1 IGH Framework 2 IGH Framework 3 Amplification Control	Target Framework 1 + Jн Framework 2 + Jн Framework 3 + Jн HLA-DQa	Units in Assay 1 x 1500 μL tube 1 x 1500 μL tube 1 x 1500 μL tube 1 x 1500 μL tube	Units in Assay MegaKit 10 x 1500 µL tubes 10 x 1500 µL tubes 10 x 1500 µL tubes 10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
1–101–0051	IGH Gene Rearrangement Assay - ABI Fluorescence Detection	30 reactions
1-101-0071	IGH Gene Rearrangement Assay MegaKit - ABI Fluorescence Detection	300 reactions

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.

IGH Gene Clonality Assays

Assay Use

Immunoglobulin heavy chain (*IGH*) clonality assays are useful for identification of B-cell clonality, studying clonal B-cell populations, and evaluation of new research and methods in malignancy studies.

Summary and Explanation of the Test

Five master mixes target conserved regions within the variable (V_H), diversity (D_H), and the joining (J_H) regions that flank the unique hypervariable, antigen-binding, complementarity determining region 3 (CDR3). Tube A contains six framework region 1 (FR1) primers and a consensus J_H region primer. Tube B contains seven framework region 2 (FR2) primers and a consensus J_H primer. Tube C contains seven framework region 3 (FR1) primers and a consensus J_H primer. Tube D contains six D_H region primers and a consensus J_H region primer. Tube E contains a D_H7 region primer and a consensus J_H primer. Positive and negative controls, as well as the Specimen Control Size Ladder Master Mix are included. PCR products can be analyzed by capillary electrophoresis or heteroduplex analysis. Clonality is indicated if any one of the master mixes generates a clonal product.

Background

The immunoglobulin heavy chain (*IGH*) gene locus on chromosome 14 (14q32.33, formerly 14q32.3) includes 46-52 functional and 30 nonfunctional variable (V_H), 27 functional diversity (D_H), and 6 functional joining (J_H) gene segments spread over 1250 kilobases.^{1,2} The most frequently used V_H gene segments in normal and malignant B cells belong to the V_H3 , V_H4 , and V_H1 family, together covering 75–95% of V_H usage. The V_H gene segments contain three framework regions (FR) and two complementarity determining regions (CDR).

The FRs are characterized by their similarity among the various V_H segments, whereas the CDRs are highly different even within the same V_H family. The CDRs represent the preferred target sequences for somatic hypermutations; however, somatic mutations can also occur in the FRs. Therefore, family-specific primers in the three different FRs were designed to increase the detection rate of clonal *IGH* B-cell populations and decrease the occurrence of false-negative results due to somatic hypermutation in primer binding sites.¹ In addition to V_H -J_H rearrangements, incomplete D_H -J_H rearrangements have been found in mature and immature B-cell malignancies. Therefore, D_H -J_H PCR analysis may be of added value for clonality assessment.²

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. M Hummel et al., Leukemia 17:2266-2272 (2003).
- 2. AW Langerak et al., Leukemia 17:2272-2275 (2003)
- 3. JJM van Dongen et al., Leukemia 17:2257-2317 (2003).

EuroClonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.



Tube E: D_H 7 Primer + J_H Consensus Primer

Figure Legend: Simple representation of the organization of a rearranged immunoglobulin heavy chain gene on chromosome 14. Black arrows represent the relative positions of primers that target the conserved framework (FR1-3) and diversity (D_H1-7) regions, and the downstream consensus J_H gene segments. The amplicon products generated from each of these regions can be differentially detected when fluorescent primer sets are used with capillary electrophoresis instruments that employ differential fluorescence detection.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0030 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0019 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0024 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0008 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg /mL	1 x 100 µL tube	$5 \times 100 \ \mu L$ tubes
Master Mixes	Taraet	Units in Assay	Units in Assay MeaaKit
	laigei	onns in Assay	onna ni Assay Megakii
IGH Tube A	Framework 1 + J _H	1 x 1500 µL tube	10 x 1500 µL tubes
IGH Tube A IGH Tube B	Framework 1 + J _H Framework 2 + J _H	1 x 1500 µL tube	10 x 1500 µL tubes 10 x 1500 µL tubes
IGH Tube A IGH Tube B IGH Tube C	Framework 1 + J _H Framework 2 + J _H Framework 3 + J _H	1 x 1500 µL tube 1 x 1500 µL tube 1 x 1500 µL tube	10 x 1500 μL tubes 10 x 1500 μL tubes 10 x 1500 μL tubes
IGH Tube A IGH Tube B IGH Tube C IGH Tube D	Framework 1 + J _H Framework 2 + J _H Framework 3 + J _H D _H 1-6 + J _H	1 x 1500 μL tube 1 x 1500 μL tube 1 x 1500 μL tube 1 x 1500 μL tube 1 x 1500 μL tube	10 x 1500 μL tubes 10 x 1500 μL tubes 10 x 1500 μL tubes 10 x 1500 μL tubes
IGH Tube A IGH Tube B IGH Tube C IGH Tube D IGH Tube E	Framework 1 + J _H Framework 2 + J _H Framework 3 + J _H D _H 1-6 + J _H D _H 7 + J _H	1 x 1500 µL tube 1 x 1500 µL tube	10 x 1500 μL tubes 10 x 1500 μL tubes

Ordering Information

Catalog #	Products	Quantity
1-101-0020	IGH Gene Clonality Assay - Gel Detection	33 reactions
1-101-0061	IGH Gene Clonality Assay - ABI Fluorescence Detection	33 reactions
1-101-0081	IGH Gene Clonality Assay MegaKit - ABI Fluorescence Detection	330 reactions

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.

IGK Gene Clonality Assays

Assay Use

Immunoglobulin kappa light chain (*IGK*) gene clonality assays are useful for the identification of B-cell clonality, studying clonal B-cell populations, and evaluation of new research and methods in malignancy studies.

Summary and Explanation of the Test

Two master mixes target conserved regions within the variable (Vk1-7) and the joining (Jk1-5) regions that flank the unique hypervariable, antigen-binding, complementarity determining region 3 (CDR3). Other primers target the K_{de} and intragenic regions.

Tube A contains six upstream primers and two Jk region primers. Tube B contains six upstream Vk region primers, an upstream intragenic primer and a downstream K_{de} primer. Positive and negative controls, as well as a Specimen Control Size Ladder Master Mix, are included. PCR products can be analyzed by capillary electrophoresis or heteroduplex analysis. Clonality is indicated if any one of the master mixes generates clonal products.

Background

The human immunoglobulin kappa (*IGK*) light chain locus on the short arm of chromosome 2 (2p12, formerly 2p11.2) spans 1820 kb. It is made up of 76 variable (Vk) gene segments belonging to 7 subgroups, 5 joining (Jk) gene segments, and one constant (Ck) gene segment. Productive assembly of the kappa gene is successful in about 60% of human B lymphocytes.¹ However, even when unsuccessful, clonal B cells generally retain the rearranged kappa genes. The Vk segments encode the first 95 N-terminal amino acids. Positions 96-108 are encoded by one of five joining (Jk) gene segments. The constant (Ck) portion of the kappa light chain (amino acids 109-214) is encoded by a single constant (Ck) region separated from the Jk region by an intron. The length of the hypervariable complementarity determining region 3 (CDR3) in kappa light chain genes is limited and rearrangements in this region display significant skewing (platykurtosis).²

Therefore, clonal CDR3 products generated from this region are most easily and reliably identified by heteroduplex analysis using standard polyacrylamide gels. Alternatively, capillary electrophoresis or gene sequencing instruments coupled with differential fluorescence detection can be used for analysis.

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. AW Langerak et al., Leukemia 17: 2275-2280 (2003).
- 2. EP Rock, PR Sibbald, MM Davis, and YH Chien. J. Exp. Med. 179(1): 323-328 (1994).
- 3. JJM van Dongen et al., Leukemia 17: 2257-2317 (2003).

EuroClonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.



IGK tube B: 6 Vk primers and INTR primer + 1 Kde primer

Figure Legend: Schematic diagram of the immunoglobulin kappa light chain gene complex on chromosome 2p11.2. Shown are the relative positions and orientations for the Vk-Jk, and K_{de} primers, which are included in the *IGK* master mix tubes.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0007 Clonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg /mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
Master Mixes IGK Tube A	Target VR-JK	Units in Assay 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes
Master Mixes IGK Tube A IGK Tube B	Target Vk-Jk Vk-Kde	Units in Assay 1 x 1500 µL tube 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes 10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
1-102-0020	IGK Gene Clonality Assay - Gel Detection	33 reactions
1-102-0021	IGK Gene Clonality Assay - ABI Fluorescence Detection	33 reactions
1-102-0031	IGK Gene Clonality Assay MegaKit - ABI Fluorescence Detection	330 reactions

These products are for Research Use Only (RUO). Not intended for diagnostic purposes.

IGL Gene Clonality Assays

Assay Use

Immunoglobulin lambda light chain (*IGL*) gene clonality assays are useful for the identification of B-cell clonality, studying clonal B-cell populations, and evaluation of new research and methods in malignancy studies.

Summary and Explanation of the Test

The *IGL* Tube master mix targets conserved regions within the variable (V λ 1-3) and the joining (J λ 1-3) regions that flank the unique, hypervariable, antigen-binding, complementarity determining region 3 (CDR3). Positive and negative controls, as well as a Specimen Control Size Ladder Master Mix, are included. PCR products can be analyzed by capillary electrophoresis or heteroduplex analysis. Clonality is indicated if the master mix generates clonal products.

Background

The human immunoglobulin lambda (*IGL*) light chain locus is located on the long arm of chromosome 22 (22q11.2) and spans 1050 kilobases. It is made up of 73–74 variable (V λ) gene segments (spread over 900 kilobases), 7–11 joining (J λ) gene segments and 7–11 constant (C λ) gene segments depending on the haplotypes. Of the 73–74 V λ region genes, only 30–33 are functional and can be grouped into 11 families and 3 clans.¹ The J λ and C λ region genes are organized in tandem with a J λ segment preceding a C λ gene. Typically there are 7 J λ –C λ segments of which four are functional and encode the four Ig lambda isotypes.

IGL gene rearrangements (V λ -J λ) rearrangements potentially represent an attractive PCR target for clonality studies to compensate for false-negative $IGH V_{H}$ -J_H PCR results mainly caused by somatic hypermutations. The limited size of the junctional region may create a challenge to distinguish polyclonal from monoclonal rearrangements when running a simple agarose or polyacrylamide gel.¹ Therefore, clonal V λ -J λ PCR products are most easily and reliably identified by heteroduplex analysis using standard polyacrylamide gels. Alternatively, capillary electrophoresis or gene sequencing instruments coupled with differential fluorescence detection can also be used for analysis.¹

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. F Davi et al., Leukemia 17:2280-2283 (2003).
- 2. JJM van Dongen et al., Leukemia 17:2257-2317 (2003).

EuroClonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.


IGL tube: 2 V λ primers + 1 J λ primer

Figure Legend: Schematic diagram of the immunoglobulin lambda light chain gene complex on chromosome 22q11.2. Shown are the relative positions and orientations for the Vλ and Jλ primers, which are included in the *IGL* master mix tube. The two Vλ primers only target Vλ1, 2, and 3 because these three V families cover approximately 70% of rearrangeable Vλ gene segments, and approximately 90% of all *IGL* gene rearrangements involve these three families. Similarly, the single Jλ primer only targets Jλ1, 2, and 3 because these three J segments are involved in 98% of all *IGL* gene rearrangements.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0010 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0029 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
<i>IGL</i> Tube	Vλ–Jλ	1 x 1500 µL tube	10 x 1500 µL tubes
Specimen Control Size Ladder	Multiple Genes	1 x 1500 ul tube	10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
1–103–0010	IGL Gene Clonality Assay - Gel Detection	33 reactions
1-103-0011	IGL Gene Clonality Assay - ABI Fluorescence Detection	33 reactions
1-103-0021	IGL Gene Clonality Assay MegaKit - ABI Fluorescence Detection	330 reactions

TCRB Gene Clonality Assays

Assay Use

T-Cell receptor beta (*TCRB*) gene clonality assays are useful for the identification of T-cell clonality, studying clonal T-cell populations, and evaluation of new research and methods in malignancy studies.

Summary and Explanation of the Test

Three multiplex master mixes target conserved regions within the variable (V β), diversity (D β), and the joining (J β) regions that flank the unique hypervariable, antigen-binding, complementarity determining region 3 (CDR3) of the T-cell receptor beta locus. Tube A contains 23 V β primers, six J β 1 primers, and three J β 2 primers. Tube B contains 23 V β and four J β 2 primers. Tube C contains two D β and 13 J β primers. Positive and negative DNA controls, as well as a Specimen Control Size Ladder Master Mix, are included. PCR products can be analyzed by capillary electrophoresis or heteroduplex analysis. Clonality is indicated, if any one of the master mixes generates clonal products.

Background

The human T-cell receptor beta (*TRB*, formerly known as *TCRB*) gene locus on chromosome 7 (7q34, formerly 7q35) includes 64-67 variable (V_β) gene segments (belonging to 30 subgroups), two diversity (D_β) gene segments, and 13 joining (J_β) gene segments, spread over 685 kilobases. The diversity of this locus has complicated PCR-based testing and extended dependence on Southern blot analysis in many testing centers. However, this standardized multiplex PCR assay detects the vast majority of clonal *TRB* gene rearrangements using only three multiplex master mixes.¹

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. M Brüggemann et al., Leukemia 17: 2283-2289 (2003).
- 2. JJM van Dongen et al., Leukemia 17: 2257-2317 (2003).

EuroClonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.



Figure Legend: Simplified diagram of a representative rearranged T-cell receptor beta gene showing the approximate placement of the upstream and downstream DNA primers. The numbers of primers and their specificity are listed for Master Mix Tubes A, B, and C.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0009 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0004 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0021 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
TCRB Tube A	Maultin La V/0 - 101/2		
	Multiple vb + Jbl/ 2	1 x 1500 µL tube	10 x 1500 µL tubes
<i>TCRB</i> Tube B	Multiple V β + J β 1/2 Multiple V β + J β 2	1 x 1500 μL tube 1 x 1500 μL tube	10 x 1500 µL tubes 10 x 1500 µL tubes
<i>TCRB</i> Tube B <i>TCRB</i> Tube C	Multiple V β + J β 1/2 Multiple V β + J β 1/2 Multiple D β + J β 1/2	1 x 1500 µL tube 1 x 1500 µL tube 1 x 1500 µL tube	10 x 1500 μL tubes 10 x 1500 μL tubes 10 x 1500 μL tubes

Ordering Information

Catalog #	Products	Quantity
1-205-0010	TCRB Gene Clonality Assay - Gel Detection	33 reactions
1-205-0020	TCRB Gene Clonality Assay MegaKit - Gel Detection	330 reactions
1-205-0011	TCRB Gene Clonality Assay - ABI Fluorescence Detection	33 reactions
1-205-0021	TCRB Gene Clonality Assay MegaKit - ABI Fluorescence Detection	330 reactions

T-Cell Receptor Gamma Gene Rearrangement Assay 2.0

Assay Use

This Research Use Only assay identifies T-cell receptor gamma (*TCRG*) chain gene rearrangements and is useful for the identification of clonal T-cell populations and evaluation of new research and methods in malignancy studies.

Summary and Explanation of the Test

This T-Cell Receptor Gamma Gene Rearrangement Assay 2.0 contains a single multiplex master mix which includes primers for all known groups of *TCR* gamma variable (Vy) region genes and joining (Jy) region genes involved in rearrangements. In addition, positive and negative controls, as well as a Specimen Control Size Ladder Master Mix are included. PCR products can be analyzed by polyacrylamide or capillary electrophoresis detection.

Background

The human T-cell receptor gamma (*TRG*, formerly known as *TCRG*) gene locus on chromosome 7 (7q14) includes 14 Vy genes belonging to four subgroups, five Jy segments, and two Cy genes spread over 200 kilobases. The diversity of this locus has historically complicated PCR-based testing. This multiplex PCR assay represents an improvement over other assays as it can detect the vast majority of *TCR* gamma gene rearrangements with a single multiplex master mix. This master mix targets all conserved regions within the variable (Vy) and joining (Jy) region genes, providing a more comprehensive analysis to include Vy and Jy regions that would not be identified with a single Vy (1-8) and Jy 1/Jy 2 primer set.

In addition, competitive amplification of all *TRG* gene rearrangements allows for identification of a quantitative threshold for a positive result and helps to avoid false positive results. The average size of the *TRG* gene rearrangement PCR amplicons is 190 nucleotides, with a normal distribution of product sizes between 159 and 207 nucleotides. This protocol leads to improved product formation from formalin-fixed, paraffin-embedded (FFPE) samples compared to other protocols that yield products of 260 nucleotides or larger.

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. TC Greiner et al., JMD 4: 137-143 (2002).
- 2. LC Lawnickie et al., JMD 5: 82-87 (2003).
- 3. Y Sandberg et al., Leukemia 21: 21 (2007).
- 4. Armand, Marine et al. HemaSphere, 2019;3:3.



This assay was developed by Invivoscribe. The performance of this assay was reviewed and validated by the EuroClonality/BIOMED-2 Group.⁴



Figure Legend: Simple representation of the organization of the T-cell receptor gamma gene on chromosome 7. Black arrows represent the relative positions of primers that target the variable region genes and the downstream joining region gene segments that are involved in rearrangements in T-cell lymphomas. The downstream primers are fluorescently labeled through the incorporation of a 6FAM fluorophore. The amplicon products generated from these rearrangements are detected by capillary electrophoresis.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
5% TCRG Positive Control DNA	50 µg/mL	1 x 50 µL tube	5 x 50 µL tube
TCRG Negative Control DNA	50 µg/mL	1 x 50 µL tube	5 x 50 µL tube
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
Master Mixes TCRG - 6FAM	Target Vy1-Vy11 + Jy1/Jy2, JyP, JyP1/JyP2	Units in Assay 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
1-207-0101	T-Cell Receptor Gamma Gene Rearrangement Assay 2.0 – ABI Fluorescence Detection	33 reactions
1-207-0111	T-Cell Receptor Gamma Gene Rearrangement Assay MegaKit 2.0 - ABI Fluorescence Detection	330 reactions
S100001	IdentiClone® RUO TCRG 2.0 Software	1 Software Package

T-Cell Receptor Gamma Gene Rearrangement Assays

Assay Use

T-Cell receptor gamma (*TCRG*) gene clonality assays are useful for the identification of T-cell clonality, studying clonal T-cell populations, and evaluation of new research and methods in malignancy studies.

Summary and Explanation of the Test

Sample genomic DNA is amplified using two master mixes that independently target conserved regions within the variable (Vy) and joining (Jy) regions that flank the unique, hypervariable, antigen-binding, complementarity determining region 3 (CDR3). This assay targets Vy1-9 and Jy gene segments. Positive and negative DNA controls, as well as an internal Amplification Control Master Mix, are included. PCR products can be analyzed by capillary electrophoresis.

Background

The T-cell receptor gamma (*TRG*, formerly known as *TCRG*) chain locus spans 160 kilobases on chromosome 7 (7p14). The locus consists of 14 variable (Vy) gene segments in six subgroups, and five joining (Jy) gene segments interspersed between two constant (Cy) gene segments. However, the repertoire of functional *TRG* molecules is limited to 4-6 functional Vy gene segments that belong to two subgroups.²

Rearrangement of the Vy and Jy gene segments of the *TRG* locus results in Vy-Jy products of unique length and sequence. Clonal *TRG* rearrangements can be most rapidly identified by analyzing the size distribution of DNA products amplified from conserved sequences that flank this Vy-Jy region.¹ DNA isolated from a normal heterogeneous population of polyclonal T-cells produces a Gaussian distribution (bell-shaped size curve) of amplified products. DNA amplified from a clonal T-cell population generates one or two product(s) of unique size that reflects proliferation of a single rearranged clone.^{1,2}

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. JE Miller, SS Wilson, DL Jaye, and M Kronenberg. J. Mol. Diag. 4: 101-117 (1999).
- 2. K Beldjord et al., Leukemia 17: 2289-2292 (2003).



Figure Legend: Simplified figure representing the organization of a rearranged T-cell receptor gamma chain gene on chromosome 7. Colored arrows represent conserved regions within the variable region gene segments targeted by primers. Primers are represented by arrows with the size range of valid products generated with each of the master mixes indicated below the figure. Colors correspond to the peak colors assigned to products when differential fluorescence detection methods are used.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0009 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
Master Mixes T-Cell Receptor Gamma Mix 1	Target Vy1-8,9 + Jy1/2	Units in Assay 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes
Master Mixes T-Cell Receptor Gamma Mix 1 T-Cell Receptor Gamma Mix 2	Target Vy1-8,9 + Jy1/2 Alt Vy + Jy1/2	Units in Assay 1 x 1500 µL tube 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes 10 x 1500 µL tubes
Master Mixes T-Cell Receptor Gamma Mix 1 T-Cell Receptor Gamma Mix 2 Amplification Control	Target Vy1-8,9 + Jy1/2 Alt Vy + Jy1/2 HLA-DQa	Units in Assay 1 x 1500 µL tube 1 x 1500 µL tube 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes 10 x 1500 µL tubes 10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
1-207-0051	T-Cell Receptor Gene Rearrangement Assay - ABI Fluorescence Detection	30 reactions

TCRD Gene Clonality Assays

Assay Use

T-Cell receptor delta (*TCRD*) gene clonality assays are useful for the identification of T-cell clonality, studying clonal T-cell populations, and evaluation of new research and methods in malignancy studies.

Summary and Explanation of the Test

The *TCRD* Tube Master Mix targets conserved regions within the variable ($V\delta$ 1-6), the diversity ($D\delta$ 2-3) and the joining ($J\delta$ 1-4) regions that flank the unique, hypervariable, antigen-binding, complementarity determining region 3 (CDR3) of the T-cell receptor delta (*TRD*, formerly known as *TCRD*). Positive and negative controls, as well as a Specimen Control Size Ladder Master Mix, are included. PCR products can be analyzed by capillary electrophoresis or heteroduplex analysis. Clonality is indicated if the master mix generates clonal products.

Background

The human T-cell receptor delta (*TRD*, formerly known as *TCRD*) gene locus is comprised of a cluster of 10 genes located on chromosome 14 (14q11.2) spread over 60 kilobases, localized between the T-cell receptor alpha (*TRA*, formerly known as *TCRA*) variable (Va) and joining (Ja) gene segments. It is made up of eight variable (Vo), three diversity (Do), and four joining (Jo) gene segments.¹ At least five of the eight Vo gene segments can also rearrange to Jo gene segments and other Vo gene segments may also be utilized in *TRD* gene rearrangements in rare cases. Although the small number of Vo, Do, and Jo gene segments available for recombination limits the potential combinatorial diversity, the complementarity determining region 3 (CDR3) or junctional diversity is extensive due to the addition of N regions, P regions, and random deletion of nucleotides by recombinases. This diversity is also extended by the recombination of up to three Do segments and therefore up to four N regions within the rearranged *TRD* locus. This limited germline diversity encoded at the *TRD* locus in conjunction with extensive junctional diversity results in a useful target for PCR analysis. *TRD* recombination events have been used most extensively as clonal markers in both T- and B-cell ALL. This standardized multiplex PCR assay detects the vast majority of clonal *TRD* gene rearrangements using a single multiplex master mix¹.

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. FL Lavender et al., Leukemia 17: 2292-2296 (2003).
- 2. JJM van Dongen et al., Leukemia 17: 2257-2317 (2003).

Second State of the State of th



TCRD tube: 6 V δ and 1 D δ 2 primers + 4 J δ and 1 D δ 3 primers

Figure Legend: Simplified diagram of a representative rearranged T-cell receptor delta gene showing the approximate placement of the upstream and downstream DNA primers. The numbers of primers and their specificity are listed for the TRD Tube Master Mix tube.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0021 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
<i>TCRD</i> Tube	Multiple Võ + Dõ + Jõ	1 x 1500 µL tube	10 x 1500 µL tubes
Specimen Control Size Ladder	Multiple Genes	1 x 1500 µL tube	10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
1-206-0010	TCRD Gene Clonality Assay - Gel Detection	33 reactions
1-206-0011	TCRD Gene Clonality Assay - ABI Fluorescence Detection	33 reactions
1-206-0021	TCRD Gene Clonality Assay MegaKit - ABI Fluorescence Detection	330 reactions

BCL1/J_H Translocation Assay

Assay Use

These assays identify BCL1/J_H t(11;14) translocations and are useful for the evaluation of new research and methods in malignancy studies.

Summary and Explanation of the Test

Two master mixes are included in this assay kit. The $BCL1/J_H$ Master Mix targets the major translocation cluster (MTC) of the *CCND1* locus (formerly known as BCL1) and the joining region (J_H) of the immunoglobulin heavy chain locus (*IGH*). The Specimen Control Size Ladder Master Mix targets multiple genes and generates a series of amplicons of approximately 100, 200, 300, 400, and 600 base pairs to ensure that the quality and quantity of input DNA is adequate to yield a valid result. Positive and negative controls are included. PCR products can be analyzed using standard gel electrophoresis with ethidium bromide staining. A *CCND1* translocation is indicated if the master mix generates product(s) within the valid size range.

Background

This aberrant $BCL1/J_{H}$ t(11;14) translocation juxtaposes genes of the immunoglobulin heavy chain (*IGH*) gene on chromosome 14q32 with the cyclin D1 gene on chromosome 11q13. The juxtaposition of *IGH*-sequences results in the transcriptional activation of cyclin D1.^{2,3} Cyclin D1 is involved in the regulation of the G1 progression and G1/S transition of the cell cycle.³ Translocation does not lead to expression of a fusion protein. In fact, oncogenesis is due to a promoter/enhancer exchange, wherein the immunoglobulin gene enhancer stimulates the expression of cyclin D1. Overexpression of cyclin D1, in turn, accelerates passage of transformed cells through the G1 phase.

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. P Wijers et al., *Leukemia* 17: 2296-2298 (2003).
- 2. JJM van Dongen et al., Leukemia 17: 2257–2317 (2003).
- 3. Shimazaki C, et. al., (1997). International Journal of Hematology. 66(1):111-5.

See Euro Clonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.



t(11;14) tube: 1 BCL1 MTC primer + 1 JH primer

Figure Legend: Schematic diagram of the *IGH-CCND1* t(11;14) translocation showing the cyclin D1 (*CCND1*) gene on the left and the Ig heavy chain (*IGH*) gene on the right. Shown are the relative positions and orientations for the *BCL1*/MTC primer and the J_H primer, which are included in the *BCL1*/J_H Master Mix tube.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0010 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0000 Polyclonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
Master Mixes BCL1/J _H Tube	Target MTC of <i>CCND1</i> + <i>IGH</i> J _H	Units in Assay 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes
Master Mixes <i>BCL1/J</i> H Tube Specimen Control Size Ladder	Target MTC of <i>CCND1</i> + <i>IGH</i> J _H Multiple Genes	Units in Assay 1 x 1500 µL tube 1 x 1500 µL tube	Units in Assay MegaKit 10 x 1500 µL tubes 10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
1-308-0010	BCL1/J _H Translocation Assay - Gel Detection	33 reactions
1-308-0020	BCL1/JH Translocation Assay MegaKit - Gel Detection	330 reactions

BCL2/J_H t(14;18) Translocation Assays

Assay Use

This Research Use Only assay identifies $BCL2/J_{H}$ t(14;18) translocations.

Summary and Explanation of the Test

Five master mixes are included in this assay kit. Two master mixes target *BCL2* major break point (MBR) translocations and two target *BCL2* minor cluster region (mcr) translocations. An Amplification Control Master Mix is also included to ensure the quality and quantity of sample DNA. Positive and negative controls are also included. This assay can be run either in a standard or nested assay format. PCR products can be analyzed by standard gel electrophoresis with ethidium bromide staining. A *BCL2* translocation is indicated if just one of the 2nd round master mixes (mixes ending in b) generates product(s) within the valid size range.

Background

BCL2 translocations are reciprocal chromosome exchanges that place the *bcl*-2 proto-oncogene, located on chromosome 18, under aberrant transcriptional control of the immunoglobulin heavy chain gene, located on chromosome 14. The *bcl*-2 protein is an antagonist to apoptosis (programmed cell death), a normal process designed to eliminate unneeded and damaged cells during hematopoiesis. Increased expression of the *bcl*-2 protein leads to an increase in the levels of B-cells in the body.

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. MS Lee et al., Science 237: 175-178 (1987).
- 2. M Crescenzi et al., Proc. Natl. Acad. Sci. USA 85: 4869-4873 (1988).



Figure Legend: Simplified view of the genomic organization of the *BCL2* and *IGH* genes on chromosomes 18 and 14, respectively. Yellow boxes represent the exon regions of the *BCL2* gene. Exons of the immunoglobulin heavy chain gene are represented in other colors. The solid black lines represents intron regions, which have been left incompletely spliced to assist in demarcation of the exon segments. MBR and mcr type t(14;18) translocations are shown in the lower portions of the figure with the relative positions of primers and the size of the amplicons generated from the positive control DNAs indicated.

Reagents

Controls

IVS-0030 Clonal Control DNA IVS-0031 Clonal Control DNA IVS-0009 Clonal Control DNA

Master Mixes

BCL2/JH t(14;18) (MBR) Mix 1b BCL2/JH t(14;18) (mcr) Mix 2b BCL2/JH t(14;18) (MBR) Mix 1a BCL2/JH t(14;18) (mcr) Mix 2a Amplification Control

trol DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
trol DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
itrol DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
	Target	Units in Assay	Units in Assay MegaKit
R) Mix 1b	Inside BCL2 MBR	1 x 1500 µL tube	10 x 1500 µL tubes
) Mix 2b	Inside BCL2 mcr	1 x 1500 µL tube	10 x 1500 µL tubes
R) Mix 1a	Outside BCL2 MBR	1 x 1500 µL tube	10 x 1500 µL tubes
) Mix 2a	Outside BCL2 mcr	1 x 1500 µL tube	10 x 1500 µL tubes
I	HLA-DQa	1 x 1500 µL tube	10 x 1500 µL tubes

Units in Assay

Ordering Information

Catalog #

Products

1-309-0010

BCL2/JH t(14;18) Translocation Assay - Gel Detection

Concentration

Quantity

30 reactions

Units in Assay MegaKit

BCL2/J_H Translocation Assay

Assay Use

This Research Use Only assay identifies *BCL2/*J_H translocations.

Summary and Explanation of the Test

Four master mixes are included in this assay. Three are used to identify translocations in the major breakpoint region (MBR) and minor cluster region (mcr) of *BCL2*. The Specimen Control Size Ladder master mix targets multiple genes and generates a series of amplicons of approximately 100, 200, 300, 400, and 600 base pairs to ensure that the quality and quantity of input DNA is adequate to yield a valid result. This assay includes negative control DNA and positive control DNAs for both the MBR and mcr. PCR products can be analyzed using standard gel electrophoresis with ethidium bromide staining. A *BCL2* translocation is indicated if any one of the master mixes generates product(s) within the valid size range.

Background

BCL2 translocations are reciprocal chromosome exchanges that place the *bcl*-2 proto-oncogene, located on chromosome 18, under aberrant transcriptional control of the immunoglobulin heavy chain gene, located on chromosome 14. The *bcl*-2 protein is an antagonist to apoptosis (programmed cell death), a normal process designed to eliminate unneeded and damaged cells during hematopoiesis. Increased expression of the *bcl*-2 protein leads to an increase in the levels of B-cells in the body.

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. PAS Evans et al., Leukemia 17: 2298-2301 (2003).
- 2. JJM van Dongen et al., Leukemia 17: 2257-2317 (2003).



t(14;18) tube A: 2 *BCL2* MBR primers + 1 J_{H} primer

t(14;18) tube B: 4 BCL2 3'MBR primers + 1 J_H primer

t(14;18) tube C: 3 *BCL2* mcr primers + 1 J_H primer

Figure Legend: Schematic diagram of the *IGH-BCL2* t(14;18) translocation showing the *BCL2* gene on the left and the Ig heavy chain (*IGH*) gene on the right. Shown are the relative positions and orientations for the major breakpoint region (MBR) primers, the minor cluster region (mcr) primers, and the JH primer, which are included in the 3 *BCL2/*JH master mix tubes.

EuroClonality This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.

Reagents

Controls

IVS-0030 Clonal Control DNA IVS-P002 Clonal Control DNA IVS-0031 Clonal Control DNA IVS-0000 Polyclonal Control DNA

Master Mixes

BCL2/Jн Tube A BCL2/Jн Tube B BCL2/Jн Tube C Specimen Control Size Ladder

	Concentration	Units in Assay	Units in Assay MegaKit
Д	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
Δ.	1600 pg/mL	1 x 100 µL tube	5 x 100 µL tubes
A Contraction of the second se	200 µg/mL	1 x 100 µL tube	$5 \times 100 \ \mu L \ tubes$
DNA	200 µg/mL	1 x 100 µL tube	$5 \times 100 \ \mu L \ tubes$
	Target	Units in Assay	Units in Assay MegaKit
	BCL2 MBR + IGH JH	1 x 1500 µL tube	10 x 1500 µL tubes
	BCL2 3' MBR + IGH JH	1 x 1500 µL tube	10 x 1500 µL tubes
	BCL2 mcr + IGH JH	1 x 1500 µL tube	10 x 1500 µL tubes
er	Multiple Genes	1 x 1500 µL tube	10 x 1500 µL tubes

Ordering Information

Catalog

1-309-0020 1-309-0040 **Products** BCL2/J_H Translocation Assay – Gel Detection BCL2/J_H Translocation Assay MegaKit – Gel Detection Quantity

33 reactions330 reactions

BCR/ABL t(9;22) Translocation Assays

Assay Use

This Research Use Only assay identifies BCR/ABL t(9;22) translocations.

Summary and Explanation of the Test

The master mixes are included in these assay kits used to amplify complementary DNA (cDNA) produced from specimen(s), and positive and negative RNA controls (included). Primers target an internal control transcript (*AbI*) and p190-, p210-, and p230-type transcripts expressed from *BCR-ABL1* translocations. Amplicon products can be analyzed by capillary electrophoresis or standard gel electrophoresis with ethidium bromide staining. A *BCR-ABL1* translocation is indicated if just one of the 2nd round master mixes (Mix 2b, Mix 2c, Mix 3b, Mix 3c, or Mix 3d) generates product(s) of the valid size. Reagents for RNA extraction and reverse transcription are not included. This assay is compatible with all standard RNA extraction and cDNA synthesis methods. This is a qualitative assay and has not been validated for quantitative use.

Background

BCR/ABL translocations are associated with a variety of hematologic malignancies. The Philadelphia chromosome (Ph1) is a specific chromosomal abnormality that results from reciprocal t(9;22)(q34;q11) chromosomal rearrangements fuse coding regions of the BCR gene, located on chromosome 22, with the *ABL* receptor-independent tyrosine kinase gene on chromosome 9. This assay detects and identifies the variety of p190-, p210- and p230-type transcripts produced from all known *BCR/ABL* translocations.

Specimen Requirements

This assay tests complementary DNA (cDNA) template.

References

- 1. R Kurzrock et al., Ann. Intern. Med. 138: 819-30 (2003).
- 2. JV Melo. Blood 88: 2375-2384 (1996).
- 3. JP Radich et al., *Blood* 85: 2632-2638 (1995).



Figure Legend: This figure shows the genomic organization of the BCR and ABL genes on chromosomes 22 and 9, respectively. Boxes represent exon regions of the ABL (red boxes) and BCR encoding exons (other colors). The solid black line represents intron regions, which have been left incompletely spliced to assist in demarcation of the exon segments. The location of exon regions targeted by labeled and unlabeled primers are indicated by arrows. A p210-type BCR-ABL1 translocation (b3a2 junction) is depicted in the lower portion of the figure along with the control ABL transcript control.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0032 Clonal Control RNA	400 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0011 Clonal Control RNA	400 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0035 Clonal Control RNA	400 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
BCR/ABL t(9;22) Mix 1a	Abl	1 x 1500 µL tube	10 x 1500 µL tubes
<i>BCR/ABL</i> t(9;22) Mix 2a	p190	1 x 1500 µL tube	10 x 1500 µL tubes
<i>BCR/ABL</i> t(9;22) Mix 3a	p210+230	1 x 1500 µL tube	10 x 1500 µL tubes
<i>BCR/ABL</i> t(9;22) Mix 1b	Abl	1 x 1500 µL tube	10 x 1500 µL tubes
<i>BCR/ABL</i> †(9;22) Mix 2b	p190	1 x 1500 µL tube	10 x 1500 µL tubes
<i>BCR/ABL</i> t(9;22) Mix 2c	p190	1 x 1500 µL tube	10 x 1500 µL tubes
<i>BCR/ABL</i> t(9;22) Mix 3b	p210+230	1 x 1500 µL tube	10 x 1500 µL tubes
<i>BCR/ABL</i> t(9;22) Mix 3c	p210+230	1 x 1500 µL tube	10 x 1500 µL tubes
<i>BCR/ABL</i> t(9;22) Mix 3d	p210+230	1 x 1500 µL tube	10 x 1500 µL tubes

Ordering Information

Catalog #	Products	Quantity
1-310-0010	BCR/ABL t(9;22) Translocation Assay - Gel Detection	30 reactions
1-310-0031	BCR/ABL t(9;22) Translocation Assay – ABI Fluorescence Detection	30 reactions

PML/RARa t(9;22) Translocation Assays

Assay Use

This Research Use Only assay identifies PML/RARa t(15;17) translocations.

Summary and Explanation of the Test

Four master mixes are included in these assay kits. Master mixes are used to amplify complementary DNA (cDNA) produced from specimen(s), as well as positive and negative RNA controls (included). Primers target an internal control transcript and the variety of Bcr1, Bcr2, and Bcr3 type transcripts expressed from *PML-RARa* translocations. Amplicon products can be analyzed by differential fluorescence detection using capillary electrophoresis. A *PML-RARa* translocation is indicated if just one of the 2nd round master mixes (Mix 2b or Mix 2c) generates product(s) of the valid size. Reagents for RNA extraction and reverse transcription are not included. This assay is compatible with all standard RNA extraction and cDNA synthesis methods. This is a qualitative assay and has not been validated for quantitative use.

Background

Three *PML/RARa* translocation patterns have been identified in samples with acute myelogenous leukemia (AML): type A is the short (S-form); the breakpoint occurs within breakpoint cluster region 3 (Bcr-3). Type B is the long (L-form); the breakpoint occurs within Bcr-1. There is a third type B variant or variable (V-form); the breakpoint is within Bcr-2. Identification of the *PML/RARa* t(15;17) rearrangements are commonly used in the study of APL because it is correlated with responsiveness to treatment. This RT-PCR method directly identifies the chimeric *PML/RARa* transcripts expressed from all three forms of *PML/RARa* translocations.

Specimen Requirements

This assay tests complementary DNA (cDNA) template.

References

- 1. H De Thé et al., Nature 347: 558-561 (1990).
- 2. H De Thé et al., Cell 66: 675-684 (1991).
- 3. A Kakizuka et al., Cell 66: 663-674 (1991).
- 4. WH Miller et al., Proc. Natl. Acad. Sci. 89: 2694-2698 (1992).



Figure Legend: This figure shows the genomic organization of the PML and RARa genes on chromosomes 15 and 17, respectively. Boxes represent exon regions of the PML (red boxes) and RARa (orange) encoding exons. The solid black line represents intron regions, which were left incompletely spliced to assist in demarcation of the exon segments. Primers are indicated by arrows, and the size of several of the products are indicated below the translocated gene segments. S-form (Bcr3) and L-form (Bcr1) PML-RARa translocations are depicted in the lower portion of the figure.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0020 Clonal Control RNA	400 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0035 Clonal Control RNA	400 μg/mL 1 x 100 μL tube 5 x 100 μL tube		5 x 100 µL tubes
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
PML/RARα t(15;17) Mix 1	RARA	1 x 1500 µL tube	10 x 1500 µL tubes
<i>PML/RAR</i> α t(15;17) Mix 2α	PML-RARA	1 x 1500 µL tube	10 x 1500 µL tubes
PML/RARa †(15;17) Mix 2b	S- and L-Forms	1 x 1500 µL tube	10 x 1500 µL tubes
<i>PML/RAR</i> α t(15;17) Mix 2c	L-Form	1 x 1500 µL tube	10 x 1500 µL tubes

Ordering Information

Products

Catalog #

PML/RARat(15;17) Translocation Assay - ABI Fluorescence Detection

Quantity

30 reactions

IGH Somatic HyperMutation Assays v2.0

Assay Use

The Research Use Only *IGH* Somatic Hypermutation Assay v2.0 is used to identify clonal rearrangements of the immunoglobulin heavy (*IGH*) chain gene and determine the somatic mutation status of the variable (V) gene sequence and is useful for the study of:

- 1. Identifying clonal rearrangements of the IGH chain gene
- 2. Assessing the extent of somatic hypermutation in the variable region of the immunoglobulin heavy chain gene
- 3. Evaluating new research and methods in malignancy studies

Summary and Explanation of the Test

These assays amplify either genomic DNA or complementary DNA (cDNA) that lies between the upstream leader (VHL) or framework 1 (FR1) regions and the downstream joining (JH) region of the *IGH* gene. The assays employ two different master mixes: Hypermutation Mix 1 and Hypermutation Mix 2. The Hypermutation Mix 1 targets sequences between the leader (VHL) and joining (JH) regions. Therefore the amplicon product(s) span the entire variable (VH) region, which contains all framework (FR) and complementarity-determining regions (CDR). The Hypermutation Mix 2 targets sequences between the framework 1 (FR1) and joining (JH) regions. The resulting amplicons include a portion of the FR1 region to the downstream JH region.

The primers that target the VHL and FR1 regions have been redesigned to include a universal sequencing tag at the 5'end. This design allows for bi-directional sequencing of clonal PCR products with just one sequencing-tag specific forward primer and one JH reverse primer, thus ensuring a more reliable and complete coverage of clonal products. Positive and negative DNA, positive RNA, as well as an amplification control are included in the assay. Clonality is indicated if any one of the master mixes generates clonal products.

Background

Rearrangements of the antigen receptor genes occur during ontogeny in B and T lymphocytes. These gene rearrangements are unique in length and sequence for each cell. Therefore, polymerase chain reaction (PCR) assays can be used to identify lymphocyte populations derived from a single cell by detecting the unique V-J gene rearrangements present within these antigen receptor loci.¹ This PCR-based assay employs multiple consensus DNA primers that target conserved genetic regions within the immunoglobulin heavy chain (*IGH*) gene. This test is used to detect and sequence the majority of clonal *IGH* rearrangements from either genomic DNA (gDNA) or complementary DNA (cDNA). Clonal products can be detected using a variety of methods, including gel and capillary electrophoresis. The presence of *IGH* somatic hypermutation (SHM) is defined as greater or equal to 2% difference from the germline variable (V) gene sequence, whereas less than 2% difference is considered evidence of no somatic hypermutation.³

Specimen Requirements

This assay tests extracted and purified genomic DNA (gDNA).

References

- 1. P Ghia et al., Leukemia 21: 1-3 (2007).
- 2. P Ghia et al., Blood 105: 1678-1685 (2005).
- 3. F Davi et al., *Leukemia* 22: 212–214 (2008).



- Universal Sequencing Tag

Figure Legend: Simple representation of the organization of a rearranged immunoglobulin heavy chain gene on chromosome 14. Black arrows represent the relative positions of primers that target the conserved Leader (L) and Framework 1 (FR1) regions, and the downstream consensus J_H gene segments.

Reagents

Controls	Concentration	Units in Assay	Units in Assay MegaKit
IVS-0013 Clonal Control DNA	200 µg/mL	1 x 100 µL tube	5 x 100 µL tubes
IVS-0013 Clonal Control RNA	400 µg/mL	1 x 100 µL tube	$5 \times 100 \ \mu L \ tubes$
IVS-0000 Polyclonal Control DNA	200 µg/mL	1 x 100 µL tube	$5 \times 100 \ \mu L \ tubes$
Master Mixes	Target	Units in Assay	Units in Assay MegaKit
Hypermutation Mix 1 v2.0	Leader + JH	1 x 1500 µL tube	10 x 1500 µL tubes
Hypermutation Mix 2 v2.0	Framework 1 + J _H	1 x 1500 µL tube	10 x 1500 µL tubes
Specimen Control Size Ladder	Multiple Genes	1 x 1500 µL tube	10 x 1500 µL tubes
Primers	Target	Units in Assay	Units in Assay MegaKit
Primer - Hypermutation	Leader + Framework 1	1 x 10 μ L tube at 100 μ M	5 x 10 µL tube
IGH J _H Primer	Јн	1 x 10 μ L tube at 100 μ M	5 x 10 µL tube

Ordering Information

Catalog #	Products	Quantity
5-101-0030	IGH Somatic Hypermutation Assay v2.0 - Gel Detection	33 reactions
5-101-0040	IGH Somatic Hypermutation Assay v2.0 MegaKit - Gel Detection	330 reactions
5-101-0031	IGH Somatic Hypermutation Assay v2.0 - ABI Fluorescence Detection	33 reactions
5-101-0041	IGH Somatic Hypermutation Assay v2.0 MegaKit - ABI Fluorescence Detection	330 reactions

FLT3 Mutation Assays

Assay Use

These Research Use Only assays identify FLT3 mutations.

Summary and Explanation of the Test

FLT3 Mutation Assays target regions of the *FLT3* gene to identify internal tandem duplication (ITD) mutations and tyrosine kinase domain (TKD) mutations, such as the D835 and I836 mutations. DNA is amplified by PCR, TKD amplicon is enzymatically digested, and *FLT3* mutations are detected via agarose gel (catalog #14120010) or capillary (catalog #14120031) electrophoresis.

Assay kits include three PCR master mixes, along with positive and negative controls. *FLT3* ITD master mix tests for internal tandem duplication mutations. *FLT3* D835 master mix tests for TKD region mutations. The Specimen Control Size Ladder master mix targets multiple genes and generates a series of amplicons of approximately 100, 200, 300, 400, and 600 base pairs to ensure that the quality and quantity of input DNA is adequate to yield a valid result.

Master mixes contain fluorophore-labeled (capillary) or unlabeled (gel) primer sets as appropriate to kit detection method.

Background

FLT3 is a receptor tyrosine kinase that is normally expressed on many cell types including hematologic stem cells. Mutation of the *FLT3* receptor, by either internal tandem duplication (ITD) of the juxtamembrane domain or point mutation in the activation loop of the tyrosine kinase domain (TKD), causes constitutive activation of the *FLT3* receptor.

Such gain-of-function mutations in the FMS related tyrosine kinase 3 (*FLT3*) gene are the subject of research studies and multiple clinical trials targeting Acute Myeloid Leukemia (AML) subjects. The most prevalent type of *FLT3* mutation is an internal tandem duplication in and around the juxtamembrane domain. The second most common mutation type in the *FLT3* gene is a TKD point mutation in aspartate (D835) or isoleucine (1836).

Specimen Requirements

High quality genomic DNA

References

https://clinicaltrials.gov

- 1. Acute Myeloid Leukemia, Clinical Practice Guidelines in Oncology, National Comprehensive Cancer Network (v.2.2014).
- 2. Lowenberg, B. et al. "Acute myeloid leukemia." N Engl J Med 341(14):1051-62 (1999).
- 3. Thiede, C. et al. "Analysis of *FLT3*-activating mutations in 979 patients with acute myelogenous leukemia: association with FAB and identification of subgroups with poor prognosis." *Blood* 99(12): 4326-35 (2002).
- 4. Nakao, M. et al. "Internal tandem duplication of the FLT3 gene found in acute myeloid leukemia." Leukemia 10(12):1911-18 (1996).
- 5. Yamamoto, Y et al. Activating mutation of D835 within the activation loop of *FLT3* in human hematologic malignancies. *Blood*, 97(8):2434–9 (2001). 6. Gilliland, DG et al. The roles of *FLT3* in hematopoiesis and leukemia. *Blood* 100(5):1532–154 (2002).

These products are not available for sale or use in the United States.

Reagents

Controls	Concentration	Units in 33 Reaction Assay
IVS-0017 Clonal Control DNA	200 µg/mL	1 x 100 µL tube
IVS-P001 Clonal Control DNA	200 µg/mL	1 x 100 µL tube
IVS-0000 Polyclonal Control DNA	200 µg/mL	1 x 100 µL tube
Master Mixes	Target	Units in 33 Reaction Assay
FLT3 ITD Master Mix	FLT3 ITD	1 x 1500 µL tube
FLT3 D835 Master Mix	FLT3 TKD	1 x 1500 µL tube
FLT3 D835 Master Mix Specimen Control Size Ladder	<i>FLT3</i> TKD Multiple Genes	1 x 1500 μL tube

Gel Electrophoresis Detection

Data was generated using the *FLT3* ITD Master Mix, and amplified products were run on a 2% TBE agarose gel alongside a 100bp DNA size ladder. Lane 1 is a *FLT3* ITD control*; lane 2 is a 10% dilution of a *FLT3* ITD control; and lane 3 is IVS-0000, which is representative of a wild type product.

*IVS-0050 performs comparable to IVS-0017 clonal control DNA, which is included in the kit as the positive control.



Capillary Electrophoresis Detection (ABI)

Data was generated using the *FLT3* ITD Master Mix and amplified products were run on an ABI 3100 instrument. Templates analyzed follow: Panel 1 is the *FLT3* ITD positive control; panel 2 is a 10% dilution of the positive control; and Panel 3 is IVS-0000, which is representative of a wild type product.

Ordering Information

Catalog #	Products
1-412-0010	FLT3 Mutation Assay - Gel Detection
1-412-0031	FLT3 Mutation Assay - ABI Fluorescence Detection

Quantity

33 reactions 33 reactions *FLT3* ITD Master Mix Lane 1 = 100% IVS-0050 Lane 2 = 10% IVS-0050 Lane 3 = 100% IVS-0000

> Mutant: 360 bp 🔶 Wild Type: 330 bp 🔶



These products are for Research Use Only (RUO). Not intended for diagnostic purposes., and additionally not available for sale or use in regions where CE-IVD products are registered.

These products are not available for sale or use in the United States. Please see the LeukoStrat CDx *FLT3* Mutation Assay product pages for IVD product availability in the United States and other global regions.

Analyte Specific Reagents

Invivoscribe[®] Analyte Specific Reagents (ASRs) target B- and T-cell antigen receptor loci, or chromosome translocations (*IGH - BCL2*, *BCR-ABL1*, *PML-RARA*).

LEARN MORE: invivoscribe.com/products/reagents

Analyte Specific Reagents

Invivoscribe Analyte Specific Reagents (ASRs) target B- and T-cell antigen receptor loci.

The ASRs are available as a single tube containing a volume of 1500 µL. To ensure the highest quality and reliability of reagents ASRs are manufactured under cGMP and ISO 13485 standards.

Per the current US FDA regulations, ASRs may only be sold to *in vitro* diagnostic manufacturers, CLIA accredited high complexity laboratories, VHA regulated clinical laboratories, and laboratories not intending to use ASRs as a part of a diagnostic test.

Analytical and performance characteristics are not established. ASRs are not available for sale or use outside of the USA.

IGH (Immunoglobulin Heavy Chain Gene Locus)

Description IGH Framework 1 – 6FAM IGH Framework 2 – 6FAM IGH Framework 3 – HEX IGH FR1 – 6FAM IGH FR2 – 6FAM IGH FR3 – HEX IGH DH1 – 6 – HEX IGH DH7 – 6FAM Catalog # A-101-0061 A-101-0091 A-101-0081 A-101-0011 A-101-0101 A-101-0031 A-101-0041 A-101-0051

TRB (T-Cell Receptor Beta Chain Gene Locus)

Description

TCRB V - J1 + 2 - 6FAM & HEX TCRB V - J2 - 6FAM TCRB D - J1 + 2 - 6FAM & HEX Catalog # A-205-0011 A-205-0021 A-205-0031

IGK (Immunoglobulin Kappa Light Chain Gene Locus)

Description *IGK* V - | - 6FAM

IGKV - K_{de} - 6FAM

Catalog # A-102-0011 A-102-0021

TRG (T-Cell Receptor Gamma Chain Gene Locus)

Description

TCRG V(2-5,8-11) J 1 + 2+P - 6FAM TCRG V(1-8,9) J - 6FAM TCRG V(1-8) J - HEX Catalog # A-207-0091 A-207-0071 A-207-0021

Controls, Reagents, and Enzymes

Invivoscribe[®] offers an extensive range of General Purpose Reagents (GPRs) and Research Use Only (RUO) nucleic acid controls.

LEARN MORE: invivoscribe.com/products/controls

For Legal Notices related to this product, visit: https://invivoscribe.com/legal-notices/

Introduction

Controls are available in DNA extracted from tissue or cell lines, or RNA extracted from cell lines. These controls can be purchased in various dilutions or as complete dilution sets and panels for several purposes, such as to help with assay validation, sensitivity or proficiency testing, or troubleshooting.

The following pages will provide an overview of available controls, along with a number of tables and reference guides, to help you decide which Invivoscribe control(s) will be suitable for your application.

Chapter Contents

DNA Controls

136 Quick Reference for DNA Controls139 LymphoTrack Low Positive Controls139 LymphoQuant Internal Controls

Control Panels

143 RNA Sensitivity Panels

Mastermix Controls

144 Specimen Control Size Ladder

RNA Controls

140 Quick Reference for RNA Controls **142** *BCR/ABL1* RNA Dilution Sets

Reagents

144 ABI Detection Reagents

Enzyme 144 FalconTaq DNA Polymerase

DNA Controls

Every laboratory needs suitable controls (positive and negative) for sensitivity and proficiency testing, as well as for troubleshooting. Since patient samples cannot serve as true controls (due to a lack of characterization and inter-sample variability), Invivoscribe offers a multitude of high quality, reliable DNA controls manufactured under cGMP conditions.

These controls can be used for most assays targeting B- and T-cell antigen receptor loci, or IGH-BCL2, BCR-ABL1, and PML-RARa chromosome translocations.

Quick Reference for DNA Controls

The vast majority of our high-quality DNA controls, including sensitivity controls and panels, are supplied in aliquots of 100 µL and are adjusted to a final concentration of 200 µg/mL in 1/10 TE (1 mM Tris- HCl (pH 8.0), 0.1 mM EDTA).

Positive for	lmr Rea	nunogla Irrange	bulin ments	Mutations	Transloc	ations	Rece Reari	T-Cell eptor Ge rangeme	ne ents
	IGH	IGK	IGL	IGHV SHM	IGH-CCND1**	IGH-BCL2	TRB	TRG	TRD
IVS-0001						•			
IVS-0004							•	•	
IVS-0007	•	٠	•			•			
IVS-0008€	0						•	•	
IVS-0009							•	•	
IVS-0010	•	•			•				
IVS-0013	•	٠							
IVS-0019	•	٠							
IVS-0021		٠					•	•	•
IVS-0024	٠	•							
IVS-0029	٠	•							
IVS-0030 ‡	•	٠		•					
IVS-0031	•	٠							
LymphoQuant B-cell Internal Control	~	~							
LymphoQuant T-cell Internal Control							~	~	
LymphoTrack B-cell Low Positive Control	~	~							
LymphoTrack T-cell Low Positive Control							~	~	

^tThese controls can be used as SHM positive controls with $\geq 2\%$ mutational rates compared to the germline sequence. ^cThis control does not contain a complete *IGH* V_H-J_H rearrangement and may only be suitable for *IGH* D_H-J_H rearrangements. ***IGH-CCND1* was previously referred to as *BCL1/J*H. Gene rearrangement
Partial IGH D_H-J_H rearrangement

∼ Recommended for NGS

Tissue DNA

Standard Concentrations

Our high-quality DNA controls are supplied in aliquots of 100 μ L and are adjusted to a final concentration of 200 μ g/mL in 1/10 TE (1 mM Tris- HCl (pH 8.0), 0.1 mM EDTA). This diluent provides sufficient buffering capacity and EDTA to protect the DNA without interfering with the Mg²⁺ concentrations required for robust amplification reactions.

IVS-0000 Polyclonal Control DNA

Tissue DNA controls are extracted from normal, disease-free tissue and are tested extensively to ensure quality and reproducibility of your test results. IVS-0000 Polyclonal Control DNA consists of genomic DNA isolated from the tissue of normal human tonsils. This control represents an excellent negative control for gene rearrangements, chromosome translocations, and mutation tests and is included in all of our PCR DNA-based assay kits. This DNA is supplied at a volume of 100 μ L and at a concentration of 200 μ g/mL.

Catalog

4-092-0010

IVS-0000 Polyclonal Control DNA

Description

*These controls are general purpose reagents (GPRs).

Cell Line DNA

Reliable Positive Controls

Cell Line DNA controls are extracted from established cell lines grown under cell culture conditions recommended by the supplier. Our controls are tested extensively to ensure quality and reproducibility of your test results. Please note, these controls are for qualitative use only.

Note: n/c is used to indicate that the control has not been fully characterized; there may be additional rearrangements, translocations or mutations associated with the control.

Standard Concentrations

Our high-quality DNA controls are supplied in aliquots of 100 µL and are provided at a final concentration of 200 µg/mL in 1/10 TE (1 mM Tris-HCl (pH 8.0), 0.1 mM EDTA). This diluent provides sufficient buffering capacity and EDTA to protect the DNA controls without interfering with the Mg²⁺ concentrations required for robust amplification reactions.

IVS-0001 Clonal Control DNA

IVS-0001 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

n/c IGH-BCL2 t(14;18) mcr n/c

Catalog # 4-088-0010

IVS-0001 Clonal Control DNA

Description

IVS-0004 Clonal Control DNA

IVS-0004 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

TRB. TRG n/c n/c

Catalog

4-088-0190

Description

IVS-0004 Clonal Control DNA

IVS-0007 Clonal Control DNA

IVS-0007 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

IGH, IGK, IGL IGH-BCL2 t(14;18) Mbr n/c

Catalog

4-088-0370

Description IVS-0007 Clonal Control DNA

IVS-0008 Clonal Control DNA[‡]

IVS-0008 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements:	IGH
Chromosome Translocations:	n/c
Mutations:	n/c

DH-H[‡], TRB, TRG

Catalog # Description

4-088-0430 IVS-0008 Clonal Control DNA

IVS-0009 Clonal Control DNA

IVS-0009 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

TRB. TRG n/c n/c

Cataloa # Description

4-088-0490

IVS-0009 Clonal Control DNA

These controls are general purpose reagents (GPRs). * This control does not contain a complete IGH V $_{\rm H}$ -J $_{\rm H}$ rearrangement and may only be suitable for $IGH D_{H}-J_{H}$ rearrangements.

Cell Line DNA

IVS-0010 Clonal Control DNA

IVS-0010 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

IGH. IGK. IGL IGH-BCL1 t(11:14) n/c

Catalog # 4-088-0550

Description

IVS-0010 Clonal Control DNA

IVS-0013 Clonal Control DNA

IVS-0013 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

IGH, IGK, IGL n/c n/c

Catalog

4-088-0730

Description IVS-0013 Clonal Control DNA

IVS-0019 Clonal Control DNA

IVS-0019 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

IGH, IGK n/c n/c

Catalog #

4-088-1090

Description IVS-0019 Clonal Control DNA

IVS-0021 Clonal Control DNA

IVS-0021 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

TRB, TRD, TRG n/c n/c

Description

Catalog # 4-088-1210

Catalog #

4-088-1390

IVS-0021 Clonal Control DNA

IVS-0024 Clonal Control DNA

IVS-0024 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

IGH, IGK n/c n/c

Description IVS-0024 Clonal Control DNA

IVS-0029 Clonal Control DNA

IVS-0029 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements:
Chromosome Translocations:
Mutations:

IGH, IGK, IGL n/c n/c

IGH

n/c

IGH SHM

Catalog

Description

4-088-1690

IVS-0029 Clonal Control DNA

IGH SHM Positive Control DNA

IGH SHM Postive Control can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

Catalog 4-088-0008 Description

IGH SHM Positive Control DNA

IVS-0030 Clonal Control DNA

IVS-0030 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

IGH. IGK IGH-BCL2 t(14;18) Mbr n/c

Catalog 4-088-1750

Description IVS-0030 Clonal Control DNA

IVS-0031 Clonal Control DNA

IVS-0031 Clonal Control DNA can be used as a positive control for:

Gene Rearrangements: Chromosome Translocations: Mutations:

IGH. IGK IGH-BCL2 t(14;18) mcr n/c

Catalog

Description

4-088-1810

IVS-0031 Clonal Control DNA

These controls are general purpose reagents (GPRs).

LymphoTrack[®] Low Positive Controls

Measurable Residual Disease (MRD) testing is a valuable tool that allows investigators to study and monitor multiple myeloma (MM), chronic lymphocytic leukemia (CLL), acute lymphoblastic leukemia (ALL), acute myelogenous leukemia (AML) and other hematologic diseases. Recent treatment advances have led to significantly increased clinical response and overall survival, but ultimately most subjects will relapse, driving the need for sensitive MRD monitoring. Sensitive and standardized testing such as NGS-based MRD may one day enable identification of those cases that will eventually relapse versus those who are potentially cured. In addition to the need for more sensitive tracking, it is clear that standardized methods are needed. Currently, MRD methods are highly subjective and recommendations are often based on consensus expert-shared knowledge and experience, not on a validated, objective method. Once specific rearrangements have been identified, LymphoTrack assays can be used with LymphoQuant and LymphoTrack Low Positive Controls to track these clonotype populations to a sensitivity as low as 10⁻⁴.

LymphoTrack[®] B-cell Low Positive Control

LymphoTrack B-cell Low Positive Control can be used as a control for:

Gene Rearrangements:
Chromosome Translocations:
Mutations:

IGH n/c n/c

LymphoTrack T-cell Low Positive Control can be used as a control for:

LymphoTrack[®] T-cell Low Positive Control

Gene Rearrangements: Chromosome Translocations: Mutations: *TRB, TRG* n/c n/c

Catalog #

4-088-0098

Description

LymphoTrack[®] B-cell Low Positive Control*

Cata	log #
------	-------

Description

4-088-0108

LymphoTrack[®] T-cell Low Positive Control*

LymphoQuant[®] Internal Controls

LymphoQuant T-cell or B-cell Internal Controls may be spiked into specimens to estimate the respective number of clonotype T-cell or *IGH* equivalents present. Addition of the LymphoQuant Internal Control to the specimen PCR facilitates clonotype tracking over time without any additional sequencing cost. Consistent use of a LymphoQuant Internal Control enables investigators to objectively monitor the disease over time with a highly standardized, sensitive method. The LymphoTrack MRD software helps researchers that use the LymphoQuant Internal Control, calculate and report an estimated number of clonotype cell equivalents and the percent clonotype in the sample, enabling researchers and pharmaceutical companies to accurately monitor hematologic disease in longitudinal studies.

LymphoQuant® B-cell Internal Control

LymphoQuant B-cell Internal Control can be used to objectively track Ig clonotypes.

IGH

n/c

n/c

Gene Rearrangements:	
Chromosome Translocations:	
Mutations:	

Catalog # Description

4-088-0118 LymphoQuant® B-cell Internal Control*

LymphoQuant[®] T-cell Internal Control

LymphoQuant T-cell Internal Control can be used to objectively track TCR clonotypes.

Gene Rearrangements: Chromosome Translocations: Mutations: *TRB, TRG* n/c n/c

Catalog

4-088-0128

Description LymphoQuant® T-cell Internal Control*

*LymphoTrack® Low Positive Controls and LymphoQuant® Internal Controls are research use only (RUO), not for diagnostic procedures.



RNA Controls

Quick Reference for RNA Controls

Reliable Assay Controls

Our RNA controls are extracted from well characterized cell lines grown under standard and carefully controlled culture conditions. The general purpose reagent (GPR) controls are tested to ensure linearity and reproducible results. Since this RNA is extracted from cell lines, these controls can be used with any of the standard housekeeping genes.

Standard Concentrations

Each RNA single control tube (as separate control tube, RNA sensitivity panel and proficiency panel) is supplied in aliquots of 100 μ L at a final concentration of 400 μ g/mL in water. Each *BCR/ABL* RNA dilution set member is supplied in aliquots of 50 μ L at a final concentration of 400 μ g/mL in water. To ensure maximum stability, store the dilutions at -85 °C to -65 °C and minimize the number of freeze-thaw cycles .

RNAs positive for chromosome translocations

Chromosome Translocation	Clonal Control RNA	Chromosome Translocation	Clonal Control RNA
<i>BCR-ABL1</i> †(9;22) p210 e13a2 (b2a2)	IVS-0003	CBFB-MYH11 inv(16)	IVS-0015
<i>BCR-ABL1</i> †(9;22) p210 e14a2 (b3a2)	IVS-0011	<i>E2A-PBX1</i> t(1;19)(q23;p13)	IVS-0002
<i>BCR-ABL1</i> t(9;22) p190 e1a2	IVS-0032	PML-RARA t(15;17)(q22;q11)	IVS-0020

RNAs negative for chromosome translocations

IVS-0035 Clonal Control RNA is negative for BCR-ABL t(9;22) and PML-RARa t(15;17) chromosome translocations.

IVS-0035 may be used as a negative control for other chromosome translocations or diluents for other chromosome translocation positive controls. Please do not hesitate to contact us at support@invivoscribe.com so we can evaluate whether this control may work for your testing needs.

Cell Line RNA

Reliable Positive and Negative Controls

Cell Line RNA controls are extracted from established cell lines grown under cell culture conditions recommended by the supplier. Our GPR controls are tested extensively to ensure quality and reproducibility of your test results. Please note, these controls are for qualitative use only.

Standard Concentrations

Our GMP-manufactured high-quality RNA controls, including sensitivity controls and proficiency panel samples, are supplied in aliquots of 100 µL and are adjusted to a final concentration of 400 µg/mL in RNase-free glass-distilled water. The pH of distilled water is slightly acidic; this protects the RNA from hydrolysis. RNA dilutions are diluted volume to volume in our negative control RNA, IVS-0035 Clonal Control RNA.

IVS-0002 Clonal Control RNA

IVS-0002 Clonal Control RNA can be used as a positive control for the chromosome translocation: *E2A-PBX1* t(1;19) (q23;p13).

Catalog # Description

4-089-0100

IVS-0002 Clonal Control RNA*

IVS-0003 Clonal Control RNA

IVS-0003 Clonal Control RNA can be used as a positive control for the chromosome translocation: *BCR-ABL1* t(9;22) p210 e13a2 (b2a2).

This control is also available as several ready-to-use dilutions into a standard negative control as listed in the table below.

Catalog #	Description
4-089-0190	IVS-0003 Clonal Control RNA*
4-089-0200	10 ⁻¹ IVS-0003 Clonal Control RNA
4-089-0210	10 ⁻² IVS-0003 Clonal Control RNA
4-089-0220	10 ⁻³ IVS-0003 Clonal Control RNA
4-089-0230	10 ⁻⁴ IVS-0003 Clonal Control RNA
4-089-0240	10 ⁻⁵ IVS-0003 Clonal Control RNA

IVS-0011 Clonal Control RNA

IVS-0011 Clonal Control RNA can be used as a positive control for the chromosome translocation: *BCR-ABL1* t(9;22) p210 e14a2 (b3a2).

This control is also available as several ready-to-use dilutions into a standard negative control as listed in the table below.

Catalog #	Description
4-089-0910	IVS-0011 Clonal Control RNA*
4-089-0920	10 ⁻¹ IVS-0011 Clonal Control RNA
4-089-0930	10 ⁻² IVS-0011 Clonal Control RNA
4-089-0940	10 ⁻³ IVS-0011 Clonal Control RNA
4-089-0950	10 ⁻⁴ IVS-0011 Clonal Control RNA
4-089-0960	10 ⁻⁵ IVS-0011 Clonal Control RNA

*These controls are general purpose reagents (GPRs). All others are research use only (RUO).

IVS-0015 Clonal Control RNA

IVS-0015 Clonal Control RNA can be used as a positive control for the chromosome translocation: *CBFB-MYH11* inv(16)

Catalog #	Description
4-089-1270	IVS-0015 Clonal Control RNA*

IVS-0020 Clonal Control RNA

IVS-0020 Clonal Control RNA can be used as a positive control for the chromosome translocation: *PML-RARA* t(15;17) L-Form.

Catalog

4-089-1720

Description

IVS-0020 Clonal Control RNA*

IVS-0032 Clonal Control RNA

IVS-0032 Clonal Control RNA can be used as a positive control for the chromosome translocation: *BCR-ABL1* t(9;22) p190 e1a2.

This control is also available as several ready-to-use dilutions into a standard negative control as listed in the table below.

Catalog #	Description
4-089-2800	IVS-0032 Clonal Control RNA*
4-089-2810	10 ⁻¹ IVS-0032 Clonal Control RNA
4-089-2820	10 ⁻² IVS-0032 Clonal Control RNA
4-089-2830	10 ⁻³ IVS-0032 Clonal Control RNA
4-089-2840	10 ⁻⁴ IVS-0032 Clonal Control RNA
4-089-2850	10 ⁻⁵ IVS-0032 Clonal Control RNA

IVS-0035 Clonal Control RNA

IVS-0035 Clonal Control RNA can be used as a negative control for *BCR-ABL1* t(9;22) and PML-RARa t(15; 17) chromosome translocations.

Catalog #	D
4-089-3070	IVS

Description IVS-0035 Clonal Control RNA*

-003-3070					

BCR/ABL RNA Dilution Sets

Our *BCR/ABL* b2a2, b3a2, and e1a2 RNA Dilution Sets consist of RNA that has been extracted from *BCR-ABL1* expressing and *BCR-ABL1* negative cell lines. Each set is composed of several dilutions (10^{-1} , 10^{-2} , 10^{-3} , 10^{-4} , 10^{-5}) of the *BCR-ABL1* positive RNA diluted (v/v) into RNA purified from a cell line that does not contain a *BCR-ABL1* translocation. Also included in these sets is a 100% *BCR-ABL1* negative RNA.

The individual *BCR/ABL* b2a2, b3a2, and e1a2 RNA Dilution Sets can be used as reference and validation materials with assays that target the main transcripts of *BCR-ABL1* t(9;22) translocations: p210 (e13a2 (b2a2), e14a2 (b3a2), and p190 (e1a2). These products may be used as the following:

- 1. Routine testing controls for cDNA synthesis, amplification and detection
- 2. Controls to establish a standard reference curve
- 3. Proficiency controls
- 4. Sensitivity controls for specific target assays

Data

Plot of Ct values (5 replicates) for the 10⁻¹, 10⁻², 10⁻³, 10⁻⁴, and 10⁻⁵ dilutions.



RNA Control Panels

RNA Sensitivity Panels

RNA sensitivity panels are 7 member panels that consist of 100% clonal RNA extracted from a positive control cell line and 10^{-1} , 10^{-2} , 10^{-3} , 10^{-4} , 10^{-5} , and 10^{-6} (1:10 – 1:1 000 000) dilutions of the positive clonal RNA diluted (v/v) into our standard negative control RNA, IVS-0035 Clonal Control RNA. Each tube contains 100 μ L of RNA at 400 μ g/mL in RNase-free glass-distilled water. The pH of distilled water is slightly acidic thereby protecting the RNA from hydrolysis. Please note, these controls are for qualitative use only.

RNA Sensitivity Panels

Catalog #	Description	Can be used as a positive control for:
4-087-0030	Sensitivity Panel for IVS-0003 Clonal Control RNA	<i>BCR-ABL1 t(9;22)</i> p210 e13a2 (b2a2)
4-087-0110	Sensitivity Panel for IVS-0011 Clonal Control RNA	<i>BCR-ABL1 t(9;22)</i> p210 e14a2 (b3a2)
4-087-0150	Sensitivity Panel for IVS-0015 Clonal Control RNA	CBFB/MYH11 inv16
4-087-0200	Sensitivity Panel for IVS-0020 Clonal Control RNA	PML-RARA t(15;17) L-form
4-087-0320	Sensitivity Panel for IVS-0032 Clonal Control RNA	<i>BCR-ABL1</i> t(9;22) p190 e1a2

Master Mix Controls

These master mixes serve as control for many of our DNA assays to ensure that sample DNA is of sufficient quality and integrity to generate a valid result.

Specimen Control Size Ladder

Our Specimen Control Size Ladder master mix targets four different housekeeping genes producing products of approximately 100, 200, 300, 400, and 600 base pair in size to ensure that the quality and quantity of the sample DNA is adequate to yield a valid result with the specific assay(s). This master mix is based on the BIOMED-2 Concerted Action BMH4-CT98-3936 from the EuroClonality Group and is available for Gel Detection (unlabeled) or ABI detection (labeled with 6FAM).

Catalog #	Description	
2-096-0020	Specimen Control Size Ladder - Unlabeled	
2-096-0021	Specimen Control Size Ladder - 6FAM	
These master mixes are general purpose reagents (GPRs).		

Reagents

ABI Detection Reagents

Invivoscribe also offers highly deionized (Hi-Di) Formamide with ROX size standards for ABI fluorescence detection with the ABI 310 or 3100 series. Hi-Di Formamide is used to stabilize single strands of denatured PCR amplicons. The ROX size standards are fluorescent labeled DNA standards which cover the 50 to 400 base pair size range. Sizes of the individual standards are: 50, 60, 90, 100, 120, 150, 160, 180, 190, 200, 220, 240, 260, 280, 290, 300, 320, 340, 360, 380, and 400 base pair.

For samples tested on an ABI 3100 series, we recommend using 10 µL of the Hi-Deionized Formamide with ROX Size Standards mixture for each microliter of PCR product. Please note that the ABI 3100 series require different concentrations of ROX size standards and the different Hi-Deionized Formamide with ROX Size Standards cannot be used interchangeably.

For samples tested on an ABI 3500 series, GeneScan[™] 600[®] LIZ dye Size Standard v2.0 can be purchased from Thermo Fisher Scientific.

Catalog #

6-098-0061

Available through Thermo Fisher Scientific®: 4408399

Description

Hi-Deionized Formamide with ROX Size | Standard (ABI 3100), 1 mL GeneScan^{**} 600 LIZ[®] dye v2.0 | Standard (ABI 3500), 800 reactions

Enzyme

FalconTaq[™] DNA Polymerase

FalconTaq[®] DNA Polymerase can be used for amplification using PCR to obtain high specificity, sensitivity, and yield. This enzyme has been proven to minimize extension of non-specifically bound primers. Generate reliable results by using FalconTaq[®] DNA Polymerase for robust performance.

Catalog

Description

6-097-0130 FalconTaq[™] DNA Polymerase 250 U, 5 U/μL

FalconTaq $^{\mbox{\tiny VI}}$ DNA Polymerase is replacing EagleTaq DNA Polymerase, please contact support@invivoscribe.com if you have questions.
Custom Products

The Invivoscribe® team of experts can help develop your ideas into customized products. Allow us to partner with you to take a basic concept through design, development, validation, regulatory approval, and release.

LEARN MORE: invivoscribe.com/clinical-lab-services



Customized Products to Meet Your Needs

The Invivoscribe team of experts can help develop your ideas into customized products. Allow us to partner with you to take a basic concept through design, development, validation, regulatory approval (if applicable), and release. For more information, please call our San Diego office at +1 858.224.6600 or send an email to inquiry@invivoscribe.com.

Custom Designed Assays

In response to the FDA announcing its intention to dramatically expand its regulatory oversight of laboratory developed tests (LDTs), Invivoscribe is partnering with laboratories worldwide to help facilitate the conversion of LDTs into FDA-cleared assays, as we know the barriers to bringing new assays online are often the availability of resources and the cost of validation.

By leveraging the power of our regulatory expertise, provided through each milestone, we can help ensure safety, efficacy and quality. Our customizable reagent manufacturing capabilities can reduce your LDT costs and lead to higher-quality testing.

To date, we have partnered with more than 40 laboratories around the world to develop, validate, and launch a variety of molecular products. A number of these partnerships have also led to the release of US and CE-marked *in vitro* diagnostic products and services. Learn how Invivoscribe can help you develop assays for new products, services, and novel applications.

Custom Controls and Validation Panels

We offer a large selection of well-characterized DNA and RNA controls that are used to define the performance characteristics of a wide variety of molecular reagents. To address your specific requirements, we can partner with you to design, validate, and provide custom controls and validation panels. If necessary, we are willing to acquire, characterize, and engineer custom controls for your specific application. We can produce DNA, RNA, or cDNA at any specified concentration, dilution, or volume. Please contact us with your requirements and we will be happy to provide controls to suit your needs.

Invivoscribe is a Comprehensive Partner for Companion Diagnostic Development

From biomarker identification through commercialization, Invivoscribe has expertise at every stage of companion diagnostics development.

- Discovery and Patient Stratification: We offer comprehensive gene panels to identify biomarkers and define patient populations, thus reducing development costs and improving the success of clinical trials.
- Clinical Trials: Our network of global laboratories accelerates sample acquisition and harmonizes testing to ensure accurate results.
- Regulatory Approval: Our in-house experts have experience seeking approval with global agencies.
- Commercialization: Our cGMP manufacturing expertise and distribution channels allow approved CDx to reach all global markets.

Invivoscribe is an ISO 13485-accredited and FDA/CDRH registered medical device manufacturer with a long record of successful partnerships. We are the industry-leading assay and software development company, providing full QSR design control and a complete range of cGMP manufactured assays, controls, reagents, and services to CLIA-accredited clinical laboratory and pharmaceutical communities.

Please contact us at **+1858.224.6600** or **inquiry@invivoscribe.com** for more details about partnering with Invivoscribe for the development and manufacturing of companion diagnostics, *in vitro* diagnostics, molecular reagents, and/or nucleic acid controls.

B- or T-cell Clonality Testing Workflow

The Invivoscribe[®] European Conformity marked in vitro diagnostics (CE-IVD) and Research Use Only (RUO) clonality assays detect clonal populations in just a few easy steps.

LEARN MORE: invivoscribe.com/products/controls



Peripheral blood, Tissue, or FFPE sample



Data Interpretation

The Invivoscribe European Conformity marked *in vitro* diagnostics (CE-IVD) and Research Use Only (RUO) clonality assays detect clonal populations in just a few easy steps. These steps include PCR amplification of the immunoglobulin or T-cell receptor genes of interest, followed by detection with non-denaturing polyacrylamide gels, capillary electrophoresis, or next-generation sequencing using an Illumina[®] MiSeq[™], Thermo Fisher Scientific[®] Ion S5[™] or Ion PGM[™] instrument. A flowchart illustrating this workflow is shown below.

*Or equivalent DNA Polymerase.

**For LymphoTrack and LymphoTrack Dx Assays run on the Ion S5 and PGM only. CE-marked *in vitro* diagnostic products are not available for sale or use within North America.

LymphoTrack[®] Dx and LymphoTrack[®] Workflow Summary

Illumina[®] MiSeq™



STEP1 Prepare DNA samples

- **STEP 2** Run PCR with a LymphoTrack[®] Assay
- **STEP 3** Purify amplicons
- **STEP 4** Quantify and prepare library
- STEP 5 Load onto MiSeq[™] flow cell and run
- STEP 6 Process raw data with bioinformatics software

STEP 7 Analyze results

Notice: The LymphoTrack Dx Assays are *in vitro* diagnostic products and are available in regions that accept CE-IVD products. *Image courtesy of Illumina, Inc.

LymphoTrack[®] Dx and LymphoTrack[®] Workflow Summary

Thermo Fisher Scientific[®] Ion S5/PGM[™]



- **STEP1** Prepare DNA samples
- STEP 2 Run PCR with a LymphoTrack[®] Assay
- **STEP 3** Purify amplicons
- STEP 4 Quantify, prepare, and enrich library
- **STEP 5** Load onto Ion S5 /PGM[™] chip and run
- STEP 6 Process raw data with bioinformatics software

STEP 7 Analyze results

Notice: The LymphoTrack Dx Assays are *in vitro* diagnostic products and are available in regions that accept CE-IVD products. *Image courtesy of Thermo Fisher Scientific

LymphoTrack[®] Flex Workflow Summary

Illumina[®] NextSeq[™]1000



- **STEP 1** Amplify targets via PCR1
- STEP 2 Purify amplicons with AMPure® XP
- **STEP 3** Add index primers via PCR2
- **STEP 4** Pool then purify single-target libraries
- **STEP 5** Quantify single-target libraries
- STEP 6 Load library pool, configure run and sequence
- **STEP 7** Analyze raw data with LymphoTrack[®] Enterprise Software

Notice: The LymphoTrack Flex Assays are for Research Use Only; not intended for diagnostic or therapeutci procedures.

Next-Generation Sequencing Menu

Invivoscribe offers LymphoTrack®, LymphoTrack® Flex and LymphoTrack Dx Assays for the analysis of B- and T-cell clonality, somatic hypermutation, and measurable residual disease studies**. Assays are designed for use on both industry standard next-generation sequencing (NGS) platforms: the Illumina® MiSeq[™] and Thermo Fisher Scientific® Ion PGM[™] and Ion S5[™] instruments.

Invivoscribe assays for the Illumina® MiSeq[™] platform offer the ability to analyze up to 46 (with FR1) and 22 samples and two controls per gene target and the multiplexing capabilities to generate a sequencing library that combines amplicons from different Invivoscribe LymphoTrack and LymphoTrack Dx Assays onto the same flow cell. Our software then sorts and assigns the correct sequences to their corresponding sample.

Invivoscribe assays for the Ion PGM^{**} and Ion S5^{**} platform offer the ability to analyze up to ten samples and two controls per gene target and the multiplexing capability to generate a sequencing library that combines amplicons from different Invivoscribe LymphoTrack and LymphoTrack Dx Assays onto the same sequencing chip, reducing per sample testing costs.

All LymphoTrack and LymphoTrack Dx Assays allow for fast and easy analysis and data visualization using the LymphoTrack bioinformatics software. The LymphoTrack software sorts and assigns the sequences to their corresponding sample and provides information such as the prevalence, gene segment usage, and the mutation rate (*IGH* Leader and *IGH* FR1 only). In addition, the LymphoTrack Measurable Residual Disease (MRD) Software allows for clonotype sequences to be tracked in subsequent samples for research applications.

The table below indicates which LymphoTrack (Research Use Only) and LymphoTrack Dx (CE-IVD Marked) Assays are currently available.

CE-Marked IVD Assays	MiSeq [™]	lon S5/PGM™
LymphoTrack® Dx IGHV Leader Somatic Hypermutation Assays	AVAILABLE	NOT AVAILABLE
LymphoTrack® Dx <i>IGH</i> FR1 Assays	AVAILABLE	AVAILABLE
LymphoTrack® Dx <i>IGH</i> FR2 Assays	AVAILABLE	AVAILABLE
LymphoTrack® Dx <i>IGH</i> FR3 Assays	AVAILABLE	AVAILABLE
LymphoTrack® Dx <i>IGH</i> FR1/2/3 Assays	AVAILABLE	AVAILABLE
LymphoTrack® Dx <i>IGK</i> Assays	AVAILABLE	AVAILABLE
LymphoTrack® Dx <i>TRG</i> Assays	AVAILABLE	AVAILABLE
LymphoTrack® Dx <i>TRB</i> Assays	AVAILABLE	NOT AVAILABLE
Research Use Only (RUO) Assays	MiSeq™	lon S5/PGM™
LymphoTrack® IGHV Somatic Hypermutation Assays	AVAILABLE	NOT AVAILABLE
LymphoTrack® <i>IGH</i> FR1 Assays	AVAILABLE	AVAILABLE
LymphoTrack® <i>IGH</i> FR2 Assays	AVAILABLE	AVAILABLE
LymphoTrack [®] <i>IGH</i> FR3 Assays	AVAILABLE	AVAILABLE
LymphoTrack® <i>IGH</i> FR1/2/3 Assays	AVAILABLE	AVAILABLE
LymphoTrack® <i>IGK</i> Assays	AVAILABLE	AVAILABLE

LymphoTrack[®] Flex Assays provide additional flexibility as compared to LymphoTrack[®] Assays, including higher-sample throughput, automationfriendly workflows, and variable input DNA volume. The resulting FASTQ or FASTQ.GZ data can be analyzed by the LymphoTrack[®] Enterprise Software, which applies additional analysis parameters and provides sequence prevalence, gene segment usage, and the mutation rate (*IGH* Leader and *IGH* FR1 only).

The table below indicates which LymphoTrack® Flex Assays (Research Use Only) are currently available.

LymphoTrack Flex [®] Assays	
LymphoTrack Flex TRG Assay	(compatible with the NextSeq™1000)
LymphoTrack Flex IGH FR1 Assay	(compatible with the NextSeq™1000)
LymphoTrack Flex IGHV Leader Assay	(compatible with the NextSeq™1000)

CE-marked assays are *in vitro* diagnostic products and are not available for sale or use within North America. **Measurable residual disease (MRD) applications are currently for research use only.

Gel and Capillary Electrophoresis Menu

Invivoscribe offers assays that can be analyzed using two conventional methods of fragment analysis: gel electrophoresis or capillary electrophoresis. Gel electrophoresis kits offer a comparatively easy and inexpensive solution for clonality, translocation, and mutational testing and are often the method of choice for laboratories new to using these methods and techniques. PCR products are analyzed using non-denaturing polyacrylamide gels (PAGE) and often require a heteroduplex step for resolution of generated amplicons.

Capillary electrophoresis kits are supplied with fluorescently labeled primers, allowing the resulting PCR products to be analyzed on Applied Biosystems (ABI) platforms e.g. 3130, 3500, 3500xL, 3500xL Dx. Fragment analysis by capillary electrophoresis offers the ability to detect fragments with a high level of accuracy and analytical sensitivity and allows for greater sample throughput compared to gel detection methods. In addition, capillary electrophoresis detection often facilitates a more objective interpretation of results than gel-based detection.

The table below summarizes which detection methods are available for our clonality, translocation and *FLT3* mutation assays either as RUO, CE-IVD, or IVD.

CE-Marked IVD Assays	Gel	ABI
IdentiClone® IGH + IGK B-Cell Clonality Assay	NOT AVAILABLE	AVAILABLE
IdentiClone® <i>IGH</i> Gene Clonality Assay	NOT AVAILABLE	AVAILABLE
IdentiClone® <i>IGK</i> Gene Clonality Assay	NOT AVAILABLE	AVAILABLE
IdentiClone® <i>IGL</i> Gene Clonality Assay	NOT AVAILABLE	AVAILABLE
IdentiClone® TCRB Gene Clonality Assay	NOT AVAILABLE	AVAILABLE
IdentiClone® T-Cell Receptor Gamma Gene Rearrangement Assay 2.0	NOT AVAILABLE	AVAILABLE
IdentiClone® TCRD Gene Clonality Assay	NOT AVAILABLE	AVAILABLE
LeukoStrat® CDx <i>FLT3</i> Mutation Assay (CE 2797 IVD)	NOT AVAILABLE	AVAILABLE
IVD Companion Diagnostic Assays	Gel	ABI
LeukoStrat® CDx <i>FLT3</i> Mutation Assay IVD (USA)	NOT AVAILABLE	AVAILABLE
LeukoStrat® CDx <i>FLT3</i> Mutation Assay (JP)	NOT AVAILABLE	AVAILABLE
LeukoStrat® CDx FLT3 Mutation Assay (CE 2797 IVD)	NOT AVAILABLE	AVAILABLE
Research Use Only (RUO) Assays	Gel	ABI
Research Use Only (RUO) Assays	Gel	ABI
Research Use Only (RUO) Assays FLT3 Mutation Assay IGH + IGK B-Cell Clonality Assay	Gel AVAILABLE AVAILABLE	ABI AVAILABLE AVAILABLE
Research Use Only (RUO) Assays FLT3 Mutation Assay IGH + IGK B-Cell Clonality Assay IGH Gene Rearrangement Assay	Gel AVAILABLE AVAILABLE NOT AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE
Research Use Only (RUO) Assays FLT3 Mutation Assay IGH + IGK B-Cell Clonality Assay IGH Gene Rearrangement Assay IGH Gene Clonality Assay	Gel AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE
Research Use Only (RUO) Assays FLT3 Mutation Assay IGH + IGK B-Cell Clonality Assay IGH Gene Rearrangement Assay IGH Gene Clonality Assay IGH Somatic Hypermutation Assay v2.0	Gel AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE AVAILABLE
Research Use Only (RUO) Assays FLT3 Mutation Assay IGH + IGK B-Cell Clonality Assay IGH Gene Rearrangement Assay IGH Gene Clonality Assay IGH Somatic Hypermutation Assay v2.0 IGL Gene Clonality Assay	Gel AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE
Research Use Only (RUO) Assays FLT3 Mutation Assay IGH + IGK B-Cell Clonality Assay IGH Gene Rearrangement Assay IGH Gene Clonality Assay IGH Somatic Hypermutation Assay v2.0 IGL Gene Clonality Assay TCRB Gene Clonality Assay	Gel AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE AVAILABLE AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE
Research Use Only (RUO) Assays FLT3 Mutation Assay IGH + IGK B-Cell Clonality Assay IGH Gene Rearrangement Assay IGH Gene Clonality Assay IGH Somatic Hypermutation Assay v2.0 IGL Gene Clonality Assay TCRB Gene Clonality Assay T-Cell Receptor Gamma Gene Rearrangement Assay	Gel AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE AVAILABLE AVAILABLE NOT AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE
Research Use Only (RUO) AssaysFLT3 Mutation AssayIGH + IGK B-Cell Clonality AssayIGH Gene Rearrangement AssayIGH Gene Clonality AssayIGH Somatic Hypermutation Assay v2.0IGL Gene Clonality AssayTCRB Gene Clonality AssayT-Cell Receptor Gamma Gene Rearrangement Assay 2.0T-Cell Receptor Gamma Gene Rearrangement Assay 2.0	Gel AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE
Research Use Only (RUO) AssaysFLT3 Mutation AssayIGH + IGK B-Cell Clonality AssayIGH Gene Rearrangement AssayIGH Gene Clonality AssayIGH Gene Clonality AssayIGH Somatic Hypermutation Assay v2.0IGL Gene Clonality AssayTCRB Gene Clonality AssayT-Cell Receptor Gamma Gene Rearrangement Assay 2.0T-Cell Receptor Gamma Gene Rearrangement Assay 2.0T-Cell Receptor Gamma Gene Rearrangement Assay 2.0TCRD Gene Clonality Assay	Gel AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE
Research Use Only (RUO) AssaysFLT3 Mutation AssayIGH + IGK B-Cell Clonality AssayIGH Gene Rearrangement AssayIGH Gene Clonality AssayIGH Somatic Hypermutation Assay v2.0IGL Gene Clonality AssayTCRB Gene Clonality AssayT-Cell Receptor Gamma Gene Rearrangement Assay 2.0T-Cell Receptor Gamma Gene Rearrangement Assay 2.0TCRD Gene Clonality AssayBCL1/JH Translocation Assay	Gel AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE AVAILABLE AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE NOT AVAILABLE
Research Use Only (RUO) AssaysFLT3 Mutation AssayIGH + IGK B-Cell Clonality AssayIGH Gene Rearrangement AssayIGH Gene Clonality AssayIGH Gene Clonality AssayIGH Somatic Hypermutation Assay v2.0IGL Gene Clonality AssayTCRB Gene Clonality AssayT-Cell Receptor Gamma Gene Rearrangement Assay 2.0T-Cell Receptor Gamma Gene Rearrangement Assay 2.0TCRD Gene Clonality AssayBCL1/JH Translocation AssayBCL2/JH Translocation Assay	Gel AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE NOT AVAILABLE NOT AVAILABLE
Research Use Only (RUO) AssaysFLT3 Mutation AssayIGH + IGK B-Cell Clonality AssayIGH Gene Rearrangement AssayIGH Gene Clonality AssayIGH Somatic Hypermutation Assay v2.0IGL Gene Clonality AssayTCRB Gene Clonality AssayT-Cell Receptor Gamma Gene Rearrangement Assay 2.0TCRD Gene Clonality AssayBCL1/JH Translocation AssayBCL2/JH Translocation AssayBCL2/JH t(14;18) Translocation Assay	Gel AVAILABLE AVAILABLE NOT AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE NOT AVAILABLE NOT AVAILABLE
Research Use Only (RUO) AssaysFLT3 Mutation AssayIGH + IGK B-Cell Clonality AssayIGH Gene Rearrangement AssayIGH Gene Clonality AssayIGH Gene Clonality AssayIGH Somatic Hypermutation Assay v2.0IGL Gene Clonality AssayTCRB Gene Clonality AssayT-Cell Receptor Gamma Gene Rearrangement Assay 2.0TCRD Gene Clonality AssayBCL1/JH Translocation AssayBCL2/JH Translocation AssayBCL2/JH t(14;18) Translocation AssayBCR/ABL t(9;22) Translocation Assay	Gel AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE	ABI AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE AVAILABLE NOT AVAILABLE NOT AVAILABLE NOT AVAILABLE AVAILABLE

CE-marked assays are *in vitro* diagnostic products and are available in regions that accept CE-IVD products. IVD assays are *in vitro* diagnostic products, and are available for sale or use in regions indicated.

Common Technical Support Questions

1. What sample types may be suitable for analysis with Invivoscribe Gel and Capillary assays?

We recommend high-quality DNA for clonality testing with our assays. This can be extracted from frozen or fresh tissue, peripheral blood, bone marrow, skin biopsies, etc.

2. When should the recommended controls be run with our assays?

The no template, positive, and negative controls should be included in every run for each target, per the product insert or instructions for use.

3. What is the purpose of the Specimen Control Size Ladder and Amplification Control master mix? What is the difference between these master mixes?

The Specimen Control Size Ladder and Amplification Control master mixes are used as troubleshooting tools that allow you to determine if the quality and quantity of your DNA sample is suitable for use with our assays. The Specimen Control Size Ladder amplifies DNA at approximately 100, 200, 300, 400, and 600 base pairs; whereas, the Amplification Control amplifies DNA at 235 bp.

4. How should the master mix and controls be stored and thawed?

The master mixes should be stored at -65 to -85 °C and should be thawed at room temperature and vortexed prior to use. If you intend to use master mixes multiple times, we recommend aliquoting the master mixes to minimize the number of freeze/ thaw cycles. For the *FLT3* CDx Mutation Assay: Opened vials of master mixes stored frozen may incur up to 4 freeze thaw cycles. Opened vials of controls stored frozen may incur up to 8 freeze thaw cycles.

5. Where can more information about the primers used in our assays be found?

Most primer information is proprietary to Invivoscribe and cannot be disclosed. We can, however, tell you the target area for the primers in each master mix, if you contact our support team by emailing support@invivoscribe.com or by calling +1 858-224-6600.

6. Which targets are recommended for the study of B-cell malignancies?

The EuroClonality/BIOMED-2 Group has shown that combined testing of *IGH* and *IGK* achieves a clinical sensitivity of 99%. If purchasing these assays separately is cost prohibitive, our *IGH* + *IGK* Gene Clonality Assay (does not include *IGH* Tubes D and E) may be a feasible alternative option (see Figure 2 and Table 1 in *Leukemia* (2007) 21, 201-206). We also offer next-generation sequencing LymphoTrack® Assays for *IGH* and *IGK* for use with MiSeq[™] or Ion S5/PGM[™] instruments. In addition, a high percentage of B-ALL patients have *TRG* rearrangements, which can be detected using our assays to detect *TRG* gene rearrangements.

7. What are the differences between our *IGH* Gene Rearrangement Assays and the *IGH* Gene Clonality Assays?

The *IGH* Gene Rearrangement Assay was designed by Invivoscribe; whereas, the *IGH* Gene Clonality Assay was designed by the EuroClonality/BIOMED-2 Group. Both assays target the conserved *IGH* framework regions, Framework 1, Framework 2, and Framework 3. The *IGH* Gene Clonality Assay also targets incomplete $D_{H} - J_{H}$ rearrangements.

The IGH Gene Clonality Assay includes 33 reactions per master mix and the IGH Gene Rearrangement Assay includes 30 reactions per master mix.

8. What do IGH Tubes D and E target do and why are they challenging to interpret?

Tubes D and E of our *IGH* Gene Clonality Assays target incomplete *IGH* D_H – J_H rearrangements. It is common to see known amplicans listed in the instructions for use in cases where a polyclonal background is absent (this is likely because these rearrangements are rare). Some of our customers are concerned by this, especially because there may be some samples that have robust germline amplification greater than the valid size range. We do not expect the germline amplification to outcompete true D_H – J_H rearrangements. PCR amplicons generated from germline templates are much larger than true D_H – J_H rearrangements. As a result, PCR products of germline amplifications are less robust when a specific target is present in samples.

9. Why does the polyclonal control produce a peak around 148 bp when amplified with *IGK* Tube A – 6FAM?

The 148 bp peak is a result of the restricted repertoire of *IGK* and this peak commonly appears flanked by several smaller peaks on each side. It is still possible to have a true clonal rearrangement at this size in samples. If you suspect that this peak is clonal in one of your samples, we recommend following up with heteroduplex analysis. Alternatively, NGS-based LymphoTrack® and LymphoTrack® Dx Assays provide an easier interpretation for *IGK* and reduces the number of master mixes to just one reaction.

10. What T-cell receptor kits would you recommend to detect T-cell clonal rearrangements?

Ideally, you should perform tests for *TRB*, *TRG*, and *TRD* to achieve the highest sensitivity. The EuroClonality/BIOMED-2 Group has shown that testing both *TRB* and *TRG* offers roughly the same sensitivity for the detection of T-cell malignancies as testing all three targets; however, they highly recommend testing all three assays in parallel to achieve optimal clinical sensitivity. *TRD* is especially useful in cases of suspected immature T-cell proliferations (see Figure 2 and Table 2 in *Leukemia* (2007) 21, 201-206). We also offer NGS kits for *TRG* for use with MiSeq[®].

11. What are the differences between the *TCRG* Gene Clonality Assay and the T-Cell Receptor Gamma Gene Rearrangement Assay 2.0?

The TCRG Gene Clonality Assay was designed by the EuroClonality/ BIOMED-2 Group and consists of two master mixes. For polyclonal populations, four Gaussian distributions are generated. The T-Cell Receptor Gamma Gene Rearrangement Assay 2.0 was designed by Invivoscribe and it's performance was subsequently reviewed and validated by the EuroClonality/BIOMED-2 Group. It targets all functional V_H - J_H rearrangements in a single master mix and produces smaller amplicons grouped under a single Gaussian distribution. This allows for easier interpretation and makes the assay more suitable for DNA extracted from FFPE tissue, which may consist of partially degraded DNA that would not amplify well with the larger valid size range of the TCRG Gene Clonality Assay.

12. What are the differences between the *IGH-BCL2* Translocation Assay and the *IGH-BCL2* t(14;18) Translocation Assay?

The *IGH-BCL2* Translocation Assay was designed by the EuroClonality/BIOMED-2 Group and is available as either a CE-IVD or research use only assay whereas; the *IGH-BCL2* (14;18) Translocation Assay was designed by Invivoscribe and is only available for research use only. Both of these assays target MBR and Mcr translocations, but the *IGH-BCL2* Translocation Assay also targets the translocations at the 3' Mbr. The *IGH-BCL2* (14;18) Translocation Assay was designed as a nested PCR allowing greater sensitivities (1 clonal cell per 10,000 normal cells) to be achieved. The limit of detection of the *IGH-BCL2* Translocation Assay is 1 clonal cell per 100 normal cells. Lastly, the *IGH-BCL2* Translocation Assay includes 33 reactions, whereas the *IGH-BCL2* 1(14;18) Translocation Assay includes 30 reactions.

13. Do you offer quantitative chromosome translocation (e.g., *BCR-ABL1*) controls?

Our controls are validated for qualitative use, although our customers do successfully use them with quantitative assays. Unfortunately, we cannot guarantee their performance with any assay that was not designed by Invivoscribe.

14. Which capillary electrophoresis instruments are currently validated for use with our assay kits?

Currently the capillary electrophoresis instruments Invivoscribe has validated include: ABI 3100 and 3130 series for all capillary electrophoresis detection assays. The ABI 310 and 3500 instrument series have also been validated for the majority of our capillary electrophoresis detection assays. We are not able to support using instruments not listed as validated in the instructions for use of our CE-IVD assays.

15. What are the recommended settings for my ABI instrument?

Instruments should be calibrated with the DS-30 matrix standards (Dye set D) for the ABI 310, 3100, or 3130 instrument series. For the ABI 3500 sequencer series, we advise that you calibrate the instrument with DS-33 matrix standards. We also recommend using either POP-4 or POP-7 depending on which ABI instrument you are using. If your equipment supports POP-7, we recommend using this polymer as it can be utilized for both fragment analysis and sequencing; whereas, POP-4 can only be utilized for fragment analysis.

16. How should peaks outside the valid size range be interpreted when using assay kits?

You should not interpret peaks outside of the valid size range; although, in theory, it is possible to have a true rearrangement fall outside this region. If you are concerned about a suspect peak, you may sequence your product for confirmation. Please note that samples should always be interpreted within the context of all available clinical information.

17. Is cell-free DNA (cfDNA) a suitable sample type for Invivoscribe LymphoTrack[®] or LymphoTrack[®] Dx Assays?

The average size of cfDNA (~170 bps) makes it a suitable sample type to run with *IGH* FR3 master mixes. The use of cfDNA with *TRG* master mixes might be possible, but expected amplicon sizes generated with this assay are near the upper limits of the fragment lengths typically found with this sample type.

18. Is DNA extracted from FFPE tissue suitable to use with Invivoscribe LymphoTrack[®] or LymphoTrack[®] Dx Assays?

To ensure DNA from challenging specimens is of sufficient quality and quantity to generate a valid result, samples may be tested with the Specimen Control Size Ladder master mix.

19. On which instruments can I use the LymphoTrack[®] and LymphoTrack[®] Dx Assays?

We have different versions of our assays for the S5/PGM" and MiSeq" instruments (LymphoTrack TRB is currently available only on MiSeq"). No other DNA sequencers (e.g. 454) are currently supported. Assays for the lon S5/PGM" and MiSeq" platforms differ slightly in terms of the total number of indices, etc., but both have similar benefits such as a one-step PCR reaction and available bioinformatics software.

20. How much DNA is needed for the LymphoTrack[®] and LymphoTrack[®] Dx Assays?

50 ng of high-quality genomic DNA is required for the lon S5/PGM[™] and MiSeq[™] LymphoTrack and LymphoTrack Dx Assays for clonality and somatic hypermutation applications.

21. Can I use a different library quantification method or kit?

We recommend using the KAPA™ kit for MiSeq[™] assays and either the 2100 Bioanalyzer[®] or the LabChip[®] GX for the Ion S5/PGM[™] assays.

22. Will the LymphoTrack[®] or LymphoTrack[®] Dx analysis software work on my computer?

The software requires Microsoft Windows 7 and Windows 10 (64-bit) and Excel 2007, 2010, or 2013 and will work with most desktop or laptop PCs. For specific requirements please refer to the software instructions for use.

23. Can I use the LymphoTrack[®] or LymphoTrack[®] Dx bioinformatics software with a different assay?

No, the software will only work with datasets obtained by our LymphoTrack and LymphoTrack Dx Assays.

24. What characters can I use when naming my samples and the file pathways? What types of files are accepted by the LymphoTrack[®] and LymphoTrack[®] Dx Software – MiSeq[™]?

Our software only recognizes file names and pathways that contain the following characters (A-Z, a-z, 0-9, . (dot), _ (underscore), - (hyphen)). In addition, spaces in the pathname for the data files or software (pathnames include file folders and file names) should be avoided. If the software encounters a character that is not listed above or extra spaces, an error message may be generated. Furthermore, the software is only compatible with adaptor-trimmed fastq.gz files that are generated by the MiSeq[®] Reporter Software when the MiSeq[®] instrument is used. An example of the naming format that the MiSeq[®] Reporter uses: SampleName_S1_L001_R1_001.fastq.gz and SampleName_S1_L001_R2_001.fastq.gz.

25. Do the Invivoscribe MiSeq[™] indices correspond to the Illumina[®] indices?

The indices included in our MiSeq[™] master mixes follow Illumina[®]'s TruSeq LT nomenclature. For instance, *IGH* FR1 MiSeq[™] 01 corresponds to A001. Information for the other indices can be found in the instructions for use on how to set up the MiSeq[™] Sample Sheet to detect the appropriate indices.

26. Why am I getting a low percent passing filter and Q30 score?

Low Q30 and percent passing filter (%PF) scores could be an indication that the flow cell is overloaded. If this is suspected, verify your amplicon and library calculations and quantifications are correct. Low run metrics can also be attributed to many additional factors including poor quality DNA, contamination, flow cell or instrument issues, etc. Please refer to your Illumina MiSeq[®] user guides and contact Illumina[®] Support.

27. Why is the same VH-JH rearrangement combination and sequence shared by two groups of reads, one of which is several bases shorter than the other when looking at the Read Summary tab of the excel document created by the LymphoTrack[®] Visualization Tool?

Our software was designed to list every unique sequence separately in order for the customer to see all of the data and make their own determination on how to interpret it. The several base pair difference can be due to a number of factors including amplification errors and sequencing errors. It could also be a result of similarities between some of the primer sequences that were designed to ensure maximum coverage. We also include a Merged Read Summary report for your reference that combines sequences that only differ by 1 or 2 basepairs.

28. Do I need to perform an adapter ligation prior to sequencing my products?

Performing an adapter ligation is not needed. The primers included in our LymphoTrack and LymphoTrack Dx master mixes already include the appropriate index barcodes and adapter sequences. After PCR amplification, you will be able to proceed with amplicon purification, amplicon quantification, library pooling, and sequencing.

29. If the LymphoTrack[®] or LymphoTrack[®] Dx software generated an error, what information should I submit to Technical Support?

Please submit the *.txt Log file that should have been created by the software in the output folder, a screenshot of the sample directory, and the Lot Number of the software you are using to <u>support@invivoscribe.com</u>.

30. Are controls provided with the kits? Can you purchase additional controls? How are they supplied?

Each kit contains the necessary positive and negative controls required to perform the assay; additional controls may also be purchased separately. Single-tube DNA controls are provided as 100 µL aliquots of 200 µg/mL in 1/10 TE Buffer, 50 µL aliquots of 50 µg/ mL in 1/10 TE Buffer, and 45 µL aliquots of 15 µg/mL in 1/10 TE Buffer. Single-tube RNA controls are provided as 100 µL aliquots of 400 µg/ml in RNAse free in glass distilled water.

31. What are the differences between dilution sets, sensitivity panels, and proficiency panels?

31a. RNA Dilution Sets

BCR/ABL b3a2 (Cat# 4-085-0210), BCR/ABL b2a2 (Cat# 4-085-0310), and BCR/ABL e1a2 (Cat# 4-085-0110). These sets contain six tubes: 100% negative control RNA and volume to volume (v/v) dilutions (10⁻¹, 10⁻², 10⁻³, 10⁻⁴, and 10⁻⁵) of the positive control RNA into the negative control RNA (IVS-0048). The RNA Dilution Sets are supplied at a concentration of 400 µg/mL, and each tube contains 50 µL. These dilution sets may be used to establish a standard reference curve, as proficiency controls, as sensitivity controls for specific target assays, and as routine testing controls for cDNA synthesis, amplification and detection.

31b. RNA Sensitivity Panels

These panels consist of seven tubes: 100% positive control RNA and v/v dilutions (10-¹, 10-², 10-³, 10-⁴, 10-⁵, and 10-⁶) of the positive control RNA into the negative control RNA (IVS-0035). The RNA Sensitivity Panels are supplied at a concentration of 400 µg/mL, and each tube contains 100 µL. The RNA Sensitivity Panels may be used as sensitivity controls for specific target assays, and as routine testing controls for cDNA synthesis, amplification and detection.

31c. DNA Sensitivity Panels

Consist of six tubes: 100% clonal DNA and v/v dilutions of the clonal DNA into negative polyclonal DNA (IVS-0000) to make 30%, 20%, 10%, 5%, and 1% dilutions. The DNA Sensitivity Panels are supplied at a concentration of 200 μ g/mL and each tube contains 100 μ L. The DNA Sensitivity Panels may be used as sensitivity controls for specific target assays.

31d. RNA Proficiency Panel

The proficiency panel for BCR-ABL1 t(9;22) can be used as a sensitivity control for specific target assays, and as routine testing controls for cDNA synthesis, amplification and detection. It consists of ten tubes: 100% positive control RNA and v/v dilutions (10- $^{\circ}$ and 10- $^{\circ}$) of IVS-0003, IVS-0011 and IVS-0032. It also includes BCR-ABL1 Negative Clonal Control RNA (IVS-0035).

Product List by Catalog Number

Companion Diagnostic Assays Capillary Analysis

K-412-0331	LeukoStrat® CDx FLT3 Mutation Assay (JP)
K-412-0341	LeukoStrat® CDx FLT3 Software (JP)
K-412-0361	LeukoStrat® CDx FLT3 Mutation Assay (USA)
K-412-0371	LeukoStrat® CDx FLT3 Software (USA)
K-412-0431	LeukoStrat® CDx FLT3 Mutation Assay (CE 2797 IVD)
K-412-0441	LeukoStrat® CDx FLT3 Software (CE 2797 IVD)

Research Use Only Assays Capillary & Gel Fragment Analysis

1-100-0010	IGH + IGK B-Cell Clonality Assay – Gel Detection
1-100-0031	IGH + IGK B-Cell Clonality Assay – ABI Fluorescence Detection
1-100-0041	IGH + IGK B-Cell Clonality Assay MegaKit – ABI Fluorescence Detection
1-101-0020	IGH Gene Clonality Assay – Gel Detection
1-101-0051	IGH Gene Rearrangement Assay – ABI Fluorescence Detection
1-101-0061	IGH Gene Clonality Assay – ABI Fluorescence Detection
1-101-0071	IGH Gene Rearrangement Assay MegaKit – ABI Fluorescence Detection
1-101-0081	IGH Gene Clonality Assay MegaKit – ABI Fluorescence Detection
1-102-0020	IGK Gene Clonality Assay – Gel Detection
1-102-0021	IGK Gene Clonality Assay – ABI Fluorescence Detection
1-102-0031	IGK Gene Clonality Assay MegaKit – ABI Fluorescence Detection
1-103-0010	IGL Gene Clonality Assay – Gel Detection
1-103-0011	IGL Gene Clonality Assay – ABI Fluorescence Detection
1-103-0021	IGL Gene Clonality Assay MegaKit – ABI Fluorescence Detection
1-205-0010	TCRB Gene Clonality Assay – Gel Detection
1-205-0011	TCRB Gene Clonality Assay – ABI Fluorescence Detection
1-205-0020	TCRB Gene Clonality Assay MegaKit – Gel Detection
1-205-0021	TCRB Gene Clonality Assay MegaKit – ABI Fluorescence Detection
1-206-0010	TCRD Gene Clonality Assay – Gel Detection
1-206-0011	TCRD Gene Clonality Assay – ABI Fluorescence Detection
1-206-0021	TCRD Gene Clonality Assay MegaKit – ABI Fluorescence Detection
1-207-0051	T-Cell Receptor Gamma Gene Rearrangement Assay – ABI Fluorescence Detection
1-207-0101	T-Cell Receptor Gamma Gene Rearrangement Assay 2.0 –
	ABI Fluorescence Detection
1-207-0111	T-Cell Receptor Gamma Gene Rearrangement Assay 2.0 MegaKit –
	ABI Fluorescence Detection
S100001	IdentiClone® RUO TCRG 2.0 Software
1-308-0010	BCL1/JH Translocation Assay – Gel Detection
1-308-0020	BCL1/JH Translocation Assay MegaKit – Gel Detection
1-309-0010	BCL2/JH t(14;18) Translocation Assay – Gel Detection
1-309-0020	BCL2/JH Translocation Assay – Gel Detection
1-309-0040	BCL2/JH Translocation Assay MegaKit – Gel Detection
1-310-0010	BCR/ABL t(9;22) Translocation Assay – Gel Detection
1-310-0031	BCR/ABL t(9;22) Translocation Assay – ABI Fluorescence Detection
1-311-0011	PML/RARat(15;17) Translocation Assay – ABI Fluorescence Detection
1-412-0010	FLT3 Mutation Assay - Gel Detection
	ELTERNA A MOLEL D. A. M.

1-412-0031 FLT3 Mutation Assay - ABI Fluorescence Detection

Master Mixes

2-096-0020 2-096-0021 2-101-0011 2-101-0031 2-101-0041 2-101-0051	Specimen Control Size Ladder – Unlabeled Specimen Control Size Ladder – 6FAM IGH Tube A – 6FAM IGH Tube C – HEX IGH Tube D – HEX IGH Tube E – 6FAM
2-101-0061	<i>IGH</i> Framework 1 (FR1) – 6FAM
2-101-0081	IGH Framework 3 (FR3) – HEX
2-101-0091	<i>IGH</i> Framework 2 (FR2) – 6FAM
2-101-0101	<i>IGH</i> Tube B – 6FAM
2-101-0180	Hypermutation Mix 2 v2.0 – Unlabeled
2-102-0011	IGK Tube A – 6FAM
2-102-0021	IGK Tube B – 6FAM
2-103-0011	<i>IGL</i> Tube – 6FAM
2-205-0011	TCRB Tube A – 6FAM & HEX
2-205-0021	TCRB Tube B – 6FAM
2-205-0031	TCRB Tube C – 6FAM & HEX
2-206-0011	TCRD Tube – 6FAM & HEX
2-207-0010	T-Cell Receptor Gamma Mix 1 – Unlabeled
2-207-0020	T-Cell Receptor Gamma Mix 2 – Unlabeled
2-207-0091	TCRG – 6FAM
2-308-0010	BCL1/JH Tube – Unlabeled
2-309-0010	BCL2/JH t(14;18) (Mbr) Mix 1b – Unlabeled
2-309-0020	BCL2/JH t(14;18) (mcr) Mix 2b – Unlabeled
2-309-0030	BCL2/JH t(14;18) (Mbr) Mix 1a – Unlabeled
2-309-0040	BCL2/JH t(14;18) (mcr) Mix 2a – Unlabeled

2-309-0050	<i>BCL2/</i> Jн Tube A – Unlabeled
2-309-0060	<i>BCL2/</i> Jн Tube B – Unlabeled
2-309-0070	<i>BCL2/</i> Jн Tube C – Unlabeled
2-310-0040	BCR/ABL t(9;22) Mix 1b – Unlabeled
2-310-0050	BCR/ABL t(9;22) Mix 2b – Unlabeled
2-310-0060	BCR/ABL t(9;22) Mix 2c – Unlabeled
2-310-0070	BCR/ABL t(9;22) Mix 3b - Unlabeled
2-310-0080	BCR/ABL t(9;22) Mix 3c – Unlabeled
2-310-0090	BCR/ABL t(9;22) Mix 3d - Unlabeled
2-311-0011	PML/RARa t(15;17) Mix 1 – HEX
2-311-0031	PML/RARa t(15;17) Mix 2b – HEX
2-311-0041	PML/RARa t(15;17) Mix 2c – HEX

BCR/ABL RNA Dilution Sets

4-085-0110	BCR/ABL e1a2 RNA Dilution Set
4-085-0210	BCR/ABL b3a2 RNA Dilution Set
4-085-0310	BCR/ABL b2a2 RNA Dilution Set

Next-Generation Sequencing RUO Assays & Software

7-121-0059	LymphoTrack® <i>IGHV</i> Somatic Hypermutation Assay Kit A – MiSeq™
7-121-0069	LymphoTrack [®] IGHV Somatic Hypermutation Assay Panel – MiSeq™
7-121-0129	LymphoTrack® <i>IGH</i> FR1/2/3 Assay Kit A – MiSeq™
7-121-0139	LymphoTrack® <i>IGH</i> FR1/2/3 Assay Panel – MiSeq™
7-121-0009	LymphoTrack® <i>IGH</i> FR1 Assay Kit A – MiSeq™
7-121-0039	LymphoTrack® <i>IGH</i> FR1 Assay Panel – MiSeq™
7-121-0089	LymphoTrack® <i>IGH</i> FR2 Assay Kit A – MiSeq™
7-121-0099	LymphoTrack® <i>IGH</i> FR2 Assay Panel – MiSeq™
7-121-0109	LymphoTrack® <i>IGH</i> FR3 Assay Kit A – MiSeq™
7-121-0119	LymphoTrack® <i>IGH</i> FR3 Assay Panel – MiSeq™
7-121-0057	LymphoTrack® <i>IGH</i> FR1/2/3 Assay – S5/PGM™
7-121-0007	LymphoTrack® <i>IGH</i> FR1 Assay – S5/PGM™
7-121-0037	LymphoTrack® <i>IGH</i> FR2 Assay – S5/PGM™
7-121-0047	LymphoTrack® <i>IGH</i> FR3 Assay – S5/PGM™
7-122-0009	LymphoTrack® <i>IGK</i> Assay Kit A – MiSeq™
7-122-0019	LymphoTrack® <i>IGK</i> Assay Panel – MiSeq™
7-122-0007	LymphoTrack® <i>IGK</i> Assay – S5/PGM™
7-225-0009	LymphoTrack® <i>TRB</i> Assay Kit A – MiSeq™
7-225-0019	LymphoTrack® <i>TRB</i> Assay Panel – MiSeq™
7-227-0019	LymphoTrack® <i>TRG</i> Assay Kit A – MiSeq™
7-227-0009	LymphoTrack® <i>TRG</i> Assay Panel – MiSeq™
7-227-0007	LymphoTrack® <i>TRG</i> Assay – S5/PGM™
7-500-0007	LymphoTrack® Software – S5/PGM™
7-500-0008	LymphoTrack® MRD Software
7-500-0009	LymphoTrack® Software – MiSeq™
7100004	LymphoTrack® Flex IGHV Leader Assay
7100003	LymphoTrack® Flex IGH FR1 Assay
7100002	LymphoTrack® Flex <i>TRG</i> Assay
S100003	LymphoTrack [®] Enterprise Software
1-412-0019	FLT3 ITD MRD Assay
S100005	FLT3 ITD MRD Software
1-416-0019	NPM1 MRD Assay
S100004	NPM1 MRD Software

Capillary Fragment Analysis CE-IVD IdentiClone® Assays

9-100-0031	IdentiClone® IGH + IGK B-Cell Clonality Assay – ABI Fluorescence Detection
9-100-0041	IdentiClone® IGH + IGK B-Cell Clonality Assay MegaKit –
	ABI Fluorescence Detection
9-101-0061	IdentiClone [®] IGH Gene Clonality Assay – ABI Fluorescence Detection
9-101-0081	IdentiClone® IGH Gene Clonality Assay MegaKit –
	ABI Fluorescence Detection
9-102-0021	IdentiClone [®] IGK Gene Clonality Assay – ABI Fluorescence Detection
9-102-0031	IdentiClone® IGK Gene Clonality Assay MegaKit –
	ABI Fluorescence Detection
9-103-0011	IdentiClone [®] IGL Gene Clonality Assay – ABI Fluorescence Detection
9-205-0011	IdentiClone® TCRB Gene Clonality Assay – ABI Fluorescence Detection
9-205-0021	IdentiClone® TCRB Gene Clonality Assay MegaKit – ABI Fluorescence Detection
9-206-0011	IdentiClone® TCRD Gene Clonality Assay – ABI Fluorescence Detection
9-206-0021	IdentiClone® TCRD Gene Clonality Assay MegaKit –
	ABI Fluorescence Detection
9-207-0101	IdentiClone® T-Cell Receptor Gamma Gene Rearrangement Assay 2.0 –
	ABI Fluorescence Detection
9-207-0111	IdentiClone® T-Cell Receptor Gamma Gene Rearrangement Assay 2.0
	MegaKit – ABI Fluorescence Detection

Analyte Specific Reagents

A-412-0071	FLT3 ITD MM - 6FAM & HEX - ASF
A-412-0081	FLT3 TKD MM - 6FAM - ASR
A-101-0011	IGH FR1 – 6FAM
A-101-0020	IGH FR2 – Unlabeled
A-101-0030	IGH FR3 – Unlabeled
A-101-0031	IGH FR3 – HEX
A-101-0041	IGH DH1 - 6 - HEX
A-101-0051	IGH DH7 – 6FAM
A-101-0061	IGH Framework 1 – 6FAM
A-101-0070	IGH Framework 2 – Unlabeled
A-101-0080	IGH Framework 3 – Unlabeled
A-101-0081	IGH Framework 3 – HEX
A-101-0091	IGH Framework 2 – 6FAM
A-101-0101	IGH FR2 – 6FAM
A-102-0010	IGKV – J – Unlabeled
A-102-0011	IGK V – J – 6FAM
A-102-0020	<i>IGK</i> V – K _{de} – Unlabeled
A-102-0021	IGKV – K _{de} – 6FAM
A-103-0011	IGL V – J – 6FAM
A-205-0010	TCRB V - J1+2 – Unlabeled
A-205-0011	TCRB V – J1+2 – 6FAM & HEX
A-205-0020	TCRBV - J2 – Unlabeled
A-205-0021	TCRB V – J2 – 6FAM
A-205-0031	<i>TCRB</i> D – J1+ 2 – 6FAM & HEX
A-207-0021	TCRG V(1-8)J – HEX
A-207-0071	TCRG V(1-8,9)J – 6FAM
A-207-0091	TCRG V(2-5,8-11)J1+2+P – 6FAM
A-309-0050	BCL2/JH Mbr – Unlabeled
A-309-0060	<i>BCL2/</i> Jн 3'Mbr – Unlabeled
A-309-0070	BCL2/JH mcr – Unlabeled

BCR/ABL RNA Dilution Sets

4-085-0110 BCR/ABL e1a2 RNA Dilution Set 4-085-0210 BCR/ABL b3a2 RNA Dilution Set 4-085-0310 BCR/ABL b2a2 RNA Dilution Set

RNA Sensitivity Panels

4-087-0030	Sensitivity Panel – IVS-0003 Clonal Control RNA
4-087-0110	Sensitivity Panel – IVS-0011 Clonal Control RNA
4-087-0150	Sensitivity Panel – IVS-0015 Clonal Control RNA
4-087-0200	Sensitivity Panel – IVS-0020 Clonal Control RNA
4-087-0320	Sensitivity Panel – IVS-0032 Clonal Control RNA

Cell Line DNA Controls

4-088-0008	IGH SHM Positive Control DNA
4-088-0010	IVS-0001 Clonal Control DNA
4-088-0098	LymphoTrack® B-cell Low Positive Control
4-088-0108	LymphoTrack® T-cell Low Positive Control
4-088-0118	LymphoQuant® B-cell Internal Control
4-088-0128	LymphoQuant® T-cell Internal Control
4-088-0190	IVS-0004 Clonal Control DNA
4-088-0370	IVS-0007 Clonal Control DNA
4-088-0430	IVS-0008 Clonal Control DNA
4-088-0490	IVS-0009 Clonal Control DNA
4-088-0550	IVS-0010 Clonal Control DNA
4-088-0730	IVS-0013 Clonal Control DNA
4-088-1090	IVS-0019 Clonal Control DNA
4-088-1210	IVS-0021 Clonal Control DNA
4-088-1390	IVS-0024 Clonal Control DNA
4-088-1690	IVS-0029 Clonal Control DNA
4-088-1750	IVS-0030 Clonal Control DNA
4-088-1810	IVS-0031 Clonal Control DNA

Cell Line RNA Controls

4-089-0100	IVS-0002 Clonal Control RNA
4-089-0190	IVS-0003 Clonal Control RNA
4-089-0200	10 ⁻¹ IVS-0003 Clonal Control RNA
4-089-0210	10 ⁻² IVS-0003 Clonal Control RNA
4-089-0220	10 ⁻³ IVS-0003 Clonal Control RNA
4-089-0230	10 ⁻⁴ IVS-0003 Clonal Control RNA
4-089-0240	10 ⁻⁵ IVS-0003 Clonal Control RNA
4-089-0910	IVS-0011 Clonal Control RNA

 4-089-0920
 10⁻¹ IVS-0011 Clonal Control RNA

 4-089-0930
 10⁻² IVS-0011 Clonal Control RNA

 4-089-0940
 10⁻³ IVS-0011 Clonal Control RNA

 4-089-0950
 10⁻⁴ IVS-0011 Clonal Control RNA

 4-089-0960
 10⁻⁵ IVS-0011 Clonal Control RNA

 4-089-0950
 10⁻⁴ IVS-0011 Clonal Control RNA

 4-089-0960
 10⁻⁵ IVS-0011 Clonal Control RNA

 4-089-1270
 IVS-0015 Clonal Control RNA

 4-089-2810
 IVS-0022 Clonal Control RNA

 4-089-2800
 IVS-0032 Clonal Control RNA

 4-089-2810
 10⁻¹ IVS-0032 Clonal Control RNA

 4-089-2820
 10⁻² IVS-0032 Clonal Control RNA

 4-089-2830
 10⁻³ IVS-0032 Clonal Control RNA

 4-089-2840
 10⁻⁴ IVS-0032 Clonal Control RNA

 4-089-2850
 10⁻⁵ IVS-0032 Clonal Control RNA

Tissue DNA Control

4-092-0010 IVS-0000 Polyclonal Control DNA

LymphoTrack[®] Low Positive Controls

4-088-0098 LymphoTrack® B-cell Low Positive Control 4-088-0108 LymphoTrack® T-cell Low Positive Control

LymphoQuant[®] Internal Controls

4-088-0118 LymphoQuant® B-cell Internal Control 4-088-0128 LymphoQuant® T-cell Internal Control

Somatic Hypermutation Sequencing Assays

5-101-0030	IGH Somatic Hypermutation Assay v2.0 – Gel Detection
5-101-0031	IGH Somatic Hypermutation Assay v2.0 – ABI Fluorescence Detection
5-101-0040	IGH Somatic Hypermutation Assay v2.0 MegaKit – Gel Detection
5-101-0041	IGH Somatic Hypermutation Assay v2.0 MegaKit – ABI Fluorescence

ABI Reagents

6-098-0061 HI-Deionized Formamide with ROX Size Standard (ABI 3100)

Taq DNA Polymerases

6-097-0130 FalconTaq DNA Polymerase

Next-Generation Sequencing CE-IVD LymphoTrack® Dx Assays

9-121-0007	LymphoTrack® Dx <i>IGH</i> FR1 Assay – S5/PGM™
9-121-0009	LymphoTrack® Dx <i>IGH</i> FR1 Assay Kit A – MiSeq™
9-121-0037	LymphoTrack® Dx <i>IGH</i> FR2 Assay – S5/PGM™
9-121-0039	LymphoTrack® Dx <i>IGH</i> FR1 Assay Panel – MiSeq™
9-121-0047	LymphoTrack® Dx <i>IGH</i> FR3 Assay – S5/PGM™
9-121-0057	LymphoTrack® Dx <i>IGH</i> FR1/2/3 Assay – S5/PGM™
9-121-0059	LymphoTrack® Dx <i>IGHV</i> Leader Somatic Hypermutation Assay Kit A – MiSeq™
9-121-0069	LymphoTrack® Dx IGHV Leader Somatic Hypermutation Assay Panel – MiSeq™
9-121-0089	LymphoTrack® Dx <i>IGH</i> FR2 Assay Kit A – MiSeq™
9-121-0099	LymphoTrack® Dx <i>IGH</i> FR2 Assay Panel – MiSeq™
9-121-0109	LymphoTrack® Dx <i>IGH</i> FR3 Assay Kit A – MiSeq™
9-121-0119	LymphoTrack® Dx <i>IGH</i> FR3 Assay Panel – MiSeq™
9-121-0129	LymphoTrack® Dx <i>IGH</i> FR1/2/3 Assay Kit A – MiSeq™
9-121-0139	LymphoTrack® Dx <i>IGH</i> FR1/2/3 Assay Panel – MiSeq™
9-122-0009	LymphoTrack® Dx <i>IGK</i> Assay Kit A – MiSeq™
9-122-0019	LymphoTrack® Dx <i>IGK</i> Assay Panel – MiSeq™
9-122-0007	LymphoTrack® Dx <i>IGK</i> Assay – S5/PGM™
9-225-0009	LymphoTrack® Dx <i>TRB</i> Assay Kit A – MiSeq™
9-225-0019	LymphoTrack® Dx <i>TRB</i> Assay Panel – MiSeq™
9-227-0007	LymphoTrack® Dx <i>TRG</i> Assay – S5/PGM™
9-227-0009	LymphoTrack® Dx <i>TRG</i> Assay Panel – MiSeq™
9-227-0019	LymphoTrack® Dx <i>TRG</i> Assay Kit A – MiSeq™
9-500-0007	LymphoTrack® Dx Software – S5/PGM™

9-500-0009 LymphoTrack® Dx Software – MiSeq