

LeukoStrat[®]

CDx *FLT3* Mutation Assay

FDA Approved | CE 2797 IVDR Certified

Enable Clinical Trial Enrollment and Patient Management with LabPMM[®]'s LeukoStrat[®] CDx *FLT3* Mutation Assay Service

Clinical Information

Comprehensive *FLT3* mutation testing plays a critical role in the management of patients with acute myeloid leukemia (AML). LabPMM's LeukoStrat[®] CDx *FLT3* Mutation Assay service detects both *FLT3* internal tandem duplications (ITD) and tyrosine kinase domain (TKD) mutations in peripheral blood or bone marrow specimens. The only regulatory approved *FLT3*-mutation analysis test in the United States and the European Union, the globally standardized test identifies patients who may benefit from *FLT3*-targeted therapies*, including VANFLTYA[®], RYDAPT[®], and XOSPATA[®].

LeukoStrat[®] CDx *FLT3* Mutation provides physicians with clinically relevant molecular insight to support personalized management of patients with *FLT3*-mutated AML.

LabPMM's Report Includes

- ✓ *FLT3*-ITD detection and size characterization**
- ✓ *FLT3*-TKD mutation result
- ✓ *FLT3*-ITD and -TKD signal ratio
- ✓ Objective, actionable result interpretation

Partner with LabPMM[®] today for trusted *FLT3* testing.

For questions connect with a scientist at inquiry@invivoscribe.com

* The availability of XOSPATA[®], VANFLTYA[®], and RYDAPT[®] varies by region and is subject to local regulatory approval. Contact Invivoscribe for details on local regulatory approvals.

**The *FLT3*-ITD base pair length is reported exclusively by the LabPMM's US location.

***LabPMM Results Portal is accessible to US customers only.



Specimen Requirements

- » 1-2 mL of peripheral blood in EDTA or Sodium Heparin
- » 0.25-0.5 mL of bone marrow in EDTA or Sodium Heparin



Turnaround Time

- » 2 to 3 business days



Shipping Conditions

- » 2-8 °C up to 72 hours; do not freeze (peripheral blood or bone marrow)



Specimen Stability

- » 2-8 °C up to 7 days **prior to testing**

Expand your level of *FLT3* characterization by ordering today. Follow our simple submission process:

1

Complete LabPMM's TRF



2

Send us your bone marrow or peripheral sample

3

Receive on interpretive report in only 2-3 days

4

Conveniently access secure test results through the LabPMM[®] Results Portal***



Send completed TRF and results portal forms to customerservice@invivoscribe.com

LabPMM[®]
an invivoscribe company

[invivoscribe.com](https://www.invivoscribe.com)

Tel +1 858.224.6666 | Fax +1 858.408.3964
support@labpmm.com