

TCRG Gene Clonality Master Mixes

For identification of clonal T-cell receptor gamma chain gene rearrangements.

RUO For RESEARCH USE ONLY. Not for use in diagnostic procedures.

Master mixes are components of complete assay/s. In general, the protocols used for all master mixes are identical, but some assays require different DNA polymerase volumes, methods of detection, and/or PCR program protocols. Therefore, it is important to refer to the specific assay product insert for the operating procedures. Product inserts are available on the web at www.invivoscribe.com. The product insert of interest will be found on the specific webpage of the respective assay.


This assay is based on the EuroClonality/BIOMED-2 Concerted Action BMH4-CT98-3936.



Invivoscribe, Inc.
10222 Barnes Canyon Road, Bldg.1
San Diego, CA 92121-2711 USA
invivoscribe.com | support@invivoscribe.com
Tel 858.224.6600 | Fax 858.224.6601

Manufactured in U.S.A.



 Storage Conditions: **-85°C to -65°C**

Catalog #	Product	Quantity
REF 22070030	TCRG Tube A - Unlabeled for Gel Detection	33 reactions
REF 22070040	TCRG Tube B – Unlabeled for Gel Detection	33 reactions
REF 22070031	TCRG Tube A – 6FAM & HEX for ABI Fluorescence Detection	33 reactions
REF 22070041	TCRG Tube B – 6FAM & HEX for ABI Fluorescence Detection	33 reactions